



Document Number: FB003289

Revision: A

Group: Fact Book

Type: External

State: Released

Latest Released: YES

Implemented Date: 11/18/2019

Stamp Date: Monday, November 18, 2019 10:27:39 AM EST



Johnson & Johnson Medical KK

Sugakawa, Japan

Service and Repair Center for the Megapower 1000 Product Code

Factbook # FB003289



Factbook Approval

We have reviewed and do approve Factbook #FB003289, Addendum to Factbook F-261 (EpiCenter document Number FB001103), for the service and repair of the MegaPower 1000 Generator at Johnson & Johnson Medical KK Sugakawa, Japan. We find that the documentation contained in this factbook meets the requirements as defined by WE001534, Rev F, Establishment and Certification of Service Depot, Field Service, and Parts Depot Centers. The signoff of this factbook represents the certification process and indicates the service center is to be considered authorized/qualified to perform service on the identified product(s).

E-Sig in EpiCenter

Shannon Gillespie
Manager, Service Center

E-Sig in EpiCenter

Date

E-Sig in EpiCenter

Robert Peters
Team Leader, Customer Quality
Worldwide Service and Repair

E-Sig in EpiCenter

Date



FACTBOOK CHECKLIST

Date: October 29, 2019

From: Jason Stivers

Re: 1000

Attention: Factbook # FB003289

Activity	Complete	Not Applicable
Service System Quality Assessment	✓	
Technical Training Program	✓	
Documentation System	✓	
Equipment Installation Qualification (EIQ)	✓	
Physical Requirements	✓	
Operating Agreement	✓	
Start-Up Activities	✓	
Supplier Approval	✓	

Indicate Activity status with a single “✓” mark.

E-Sig in EpiCenter

Jason Stivers
Service Engineer, EES – Service Staff Engineer



FACTBOOK STRATEGY

Date: October 29, 2019

From: Jason Stivers

Re: 1000

Attention: Factbook # FB003289

The support data establishing Johnson & Johnson Medical KK Sugakawa, Japan as an Ethicon Endo-Surgery authorized service center for the 1000 Megadyne capital equipment is contained in this Factbook.

The process follows WE001534 Rev F, Establishment and Certification of Service Depot, Field Service, and Parts Depot Centers. A Service and Repair Facility Qualification Record has been created which serves as the index for this factbook. Approval signatures will appear on each memorandum page as outlined by the Service and Repair Facility Qualification Record.

This information is an Addendum to the original Factbook F-261 that qualified Johnson & Johnson Medical KK Sugakawa, Japan as an authorized service center for the GEN04, SCM12, SCM23, and STHC1. This Factbook documents their training and qualification to repair the product codes 1000 generator. The original product codes SCM12, SCM23, and STHC1 were sold to another company and are not part of this factbook process. Johnson & Johnson Medical KK also services the RF60 generator approved in factbook number FB002269 and GEN11 approved in factbook number FB002911. Those factbooks and associated service of products are not affected by this factbook.

E-Sig in EpiCenter

Jason Stivers
Service Engineer, EES – Service Staff Engineer



Service & Repair Facility Qualification Record

The table of contents of this Factbook is listed below. The documentation contained meets the intent of WE001534 Rev F, Establishment and Certification of Service Depot, Field Service, and Parts Depot Centers.

Section	Pg.
1.0 Service System Quality Assessment	7
1.1 Reference Memorandum	8
2.0 Technical Training Results	9
2.1 Training Results Memorandum	10-11
2.1.1 Training Record Form	12-13
2.1.2 Training Agenda	14
2.1.3 The Science of Electrosurgery – Training Presentation	15-59
2.1.4 Service and Repair Forms	60-70
2.1.5 Test Results	71-72
2.1.6 Training Certificate	73
3.0 Documentation System	74
3.1 Documentation System Memorandum	75
4.0 Equipment Installation Qualification (EIQ) RESULTS	76
4.1 EIQ Results Memorandum	77
4.1.1 Mega Power 1000 Installation Qualification Test Protocol ENG-PRT-473	78-79
4.1.2 Mega Power 1000 Installation Qualification ENG-PRT-473 Appendix A	80-89
4.1.3 Mega Power 1000 Performance/Operational Qualification Test Protocol ENG-PRT-502	90-91
4.1.4 Mega Power 1000 Performance/Operational Qualification Test Protocol Results – CS-FRM-034	92-102
4.1.5 Mega Power 1000 Performance/Operational Qualification Test Protocol Results – Master Error Log	103
4.1.6 Mega Power 1000 Performance/Operational Qualification Test Protocol Results – ESU Test Results	104
4.1.7 Specification sheet – Fluke QA-ESII Electrosurgery Analyzer	105-108
4.1.8 Calibration sheet – Fluke QA-ESII Electrosurgery Analyzer	109-111
4.1.9 Specification sheet – Yokogawa Digital Multimeter 733	112
4.1.10 Calibration sheet - Yokogawa Digital Multimeter 733	113-116
4.1.11 Specification sheet – Kikusui Ground Bond Tester TOS6100	117
4.1.12 Calibration sheet - Kikusui Ground Bond Tester TOS6100	118-120
4.1.13 Specification sheet – Hioki Leak Current HiTester 3156	121
4.1.14 Calibration sheet - Hioki Leak Current HiTester 3156	122-124
5.0 Physical Requirements	125
5.1 Reference Memorandum	126
5.1.1 J&J Medical KK Service Center Floor Plan	127
6.0 Operating Agreement	128
6.1 Reference Memorandum	129
6.1.1 Operating Agreement	130-146
7.0 Start-Up Activities	147

7.1	Reference Memorandum Start up Activities	148
8.0	Supplier Approval	149
8.1	Reference Memorandum Supplier Approval	150

1.0 SERVICE SYSTEM QUALITY ASSESSMENT



SERVICE SYSTEM QUALITY ASSESSMENT

Date: October 29, 2019

From: Jason Stivers

Attention: Factbook # FB003289

Re: 1000

Refer to section 1.0 of Factbook F-261 (EPI Center document number FB001103) for the support documentation that a servicing quality system assessment was conducted according to established procedures. Ongoing quality assessments of this facility will be maintained on file in EtQ per schedule.

E-Sig in EpiCenter

Jason Stivers
Service Engineer, EES – Service Staff Engineer

2.0 TECHNICAL TRAINING RESULTS



TECHNICAL TRAINING RESULTS

Date: October 29, 2019

From: Jason Stivers

Re: 1000

Attention: Factbook # FB003289

On February 26-28, 2019, one service repair technician (Service Manager) from the J&J K.K. Medical Company Ltd. Japan Service and Repair Depot was provided training at Megadyne Medical Products, Draper, Utah for the following products:

- 1000, Mega Power Electrosurgical Generator

The training was conducted by Megadyne Service and Repair Trainers John Minuth and Tyler Skinner with assistance from Megadyne Repair Technician, Bruce Hevelone. Training began with a basic introduction to the product application and use, theory of operation of the system and a general product description and functional description. The attached training log contains the list of Megadyne forms, work instructions, software, and service bulletins required to be covered during the training for product codes. Additionally, procedures/manuals that included disassembly, reassembly, repair, testing, quality inspection, and product release were covered in the training. Troubleshooting information was also covered, which included identifying common causes of failure, hardware troubleshooting, and service testing. Afterward, standard service center processes, such as bench tests, electrical safety tests, and product release tests were demonstrated. To demonstrate the ability to repair the products the trainee was provided and passed a written test post training.

Finally, complaint awareness training is conducted on an annual basis and thus was not a needed deliverable for this specific training. Training records for complaint awareness are maintained within the training management system at the J&J KK Medical.

With this successful completion of the activities referenced above, the following individual(s) should now be considered trained as an authorized EES representative capable of the analysis, service, and repair of the Megadyne products listed above, and as a qualified and authorized trainer for the product.

1. Kohei Seki, Service Manager

Additionally, with this successful completion of the activities referenced above, the following individual(s) should now be considered trained as an authorized EES representative capable of the final release of product to inventory and authorized to train the quality release person(s) within their center.

1. Kohei Seki, Service Manager

Attached evidence of completion of these activities is:

- Franchise Qualification and Training Record Form (Shared)
- Training Agenda
- The Science of Electrosurgery – training presentation
- Megadyne Mega Power Service Center Repair Form, New Faceplate – objective evidence of completed device testing
- Mega Power Training Exam Results

- Certificate of Training

Note: Training Agenda and Training Record Form, include references to the Megadyne Mega Soft Patient Return Electrodes, and Mega Vac and Mega Vac Plus smoke evacuators. These references are not applicable for this factbook. The smoke evacuator products are not sold or serviced in the Japan market but were included in Megadyne training session. The Mega Soft products will be included in a future factbook.

E-Sig in EpiCenter

Jason Stivers
Service Engineer, EES – Staff Service Engineer

Form Non-PPE
Quality System
Franchise Qualification and Training Record Form(Shared)

FM-0000809 / Rev: 11
CO: 100648063

FRANCHISE QUALIFICATION AND TRAINING RECORD FORM (Shared)

Page 1 of 2

NOTE: COMPLETE THIS FORM IN BLUE OR BLACK INK ONLY.

Franchise/Site Impacted: <input type="checkbox"/> Acclarent <input type="checkbox"/> ASP <input type="checkbox"/> Biosense Webster <input type="checkbox"/> Ethicon <input checked="" type="checkbox"/> Ethicon-Endo <input type="checkbox"/> Mentor <input type="checkbox"/> NeuWave <input type="checkbox"/> Sterilmed					
Method of Training: <input checked="" type="checkbox"/> Instructor Led <input type="checkbox"/> Self Training <input type="checkbox"/> External Training <input type="checkbox"/> Awareness					
Instructor's Name/Signature/Date (Print {N/A if Self or External training}): John Minuth <i>John Minuth 2/28/19</i>					
Course, NR/CAPA, or Document Description: (Description of content or course agenda/outline. Course description may include, but not limited to, the complete document title, version and sections/topic covered or brief description of reason for training)					
Megadyne Service & Repair Training – Mega Power 1000 Service & Repair, Mega Soft Pads Evaluation, Mega Vac Planned Maintenance, and Trainer Training					
Training ID: Course Code, Title, NR/CAPA Number, or Document Number	Document Revision Number	Training Start Date	Training End Date	Training Duration	For Training Dept Use Only (or designee)
					Class Code (Optional)
CS-FRM-034	004	2/26/19	2/28/19	N/A	
CS-FRM-035	003	2/26/19	2/28/19	N/A	
ENG-FRM-013	003	2/26/19	2/28/19	N/A	
CS-FRM-027	001	2/26/19	2/28/19	N/A	
CS-FRM-028	001	2/26/19	2/28/19	N/A	
ENG-WI-037	001	2/26/19	2/28/19	N/A	
ENG-WI-036	001	2/26/19	2/28/19	N/A	
ENG-WI-035	002	2/26/19	2/28/19	N/A	
ENG-WI-053	004	2/26/19	2/28/19	N/A	
ICM-470-9024	H	2/26/19	2/28/19	N/A	
MKT-LBL-062	001	2/26/19	2/28/19	N/A	
MKT-LBL-063	002	2/26/19	2/28/19	N/A	
MKT-LBL-081	002	2/26/19	2/28/19	N/A	
MKT-LBL-101	002	2/26/19	2/28/19	N/A	
<i>For Training Dept / Designee Use Only (Optional)</i>					

Form Non-PPE
 Quality System
 Franchise Qualification and Training Record Form(Shared)

FM-0000809 / Rev: 11
 CO: 100648063

Printed Name of Data Entry Person: _____ Signature: _____ Date Entered: _____

The signature indicates that the training records are entered into ComplianceWire by the indicated person(s).

FRANCHISE QUALIFICATION AND TRAINING RECORD FORM (Shared)

KHS 9/18/19

Page 2 of 2

GR 9/18/19

First Training ID: N/A	<i>CS-FRM-034</i>	End Date: N/A	<i>2/28/19</i>
My signature below indicates I have self-trained / reviewed or participated in Instructor Led Training as described on page 1 and I understand the information or concepts covered.			
TYPE / WRITE NAME	SIGNATURE	DATE	WWID #
1. Kohei Seki	<i>K. Seki</i>	<i>02/28/2019</i>	152816107

CONFIDENTIAL use pursuant to Company Procedures

Page 4 of 4

Megadyne Medical Products
11506 S State St | Draper, UT 84020



Megadyne, Mega Power 1000, Training Schedule, Salt Lake City: Feb 26 to Feb 28, 2018

Tuesday, Feb 26, 2019


9:00 am – 10:30 am	Principles to Electrosurgery
10:45 am – 12:00 pm	Anatomy and Troubleshooting, Mega Power 1000
12:00 pm – 1:00 pm	Lunch
1:00 pm – 2:00 pm	Disassembly of Mega Power 1000
2:00 pm – 4:00 pm	Assembly of Mega Power 1000

Wednesday, Feb 27, 2019

9:00 am – 10:00 am	Calibration, Mega Power 1000
10:00 am – 12:00 pm	Final Test
12:00 pm – 1:00 pm	Lunch
1:00 pm – 2:00 pm	Paperwork Overview

Thursday, Feb 28, 2019


9:00 am – 12:00 am	Preventive Maintenance Smoke Evacuator
12:00 am – 1:00 pm	Lunch
1:00 pm - 4:00pm	Pad Loaners



MEGADYNE
The Electrosurgical Authority

Electricity

Electricity flows in a closed path



Confidential

Confidential

Confidential


2

1



The Transformer

For Electrosurgery Current goes from Pin 3 back to Pin 4



Confidential


Confidential

Confidential

3



The Electrosurgical “Circuit”



Confidential

Confidential

Confidential

4

2



Basic Terms

- Voltage – Potential difference (Volts)
- Current – The time rate of flow of an electrical charge. (Amps)
- Power – The Rate that work is done. (Watts or Horsepower)

1 Horsepower = 746 Watts

Confidential

Confidential


Confidential

5



Basic Equations

- Ohm's Law: $V = I \cdot R$
- Power = $V \cdot I$ or V^2/R or $I^2 \cdot R$



Confidential

Confidential

Confidential

6



Generator Feedback Loop

- Generate a signal
- Measure Voltage and Current
- Calculate Power
- Change the voltage of the signal
- Start over.

Confidential

Confidential


Confidential

7



Electrocuted?

- Why doesn't electrosurgery shock people?



Confidential

Confidential

Confidential


8

MEGADYNE
The Electrosurgical Authority

The Answer: Radio Freq AC

Alternating Current

AC



Confidential

Confidential

Confidential


9

MEGADYNE
The Electrosurgical Authority

The Nervous System

Neurons (the copper wires of the body)

Structure of a Typical Neuron



Confidential

Confidential


Confidential

10

5

MEGADYNE
The Electrosurgical Authority

Sodium / Potassium Pump



Copyright © 2001 Benjamin Cummings, an imprint of Addison Wesley Longman, Inc.

Confidential


Confidential

Confidential

11

MEGADYNE
The Electrosurgical Authority

The Pump is Slow



Generator: 400 kHz
Each cycle = 2.5µS

Wall Plug: 60 Hz
Each cycle = 17mS

Neuron Cycle = 4mS

Confidential

Confidential



Confidential

12



Cutting at the Surgical Site Only

- If current is travelling through the entire body, why does it only cut at the active electrode?



Confidential

Confidential


Confidential

13



Current Density

- The current enters the body at one point, but travels out through a large pad.
- This is similar to a pressure washer.



Confidential

Confidential


Confidential

14



Current Density

- Increase the current density and you increase the surgical effect.



Confidential

Confidential

Confidential

15



What are the different modes?

- Ace
- Cut
- Coag
- Blend
- Bipolar



Confidential

Confidential

Confidential

16



Cut

- The purpose of cut is to heat up cells so that their cytoplasm boils and the cell ruptures.



Confidential

Confidential


Confidential

17



Cut Waveform

- The Cut Wave form is a sine wave. It continually applies power and ruptures cells.



Confidential

Confidential

Confidential

18



Cut Power Curves




Figure 1 - Manufacturer Published Power Output versus Impedance Graphs

Confidential

Confidential


Confidential

19



Coag

- The purpose of Coag is to heat a large area and dry out blood.



Confidential

Confidential

Confidential

20



Coag Waveform

Coag bursts pulses for two reasons.

1. To let tissue cool and prevent penetration into the tissue.
2. To get larger voltages to encourage sparking



Confidential

Confidential

Confidential

21



Coag Power Curves

Power Output vs Impedance for 60 Watt Coag Setting




Figure 3 - Example of Actual Power Output versus Impedance Graphs

Confidential

Confidential


Confidential

22



Cut Vs Coag

all @ 40W



Confidential

Confidential


Confidential

23



Blend

- Blend combines Cut and Coag.



Confidential

Confidential



Confidential

24



Bipolar

- Bipolar returns the current at the surgical site.
- Bipolar reduces power as the tissue dries.



Confidential

Confidential


Confidential

25



ACE (A better cut)

- Two requirements for ACE
 1. A special blade
 2. Special power algorithms



Confidential

Confidential

Confidential

26



The ACE Blade

ACE blades have a patented geometry that naturally funnels (or focuses) the electrosurgical current



Confidential

Confidential


Confidential

27



ACE Power Curves

- Power goes down once the tissue has been cut.



Confidential


Confidential

Confidential

28



Power Curves (Cut and ACE)



Confidential

Confidential

Confidential

29



Generator Specifications

Output Power Characteristics

Mode	Power (Watts)	Output Tolerance * (Rated Load)	Rated Load (Ohms)	Maximum Open Circuit Voltage (Vp-p)	Operating Frequency (Rated Load)	Crest Factor Nominal @ (Rated Load)
Monopolar CUT						
ACE Cut	150	20%	200	1500	400kHz	1.6
Pure Cut	300	20%	300	3000	400kHz	1.6
Blend	200	20%	300	4000	400kHz	3.0
Monopolar COAG						
COAG 1	120	20%	500	5000	2.5μs Pulse @ 32kHz	6.9
COAG 2	120	20%	500	5000	2.5μs Pulse @ 30kHz	7.1
Spray	120	20%	500	6000	2.5μs Pulse @ 22kHz	8.0
Bipolar						
Micro	80	20%	100	360	400kHz	1.6
Macro	80	20%	100	760	400kHz	1.6

Confidential

Confidential

Confidential


30



Return Pads

Three Basic Kinds

1. Single Plate Pad Electrode
2. Dual Plate Pad Electrode
3. Mega Soft



Confidential

Confidential


Confidential

31



Dual Pads

- Dual Plate Pads are used to monitor the pad connection to the patient.



Confidential

Confidential

Confidential




32

16



Pads and Edge Effects

When more current comes into a particular edge this is called an edge effect. This can lead to unbalanced heating and burns.



Confidential


Confidential

Confidential

33



The Mega Soft



Confidential

Confidential

Confidential

34



Capacitively Coupled

- The Mega Soft connects the patient to the return through a capacitance.
- The capacitance is proportional to the area of the patient in contact with the pad.
- Prevents burns due to lack of contact.
- Prevents edge effects.

Confidential


Confidential

Confidential

35



Mega Soft Setup



Confidential


Confidential

Confidential

36



Mega Soft Positions



Confidential

Confidential

Confidential

37



Questions

- Any Questions?



Confidential

Confidential

Confidential

38

19



More Terms


* Voltage

V_{pp} = peak-to-peak voltage – lowest negative to highest voltage

V_p = peak voltage – highest voltage away from zero (V_{pp} estimate = $2 \times V_p$)

$V = V_{rms} =$ root mean square voltage – effective average

$$CF = V_p / V_{rms}$$



Confidential

Confidential



The Electrosurgical Authority™


Electricity

* Voltage

(Example: home voltage)

$$120V = 120V_{rms} = 170V_p = 340V_{pp}$$

$$CF = 170/120 = 1.4$$



Confidential


Confidential

Confidential


40



Electrosurgery - Overview




Waveforms





Waveforms

a@ 40W



Confidential

Confidential

Confidential

43



Power Curves

ESU Output

Most ESUs (including the Mega Power) have modes that regulate power output:

- 1) Constant Power Output vs Resistive Load (e.g., Pure Cut or Standard Coag), or
- 2) Constant Voltage Output vs Resistive Load (e.g., ACE mode, some Erbe modes).

Confidential

Confidential

Confidential

44


22

MEGADYNE


The Electrosurgical Authority

Power Curves


CONSTANT POWER



Current Decreases as R Increases




Voltage Increase as R Increases

**MEGADYNE**


The Electrosurgical Authority

Power Curves


CONSTANT VOLTAGE



Current Decreases as R Increases



Constant Voltage





Power Curves

Power Modes

As a safety feature to prevent unexpected power delivery spikes, simultaneous activation of multiple instruments is not possible on the ForceTract™ system.

Monopolar Mode

The ForceTract™ energy platform produces five different modes of power output:

Cut Modes

Pure cut provides a clean, precise cut to my tissue with little or no hemostasis. Blend cut is a conventional blend mode of cut that provides slower cutting and additional hemostasis.

Vapor/Blow Mode

Vapor/Blow mode is a unique combination of hemostasis and desiccation and allows the user to dry down for more hemostasis and speed up for faster desiccation. Thermal spread is equal or superior to Cut or Blend modes.

Cog Mode

It generates coagulation tissue by sparking from the active electrode through air to the patient tissue. Since sparks may spray unpredictably from the electrode during fulguration, using fulguration for delicate tissue or in confined areas can complicate surgery. Accidental sparking to adjacent areas can occur as tissue at the surgical site dries and becomes more resistive to current flow.


1-6

ForceTract™ User's Guide



ESU Technical Background

Mega Power Coagulation Standard Output Power Graphs
(from User's Manual)



Confidential

Confidential


Confidential

48



ESU Technical Background

Competitive Generator - Coagulation Standard Output Power Graphs
(from their User's Manual – Force FX C)



Confidential

Confidential


Confidential

49



ESU Technical Background

ACE Power Graphs (from User's Manual)





Confidential

Confidential


Confidential

50





ESU Technical Background



Confidential


Confidential

Confidential

53



Mega Power Specifications



Confidential

Confidential

Confidential

54



Mega Power Specifications

Output Power Characteristics

Mode	Power (Watts)	Output Tolerance * (Rated Load)	Rated Load (Ohms)	Maximum Open Circuit Voltage (Vp-p)	Operating Frequency (Rated Load)	Crest Factor Nominal @ (Rated Load)
Monopolar CUT						
ACE Cut	150	20%	200	1500	400kHz	1.6
Pure Cut	300	20%	300	3000	400kHz	1.6
Blend	200	20%	300	4000	400kHz	3.0
Monopolar COAG						
COAG 1	120	20%	500	5000	2.5µs Pulse @ 32kHz	6.9
COAG 2	120	20%	500	5000	2.5µs Pulse @ 30kHz	7.1
Spray	120	20%	500	6000	2.5µs Pulse @ 22kHz	8.0
Bipolar						
Micro	80	20%	100	360	400kHz	1.6
Macro	80	20%	100	760	400kHz	1.6

Confidential

Confidential

Confidential

55



Mega Power Specifications

Power Output vs Impedance for 160 Watt Cut Setting




Figure 1 - Manufacturer Published Power Output versus Impedance Graphs

Confidential

Confidential

Confidential

56

MEGADYNE
The Electrosurgical Authority

Mega Power Specifications

Power Output vs Impedance for 60 Watt Coag Setting




Figure 2 - Manufacturer Published Power Output versus Impedance Graphs

Confidential

Confidential

Confidential

57

MEGADYNE
The Electrosurgical Authority

Mega Power Specifications

Power Output vs Impedance for 60 Watt Coag Setting




Figure 3 - Example of Actual Power Output versus Impedance Graphs

Confidential

Confidential


Confidential

58




Mega Power Specifications

MEGAPOWER
Electrosurgical Generator



OPERATOR'S MANUAL

MEGAPOWER
Electrosurgical Generator



Service Manual

Confidential

Confidential

Confidential

59



Let's Take a Peek Inside

Confidential

Confidential

Confidential

60

30



Mega Power – what is inside?

1) General Purpose Generator

- CUT Modes
- COAG Modes
- BIPOLAR Modes

2) Advanced Cutting Effect

- Smart Cutting Mode

Confidential

Confidential


Confidential

61



Advanced Cutting Effect

ACE BLADE



Special Mode



Confidential


Confidential

Confidential

62



Advanced Cutting Effect



Confidential

Confidential

Confidential

63



Advanced Cutting Effect

How Does ACE Work?

Confidential

Confidential

Confidential

64

32



Advanced Cutting Effect

How Does ACE Work?

Two Keys

1) ACE Mode (generator)

+

2) ACE Blade (blade shape)

Confidential

Confidential

Confidential

65



Advanced Cutting Effect

Let's Look at ACE Mode

Background

Most Modern Day Generators
try and maintain Constant Power
Output

Confidential

Confidential




Confidential

66



ESU Technical Background

CONSTANT POWER



Advanced Cutting Effect

ACE Mode
is
DIFFERENT!

Confidential

Confidential



Confidential

68

MEGADYNE
The Electrosurgical Authority

ESU Technical Background

CONSTANT VOLTAGE




MEGADYNE
The Electrosurgical Authority

ESU Technical Background

ACE Power Graphs (from User's Manual)

ACE Cut- Output Power Vs Load Resistance



Confidential

Confidential

Confidential

70



Advanced Cutting Effect

Thought Experiment

How do we measure Power?

(what units are used)

Confidential

Confidential

Confidential

71



Advanced Cutting Effect

For Electrosurgical Generators

Power is measured in WATTS

(metric units)

Confidential

Confidential

Confidential

72



Advanced Cutting Effect

For Electrosurgical Generators

Power is measured in WATTS

(metric units)

How about Power in
automobiles – what units?

Confidential

Confidential

Confidential

73



Advanced Cutting Effect

Power is measured in

Horse Power

1 Horse Power = 746 Watts

Confidential

Confidential

Confidential


74



Advanced Cutting Effect

Thought Experiment

What if we have an auto putting out
100 HP – how fast will it go?



Confidential

Confidential


Confidential

75



Advanced Cutting Effect

still 100 HP - How about now?



Constant Power = Different Effect

Confidential

Confidential

Confidential

76



Advanced Cutting Effect

for ACE

the only effect we want is

*** CLEAN CUTTING ***

(no thermal damage, no drag)

regardless of Load!

Confidential

Confidential


Confidential

77



Advanced Cutting Effect

Back to Thought Experiment:



Using a Constant Power Mode, would one expect the same surgical effect if you changed blade depth or tissue type?

Confidential

Confidential


Confidential

78



Advanced Cutting Effect

With Constant Power, Voltage Increases as Resistance Increases



Confidential

Confidential

Confidential

79



Advanced Cutting Effect

Only
EFFECT
we
want
with
ACE

by
yes,
this
1 GI

iod
the
the
final
) + 0
.12-1

hemostasis.
Electrosurgical generators may deliver monopolar current to produce electrosurgical cutting and/or coagulation. Electrosurgical cutting is achieved by a high voltage (>200 V) continuous current. This cutting effect begins with micro-electric arcs that are concentrated at the area of application via an active electrode (accessory device). This produces a rapid rise in the internal temperature of cells, causing intracellular water to boil and the cells to burst. This path of disrupted cells creates an electrosurgical cut. Dissipation of heat to the tissue at the margins of the cut produces a zone of coagulation along the cut edge.

VOLUME 58, NO. 5, 2003

Confidential

Confidential


Confidential

80

40



Advanced Cutting Effect



Confidential

Confidential

Confidential

81



Advanced Cutting Effect

- ACE mode is actually not as simple as Constant Voltage.
- As the resistive load increases, there are other losses in the electrosurgical circuit that rob voltage from the surgical site.
- The ACE mode uses a proprietary feedback algorithm to adjust the Mega Power to maintain constant voltage at the surgical site, thus compensating the output for other voltage losses in the system as impedance dynamically changes.

Confidential

Confidential


Confidential

82



Advanced Cutting Effect

Ace Power Curve Graphs



Confidential

Confidential

Confidential

83



Advanced Cutting Effect

If you do not get it right –

- Too much Voltage = thermal damage
- Not enough Voltage = drag
- Right amount of Voltage = ACE effect

Confidential

Confidential

Confidential

84

42



Advanced Cutting Effect

Two Keys to ACE

The other Half of the

ACE Effect

comes from the blade geometry

Confidential



Confidential

Confidential

85



Advanced Cutting Effect



Confidential

Confidential

Confidential

86



Advanced Cutting Effect

ACE blades have a patented geometry that naturally funnels (or focuses) the electrosurgical current



Confidential

Confidential


Confidential

87



Advanced Cutting Effect

Competitive Advanced Surgical Cutters



Confidential

Confidential

Confidential



88

44



Other RF-based Generators

Hyfrecators



Confidential


Confidential

Confidential

89



QUESTIONS?



TRIAL LAKE, UTAH AUG 2012

Confidential

Confidential

Confidential

90

The user must ensure that they are using the correct/current revision of this document.
Document: CS-FRM-034 Rev: 004 Effective: 15 Dec 2017

Mega Power Service Center Repair Form, New Faceplate

MEGADYNE
The Electrosurgical Authority



Mega Power Information	
RGA	16180
Serial Number	175027001
RoHS	<input type="checkbox"/> Leaded <input checked="" type="checkbox"/> RoHS

PART A: Equipment Calibration

Record calibration data:

Table A1: Calibration Data

Equipment	Brand	Model Number	Serial Number	Cal. Due Date
ESU Analyzer	FLUKE	DA-ESTI	202407	18 Jun 19
Digital Multi-meter	FLUKE	179	37220327	1/31/2020
Safety Analyzer	FLUKE	ESA 620	2543060	11/30/2019

*All testing must be performed with calibrated equipment. Do not proceed with final testing if any piece of equipment is not current with calibration schedule.

Part B: Investigation and Repair

Section 1: Before the Mega Power is disassembled

- With the ESU analyzer set to 300 Ohms, measure the output power with the Mega Power set to PURE CUT and the power set to 30 W, 75 W, and 300 W. Record the results in Table B1.
- With the ESU analyzer set to 500 Ohms, measure the output power with the Mega Power set to COAG 1 and the power set to 30 W, 50 W, and 120 W. Record the results in Table B2.
- Using Field Calibration Software, 1410010-01, download, print, and attach the error logs from the master processor and watch dog processor.

Table B1: PURE CUT

Setting	Power Output
300 Ω	30 W
300 Ω	75 W
300 Ω	300 W

Table B2: COAG 1

Setting	Power Output
500 Ω	30 W
500 Ω	50 W
500 Ω	120 W

Section 2: Investigation

Table B3: Investigation

The Unit is good working. Loauer return

Initial/Date: KS 02/27/19

Document Number	Revision	Previous Number	Page
CS-FRM-034	004	N/A	1 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 27 Feb 2019, 10:18:31 am; Printed by: MSCHROEDER.

The user must ensure that they are using the correct/current revision of this document.
Document: CS-FRM-034 Rev: 004 Effective: 15 Dec 2017

Mega Power Service Center Repair Form, New Faceplate



Section 3: Repair Actions

Table B4: Repair Actions

Table B5: Replacement Parts (Confirm all replacement parts with CS-FRM-038)

Ship Printed Circuit Boards to Megadyne

Ship replaced parts
to Megadyne for evaluation.

Megadyne Medical Product, Inc.
Attn: Service and Repair
11506 South State St.
Draper, Utah 84020
USA

Table B6: Shipped to Megadyne

Part Number	Description	PCB Number (Where Applicable)
		Ref ID: U212A12019
Shipped By:	L&G	Date:
Tracking Number:	G4U912019	

Document Number CS-FRM-034	Revision 004	Previous Number N/A	Page 2 of 11
-------------------------------	-----------------	------------------------	-----------------

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 27 Feb 2019, 10:18:31 am; Printed by: MSCHROEDER.

The user must ensure that they are using the correct/current revision of this document.
Document: CS-FRM-034 Rev: 004 Effective: 15 Dec 2017

Mega Power Service Center Repair Form, New Faceplate



Once the repair work has been completed and the Mega Power has been reassembled, test per Part C through Part G.

PART C: Display Test

Section 1: Power-Up Test

Turn on the Mega Power. Once the Mega Power has completed the POST test check the appropriate boxes in Table C1.

Table C1: Power-Up Test

Pass	Fail	Criteria
<input type="checkbox"/>	<input type="checkbox"/>	Power Switch is illuminated
<input checked="" type="checkbox"/>	<input type="checkbox"/>	PURE CUT LED is illuminated
<input checked="" type="checkbox"/>	<input type="checkbox"/>	COAG 1 LED is illuminated
<input checked="" type="checkbox"/>	<input type="checkbox"/>	MICRO LED is illuminated
<input checked="" type="checkbox"/>	<input type="checkbox"/>	No errors are displayed
<input checked="" type="checkbox"/>	<input type="checkbox"/>	All display segments show dashes.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	CQM Alarm LED is red and alarm tone sounded (Without a return pad connected to the Mega Power)

Section 2: Display Functionality

Table C2: Display Functionality

Pass	Fail	Criteria
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Press each mode button and verify that the appropriate LED button lights up.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Press the up arrow button on each of the three channels and verify power goes up.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Press the down arrow button on each of the three channels and verify power goes down.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Press the Bipolar Tone button and verify the LED lights up.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Plug in the Monopolar Return Test Cable and verify CQM alarm LED turns green.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Set all power settings to 1. Turn power off. Turn power on. Press recall and verify power setting of 1 is displayed for each mode.

Document Number	Revision	Previous Number	Page
CS-FRM-034	004	N/A	Page 3 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 27 Feb 2019, 10:18:31 am; Printed by: MSCHROEDER.

The user must ensure that they are using the correct/current revision of this document.
Document: CS-FRM-034 Rev: 004 Effective: 15 Dec 2017

Mega Power Service Center Repair Form, New Faceplate



PART D: Monopolar Output Tests

Section 1: Monopolar A, Power Outputs and Crest Factors

Procedure:

- 1) Set the ESU analyzer to the loads called out in Table D1.
- 2) Set the Mega Power to the power settings (PSET) and modes called out in Table D1.
- 3) Measure and record the appropriate output wattage and crest factors where indicated.
- 4) Verify Mega Power Tones.

Table D1: Monopolar A, Power Outputs and Crest Factors

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)	Crest Factor	Criteria
300	Pure	1	2	0 - 5		1.3 – 1.9
300	Pure	10	9	5 - 15		
300	Pure	150	149	120 - 180		
300	Pure	300	299	240 - 360	1.5	
300	Blend	1	2	0 - 5		2.5 – 3.5
300	Blend	10	10	5 - 15		
300	Blend	100	103	80 - 120		
300	Blend	200	209	160 - 240	2.8	
300	ACE		131	108 - 162	1.5	1.3 – 1.9
Verify Tone for CUT: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			Verify volume can be adjusted: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			
500	Coag 1	1	2	0 - 5		6.4 – 7.4
500	Coag 1	10	10	5 - 15		
500	Coag 1	60	60	48 - 72		
500	Coag 1	120	120	96 - 144	6.5	
500	Coag 2	1	1	0 - 5		6.6 – 7.6
500	Coag 2	10	9	5 - 15		
500	Coag 2	60	60	48 - 72		
500	Coag 2	120	121	96 - 144	6.7	
500	Spray	1	1	0 - 5		7.7 – 8.7
500	Spray	10	9	5 - 15		
500	Spray	60	60	48 - 72		
500	Spray	120	121	96 - 144	7.8	
Verify Tone for COAG: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			L.S. 02/21/19 L.S. 02/21/19			

Confidential - Controlled copy on date printed

Document Number	Revision	Previous Number	Page
CS-FRM-034	004	N/A	4 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 27 Feb 2019, 10:18:31 am; Printed by: MSCHROEDER.

The user must ensure that they are using the correct/current revision of this document.
Document: CS-FRM-034 Rev: 004 Effective: 15 Dec 2017

Mega Power Service Center Repair Form, New Faceplate



Section 2: Monopolar B, Power Outputs and Crest Factors

Procedure:

- 1) Set the Mega Power up so that Monopolar B is connected to the ESU Analyzer
- 2) Ensure that the activation button activates the Monopolar B electrode.
- 3) Set the ESU analyzer to the loads called out in Table D2.
- 4) Set the Mega Power to the power settings (PSET) and modes called out in Table D2.
- 5) Measure and record the appropriate output wattage and crest factors where indicated.
- 6) Verify Mega Power tones.

Table D2: Monopolar B, Power Outputs and Crest Factors

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)	Crest Factor	Criteria
300	Pure	1	2	0 - 5		1.3 – 1.9
300	Pure	10	9	5 - 15		
300	Pure	150	150	120 - 180		
300	Pure	300	300	240 - 360		
300	Blend	1	2	0 - 5		2.5 – 3.5
300	Blend	10	10	5 - 15		
300	Blend	100	104	80 - 120		
300	Blend	200	210	160 - 240		
300	Ace		131	108 - 162	1.5	1.3 – 1.9
		Verify Tone for Cut:		<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail		
500	Coag 1	1	2	0 - 5		6.4 – 7.4
500	Coag 1	10	10	5 - 15		
500	Coag 1	60	56	48 - 72		
500	Coag 1	120	120	96 - 144		
500	Coag 2	1	1	0 - 5		6.6 – 7.6
500	Coag 2	10	9	5 - 15		
500	Coag 2	60	57	48 - 72		
500	Coag 2	120	121	96 - 144		
500	Spray	1	1	0 - 5		7.7 – 8.7
500	Spray	10	9	5 - 15		
500	Spray	60	60	48 - 72		
500	Spray	120	117	96 - 144		
		Verify Tone for COAG:		<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail		

Confidential - Controlled copy on date printed

Document Number	Revision	Previous Number	Page
CS-FRM-034	004	N/A	5 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 27 Feb 2019, 10:18:31 am; Printed by: MSCHROEDER.

The user must ensure that they are using the correct/current revision of this document.
Document: CS-FRM-034 Rev: 004 Effective: 15 Dec 2017

Mega Power Service Center Repair Form, New Faceplate



Section 3: Monopolar B Power Curves

Procedure:

- 1) Set the ESU analyzer to the loads called out in Table D3 and Table D4.
- 2) Set the Mega Power to the power settings and modes called out in Table D3 and Table D4.
- 3) Measure and record the appropriate output wattage where indicated.

Table D3: CUT Power Curves

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)
100	ACE		159	120 – 180
200	ACE		151	120 – 180
300	ACE		See Table D2	108 – 162
500	ACE		92	76 – 114
1,000	ACE		50	42 – 64
2,000	ACE		22	18 – 28
100	Pure	300	155	105 – 175
200	Pure	300	269	240 – 360
300	Pure	300	See Table D2	240 – 360
500	Pure	300	245	240 – 360
1,000	Pure	300	286	240 – 360
2,000	Pure	300	195	160 – 240
100	Blend	200	96	75 – 125
200	Blend	200	264	160 – 240
300	Blend	200	See Table D2	160 – 240
500	Blend	200	203	160 – 240
1,000	Blend	200	190	160 – 240
2,000	Blend	200	166	130 – 190

Table D4: COAG Power Curves

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)
100	COAG 1	120	117	96 – 144
200	COAG 1	120	129	96 – 144
500	COAG 1	120	See Table D2	96 – 144
1,000	COAG 1	120	117	92 – 140
2,000	COAG 1	120	118	90 – 135
100	COAG 2	120	121	96 – 144
200	COAG 2	120	129	96 – 144
500	COAG 2	120	See Table D2	96 – 144
1,000	COAG 2	120	113	92 – 138
2,000	COAG 2	120	110	86 – 129
100	Spray	120	123	96 – 144
200	Spray	120	128	96 – 144
500	Spray	120	See Table D2	96 – 144
1,000	Spray	120	103	90 – 135
2,000	Spray	120	101	80 – 120

Document Number CS-FRM-034	Revision 004	Previous Number N/A	Page 6 of 11
-------------------------------	-----------------	------------------------	-----------------

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 27 Feb 2019, 10:18:31 am; Printed by: MSCHROEDER.

The user must ensure that they are using the correct/current revision of this document.
Document: CS-FRM-034 Rev: 004 Effective: 15 Dec 2017

Mega Power Service Center Repair Form, New Faceplate



PART E: Bipolar Output Tests

Section 1: Bipolar Power Outputs and Crest Factors

Procedure:

- 1) Set the Mega Power up so that Bipolar is connected to the ESU Analyzer
- 2) Ensure that footswitch activates the Bipolar electrodes.
- 3) Set the ESU analyzer to the loads called out in Table E1.
- 4) Set the Mega Power to the power settings (PSET) and modes called out in Table E1.
- 5) Measure and record the appropriate output wattage and crest factors where indicated.
- 6) Verify Mega Power Tones.

Table E1: Bipolar, Power Outputs and Crest Factors

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)	Crest Factor	Criteria
100	Micro	1	1	0 - 5		1.4 – 2.0
100	Micro	10	10	5 - 15		
100	Micro	40	39	32 - 48		
100	Micro	80	82	64 - 96		
100	Macro	1	1	0 - 5		1.4 – 2.0
100	Macro	10	9	5 - 15		
100	Macro	40	39	32 - 48		
100	Macro	80	81	64 - 96		
Verify normal tone*: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			Verify pulsing tone*: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			

* When bipolar is activated the Mega Power displays the estimated current flow in the current flow indicator window. In addition to the current flow indicator, the Mega Power can be set to indicate the current with an audio alarm. The current is proportional to the pulsing rate of the audio tone. By pressing the bipolar tone button (the button to the right of the current flow window) the bipolar alarm will toggle between a continuous tone and a pulsing tone. Verify that both the continuous and the pulsing tones are operational in the Macro Mode with a 100 Ω load, and the Mega Power set to 80 W.

Document Number	Revision	Previous Number	Page
CS-FRM-034	004	N/A	7 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 27 Feb 2019, 10:18:31 am; Printed by: MSCHROEDER.

The user must ensure that they are using the correct/current revision of this document.
Document: CS-FRM-034 Rev: 004 Effective: 15 Dec 2017

Mega Power Service Center Repair Form, New Faceplate



Section 2: Bipolar Power Curves

Procedure:

- 1) Set the ESU analyzer to the loads called out in Table E2.
- 2) Set the Mega Power to the power settings (PSET) and modes called out in Table E2.
- 3) Measure and record the appropriate output wattage where indicated.

Table E2: BIPOLAR Power Curves

Load (Ω)	Mode	PSET (watts)	Output (Watts)	Criteria (Watts)
50	Micro	80	66	46 – 68
75	Micro	80	81	60 – 90
100	Micro	80	See Table E1	64 – 96
200	Micro	80	76	64 – 96
650	Micro	80	7	5 – 15
1000	Micro	80	5	0 – 10
50	Macro	80	66	46 – 68
75	Macro	80	81	60 – 90
100	Macro	80	See Table E1	64 – 96
200	Macro	80	76	64 – 96
500	Macro	80	42	38 – 58
1000	Macro	80	28	20 – 32

Confidential - Controlled copy on date printed

Document Number CS-FRM-034	Revision 004	Previous Number N/A	Page 8 of 11
-------------------------------	-----------------	------------------------	-----------------

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 27 Feb 2019, 10:18:31 am; Printed by: MSCHROEDER.

The user must ensure that they are using the correct/current revision of this document.
Document: CS-FRM-034 Rev: 004 Effective: 15 Dec 2017

Mega Power Service Center Repair Form, New Faceplate



PART F: Isolation Voltages and Contact Quality Monitor (CQM) Test

Section 1: Isolation Voltages

Table F1: Isolation Voltages

Channel	Directions	Voltage (V)	Criteria (V)
Monopolar A	With a digital multi-meter measure the voltage across the active electrode and one of the hand switching inputs.	11.5	> 10.5
Monopolar B	With a digital multi-meter measure the voltage across the active electrode and one of the hand switching inputs.	11.5	> 10.5
Bipolar	With a digital multi-meter measure the voltage across the active electrode and the ground output.	12.05	> 10.5

Section 2: Verify Potentiometer, PN: 2010127-02, Calibration

Procedure:

- 1) Set the potentiometer to resistances called out in Table F2.
- 2) With a digital multi-meter measure the voltages across the potentiometer.

Table F2: Verify Potentiometer Calibration

Potentiometer Setting	Measured Resistance (Ω)	Criteria (Ω)
10 Ω	10.1	9 - 11
50 Ω	49.9	49 - 51
200 Ω	197.1	196 - 204

Section 3: CQM Test

Table F3: CQM Test

Start up	
With the Mega Power off, connect and set the potentiometer to 50 Ω .	CQM is green when power up is complete: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail
Lower CQM Alarm Trip Point	
Set the potentiometer to 0 Ω and connect to Mega Power.	CQM turns green: <input type="checkbox"/> Pass <input type="checkbox"/> Fail
Raise the resistance until CQM turns red.	Resistance: 16 < 20 Ω
Midlevel CQM Trip Point	
Unplug the potentiometer. Set the potentiometer to 50 Ω and reconnect. Raise the resistance until CQM turns red.	CQM turns green: <input type="checkbox"/> Pass <input checked="" type="checkbox"/> Fail
	Resistance: 65 58 - 72 Ω
Upper CQM Alarm Trip Point	
Unplug the potentiometer. Set the potentiometer to 200 Ω and reconnect. Lower the resistance until CQM turns green.	CQM stays red: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail
	Resistance: 136 122 - 148 Ω

Document Number	Revision	Previous Number	Page
CS-FRM-034	004	N/A	9 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 27 Feb 2019, 10:18:31 am; Printed by: MSCHROEDER.

The user must ensure that they are using the correct/current revision of this document.
Document: CS-FRM-034 Rev: 004 Effective: 15 Dec 2017

Mega Power Service Center Repair Form, New Faceplate



Mega Power Information	
RGA	16780
Serial Number	1678 17502700 1678 02/2019

PART G: Electrical Safety Test

Table G1: Protective Earth Connection to Chassis

Protective Earth ground to outlet.

Current Applied: 25 A	Impedance: 54	Criteria: < 200m Ω
-----------------------	---------------	--------------------

Table G2: Enclosure Leakage Test 1

Protective earth ground to power entry: ground open.

Voltage applied: Mains	Current: 101	Criteria: 50 – 300 μA
------------------------	--------------	-----------------------

Table G3: Enclosure Leakage Test 2

Protective earth ground to power entry: neutral open, and ground open.

Voltage applied: Mains	Current: 190.8	Criteria: < 1000 μA
------------------------	----------------	---------------------

Table G4: Enclosure Leakage Test 3

Protective earth ground to power entry: neutral open, reverse polarity, and ground open.

Voltage applied: Mains	Current: 190.9	Criteria: < 1000 μA
------------------------	----------------	---------------------

Table G5: Enclosure Leakage Test 4

Protective earth ground to power entry: reverse polarity, and ground open.

Voltage applied: Mains	Current: 101.5	Criteria: < 1000 μA
------------------------	----------------	---------------------

Table G6: Patient Leakage

All patient connections to power entry.

Voltage applied: Mains	Current: 4.3	Criteria: < 10 μA
------------------------	--------------	-------------------

Technician Performing Electrical Safety Test

The Mega Power has passed all tests stated in Sections G.

Print: KOHÉI SEKI

Signature: K. seki

Date: 02/27/2019

Document Number	Revision	Previous Number	Page
CS-FRM-034	004	N/A	10 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 27 Feb 2019, 10:18:31 am; Printed by: MSCHROEDER.

The user must ensure that they are using the correct/current revision of this document.
Document: CS-FRM-034 Rev: 004 Effective: 15 Dec 2017

Mega Power Service Center Repair Form, New Faceplate



PART H: Signatures

Technician Performing Tests

All equipment used to perform final test was calibrated. The Mega Power has passed all tests stated in Sections C through Section G and packaged per ENG-WI-035.

Print: L. Sebel

Signature: L. Sebel

Date: 02/29/2019

Megadyne Review

The document has been reviewed. The test data is within the stated criteria.

Print: Bruce Revelone

Signature: Bruce Revelone

Date: 2-27-19

Document Number	Revision	Previous Number	Page
CS-FRM-034	004	N/A	11 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 27 Feb 2019, 10:18:31 am; Printed by: MSCHROEDER.

Megadyne Medical Products
11506 S. State St | Draper, UT 84020


Score: 9.10
Passing: 8/10 or better



Mega Power 1000 Service Center Test


- 1) Rated load is the resistance at which ESU is tested. What is the rated load for:

- a. Pure Cut: 300
- b. Coag: 500
- c. Bipolar: 100



- 2) Mega Power Preventive Maintenance is performed with what form numbers:

- a. New Face Plate: CS-FRM-034
- b. Original Face Plate: CS-FRM-035



- 3) The CQM measures the Resistance between the two pads of a dual plate return electrode.



- 4) The power curves test measures power at different load (s) by means of an ESU analyzer



- 5) Crest Factor is calculated by the ESU analyzer by dividing the Peak Voltage by the RMS Voltage.




- 6) What two modes are found on the new face plate ESU, but not the old face plate ESU? MICRO and MACRO



Legal Entity
123 Street Address | Town, State, Zip

- ✓) When connecting to the ESU Analyzer, which Monopolar Pin do you plug the red lead into (Left, Middle, or Right)? RIGHT



- ✓) True or False: The Potentiometer 2010127-02 calibration needs to be verified before each use. YES

- ✓) The ESU needs to be electrical safety tested to which test standard?

IEC 60601-1

- ✓) True or False: Fields can be left blank on the form if the information is unknown or does not apply. False

Print Technician Name: KOTERI SEKI

Signature: K. Seki Date: 02/27/2019

Megadyne Review: Joe Date: 2/27/2019

Certificate of Training

This certificate is presented to:

Kohei Seki
Johnson and Johnson, Service Center

In recognition for completing service training on Mega Power 1000



Trainer (Tyler Skinner)

Aug. 19, 2019
Date



Trainer (John Minuth)

Aug. 19, 2019
Date

MEGADYNE
PART OF THE JOHNSON & JOHNSON FAMILY OF COMPANIES

ETHICON
PART OF THE JOHNSON & JOHNSON FAMILY OF COMPANIES

3.0 DOCUMENTATION SYSTEM



DOCUMENTATION SYSTEM

Date: October 29, 2019

From: Jason Stivers

Re: 1000

Attention: Factbook # FB003289

Procedures, work instructions, and items that are non-product specific to 1000 already reside in the J&J Medical KK documentation system as they were previously qualified as a service center. Refer to section 3.0 of Factbook F-261 (EpiCenter document Number FB001103). Procedures, work instructions, and items that are product specific to the 1000 have been created and implemented into the J&J Medical KK service process and are listed below.

Megadyne Reference Documents

- 3000144-01 Megadyne Mega Power Field Calibration Manual
- MKT-LBL-063 Megadyne Mega Power Trouble Shooting Guide
- 3000158-01 Megadyne Mega Power Electrosurgical Generator Operators Manual ~ www.e-ifu.com
- 3000159-01 Megadyne Mega Power Electrosurgical Generator Service Manual ~ www.e-ifu.com
- ENG-WI-035 Mega Power 1000 Packaging Instructions, Service and Repair
- ENG-WI-036 Mega Power 1000 Disassembly Instructions, Service and Repair
- ENG-WI-037 Mega Power 1000 Assembly Instructions, Service and Repair
- CS-FRM-034, Mega Power Service Center Repair Form, New Faceplate

E-Sig in EpiCenter

Jason Stivers
Service Engineer, EES – Service Staff Engineer

4.0 EQUIPMENT INSTALLATION QUALIFICATION RESULTS



EIQ RESULTS

Date: October 29, 2019

From: Jason Stivers

Re: 1000

Attention: Factbook # FB003289

All equipment necessary to perform service and repair activities on the Megadyne Mega Power 1000 Electrosurgical Generator have been installed at the J&J K.K. Medical Products Service and Repair Depot, Sukagawa, Japan location. Equipment Installation was conducted per Megadyne Protocol ENG-PRT-437 revision 001 and ENG-PRT-502 revision 001 (see attached).

As described in the service procedures, the product specific test equipment used for the repair and testing of the Megadyne Mega Power 1000 Electrosurgical Generator are as follows:

Custom

Mega Power 1900 Software Calibration Kit
Monopolar Footswitch (1400 or 1400J)
Bipolar Footswitch (1450 or 1450J)
Potentiometer (PN: 2010127-02)

Standard

Fluke QA-ESII Electrosurgery Analyzer (Reference F-124)
Yokogawa Digital Multimeter 733 (Reference F-07)
Kikusui Ground Bond Tester TOS6100 (Reference F-09)
Hioki Leak Current HiTester 3156 (Reference F-56)

The J&J K.K. Medical Products Service and Repair Depot personnel were trained and qualified to use the above-mentioned equipment as referenced in the training section of this Factbook.

The equipment can be found within the J&J K.K. Medical Products Service and Repair Depot calibration system where appropriate and custom tooling or equipment identified has been documented as installed for use. Installation Qualification (IQ) and Performance Qualification (PQ) of equipment was performed. The IQ/PQ stored in local J&J K.K. document control system represents the successful execution of test procedures which qualifies the test equipment being used as properly installed and in good working order.

Specification Sheets and calibration certificates for Standard Equipment listed above are attached to this Factbook.

Part of product testing includes performing electrical safety testing according to appropriate standards as outlined in the service documentation. The equipment used at J&J K.K. meets all applicable standards.

E-Sig in EpiCenter

Jason Stivers
Service Engineer, EES – Service Staff Engineer

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number ENG-PRT-473
	Installation Qualification for Mega Power 1000 Controlled by Service Centers	Revision: 001
		Page 1 of 2

Authored By: John Minuth

Table of Contents

1.0 References	1
2.0 Attachments.....	1
3.0 Scope.....	1
4.0 Purpose	1
5.0 Equipment.....	1
6.0 Passing Criteria.....	1
7.0 Procedure.....	2
8.0 Results and Conclusions.....	2

1.0 References

1000 Mega Power Assembly

2.0 Attachments

Appendix A Protocol Results

3.0 Scope

This protocol is intended to be used with Mega Power 1000's used by service centers as a known good unit.

4.0 Purpose

Service Centers will receive a working Mega Power 1000 that can be used during the investigation of a failure. This installation protocol will be performed upon initial receipt of the Mega Power. It will also be performed if the Mega Power 1000 is disassembled and then reassembled.

5.0 Equipment

ESU Analyzer	Fluke, QA-ESII or equivalent
Digital Multi-meter	Fluke 289 or equivalent
Safety Analyzer	Fluke ESA620 or equivalent

6.0 Passing Criteria

All data acquired from Appendix A is required to be within the listed specification

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number ENG-PRT-473
	Installation Qualification for Mega Power 1000 Controlled by Service Centers	Revision: 001
		Page 2 of 2

7.0 Procedure

See Appendix A for Procedure

8.0 Results and Conclusions

Once Appendix A is completed, the document will be scanned and sent to
SRINTL@Megadyne.com. They will be entered into the Service Center Fact Book.

Megadyne Medical Products, Inc.	TEST PROTOCOL Appendix A	Document Number ENG-PRT-473 Appendix A
	Installation Qualification for Mega Power 1000 Controlled by Service Centers	Revision: 001
		Page 1 of 10

ENG-PRT-473: Appendix A

PART A: Equipment Calibration

Record calibration data:

Table A1: Calibration Data

Equipment	Brand	Model Number	Serial Number	Cal. Due Date
ESU Analyzer	Fluke	QA-ESTI	300091	04/30/2020
Digital Multi-meter	Yokogawa	73302	1501471	03/31/2020
Safety Analyzer	Hioki	3156	051238345	11/30/2019

*All testing must be performed with calibrated equipment. Do not proceed with final testing if any piece of equipment is not current with calibration schedule.

Megadyne Medical Products, Inc.	TEST PROTOCOL Appendix A	Document Number ENG-PRT-473 Appendix A Revision: 001
	Installation Qualification for Mega Power 1000 Controlled by Service Centers	Page 2 of 10

PART A: Display Test**Section 1: Power-Up Test**

Turn on the Mega Power. Once the Mega Power has completed the POST test check the appropriate boxes in Table A1.

Table A1: Power-Up Test

Pass	Fail	Criteria
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Power Switch is illuminated
<input type="checkbox"/>	<input type="checkbox"/>	PURE CUT LED is illuminated
<input checked="" type="checkbox"/>	<input type="checkbox"/>	COAG 1 LED is illuminated
<input type="checkbox"/>	<input type="checkbox"/>	MICRO LED is illuminated
<input type="checkbox"/>	<input type="checkbox"/>	No errors are displayed
<input type="checkbox"/>	<input type="checkbox"/>	All display segments show dashes.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	CQM Alarm LED is red and alarm tone sounded (Without a return pad connected to the Mega Power)

Section 2: Display Functionality**Table A2: Display Functionality**

Pass	Fail	Criteria
<input type="checkbox"/>	<input type="checkbox"/>	Press each mode button and verify that the appropriate LED button lights up.
<input type="checkbox"/>	<input type="checkbox"/>	Press the up arrow button on each of the three channels and verify power goes up.
<input type="checkbox"/>	<input type="checkbox"/>	Press the down arrow button on each of the three channels and verify power goes down.
<input type="checkbox"/>	<input type="checkbox"/>	Press the Bipolar Tone button and verify the LED lights up.
<input type="checkbox"/>	<input type="checkbox"/>	Plug in the Monopolar Return Test Cable and verify CQM alarm LED turns green.
<input type="checkbox"/>	<input type="checkbox"/>	Set all power settings to 1. Turn power off. Turn power on. Press recall and verify power setting of 1 is displayed for each mode.

Megadyne Medical Products, Inc.	TEST PROTOCOL Appendix A	Document Number ENG-PRT-473 Appendix A
	Installation Qualification for Mega Power 1000 Controlled by Service Centers	Revision: 001
		Page 3 of 10

PART B: Monopolar Output Tests**Section 1: Monopolar A, Power Outputs and Crest Factors**

Procedure:

- 1) Set the ESU analyzer to the loads called out in Table B1.
- 2) Set the Mega Power to the power settings (PSET) and modes called out in Table B1.
- 3) Measure and record the appropriate output wattage and crest factors where indicated.
- 4) Verify Mega Power Tones.
- 5) Verify that Error 219 does not get displayed.

Table B1: Monopolar A, Power Outputs and Crest Factors

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)	Crest Factor	Criteria
300	Pure	1	2	0 - 5		1.6 1.3 - 1.9
300	Pure	10	10	5 - 15		
300	Pure	150	152	120 - 180		
300	Pure	300	306	240 - 360		
300	Blend	1	2	0 - 5		1.6 1.3 - 1.9
300	Blend	10	10	5 - 15		
300	Blend	100	104	80 - 120		
300	Blend	200	212	160 - 240		
300	ACE		133	108 - 162	1.6	1.3 - 1.9
Verify Tone for CUT: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			Verify volume can be adjusted: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			
500	Coag 1	1	2	0 - 5		6.8 6.4 - 7.4
500	Coag 1	10	10	5 - 15		
500	Coag 1	60	62	48 - 72		
500	Coag 1	120	124	96 - 144		
500	Coag 2	1	2	0 - 5		6.8 6.4 - 7.4
500	Coag 2	10	10	5 - 15		
500	Coag 2	60	62	48 - 72		
500	Coag 2	120	125	96 - 144	7.0	6.6 - 7.6
500	Spray	1	1	0 - 5		8.1 7.7 - 8.7
500	Spray	10	9	5 - 15		
500	Spray	60	61	48 - 72		
500	Spray	120	123	96 - 144	8.1	7.7 - 8.7
Verify Tone for COAG: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail						
Error 219: Set generator to Coag 1, 30 Watts. With the Generator unloaded (Disconnected from the ESU Analyzer), activate the generator for 10 to 20 seconds. Verify that an error 219 was not displayed.					<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail	

Megadyne Medical Products, Inc.	TEST PROTOCOL Appendix A	Document Number ENG-PRT-473 Appendix A Revision: 001
	Installation Qualification for Mega Power 1000 Controlled by Service Centers	Page 4 of 10

Section 2: Monopolar B, Power Outputs and Crest Factors

Procedure:

- 1) Set the Mega Power up so that Monopolar B is connected to the ESU Analyzer
- 2) Ensure that the activation button activates the Monopolar B electrode.
- 3) Set the ESU analyzer to the loads called out in Table B2.
- 4) Set the Mega Power to the power settings (PSET) and modes called out in Table B2.
- 5) Measure and record the appropriate output wattage and crest factors where indicated.
- 6) Verify Mega Power tones.
- 7) Verify Error 219 does not get displayed

Table B2: Monopolar B, Power Outputs and Crest Factors

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)	Crest Factor	Criteria
300	Pure	1	2	0 - 5		
300	Pure	10	10	5 - 15		
300	Pure	150	155	120 - 180		
300	Pure	300	309	240 - 360	1.6	1.3 - 1.9
300	Blend	1	2	0 - 5		
300	Blend	10	10	5 - 15		
300	Blend	100	105	80 - 120		
300	Blend	200	213	160 - 240	2.9	2.5 - 3.5
300	Ace		155	108 - 162	1.6	1.3 - 1.9
		Verify Tone for Cut:		<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail		
500	Coag 1	1	2	0 - 5		
500	Coag 1	10	10	5 - 15		
500	Coag 1	60	62	48 - 72		
500	Coag 1	120	124	96 - 144	6.8	6.4 - 7.4
500	Coag 2	1	2	0 - 5		
500	Coag 2	10	10	5 - 15		
500	Coag 2	60	62	48 - 72		
500	Coag 2	120	125	96 - 144	7.0	6.6 - 7.6
500	Spray	1	1	0 - 5		
500	Spray	10	9	5 - 15		
500	Spray	60	61	48 - 72		
500	Spray	120	123	96 - 144	8.1	7.7 - 8.7
		Verify Tone for COAG:		<input type="checkbox"/> Pass <input checked="" type="checkbox"/> Fail		
Error 219: Set generator to Coag 1, 20 Watts. With the Generator unloaded (Disconnected from the ESU Analyzer), activate the generator for 10 to 20 seconds. Verify that an error 219 was not displayed.					<input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail	

Megadyne Medical Products, Inc.	TEST PROTOCOL Appendix A				Document Number ENG-PRT-473 Appendix A
	Installation Qualification for Mega Power 1000 Controlled by Service Centers				Revision: 001
					Page 5 of 10

Section 3: Monopolar B Power Curves

Procedure:

- 1) Set the ESU analyzer to the loads called out in Table D3 and Table B4.
- 2) Set the Mega Power to the power settings and modes called out in Table D3 and Table B4.
- 3) Measure and record the appropriate output wattage where indicated.

Table B3: CUT Power Curves

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)
100	ACE		163	120 – 180
200	ACE		155	120 – 180
300	ACE		See Table B2	108 – 162
500	ACE		93	76 – 114
1,000	ACE		93	42 – 64
2,000	ACE		23	18 – 28
100	Pure	300	168	105 – 175
200	Pure	300	281	240 – 360
300	Pure	300	See Table B2	240 – 360
500	Pure	300	303	240 – 360
1,000	Pure	300	295	240 – 360
2,000	Pure	300	292	160 – 240
100	Blend	200	99	75 – 125
200	Blend	200	206	160 – 240
300	Blend	200	See Table B2	160 – 240
500	Blend	200	206	160 – 240
1,000	Blend	200	194	160 – 240
2,000	Blend	200	169	130 – 190

Table B4: COAG Power Curves

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)
100	COAG 1	120	120	96 – 144
200	COAG 1	120	129	96 – 144
500	COAG 1	120	See Table B2	96 – 144
1,000	COAG 1	120	122	92 – 140
2,000	COAG 1	120	120	90 – 135
100	COAG 2	120	124	96 – 144
200	COAG 2	120	131	96 – 144
500	COAG 2	120	See Table B2	96 – 144
1,000	COAG 2	120	119	92 – 138
2,000	COAG 2	120	113	86 – 129
100	Spray	120	126	96 – 144
200	Spray	120	129	96 – 144
500	Spray	120	See Table B2	96 – 144
1,000	Spray	120	114	90 – 135
2,000	Spray	120	109	80 – 120

Megadyne Medical Products, Inc.	TEST PROTOCOL Appendix A	Document Number ENG-PRT-473 Appendix A Revision: 001
	Installation Qualification for Mega Power 1000 Controlled by Service Centers	Page 6 of 10

PART C: Bipolar Output Tests**Section 1: Bipolar Power Outputs and Crest Factors**

Procedure:

- 1) Set the Mega Power up so that Bipolar is connected to the ESU Analyzer
- 2) Ensure that footswitch activates the Bipolar electrodes.
- 3) Set the ESU analyzer to the loads called out in Table C1.
- 4) Set the Mega Power to the power settings (PSET) and modes called out in Table C1.
- 5) Measure and record the appropriate output wattage and crest factors where indicated.
- 6) Verify Mega Power Tones.

Table C1: Bipolar, Power Outputs and Crest Factors

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)	Crest Factor	Criteria
100	Micro	1	1	0 - 5		
100	Micro	10	10	5 - 15		
100	Micro	40	38	32 - 48		
100	Micro	80	79	64 - 96	1.6	1.4 - 2.0
100	Macro	1	1	0 - 5		
100	Macro	10	8	5 - 15		
100	Macro	40	38	32 - 48		
100	Macro	80	76	64 - 96	1.6	1.4 - 2.0
Verify normal tone*: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			Verify pulsing tone*: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			

* When bipolar is activated the Mega Power displays the estimated current flow in the current flow indicator window. In addition to the current flow indicator, the Mega Power can be set to indicate the current with an audio alarm. The current is proportional to the pulsing rate of the audio tone. By pressing the bipolar tone button (the button to the right of the current flow window) the bipolar alarm will toggle between a continuous tone and a pulsing tone. Verify that both the continuous and the pulsing tones are operational in the Macro Mode with a 100 Ω load, and the Mega Power set to 80 W.

Megadyne Medical Products, Inc.	TEST PROTOCOL Appendix A	Document Number ENG-PRT-473 Appendix A
	Installation Qualification for Mega Power 1000 Controlled by Service Centers	Revision: 001
		Page 7 of 10

Section 2: Bipolar Power Curves

Procedure:

- 1) Set the ESU analyzer to the loads called out in Table C2.
- 2) Set the Mega Power to the power settings (PSET) and modes called out in Table C2.
- 3) Measure and record the appropriate output wattage where indicated.

Table C2: BIPOLAR Power Curves

Load (Ω)	Mode	PSET (watts)	Output (Watts)	Criteria (Watts)
50	Micro	80	62	46 – 68
75	Micro	80	78	60 – 90
100	Micro	80	See Table C1	64 – 96
200	Micro	80	75	64 – 96
650	Micro	80	8	5 – 15
1000	Micro	80	6	0 – 10
50	Macro	80	62	46 – 68
75	Macro	80	78	60 – 90
100	Macro	80	See Table C1	64 – 96
200	Macro	80	75	64 – 96
500	Macro	80	54	38 – 58
1000	Macro	80	29	20 – 32

Megadyne Medical Products, Inc.	TEST PROTOCOL Appendix A	Document Number ENG-PRT-473 Appendix A Revision: 001
	Installation Qualification for Mega Power 1000 Controlled by Service Centers	Page 8 of 10

PART D: Isolation Voltages and Contact Quality Monitor (CQM) TestSection 1: Isolation Voltages**Table D1: Isolation Voltages**

Channel	Directions	Voltage (V)	Criteria (V)
Monopolar A	With a digital multi-meter measure the voltage across the active electrode and one of the hand switching inputs.	12.61	> 10.5
Monopolar B	With a digital multi-meter measure the voltage across the active electrode and one of the hand switching inputs.	12.29	> 10.5
Bipolar	With a digital multi-meter measure the voltage across the active electrode and the ground output.	15.02	> 10.5

Section 2: Verify Potentiometer, PN: 2010127-02, Calibration

Procedure:

- 1) Set the potentiometer to resistances called out in Table D2.
- 2) With a digital multi-meter measure the voltages across the potentiometer.

Table D2: Verify Potentiometer Calibration

Potentiometer Setting	Measured Resistance (Ω)	Criteria (Ω)
10 Ω	9.2	9 - 11
50 Ω	50.3	49 - 51
200 Ω	201.6	196 - 204

Section 3: CQM Test**Table D3: CQM Test**

Start up With the Mega Power off, connect and set the potentiometer to 50 Ω .		CQM is green when power up is complete: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail
Lower CQM Alarm Trip Point Set the potentiometer to 0 Ω and connect to Mega Power. Raise the resistance until CQM turns red.		CQM turns green: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail Resistance: 17 < 20 Ω
Midlevel CQM Trip Point Unplug the potentiometer. Set the potentiometer to 50 Ω and reconnect. Raise the resistance until CQM turns red.		CQM turns green: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail Resistance: 69 58 - 72 Ω
Upper CQM Alarm Trip Point Unplug the potentiometer. Set the potentiometer to 200 Ω and reconnect. Lower the resistance until CQM turns green.		CQM stays red: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail Resistance: 143 122 - 148 Ω

Megadyne Medical Products, Inc.	TEST PROTOCOL Appendix A	Document Number ENG-PRT-473 Appendix A
	Installation Qualification for Mega Power 1000 Controlled by Service Centers	Revision: 001
		Page 9 of 10

PART E: Electrical Safety Test**Table E1: Protective Earth Connection to Chassis**

Protective Earth ground to outlet.

Current Applied: 25 A	Impedance: $1.0\ \Omega$	Criteria: < 200m Ω
-----------------------	--------------------------	---------------------------

Table E2: Enclosure Leakage Test 1

Protective earth ground to power entry: ground open.

Voltage applied: Mains	Current: $8.8\ \mu A$	Criteria: 50 – 300 μA
------------------------	-----------------------	----------------------------

Table E3: Enclosure Leakage Test 2

Protective earth ground to power entry: neutral open, and ground open.

Voltage applied: Mains	Current: $7.63\ \mu A$	Criteria: < 1000 μA
------------------------	------------------------	--------------------------

Table E4: Enclosure Leakage Test 3

Protective earth ground to power entry: neutral open, reverse polarity, and ground open.

Voltage applied: Mains	Current: $7.84\ \mu A$	Criteria: < 1000 μA
------------------------	------------------------	--------------------------

Table E5: Enclosure Leakage Test 4

Protective earth ground to power entry: reverse polarity, and ground open.

Voltage applied: Mains	Current: $8.7\ \mu A$	Criteria: < 1000 μA
------------------------	-----------------------	--------------------------

Table E6: Patient Leakage

All patient connections to power entry.

Voltage applied: Mains	Current: $3\ \mu A$	Criteria: < 10 μA
------------------------	---------------------	------------------------

Megadyne Medical Products, Inc.	TEST PROTOCOL Appendix A	Document Number ENG-PRT-473 Appendix A Revision: 001
	Installation Qualification for Mega Power 1000 Controlled by Service Centers	Page 10 of 10

PART F: Signatures**Technician Performing Tests**

All equipment used to perform final test was calibrated. The Mega Power has passed all tests stated in Sections B through Section F.

Print: KOHEI SELCI


Signature: K. selci

Date: 10/09/2019

Megadyne Review

The document has been reviewed. The test data is within the stated criteria.

Print: JASON STIVERS

Signature: 

Date: 10/25/2019

SIGNING AS SERVICE
ENGINEER FOR MEGA POWER 1000
AT ETHICON ENDO-SURGERY

The user must ensure that they are using the correct/current revision of this document.
Document: ENG-PRT-502 Rev: 001 Effective: 13 Jun 2018

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number ENG-PRT-502
	Performance/ Operational Qualification: Service Center for Mega Power 1000	Revision: 001
		Page 1 of 2

Authored By: John Minuth

1.0 Purpose

This document is intended to qualify a service center to repair the Mega Power 1000. The service center will have to be able to disassemble, reassemble, and final test a generator.

2.0 References

1000	Mega Power 1000
CS-SOP-001	Service and Repair SOP
ENG-WI-036	Mega Power 1000 Disassembly Instructions, Service and Repair
ENG-WI-037	Mega Power 1000 Assembly Instructions, Service and Repair
CS-FRM-034	Mega Power Service Center Repair Form, New Faceplate

3.0 Attachments

None

4.0 Scope

The scope of this protocol is any J&J Service Center that will repair the Mega Power 1000 both old and new faceplates.

5.0 Strategies/Assumption

- 5.1 A technician from each service center will follow the instructions in section 6 for one (1) unit.
- 5.2 Each service center to be qualified will conduct this protocol separately.
- 5.3 Although the old faceplate and new faceplates are different and there is a different Mega Power Service Center Repair form for the old faceplate, the process for disassembly, re-assembly, packaging and repair of the Mega Power is the same so the qualification will only be performed on a Mega Power with new faceplate.
- 5.4 Each of the service centers receiving this protocol have been trained.
- 5.5 Although no repair is planned to be performed as part of this qualification it is assumed that a full disassembly and re-assembly of the unit per the referenced work instructions is representative of the repair process.

The user must ensure that they are using the correct/current revision of this document.
Document: ENG-PRT-502 Rev: 001 Effective: 13 Jun 2018

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number ENG-PRT-502
	Performance/ Operational Qualification: Service Center for Mega Power 1000	Revision: 001
		Page 2 of 2

6.0 Procedural Steps

- 6.1 Steps to be conducted by J&J Service Center
 - 6.1.1 Complete Part A on CS-FRM-034
 - 6.1.2 Part B of CS-FRM-034 is not applicable to this protocol since no repair is planned. Mark each table with an N/A.
 - 6.1.3 Disassemble a Mega Power 1000 per ENG-WI-036
 - 6.1.4 Reassemble a Mega Power 1000 per ENG-WI-037
 - 6.1.5 Test the Mega Power 1000 per CS-FRM-034, Part C through H.
Complete all required fields except for the Megadyne Review of Section H.
 - 6.1.6 Scan and send completed CS-FRM-034 to Megadyne.
SRINTL@megadyne.com
Attn: Carly Rasband
- 6.2 Steps to be conducted by Megadyne
 - 6.2.1 Review completed form CS-FRM-034 sent by the service center and complete the Megadyne Review of Section H.
 - 6.2.2 Document results of each J&J Service Center qualification within a report.
 - 6.2.3 Send completed form CS-FRM-034 to Ethicon for inclusion in the factbook
Attn: Jamie Best
4545 Creek Road
Mail Location 120A
Cincinnati, Ohio 45242

7.0 Acceptance Criteria

The acceptance criterion is that after disassembly and re-assembly, the unit under test passes all the test requirements called out in CS-FRM-034.

Mega Power Service Center Repair Form, New Faceplate



Mega Power Information	
RGA	Kuseki 11/09/2019
Serial Number	182955001
RoHS	<input type="checkbox"/> Leaded <input checked="" type="checkbox"/> RoHS

PART A: Equipment Calibration

Record calibration data:

Table A1: Calibration Data

Equipment	Brand	Model Number	Serial Number	Cal. Due Date
ESU Analyzer	Fluke	8A-F5II	300091	09/30/2020
Digital Multi-meter	Yokogawa	73302	1501471	03/31/2020
Safety Analyzer	Hirochi	3156	051238345	11/30/2019

*All testing must be performed with calibrated equipment. Do not proceed with final testing if any piece of equipment is not current with calibration schedule.

Part B: Investigation and Repair**Section 1: Before the Mega Power is disassembled**

- With the ESU analyzer set to 300 Ohms, measure the output power with the Mega Power set to PURE CUT and the power set to 30 W, 75 W, and 300 W. Record the results in Table B1.
- With the ESU analyzer set to 500 Ohms, measure the output power with the Mega Power set to COAG 1 and the power set to 30 W, 50 W, and 120 W. Record the results in Table B2.
- Using Field Calibration Software, 1410010-01, download, print, and attach the error logs from the master processor and watch dog processor.

Table B1: PURE CUT

Setting	Power Output
300 Ω	30 W
300 Ω	75 W
300 Ω	300 W

Table B2: COAG 1

Setting	Power Output
500 Ω	30 W
500 Ω	50 W
500 Ω	120 W

Section 2: Investigation**Table B3: Investigation**

The Unit is good working, K.G.V

K.G.V 11/09/2019

Initial/Date: K.S 11/09/2019

Document Number
CS-FRM-034

Revision
004

Previous Number
N/A

Page
1 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Mega Power Service Center Repair Form, New Faceplate



Section 3: Repair Actions

Table B4: Repair Actions

huhe

Verdict 10/08/2011

Initial/Date:
E.S 10/08/2019

Table B5: Replacement Parts (Confirm all replacement parts with CS-FRM-038)

Ship Printed Circuit Boards to Megadyne

Ship replaced parts
to Megadyne for evaluation.

Megadyne Medical Product, Inc.
Attn: Service and Repair
11506 South State St.
Draper, Utah 84020
USA

Table B6: Shipped to Megadyne

Part Number	Description	PCB Number (Where Applicable)
	K-5025 10/08/2019	
Shipped By:	K-5025 10/08/2019	Date:
Tracking Number:		

Document Number CS-FRM-034	Revision 004	Previous Number N/A	Page 2 of 11
-------------------------------	-----------------	------------------------	-----------------

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Mega Power Service Center Repair Form, New Faceplate



Once the repair work has been completed and the Mega Power has been reassembled, test per Part C through Part G.

PART C: Display Test

Section 1: Power-Up Test

Turn on the Mega Power. Once the Mega Power has completed the POST test check the appropriate boxes in Table C1.

Table C1: Power-Up Test

Pass	Fail	Criteria
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Power Switch is illuminated
<input checked="" type="checkbox"/>	<input type="checkbox"/>	PURE CUT LED is illuminated
<input checked="" type="checkbox"/>	<input type="checkbox"/>	COAG 1 LED is illuminated
<input checked="" type="checkbox"/>	<input type="checkbox"/>	MICRO LED is illuminated
<input checked="" type="checkbox"/>	<input type="checkbox"/>	No errors are displayed
<input checked="" type="checkbox"/>	<input type="checkbox"/>	All display segments show dashes.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	CQM Alarm LED is red and alarm tone sounded (Without a return pad connected to the Mega Power)

Section 2: Display Functionality

Table C2: Display Functionality

Pass	Fail	Criteria
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Press each mode button and verify that the appropriate LED button lights up.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Press the up arrow button on each of the three channels and verify power goes up.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Press the down arrow button on each of the three channels and verify power goes down.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Press the Bipolar Tone button and verify the LED lights up.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Plug in the Monopolar Return Test Cable and verify CQM alarm LED turns green.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Set all power settings to 1. Turn power off. Turn power on. Press recall and verify power setting of 1 is displayed for each mode.

Document Number CS-FRM-034	Revision 004	Previous Number N/A	Page 3 of 11
-------------------------------	-----------------	------------------------	-----------------

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Mega Power Service Center Repair Form, New Faceplate**PART D: Monopolar Output Tests****Section 1: Monopolar A, Power Outputs and Crest Factors****Procedure:**

- 1) Set the ESU analyzer to the loads called out in Table D1.
- 2) Set the Mega Power to the power settings (PSET) and modes called out in Table D1.
- 3) Measure and record the appropriate output wattage and crest factors where indicated.
- 4) Verify Mega Power Tones.

Table D1: Monopolar A, Power Outputs and Crest Factors

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)	Crest Factor	Criteria
300	Pure	1	2	0 - 5		1.6
300	Pure	10	10	5 - 15		
300	Pure	150	152	120 - 180		
300	Pure	300	306	240 - 360		
300	Blend	1	2	0 - 5		2.5 - 3.5
300	Blend	10	10	5 - 15		
300	Blend	100	104	80 - 120		
300	Blend	200	212	160 - 240		
300	ACE		133	108 - 162	1.6	1.3 - 1.9
Verify Tone for CUT: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			Verify volume can be adjusted: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			
500	Coag 1	1	2	0 - 5		6.8
500	Coag 1	10	10	5 - 15		
500	Coag 1	60	62	48 - 72		
500	Coag 1	120	124	96 - 144		
500	Coag 2	1	2	0 - 5		6.6 - 7.6
500	Coag 2	10	10	5 - 15		
500	Coag 2	60	62	48 - 72		
500	Coag 2	120	125	96 - 144		
500	Spray	1	1	0 - 5		7.7 - 8.7
500	Spray	10	9	5 - 15		
500	Spray	60	61	48 - 72		
500	Spray	120	123	96 - 144		
Verify Tone for COAG: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail						

Document Number
CS-FRM-034Revision
004Previous Number
N/APage
4 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Mega Power Service Center Repair Form, New Faceplate**Section 2: Monopolar B, Power Outputs and Crest Factors**

Procedure:

- 1) Set the Mega Power up so that Monopolar B is connected to the ESU Analyzer
- 2) Ensure that the activation button activates the Monopolar B electrode.
- 3) Set the ESU analyzer to the loads called out in Table D2.
- 4) Set the Mega Power to the power settings (PSET) and modes called out in Table D2.
- 5) Measure and record the appropriate output wattage and crest factors where indicated.
- 6) Verify Mega Power tones.

Table D2: Monopolar B, Power Outputs and Crest Factors

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)	Crest Factor	Criteria
300	Pure	1	2	0 - 5		1.3 - 1.9
300	Pure	10	10	5 - 15		
300	Pure	150	155	120 - 180		
300	Pure	300	301	240 - 360		
300	Blend	1	2	0 - 5		2.5 - 3.5
300	Blend	10	10	5 - 15		
300	Blend	100	105	80 - 120		
300	Blend	200	213	160 - 240		
300	Ace		135	108 - 162	1.6	1.3 - 1.9
		Verify Tone for Cut: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail				
500	Coag 1	1	2	0 - 5		6.4 - 7.4
500	Coag 1	10	10	5 - 15		
500	Coag 1	60	62	48 - 72		
500	Coag 1	120	124	96 - 144		
500	Coag 2	1	2	0 - 5		6.6 - 7.6
500	Coag 2	10	10	5 - 15		
500	Coag 2	60	62	48 - 72		
500	Coag 2	120	125	96 - 144		
500	Spray	1	1	0 - 5		7.7 - 8.7
500	Spray	10	9	5 - 15		
500	Spray	60	61	48 - 72		
500	Spray	120	123	96 - 144		
		Verify Tone for COAG: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail				

Document Number CS-FRM-034	Revision 004	Previous Number N/A	Page 5 of 11
-------------------------------	-----------------	------------------------	-----------------

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Mega Power Service Center Repair Form, New Faceplate**Section 3: Monopolar B Power Curves****Procedure:**

- 1) Set the ESU analyzer to the loads called out in Table D3 and Table D4.
- 2) Set the Mega Power to the power settings and modes called out in Table D3 and Table D4.
- 3) Measure and record the appropriate output wattage where indicated.

Table D3: CUT Power Curves

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)
100	ACE		163	120 – 180
200	ACE		155	120 – 180
300	ACE		See Table D2	108 – 162
500	ACE		93	76 – 114
1,000	ACE		53	42 – 64
2,000	ACE		23	18 – 28
100	Pure	300	168	105 – 175
200	Pure	300	281	240 – 360
300	Pure	300	See Table D2	240 – 360
500	Pure	300	323	240 – 360
1,000	Pure	300	295	240 – 360
2,000	Pure	300	202	160 – 240
100	Blend	200	99	75 – 125
200	Blend	200	206	160 – 240
300	Blend	200	See Table D2	160 – 240
500	Blend	200	206	160 – 240
1,000	Blend	200	194	160 – 240
2,000	Blend	200	169	130 – 190

Table D4: COAG Power Curves

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)
100	COAG 1	120	120	96 – 144
200	COAG 1	120	129	96 – 144
500	COAG 1	120	See Table D2	96 – 144
1,000	COAG 1	120	122	92 – 140
2,000	COAG 1	120	120	90 – 135
100	COAG 2	120	124	96 – 144
200	COAG 2	120	131	96 – 144
500	COAG 2	120	See Table D2	96 – 144
1,000	COAG 2	120	119	92 – 138
2,000	COAG 2	120	113	86 – 129
100	Spray	120	126	96 – 144
200	Spray	120	129	96 – 144
500	Spray	120	See Table D2	96 – 144
1,000	Spray	120	114	90 – 135
2,000	Spray	120	104	80 – 120

Document Number CS-FRM-034	Revision 004	Previous Number N/A	Page 6 of 11
-------------------------------	-----------------	------------------------	-----------------

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Mega Power Service Center Repair Form, New Faceplate



PART E: Bipolar Output Tests

Section 1: Bipolar Power Outputs and Crest Factors

Procedure:

- 1) Set the Mega Power up so that Bipolar is connected to the ESU Analyzer
- 2) Ensure that footswitch activates the Bipolar electrodes.
- 3) Set the ESU analyzer to the loads called out in Table E1.
- 4) Set the Mega Power to the power settings (PSET) and modes called out in Table E1.
- 5) Measure and record the appropriate output wattage and crest factors where indicated.
- 6) Verify Mega Power Tones.

Table E1: Bipolar, Power Outputs and Crest Factors

Load (Ω)	Mode	PSET (Watts)	Output (Watts)	Criteria (Watts)	Crest Factor	Criteria
100	Micro	1	1	0 - 5	1.6	1.4 – 2.0
100	Micro	10	10	5 - 15		
100	Micro	40	38	32 - 48		
100	Micro	80	79	64 - 96		
100	Macro	1	1	0 - 5	1.6	1.4 – 2.0
100	Macro	10	8	5 - 15		
100	Macro	40	38	32 - 48		
100	Macro	80	76	64 - 96		
Verify normal tone*: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			Verify pulsing tone*: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail			

* When bipolar is activated the Mega Power displays the estimated current flow in the current flow indicator window. In addition to the current flow indicator, the Mega Power can be set to indicate the current with an audio alarm. The current is proportional to the pulsing rate of the audio tone. By pressing the bipolar tone button (the button to the right of the current flow window) the bipolar alarm will toggle between a continuous tone and a pulsing tone. Verify that both the continuous and the pulsing tones are operational in the Macro Mode with a 100 Ω load, and the Mega Power set to 80 W.

Document Number CS-FRM-034	Revision 004	Previous Number N/A	Page 7 of 11
-------------------------------	-----------------	------------------------	-----------------

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Mega Power Service Center Repair Form, New Faceplate**Section 2: Bipolar Power Curves****Procedure:**

- 1) Set the ESU analyzer to the loads called out in Table E2.
- 2) Set the Mega Power to the power settings (PSET) and modes called out in Table E2.
- 3) Measure and record the appropriate output wattage where indicated.

Table E2: BIPOLAR Power Curves

Load (Ω)	Mode	PSET (watts)	Output (Watts)	Criteria (Watts)
50	Micro	80	62	46 – 68
75	Micro	80	78	60 – 90
100	Micro	80	See Table E1	64 – 96
200	Micro	80	75	64 – 96
650	Micro	80	8	5 – 15
1000	Micro	80	6	0 – 10
50	Macro	80	62	46 – 68
75	Macro	80	78	60 – 90
100	Macro	80	See Table E1	64 – 96
200	Macro	80	75	64 – 96
500	Macro	80	54	38 – 58
1000	Macro	80	29	20 – 32

Document Number
CS-FRM-034Revision
004Previous Number
N/APage
8 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Mega Power Service Center Repair Form, New Faceplate**PART F: Isolation Voltages and Contact Quality Monitor (CQM) Test****Section 1: Isolation Voltages****Table F1: Isolation Voltages**

Channel	Directions	Voltage (V)	Criteria (V)
Monopolar A	With a digital multi-meter measure the voltage across the active electrode and one of the hand switching inputs.	12.6	> 10.5
Monopolar B	With a digital multi-meter measure the voltage across the active electrode and one of the hand switching inputs.	12.29	> 10.5
Bipolar	With a digital multi-meter measure the voltage across the active electrode and the ground output.	13.02	> 10.5

Section 2: Verify Potentiometer, PN: 2010127-02, Calibration**Procedure:**

- 1) Set the potentiometer to resistances called out in Table F2.
- 2) With a digital multi-meter measure the voltages across the potentiometer.

Table F2: Verify Potentiometer Calibration

Potentiometer Setting	Measured Resistance (Ω)	Criteria (Ω)
10 Ω	12.2	9 - 11
50 Ω	50.3	49 - 51
200 Ω	201.6	196 - 204

Section 3: CQM Test**Table F3: CQM Test**

Start up	
With the Mega Power off, connect and set the potentiometer to 50 Ω .	CQM is green when power up is complete: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail
Lower CQM Alarm Trip Point	
Set the potentiometer to 0 Ω and connect to Mega Power. Raise the resistance until CQM turns red.	CQM turns green: <input type="checkbox"/> Pass <input checked="" type="checkbox"/> Fail Resistance: 17 < 20 Ω
Midlevel CQM Trip Point	
Unplug the potentiometer. Set the potentiometer to 50 Ω and reconnect. Raise the resistance until CQM turns red.	CQM turns green: <input type="checkbox"/> Pass <input checked="" type="checkbox"/> Fail Resistance: 69 58 - 72 Ω
Upper CQM Alarm Trip Point	
Unplug the potentiometer. Set the potentiometer to 200 Ω and reconnect. Lower the resistance until CQM turns green.	CQM stays red: <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail Resistance: 143 122 - 148 Ω

Document Number	Revision	Previous Number	Page
CS-FRM-034	004	N/A	9 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Mega Power Service Center Repair Form, New Faceplate

Mega Power Information	
RGA	Kohsei 11/09/2019
Serial Number	18295500

PART G: Electrical Safety Test**Table G1: Protective Earth Connection to Chassis**

Protective Earth ground to outlet.

Current Applied: 25 A	Impedance: 10	Criteria: < 200m Ω
-----------------------	---------------	--------------------

Table G2: Enclosure Leakage Test 1

Protective earth ground to power entry: ground open.

Voltage applied: Mains	Current: 88	Criteria: 50 – 300 μA
------------------------	-------------	-----------------------

Table G3: Enclosure Leakage Test 2

Protective earth ground to power entry: neutral open, and ground open.

Voltage applied: Mains	Current: 7.63	Criteria: < 1000 μA
------------------------	---------------	---------------------

Table G4: Enclosure Leakage Test 3

Protective earth ground to power entry: neutral open, reverse polarity, and ground open.

Voltage applied: Mains	Current: 7.84	Criteria: < 1000 μA
------------------------	---------------	---------------------

Table G5: Enclosure Leakage Test 4

Protective earth ground to power entry: reverse polarity, and ground open.

Voltage applied: Mains	Current: 89	Criteria: < 1000 μA
------------------------	-------------	---------------------

Table G6: Patient Leakage

All patient connections to power entry.

Voltage applied: Mains	Current: 3	Criteria: < 10 μA
------------------------	------------	-------------------

Technician Performing Electrical Safety Test

The Mega Power has passed all tests stated in Sections G.

Print: Kohsei SEKI

Signature:

Date: 11/09/2019

Document Number	Revision	Previous Number	Page
CS-FRM-034	004	N/A	10 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Mega Power Service Center Repair Form, New Faceplate

PART H: Signatures**Technician Performing Tests**

All equipment used to perform final test was calibrated. The Mega Power has passed all tests stated in Sections C through Section G and packaged per ENG-WI-035.

Print: KOHÉI SEKI

Signature:

Date: 10/08/2019

Megadyne Review

The document has been reviewed. The test data is within the stated criteria.

Print: JASON STINERS

Signature:

Date: 10/25/2019

SIGNING AS
SERVICE ENGINEER
FOR MEGA POWER 1000
AT ETHICON END-SURGERY

Document Number	Revision	Previous Number	Page
CS-FRM-034	004	N/A	11 of 11

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Title : Error Log

Kureki 08/01/2019

Page 1 Thursday, May 30, 2019

Master Error Log

Start Index = 0

#	Err	Mode	VSI	VSU	MP	Pset	DAC	Vout	Iout	Z	Pmeas
00	255	000	000	000	000	000	0.000	0.000	0	0
41	255	255	255	255	255	255	0.000	0.000	0	0
40	255	255	255	255	255	255	0.000	0.000	0	0
39	255	255	255	255	255	255	0.000	0.000	0	0
38	255	255	255	255	255	255	0.000	0.000	0	0
37	255	255	255	255	255	255	0.000	0.000	0	0
36	255	255	255	255	255	255	0.000	0.000	0	0
35	255	255	255	255	255	255	0.000	0.000	0	0
34	255	255	255	255	255	255	0.000	0.000	0	0
33	255	255	255	255	255	255	0.000	0.000	0	0
32	255	255	255	255	255	255	0.000	0.000	0	0
31	255	255	255	255	255	255	0.000	0.000	0	0
30	255	255	255	255	255	255	0.000	0.000	0	0
29	255	255	255	255	255	255	0.000	0.000	0	0
28	255	255	255	255	255	255	0.000	0.000	0	0
27	255	255	255	255	255	255	0.000	0.000	0	0
26	255	255	255	255	255	255	0.000	0.000	0	0
25	255	255	255	255	255	255	0.000	0.000	0	0
24	255	255	255	255	255	255	0.000	0.000	0	0
23	255	255	255	255	255	255	0.000	0.000	0	0
22	255	255	255	255	255	255	0.000	0.000	0	0
21	255	255	255	255	255	255	0.000	0.000	0	0
20	255	255	255	255	255	255	0.000	0.000	0	0
19	255	255	255	255	255	255	0.000	0.000	0	0
18	255	255	255	255	255	255	0.000	0.000	0	0
17	255	255	255	255	255	255	0.000	0.000	0	0
16	255	255	255	255	255	255	0.000	0.000	0	0
15	255	255	255	255	255	255	0.000	0.000	0	0
14	255	255	255	255	255	255	0.000	0.000	0	0
13	255	255	255	255	255	255	0.000	0.000	0	0
12	255	255	255	255	255	255	0.000	0.000	0	0
11	255	255	255	255	255	255	0.000	0.000	0	0
10	255	255	255	255	255	255	0.000	0.000	0	0
09	255	255	255	255	255	255	0.000	0.000	0	0
08	255	255	255	255	255	255	0.000	0.000	0	0
07	255	255	255	255	255	255	0.000	0.000	0	0
06	255	255	255	255	255	255	0.000	0.000	0	0
05	255	255	255	255	255	255	0.000	0.000	0	0
04	255	255	255	255	255	255	0.000	0.000	0	0
03	255	255	255	255	255	255	0.000	0.000	0	0
02	255	255	255	255	255	255	0.000	0.000	0	0
01	255	255	255	255	255	255	0.000	0.000	0	0

SN: 18295500 |

V/L

Title : ESU Test Results from lock Staffing

Thursday, May 30, 2019

k.selsi 08/01/2019

Date = 05-30-2019 Time = 23:09:41

User ID: k.s
SN: 18295001

Monopolar A measurements:

```
132.58 ;ACE cut power measurement PASS
1.60 ;ACE Crest Factor PASS
307.26 ;PURE cut 300 Watts measurement PASS
1.60 ;PURE Cut Crest Factor PASS
210.72 ;BLEND cut 200 Watts measurement PASS
2.90 ;BLEND Cut Crest Factor PASS
124.10 ;Standard Coag 120 Watts measurement PASS
6.80 ;COAG Standard Crest Factor PASS
123.60 ;Spray Coag 120 Watts measurement PASS
8.20 ;COAG Spray Crest Factor PASS
124.70 ;CWC Coag 120 Watts measurement PASS
7.00 ;COAG CWC Crest Factor PASS
```

Monopolar B ACE measurements:

```
134.50 ;ACE cut power measurement PASS
1.60 ;ACE Crest Factor PASS
308.88 ;PURE cut 300 Watts measurement PASS
1.60 ;PURE Cut Crest Factor PASS
212.60 ;BLEND cut 200 Watts measurement PASS
2.90 ;BLEND Cut Crest Factor PASS
123.86 ;Standard Coag 120 Watts measurement PASS
6.78 ;COAG Standard Crest Factor PASS
123.30 ;Spray Coag 120 Watts measurement PASS
8.16 ;COAG Spray Crest Factor PASS
124.42 ;CWC Coag 120 Watts measurement PASS
7.00 ;COAG CWC Crest Factor PASS

78.70 ;Macro Bipolar 80 Watts measurement PASS
1.60 ;Bipolar Macro Crest Factor PASS

78.78 ;Standard Bipolar 80 Watts measurement PASS
1.60 ;Bipolar Standard Crest Factor PASS
```

51 ;CQM Reading (ohms) PASS

V1

FLUKE®

Biomedical

QA-ES Series II

Electrosurgery Analyzer

Technical Data



QA-ES Series II analyzes electrosurgical units quickly and accurately.

A wide load-resistance range provides 128 user-selectable loads, including very low loads for testing many of today's ESUs.

An accuracy of $\pm 2\%$ of reading down to 20 mA guarantees reliable high-frequency leakage results. With capability to run an automatic-power-distribution test in as little as 1 minute, the QA-ES works fast so technicians save time.

An Ansur QA-ES software plug-in allows users to create and automatically run tests, capture data, and produce easy-to-read reports with a PC.

Key features

- Automatic power distribution measurement, including power, current, peak-to-peak voltage (closed load only), and crest factor
- Oscilloscope output
- High-frequency leakage measurements with accuracy of $\pm 2\%$ of reading
- 128 internal user-selectable test loads from 10Ω to 5200Ω
- Foot-switch output for triggering the ESU under test
- Ansur QA-ES software plug-in for automated test protocols
- Large display
- RS-232 and Centronic-Printer interface

Technical specifications

Generator output

Continuous operation

Continuous measurement of power, current, peak-to-peak voltage (closed load only), and crest factor

Single operation

Single measurement after the set delay time of the ESU output of power, current, peak-to-peak voltage (closed load only), and crest factor

Power distribution

Automatic measurement of power, current, peak-to-peak voltage (closed load only), and crest factor through a user-selectable load range

RF leakage current

Provides connections and load configurations to measure HF leakage from both grounded and isolated equipment

RECQM

Test the "return electrode control quality monitoring" using the QA-ES internal loads

Modes of operation

Manual or remote controlled (via Ansur)

Measurement

True-rms value of applied waveform

RMS bandwidth

30 Hz to 10 MHz (-3 dB) for instrumentation only
30 Hz to 2.5 MHz (-3 dB) with loads

Low frequency filter

100 Hz filter to avoid low-frequency disturbance or interference

Current

20 mA to 2200 mA

Current accuracy

20 mA to 2200 mA \pm 2 % of reading

Load resistance

10 Ω to 2500 Ω in step of 25 Ω (@ dc)
2500 Ω to 5200 Ω in step of 100 Ω (@ dc)



Additional fixed load

200 Ω 400 W for 30 s; max 15 % duty cycle

Crest factor

The higher of the two peak voltage measurements is used for computation

Range

1.4 to 16 (V peak/V rms)

Foot-switch output

The foot switch output can be used to trigger the electrosurgical unit

Peak-to-peak voltage

0 kV to 10 kV (closed load only)
accuracy: \pm 10 %

Note: Measurement is taken between the active and dispersive electrodes with closed load only.

Oscilloscope output

5 V/A uncalibrated, 100 mA RF current minimum input

Ansur QA-ES Plug-In

Remote control

All functions and tests in QA-ES may be performed from the PC

User-programmable test sequences

Allows unlimited numbers of test sequences with user-programmable templates and test limits. These tests include power distribution test, output test, HF leakage, and RECQM verification

Storage and recall

Protocol formats and data may be stored, recalled, printed out, or transferred

**Temperature****Operating**

15 °C to 35 °C (59 °F to 95 °F)

Storage

0 °C to 50 °C (32 °F to 122 °F)

Display

LCD graphic display

Alphanumeric format

8 lines x 40 characters

Graphic mode

240 x 64 pixel matrix

Display control

5 F-keys, enter, cancel, control knob, and an encoder

Data input/outputs

Parallel printer port and bidirectional RS-232

Power

115/230 V ac; 48 Hz to 66 Hz, 35 VA

Housing

Metal case

Dimensions (WxDxH)39.5 cm x 34.2 cm x 13.2 cm
(15.6 in x 13.5 in x 5.2 in)**Weight**

9.8 kg (21.6 lb)

Ordering Information**Models**

- 2649769** QA-ES Series II 115 V Electrosurgery Analyzer (US)
- 2651725** QA-ES Series II 230 V Electrosurgery Analyzer (Schuko)
- 2770445** QA-ES Series II 230 V Electrosurgery Analyzer (UK)
- 2770450** QA-ES Series II 230 V Electrosurgery Analyzer (Australia)
- 3096390** QA-ES Series II 100 V Electrosurgery Analyzer (Japan)

Standard accessories

- 2716044** Manual on CD
- 2716032** Manual Hard Copy
- 2772171** ESU-Dispersive Safety Lead
- 2772180** ESU-CQM Safety Lead
- 2772209** ESU-Jumper Safety Lead
- Power Cord (country specific)
- 2826194** Test Lead with stackable plugs
- 1903307** Test Lead Set with retractable sheaths
- 1610159** Sure-Grip Large Alligator Clip Set

Optional accessories

- 2461794** Carrying Case
- 2461802** Ansur Test Software, QA-ES Plug-in License
- 2461993** Data Transfer Cable, RS-232
- 2716059** Calibration Manual
- 2523266** Clamp, crocodile style, grip C, black
- 2523275** Clamp, crocodile style, grip C, red



About Fluke Biomedical

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance.

Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, while performing their work faster and more efficiently than ever. Fluke Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

Fluke Biomedical Regulatory Commitment

As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 certified and our products are:

- CE Certified, where required
- NIST Traceable and Calibrated
- UL, CSA, ETL Certified, where required
- NRC Compliant, where required

Fluke Biomedical.

Better products. More choices. One company.

Fluke Biomedical
6045 Cochran Road
Cleveland, OH 44139-3303 U.S.A.

Fluke Biomedical Europe
Science Park Eindhoven 5110,
5692EC Son, The Netherlands

For more information, contact us:
In the U.S.A. (800) 850-4608 or
Fax (440) 349-2307
In Europe/M-East/Africa +31 40 267 5200 or
Fax +31 40 267 5436
From other countries +1 (440) 248-9300 or
Fax +1 (440) 349-2307
Email: sales@flukebiomedical.com
Web access: www.flukebiomedical.com

©2007, 2008 Fluke Biomedical. Specifications
subject to change without notice. Printed in U.S.A.
1/2017 2839598 D-EN-N Rev D

SG032(3) 様式4

計測装置社外校正記録書

管理番号 Control Number	F-124	管理部門 Division	リペアサービス Repair Service
計測装置名称 Product	Electrosurgical Analyzer QA-ES II 300091		
依頼先供給者名 Supplier	大正医科器械 株式会社 Taisho Biomed Instruments Co., Ltd.		
校正実施日/期間 Date of Calibration	2018年9月28日 09/28/2018		
校正証明書の文書番号 Calibration Certificate Document Number	<ul style="list-style-type: none"> ・サービス報告書 報告書番号: 180928 Y2 ・Calibration Report Found -Left-300091-9/28/2018 ・2018年4月11日発行 校正体系図 (電気メステストレーサビリティ体系図) 		
校正結果の確認 Verification of Calibration Result	①校正担当者は校正結果の校正項目、公称値、測定値をレビューのうえ、許容範囲(合格範囲)内であることを確認し、校正結果にレ点チェックにより記録する。 (レ点チェックは、赤、青、緑などの色つきのペンを使用すること。)		
	②校正担当者は①の確認の結果、問題がなければ下記にチェックを記載する。 <input type="checkbox"/> 校正項目、公称値、測定値をレビューし、許容範囲(合格範囲)内であることを確認しました。		
合否判定 Result	PASS	合格	fail
次回校正実施期限*1 Due Date	2019年9月 09/2019		

*1 次回校正実施期限は、前回校正実施年月に校正周期を加えた年月の月末までとする。

【作成】校正担当者 署名:

Yuki Shimma 10/15/2018

【レビュー】計測装置管理責任者 署名:

20 10/16/2018

【承認】QO 署名:

O. Ota 10/17/2018

**Calibration Report**

Found-Left - 300091 - 9/28/2018

Model # QA-ES-J

Mfr. Fluke Biomedical

S/N: 300091

Customer Name: ジョンソン・エンド・ジョンソン株式会社

Customer P.O. #: N/A

Date Received: 2018/9/28

Condition as Received: Operable

WO #: N/A

Calibration Date: 2018/9/28

Calibration Data: Found-Left

Temperature (C): 24

Pressure (mm Hg): 761.4

Relative Humidity: 55

Notes: N/A

Calibration Standards

2000 or Equiv DMM IRM#: 1586027
 2000 or Equiv DMM Cal Due Date: 2019/4/16
 Oscilloscope IRM#: DM8480044
 Oscilloscope Cal Due Date: 2019/4/16
 RF Voltmeter #: 16856
 RF Voltmeter Cal Due Date: 2018/12/28
 Pearson 411 IRM#: 167612
 Pearson 411 Cal Due Date: 2019/3/15

Calibration Description

- * The calibration of this instrument is traceable to the National Institute of Standards and Technology.
- * The uncertainty ratio for this calibration is better than 4:1 except the Current Tests which are 2.3:1.
- * The results in this report relate only to the item being tested.
- * This calibration report cannot be reproduced except in full, without written approval from the laboratory.

Calibration Data

Preliminary Checks	Version	Pass/Fail	
Power Up Routine Works Properly		Pass	
Firmware Version	2.26		
User Interface and Display Operates Normally		Pass	
Cooling Fan Operates Normally		Pass	
Serial and Parallel Interface Operates Normally		Pass	
Foot Switch Operates Normally		Pass	
Housing and Measuring Inputs Isolated		Pass	
Load Tests	Spec (Ohms)	Reading (Ohms)	Pass/Fail
Fixed Load	180 - 220	200.18	Pass ✓
10 Ohm Load	9.8	9.97	+/-5% Pass ✓
25 Ohm Load	24.8	25.1	+/-5% Pass ✓
50 Ohm Load	49	49.82	+/-5% Pass ✓
75 Ohm Load	71.8	72.57	+/-5% Pass ✓
100 Ohm Load	102	102.61	+/-5% Pass ✓
200 Ohm Load	202	201.89	+/-5% Pass ✓
300 Ohm Load	302	302.11	+/-5% Pass ✓
400 Ohm Load	402	405.37	+/-5% Pass ✓
500 Ohm Load	502	505.6	+/-5% Pass ✓
700 Ohm Load	702	705.12	+/-5% Pass ✓
800 Ohm Load	802	798.68	+/-5% Pass ✓
1500 Ohm Load	1502	1502.88	+/-5% Pass ✓
1575 Ohm Load	1571.8	1572.7	+/-5% Pass ✓
1600 Ohm Load	1622	1619.6	+/-5% Pass ✓
3000 Ohm Load	3002	2992.3	+/-5% Pass ✓
5200 Ohm Load	5222	5210.3	+/-5% Pass ✓

大正医科器械株式会社
大阪市大正区三軒家東 2-2-22

TEL.
06-6553-9666

FAX.
06-6553-9305

Internet
<http://www.taishobiomed.com/>



Current Test	Measured	Tolerance	UUT	%Error	Pass/Fail
	mA		mA		
Monopolar Cut					
500 Ohm load @ 50W Cut	316	+/-5%	317	0.41	Pass ✓
300 Ohm load @ 50W Cut	404.9	+/-5%	407	0.52	Pass ✓
300 Ohm load @ 100W Cut	572.6	+/-5%	576	0.59	Pass ✓
1500 Ohm load @ 100W Cut	202	+/-5%	198	-1.98	Pass ✓
300 Ohm load @ 150W Cut	702	+/-5%	705	0.43	Pass ✓
300 Ohm load @ 200W Cut	810.8	+/-5%	812	0.15	Pass ✓
300 Ohm load @ 300W Cut	911.7	+/-5%	913	0.14	Pass ✓
700 Ohm load @ 300 W Cut	650.5	+/-5%	649	-0.23	Pass ✓
200 Ohm load @ 300 W Cut	949	+/-5%	948	-0.05	Pass ✓
100 Ohm load @ 300 W Cut	987	+/-5%	986	-0.09	Pass ✓

Power Tests	Calculated Power	UUT Power	%Error
	Watts	Watts	
500 Ohm load @ 50W	50.4	51	1.21
300 Ohm load @ 50W	49.5	50	0.95
300 Ohm load @ 100W	99.1	100	0.96
1500 Ohm load @ 100W	61.3	59	-3.79
300 Ohm load @ 150W	148.9	150	0.75
300 Ohm load @ 200W	198.6	199	0.20
300 Ohm load @ 300W	251.1	252	0.35
700 Ohm load @ 300 W	298.4	298	-0.12
200 Ohm load @ 300 W	181.6	181	-0.35
100 Ohm load @ 300 W	99.9	100	0.06

Output Waveform Test Pass ✓

Current Test	Measured	Tolerance	UUT	Pass/Fail	
	mA	Percent	Offset (mA)	mA	
Monopolar Coag					
2000 Ohm @ 1W Coag	22.0	5%	2	20	Pass ✓
1000 Ohm @ 1W Coag	30.0	5%	2	29	Pass ✓
600 Ohm @ 1W Coag	37.2	5%	2	36	Pass ✓
400 Ohm @ 1W Coag	45.1	5%	2	45	Pass ✓
300 Ohm @ 3W Coag	96.0	5%	2	97	Pass ✓
300 Ohm @ 7W Coag	148.5	5%	2	150	Pass ✓
300 Ohm @ 12W Coag	194.6	5%	2	196	Pass ✓
300 Ohm @ 27W Coag	292.0	5%	2	294	Pass ✓
300 Ohm @ 75W Coag	497.2	5%	2	490	Pass ✓
200 Ohm @ 110W Coag	616.7	5%	2	604	Pass ✓
75 Ohm @ 110W Coag	638.8	5%	2	636	Pass ✓
50 Ohm @ 110W Coag	642.6	5%	2	638	Pass ✓

Span Tests	Measured	Tolerance	UUT	Pass/Fail	
	mA	Percent	Offset (mA)	mA	
Bipolar Coag					
600 Ohm @ 1W Coag	35.6	5%	2	34	Pass ✓
400 Ohm @ 1W Coag	43.7	5%	2	43	Pass ✓
300 Ohm @ 3W Coag	97.4	5%	2	98	Pass ✓
300 Ohm @ 7W Coag	149.9	5%	2	150	Pass ✓
300 Ohm @ 12W Coag	197.4	5%	2	198	Pass ✓
300 Ohm @ 27W Coag	296.9	5%	2	299	Pass ✓
275 Ohm @ 70W Coag	503.5	5%	2	507	Pass ✓
125 Ohm @ 70W Coag	749.1	5%	2	753	Pass ✓
50 Ohm @ 70W Coag	1153.0	5%	2	1150	Pass ✓
25 Ohm @ 70W Coag	1539.0	5%	2	1523	Pass ✓

Calibration Performed by : 天正辰治
 Technical Check by : 橋本敏仁

Date: 2018/9/28

Date: 2018/10/1

The suggested re-calibration date is:

2019/9/28



99 Washington Street
Melrose, MA 02176
Fax 781-665-0780
TestEquipmentDepot.com

General Specifications of

73301 / 02 / 03

Additional Functions	Relative and percentage value computation, data / auto hold, overvoltage warning
Display	Digital display: 4,000-count digital reading; 40-segment bar graph
Measuring Rate	Digital display: 2.3 times/sec Bar graph display: 23 times/sec
Operating Temp. and Humidity	-10°C to 50°C: 80% RH or less at -10°C to 40°C, or 70% RH or less at 40°C to 50°C (no condensation)
Storage Temp. and Humidity	-20°C to 60°C, 70% RH or less (no condensation)
Temperature Coefficient	Add the accuracy (0.1°C) to the basic accuracy at a temperature within -10°C to 18°C and 28°C to 50°C
Withstanding Voltage	5.55 kV AC for 1 minute (between input terminals and casing)
Power Supply	Two AA (R6) dry cells
Battery Life	73301: Approx. 1000 hours (for continuous DC voltage) 73302/03: Approx. 350 hours measurement with alkaline cells)
Auto Power Off	The power is automatically turned off when no operation is made for approx. 20 minutes (can be disabled).
Dimensions	85 (W) x 191 (H) x 40 (D) mm
Weight	Approximately 450 g (including batteries)
Compliance with Standards	Safety EN61010-1 (1995); EN61010-2-031 (1995) (AC/DC 1000 V, CAT II; AC/DC 600 V, CAT III) EMC: EN55011 (1998); EN61326-1 (1998) + A1 (Class B, Group 1) EMS: EN50082-1 (1997); EN61326-1 (1998) + A1
Standard Accessories	Instruction manual:1, Test lead set (RD031):1, AA (R6) dry cells(built in):2

Optional Accessories

Name	Model	Specification
Fuse	A1518EF	500 mA/600 V
	A1519EF	15 A/600 V
Test leads	RD031	Red / black (1 set)
Thermistor probe	234901	-50°C to 150°C
Carrying case (hard)	B9646HH	Houses the DMM and test leads



Performance

Test conditions: Temperature and humidity = 23°C ± 5°C, 80% RH or less; Accuracy = ±(% rdg + dg).t.

Note: A response time is the time required for achieving the accuracy specified for the corresponding range.

• DC Voltage Measurement (⎓ V)

Range	Accuracy		Input Resistance	Maximum Input Voltage
	73301	73302/03		
400.0 mV fixed			10 MΩ	
4.000 V			11 MΩ	1000 Vrms AC, 1000 V DC
40.00 V	0.3% + 1	0.2% + 1	10 MΩ	
400.0 V				
1000 V				

Response time: 1 second or less

• AC Voltage Measurement (⎓ V)

Range	Accuracy			Input Impedance	Maximum Input Voltage
	50/60 Hz	50 – 500 Hz	500 Hz – 1 kHz		
400.0 mV fixed				10 MΩ, <50 pF	
4.000 V				11 MΩ, <50 pF	
40.00 V	0.5% + 2	1% + 2	1.5% + 4		1000 Vrms AC, 1000 V DC
400.0 V				10 MΩ, <50 pF	
1000 V					

Response time: 2 seconds or less

Range	Accuracy			Input Impedance	Maximum Input Voltage
	50/60 Hz	50 – 500 Hz	500 Hz – 1 kHz		
400.0 mV fixed				10 MΩ, <50 pF	
4.000 V				11 MΩ, <50 pF	
40.00 V	0.5% + 2	0.75% + 2	1.5% + 4		1000 Vrms AC, 1000 V DC
400.0 V					
1000 V					

Response time: 2 seconds or less

Range	Accuracy			Input Impedance	Maximum Input Voltage
	50/60 Hz	50 – 500 Hz	500 Hz – 1 kHz		
400.0 mV fixed				10 MΩ, <50 pF	
4.000 V				11 MΩ, <50 pF	
40.00 V	0.5% + 2	0.75% + 2	1.5% + 4		1000 Vrms AC, 1000 V DC
400.0 V					
1000 V					

Response time: 2 seconds or less

*1: 5 to 100% of F.S., or 200 to 1000 V for 1000 V range:

Range	Accuracy			Input Impedance	Maximum Input Voltage
	50/60 Hz	50 – 500 Hz	500 Hz – 1 kHz		
400.0 mV fixed				10 MΩ, <50 pF	
4.000 V				11 MΩ, <50 pF	
40.00 V	0.5% + 5	1% + 5	1.5% + 5		1000 Vrms AC, 1000 V DC
400.0 V					
1000 V					

Response time: 2 seconds or less; crest factor: <3

Range	Accuracy			Voltage Drop	Maximum Input Current
	73301	73302/03			
400.0 μA				<0.11 mV/μA	
4000 μA				400 mA (500 mA/600 V fuse-protected)	
40.00 mA	1% + 2	0.5% + 2		<2.5 mV/mA	
400.0 mA					
4.000 A	1.2% + 2			<0.1 V/A	10 A (15 A/600 V fuse-protected)
10.00 A					

Response time: 1 second or less

Range	Accuracy			Voltage Drop	Maximum Input Current
	50/60 Hz	40 Hz – 1 kHz			
400.0 μA				<0.11 mV/μA	
4000 μA				400 mA (500 mA/600 V fuse-protected)	
40.00 mA	1% + 5	1.5% + 5		<2.5 mV/mA	
400.0 mA					
4.000 A	1.2% + 5			<0.1 V/A	10 A (15 A/600 V fuse-protected)
10.00 A					

Response time: 2 seconds or less

• Capacitor Check (-||-)

Models 73302/03 (function not available with 73301)

Range	Accuracy		Input Protection Voltage
	73301	73302/03	
10.00 nF		2% + 10 (after zero calibration)	
100.0 nF			
1000 nF		2% + 5	
10.00 μF			
100.0 μF			
1000 μF		3% + 5	

Accuracy when used in combination with optional thermistor probe (234901)

• Frequency Measurement (Hz)

Models 73302/03 (function not available with 73301)

Range	Accuracy		Input Voltage Range	Maximum Input Voltage
	73302/03			
10.00 – 99.99 Hz			0.2 – 400 VRMS	
90.0 – 999.9 Hz		0.02% + 1		
900 – 9999 Hz			0.4 – 400 VRMS	
9.00 – 99.99 kHz			0.8 – 100 VRMS	

Coupling type: AC coupling

Test Equipment Depot - 800.517.8431 - 99 Washington Street Melrose, MA 02176

FAX 781.665.0780 - TestEquipmentDepot.com

SG032(3) 様式4

計測装置社外校正記録書

管理番号 Control Number	F-07	管理部門 Division	リペアサービス Repair Service
計測装置名称 Product	YOKOGAWA デジタルマルチメーター Digital Multi Meter		
依頼先供給者名 Supplier	日本電計株式会社 Denki		
校正実施日/期間 Date of Calibration	2018年3月14日 03/14/2018		
校正証明書の文書番号 Document Number	校正証明書 発行番号: UI-4428409 試験成績書 YEW-5251-01-S206-01 / 100209B02 トレーサビリティ体系図 UTM-006-1-068-18		
※受取った校正証明書の全ての文書番号を転記			
校正結果の確認 Verification Result	①校正担当者は校正結果の校正項目、公称値、測定値をレビューのうえ、許容範囲(合格範囲)内であることを確認し、校正結果にレ点チェックにより記録する。 (レ点チェックは、赤、青、緑などの色つきのペンを使用すること。)		
	②校正担当者は①の確認の結果、問題がなければ下記にチェックを記載する。 <input checked="" type="checkbox"/> 校正項目、公称値、測定値をレビューし、許容範囲(合格範囲)内であることを確認しました。		
合否判定 Result	pass <input checked="" type="radio"/> 合格 / 不合格 fail		
次回校正実施期限*1 Due Date	2020年3月 03/2020		

*1 次回校正実施期限は、前回校正実施年月に校正周期を加えた年月の月末までとする。

【作成】校正担当者 署名:

Kazuhiko Ishikawa 3/26/2018

【レビュー】計測装置管理責任者 署名:

Yuki Ishikawa 3/26/2018

【承認】QO 署名:

E-Ola 03/26/2018

Date of Issue 03/14/2018
 発行日: 2018年3月14日
 Control Number
 発行番号: UI-4428409

校正証明書

Calibration
 Certification

顧客名: ジョンソン・エンド・ジョンソン株式会社 御中
 Customer Johnson and Johnson Co., Ltd.

事業所名:
 Plant

部門名:
 Division

製品名: デジタルマルチメータ
 Product Digital Multi Meter

型番: 733 02
 Product Code

製造者名: 横河電機
 Manufactured Yokogawa

製造番号: 1501471
 Manufacturing Number

管理番号: F-07
 Control Number

校正日: 2018年3月14日
 Date of Calibration 03/14/2018

上記の製品は、当社の管理規定にもとづき校正されています。

使用した基準器は国際度量衡委員会(CIPM)/国際度量衡局(BIPM)に加盟する、国立研究開発法人産業技術総合研究所(AIST)等の国家、国際標準にトレーサブルである事を証明します。

ISO / IEC 17025 (RCL00030) 認定校正機関
 (直流電圧・電流 / 交流電圧・電流 / 直流抵抗)
 ISO9001 (JQA-QMA15393), ISO14001 (JQA-EM7241) 認証取得



ユウアイ電子株式会社

埼玉県川越市かし野台2-16-1
 TEL:049-243-4703 FAX:049-243-4703
 E-mail:calroom@yuai.co.jp



試験成績書

(1/2)

顧客名 ジョンソン・エンド・ジョンソン株式会社 御中
 Customer Johnson and Johnson k.k.

校正日 2018年3月14日
 Date of Calibration 03/14/2018
 溫度 21 °C 湿度 41 %
 Temperature humidity

製品名 デジタルマルチメータ
 Product Digital Multi Meter
 製造者名 横河電機
 Manufacturer Yokogawa
 管理番号 F-07

型番 733 02
 Product Code
 製造番号 1501471
 Manufacturing Number

校正者 矢内和宏
 Tester

承認者 金子
 Approved

判定
 Result
 埼玉県川越市かし野台2-13-2
 ユウアイ電子株式会社

上記の製品は、当社の管理規定にもとづき校正されています。使用した基準器は国際度量衡委員会(CIPM)/国際度量衡局(BIPM)に加盟する、国立研究開発法人産業技術総合研究所(AIST)等の国家、国際標準にトレーサブルです。

試験項目 Test Contents

直流電圧測定 DC Current

レンジ Range	分解能 Resolution	試験点 Test Point	確度 Spec	下限値 Lower	測定値 measured	上限値 Upper	結果 result
400 mV	0.1	100 mV	0.2 + 1	99.7 mV	100.0 mV	100.3 mV	合格 ✓ pass
4 V	0.001	1 V	0.2 + 1	0.997 V	1.000 V	1.003 V	合格 ✓ pass
40 V	0.01	10 V	0.2 + 1	9.97 V	10.00 V	10.03 V	合格 ✓ pass
400 V	0.1	100 V	0.2 + 1	99.7 V	100.0 V	100.3 V	合格 ✓ pass
1000 V	1	1000 V	0.2 + 1	997 V	999 V	1003 V	合格 ✓ pass

確度 $\pm (\% \text{ rdg} + \text{dgt})$
 Spec

交流電圧測定 AC Current f=400Hz

レンジ Range	分解能 Resolution	試験点 Test Point	確度 Spec	下限値 Lower	測定値 measured	上限値 Upper	結果 result
400 mV	0.1	100 mV	0.75 + 2	99.1 mV	100.0 mV	100.9 mV	合格 ✓ pass
4 V	0.001	1 V	0.75 + 2	0.991 V	1.000 V	1.009 V	合格 ✓ pass
40 V	0.01	10 V	0.75 + 2	9.91 V	10.00 V	10.09 V	合格 ✓ pass
400 V	0.1	100 V	0.75 + 2	99.1 V	100.0 V	100.9 V	合格 ✓ pass
1000 V	1	1000 V	0.75 + 2	991 V	1000 V	1009 V	合格 ✓ pass

確度 $\pm (\% \text{ rdg} + \text{dgt})$

直流電流測定 DC Current

レンジ Range	分解能 Resolution	試験点 Test Point	確度 Spec	下限値 Lower	測定値 measured	上限値 Upper	結果 result
400 μA	0.1	100 μA	0.5 + 2	99.3 μA	99.6 μA	100.7 μA	合格 ✓ pass
4000 μA	1	1000 μA	0.5 + 2	993 μA	999 μA	1007 μA	合格 ✓ pass
40 mA	0.01	10 mA	0.5 + 2	9.93 mA	9.96 mA	10.07 mA	合格 ✓ pass
400 mA	0.1	100 mA	0.5 + 2	99.3 mA	99.9 mA	100.7 mA	合格 ✓ pass
4 A	0.001	1 A	0.5 + 2	0.993 A	0.994 A	1.007 A	合格 ✓ pass
10 A	0.01	10 A	0.5 + 2	9.93 A	9.98 A	10.07 A	合格 ✓ pass

確度 $\pm (\% \text{ rdg} + \text{dgt})$

試験成績書

(2/2)

製品名 デジタルマルチメータ
 Product Digital Multi Meter
 製造者名 横河電機
 Product Yokogawa
 管理番号 F-07
 Control Number

型番 733 02
 Product Code
 製造番号 1501471
 Manufacturing Number

交流電流測定 DC Current $f=400Hz$

レンジ Range	分解能 Resolution	試験点 Test Point	確度 Spec.	下限値 Lower	測定値 measured	上限値 Upper	結果 Result
400 μA	0.1	100 μA	1.5 + 5	98.0 μA	99.9 μA	102.0 μA	合格 V pass
4000 μA	1	1000 μA	1.5 + 5	980 μA	999 μA	1020 μA	合格 V pass
40 mA	0.01	10 mA	1.5 + 5	9.80 mA	9.98 mA	10.20 mA	合格 V pass
400 mA	0.1	100 mA	1.5 + 5	98.0 mA	99.7 mA	102.0 mA	合格 V pass
4 A	0.001	1 A	1.5 + 5	0.980 A	0.996 A	1.020 A	合格 V pass
10 A	0.01	10 A	1.5 + 5	9.80 A	9.97 A	10.20 A	合格 V pass

確度 $\pm (\% \text{rdg} + \text{dgt})$
Spec

抵抗測定 Resistance

レンジ Range	分解能 Resolution	試験点 Test Point	確度 Spec.	下限値 Lower	測定値 measured	上限値 Upper	結果 Result
400 Ω	0.1	100 Ω	0.4 + 1	99.5 Ω	100.0 Ω	100.5 Ω	合格 V pass
4 k Ω	0.001	1 k Ω	0.4 + 1	0.995 k Ω	1.000 k Ω	1.005 k Ω	合格 V pass
40 k Ω	0.01	10 k Ω	0.4 + 1	9.95 k Ω	9.99 k Ω	10.05 k Ω	合格 V pass
400 k Ω	0.1	100 k Ω	0.4 + 1	99.5 k Ω	99.8 k Ω	100.5 k Ω	合格 V pass
4 M Ω	0.001	1 M Ω	0.5 + 1	0.994 M Ω	1.000 M Ω	1.006 M Ω	合格 V pass
40 M Ω	0.01	10 M Ω	1.0 + 2	9.88 M Ω	9.99 M Ω	10.12 M Ω	合格 V pass

確度 $\pm (\% \text{rdg} + \text{dgt})$

data end

特記 : 顧客仕様

使用基準器 Reference Equipment

No	基準器名 Name of Equipment	製造番号, Manufacturing Number	管理番号 Control Number	基準器の有効期限 Due Date
1	マルチプロダクトキャリブレータ Multi Product Calibrator 5522A/SC600	2625904	U-1-068	19年02月

全ての測定値が許容範囲内であることを確認いたしました。

担当者	部門長
Zohir Shumra 3/26/2018	Yoko Shumra 3/26/2018

計測装置管理
責任者

TOS6100

AC LOW OHM TESTER

Specifications

■ Resistance Measured Values	0 to 0.12Ω/0 to 0.6Ω, 2 ranges	■ Signal Outputs
Accuracy	±10% of full scale (3 to 5 A) ±5% of full scale (5 to 30 A)	
■ Test Current Output	3 to 30 A AC, adjusted by dial on panel Max. current 30 A AC, max. voltage: 8 V (when line voltage is at the central value of each input voltage range) → These values cannot be obtained at the same time.	TEST Testing in progress PASS Approx. 50 ms when PASS judgement is made FAIL Continuous when FAIL judgement is made WARNING See test current monitoring function READY During standby MONITOR Normal
■ Output Ammeter Scale	3 to 30 A AC ±5% of full scale	Make contact signal, amp Make contact signal, amp, buzzer Make contact signal, amp, buzzer Make contact signal, amp Make contact signal DC 0 to 10 V
Accuracy	Mean value response/rms value scale graduation	
■ Pass/fail Judgement Judgement Method	<ul style="list-style-type: none"> ● A FAIL judgement is made and output is cut off if the measured value is greater than the reference value. ● A PASS judgement is made and a PASS signal is output if there are no abnormalities after the set time has elapsed. 	<ul style="list-style-type: none"> →(1)The contact rating of the contact signal is 100 V AC, 1 A/30V DC, 1 → →(2)Buzzer volume can be adjusted with a single dial for the PASS signal and FAIL alarm. →(3)The MONITOR output consists of an output of direct current voltage proportional to the angle of deflection of the pointer on the ohmmeter. The scale is as indicated below.
Limit Value Setting Range	Set as desired between 5 to 100% of the measuring range	● 0.5 Ω range: 10V/0.5Ω ● 0.1 Ω range: 10V/0.1Ω The absolute value of output voltage error is either 5% of the output value or 50 mV, whichever is larger.
Judgement Accuracy	±15% of full scale for reference value (3 to 5 A) ±10% of full scale for reference value (5 to 30 A)	
■ Subtraction Function	<ul style="list-style-type: none"> ● A predetermined value can be subtracted from the measured value and the result of subtraction can be displayed ● The result of subtraction can be compared with the judgement reference value and the result of comparison can be used for PASS/FAIL judgement 	■ Environment Specification temperature and humidity ranges Operating temperature and humidity ranges Storage temperature ■ Power Requirements Allowable Line Voltage Frequency: 50/60 Hz Power Consumption Max. 20 VA under no-load conditions (READY state) Approx. 280 VA during 30 A output (RM = 0.22Ω) Insulation Resistance Withstanding Voltage ■ Dimensions (MAX) 430W × 149H × 370Dmm (430W × 164H × 435Dmm) ■ Weight ■ Accessories Approx. 16 kg ● Short bars (to be attached to the tester): 2 ● Power cable set: 1 ● Operation manual: 1 ● Fuse 3 A (S.B): 1 ● Fuse 1.6 A (S.B): 2 ● Model RC01-TOS/RC02-TOS remote control box ● Model LTP-2 low resistance test probe ● Model BZ01-TOS buzzer unit ● Model PL01-TOS warning light unit ● Model BH3M-TOS rack mount bracket for JIS ● Model BH4-TOS rack mount bracket for EIA
Subtraction Range	0 to 0.1Ω	
Subtraction Error	Max. ±5 of full scale (added to measuring accuracy or PASS/FAIL judgement accuracy)	
■ Test Current Monitoring Function	<ul style="list-style-type: none"> ● Test current can be monitored during testing ● A warning alarm is generated when the test current value deviates from the range of approximately 10% of the monitoring reference value ● Reference values can be set as desired within a range of 3 to 30 A ● Continuation or stopping of testing when a warning has occurred can be selected 	5 to 35°C/20 to 85% RH 0 to 40°C/20 to 90% RH -20 to 70°C/Max. 90% RH A: 90 to 110V, B: 104 to 125V, C: 194 to 236V, D: 207 to 250V Frequency: 50/60 Hz Max. 20 VA under no-load conditions (READY state) Approx. 280 VA during 30 A output (RM = 0.22Ω) Min. 30 MΩ at 500V DC 1000V AC for 1 minute 430W × 149H × 370Dmm (430W × 164H × 435Dmm) Approx. 16 kg ● Short bars (to be attached to the tester): 2 ● Power cable set: 1 ● Operation manual: 1 ● Fuse 3 A (S.B): 1 ● Fuse 1.6 A (S.B): 2 ● Model RC01-TOS/RC02-TOS remote control box ● Model LTP-2 low resistance test probe ● Model BZ01-TOS buzzer unit ● Model PL01-TOS warning light unit ● Model BH3M-TOS rack mount bracket for JIS ● Model BH4-TOS rack mount bracket for EIA
Test Time	0.5 seconds to 10 minutes (with 4-range timer)	
■ Remote Control Start/Stop Operation	<ul style="list-style-type: none"> ● Low active control ● Input conditions: <ul style="list-style-type: none"> ♦ High level input voltage: 11 to 15 V ♦ Low level input voltage: 0 to 4 V ♦ Low level sweepout current: Max. 2 mA ♦ Input signal time width: Min. 20 ms ● Since the input terminal is pulled up to a +15V power supply by a resistor, input becomes equal to high level input if the input terminal is open. 	

SG032(2) 様式4

計測装置社外校正記録書

管理番号 <i>Control Number</i>	F-09	管理部門 <i>Division</i>	リペアサービス <i>Repair Service</i>
計測装置名称 <i>Product</i>	交流低抵抗試験器 TOS6100 S/N:29070206 <i>Ground Bond Tester</i>		
依頼先供給者名 <i>Supplier</i>	日本電計株式会社 <i>Denki</i>		
校正実施日/期間 <i>Date of Calibration</i>	08/03/2017		
校正証明書の文書番号 <i>Certificate of Calibration Document Number</i>	校正証明書 発行番号: UI-4410258 試験成績書 KIK-0008-01 040422T07 トレーサビリティ体系図 UTM-006-1-0163-03 トレーサビリティ体系図 UTM-006-1-197-19		
校正結果の確認 <i>Calibration Result Confirmation</i>	①校正担当者は校正結果の校正項目、公称値、測定値をレビューのうえ、許容範囲(合格範囲)内であることを確認し、校正結果にレ点チェックにより記録する。 (レ点チェックは、赤、青、緑などの色つきのペンを使用すること。) ②校正担当者は①の確認の結果、問題がなければ下記にチェックを記載する。 <input checked="" type="checkbox"/> 校正項目、公称値、測定値をレビューし、許容範囲(合格範囲)内であることを確認しました。		
合否判定 <i>Result</i>	PASS	合格	FAIL
次回校正実施期限*1 <i>Due Date</i>	8/2019		

*1 次回校正実施期限は、前回校正実施年月に校正周期を加えた年月の月末までとする。

【作成】校正担当者 署名: *鈴木和昭 8/4/2017*

【レビュー】計測装置管理責任者 署名: *Yuki Shimura 8/4/2017*

【承認】QRC 署名: *Elini 08/07/2017*

Date of Issue 08/03/2017
 発行日: 2017年8月3日
 Control Number
 発行番号: UI-4410258

校正証明書

Calibration
 Certification

顧客名: ジョンソン・エンド・ジョンソン株式会社 御中
 Customer: Johnson and Johnson K.K.

事業所名:
 Plant

部門名:
 Division

製品名: 交流低抵抗試験器
 Product: Ground Bond Tester

型番: TOS6100
 Product Code: TOS6100

製造者名: 菊水電子工業
 Manufactured: Kikusui Electronic Corp.

製造番号: 29070206
 Manufacturing Number: 29070206

管理番号: F-09
 Control Number: F-09

校正日: 2017年8月3日
 Date of Calibration: 08/03/2017

上記の製品は、当社の管理規定にもとづき校正されています。

使用した基準器は国際度量衡委員会(CIPM)/国際度量衡局(BIPM)に加盟する、国立研究開発法人産業技術総合研究所(AIST)等の国家、国際標準にトレーサブルである事を証明します。

ISO / IEC 17025 (RCL00030) 認定校正機関
 (直流電圧・電流/交流電圧・電流/直流抵抗)
 ISO9001 (JQA-QMA15393)、ISO14001 (JQA-EM7241)認証取得



ユウアイ電子株式会社

埼玉県川越市かし野台2-16
 TEL:049-243-4703 FAX:049-243-4704
 E-mail:calroom@yuai.co.jp



試験成績書 Test Report

(1/1)

顧客名 ジョンソン・エンド・ジョンソン株式会社 御中
Customer: Johnson and Johnson KK.

校正日 2017年8月3日
Date of Calibration 08/03/2017
温度 21 °C 湿度 42 %
temperature 21 °C humidity 42 %

製品名 交流低抵抗試験器
Product: Ground Bond Tester
製造者名 菊水電子工業
Manufacturer: Kikusui Electronic Corp
管理番号 F-09
Control Number

型番 TOS6100
Product Code: TOS6100
製造番号 29070206
Manufacturing Number

校正者 矢内和宏
Tester: Kazuhiko Yachi

承認者 鈴木
Approved: Suzuki

判定 合格 PASS
Result: Pass
埼玉県川越市かし野台2-13-2
ユウアイ電子株式会社

上記の製品は、当社の管理規定にもとづき校正されています。使用した基準器は国際度量衡委員会(CIPM)/国際度量衡局(BIPM)に加盟する、国立研究開発法人産業技術総合研究所(AIST)等の国家、国際標準にトレーサブルです。

試験項目 Test Contexts

抵抗測定 Resistance 10A

レンジ Range	試験点 Test point	精度 Spec	下限値 Lower	測定値 measured	上限値 Upper	結果 Result
0.10 Ω	0.02 Ω	5	0.015 Ω	0.019 Ω	0.025 Ω	合格 ✓ Vpa
0.10 Ω	0.06 Ω	5	0.055 Ω	0.060 Ω	0.065 Ω	合格 ✓ Vpa
0.10 Ω	0.10 Ω	5	0.095 Ω	0.100 Ω	0.105 Ω	合格 ✓ Vpa
0.50 Ω	0.10 Ω	5	0.075 Ω	0.100 Ω	0.125 Ω	合格 ✓ Vpa
0.50 Ω	0.30 Ω	5	0.275 Ω	0.296 Ω	0.325 Ω	合格 ✓ Vpa
0.50 Ω	0.50 Ω	5	0.475 Ω	0.500 Ω	0.525 Ω	合格 ✓ Vpa

精度 ±(% fs)
Spec

電流計 Current

試験点	精度	下限値	測定値	上限値	結果
10 A	5	8.50 A	9.95 A	11.50 A	合格 ✓ Vpa
20 A	5	18.50 A	19.88 A	21.50 A	合格 ✓ Vpa
30 A	5	28.50 A	29.76 A	31.50 A	合格 ✓ Vpa
fs= 30 A					

精度 ±(% fs)
Spec

data end

使用基準器 Reference Equipment

No	基準器名 Name of Equipment	製造番号 Manufacturing Number	管理番号 Control Number	基準器の有効期限 Due Date
1	マルチプロダクトキャリブレータ Multi Product Calibrator 5502A	2927802	U-1-163	17年10月 ✓ V
2	デジタルマルチメータ DMM 34410A	MY45001069	U-1-197	18年04月 ✓ V

全ての値が良好であることを確認しました。	
担当者	部門長
8/4/2017	8/4/2017
計測装置管理責任者	計測装置管理責任者
8/4/2017	8/4/2017

8/4/2017 ¥3

KIK-0008-01

ユウアイ電子株式会社

040422T07

3156 Specifications

■ Leak current measurement

Measurement mode	: Earth leakage current / Leakage current between enclosure and earth / Leakage current between enclosure and enclosure / Leakage current between enclosure and line / Patient leakage current I / Patient leakage current II / Patient leakage current III / Patient auxiliary current	Measurement system	: Indication of a current value calculated based on the measured drop in voltage caused by simulated resistance of the human body. Measurement of true effective value. The measurement section: chassis-grounded and floating.
Target current	: DC, AC, AC+DC, AC peak	A/D conversion system	: $\Delta \Sigma$ system (20 bits)
Allowable value measurement current	: 25 mA max. (DC/AC/AC+DC mode) 75 mA max. (ACpeak mode)	Input resistance	: $1 \text{ M}\Omega \pm 1\%$ (single-ended input) Excluding voltmeter section, simulated resistance of the human body (current detection circuit)
Measurement range	: DC/AC/AC+DC mode: 50 μA /500 μA /5 mA/25 mA AC peak mode: 500 μA /1 mA/10 mA/75 mA	Input capacity	: 200 pF or lower (between terminals T1 and T2) ($f = 100 \text{ kHz}$, with network circuit isolated)
Range switch	: Auto range/Hold range	Grounding capacity	: 200 pF or lower (between terminals T1/T2 and chassis)
Trigger system	: Manual: Automatic generation of internal trigger, and free-run measurement Automatic: measurement started by external start signal 1. Started by pressing Start key on the operation screen. 2. Asserted by the START terminal on the EXT I/O connector. 3. Started by the :START interface command.	CMRR (between terminals T1/T2 and chassis)	: 60 dB or higher, at 60 Hz, 10 kHz 40 dB or higher, at 100 kHz, 1 MHz
Measurement speed	: Trigger-system manual measurement: $100 \pm 1 \text{ ms}$ (indication of moving average of 16 measurements) Trigger-system automatic measurement: $100 \pm 1 \text{ ms}$ (min.) - measurement setting time		
Measurement terminals	: Terminal T1, terminal T2 (with built-in fuse holder) Terminal T3 (110% voltage application terminal)		

■ Functions

110% voltage application function	: Equipped with a voltage output terminal (T3) that applies 110% power supply voltage between the functionally insulated signal input/output section (or Type F applied part) and ground. Output ON/OFF selection Output impedance: $22.5 \pm 1 \text{ k}\Omega$ •Applied when positive phase (to input power supply voltage) •Applied when negative phase (to input power supply voltage) Automatic switching function (automatic measurement function)	Allowable value judgement function	: Allowable value: Sets the upper-limit current value Judgement: PASS measurement value \leq upper-limit value FAIL measurement value $>$ upper-limit value Processing: Indication, buzzer, judgement output from EXT I/O
Wiring check function	: Polarity check/VA check	Mode selection function	: •Current measurement function Unit of current measurement: Auto/ mA, fixed •Voltage measurement function Isolates the internal network for using the product as a voltmeter between terminals T1 and T2. Maximum measurement voltage: 25 V
Automatic measurement function	: Setting of measuring time. Setting of delay (wait) time for changing setting conditions. The power supply polarity and equipment status are automatically switched during measurement.	Beep sound setting	: •Allowable value judgement •Key input •T3 (110% voltage application terminal) output •Line voltage output from T2
Application line selection function	: Use of T2 and internal contact/Use of T1 and T2	Save/load function	: 30 panels for saving the following setting data (measurement mode, network, equipment name, control number, grounding class, applied part, measurement range, filter, target current, allowable value setting, malfunction condition setting, power supply polarity switching, automatic measurement items, automatic measuring time, measurement delay time)
Ground fault prevention function	: Pre-check of current value between connection terminals to prevent a ground fault. Only effective for leakage current measurement between enclosure and line.	Data save function	: Saved content: Sample equipment information (equipment name, control number), measurement data, date Memory capacity: Data up to 100 units
Setting of single-fault condition	: •Setting of malfunction mode for power line for sample equipment. 1. Disconnection of one wire in power line (neutral side) 2. Disconnection of protective earth conductor •Application of 110% voltage for simulated connection of malfunctioning equipment. Positive phase/negative phase •Selection of application line for leakage current measurement between enclosure and line.	Clock function	: Auto calendar, automatic leap-year adjustment, 24-hour clock Clock accuracy: Deviation of about 4 minutes a month
Power line for sample equipment (switching power supply polarity)	: Positive/negative phase (Automatic switching possible when using automatic measurement function)	Data backup function	: SRAM (test condition data), RTC Backup battery life: 4 years (reference value at 25°C)
Setting of measuring time	: Setting range: 1 sec. to 5 min., in 1-sec. increments Effective only in automatic measurement	Backlight automatic OFF function	: Constant ON/Auto OFF 1 min. to 30 min., in 1-min. increments
Measurement delay (setting) function	: Setting range: 1 sec. to 30 min., in 1-sec. increments •Wait time from the completion of measurement to power supply disconnection •Wait time from switching power supply polarity to the start of measurement •Wait time for operations other than switching polarity	Self-test function	: MEM (internal RAM)/KEY (6x6 matrix touch panel)/LCD (front LCD panel)/LED /Buzzer
Maximum value hold function	: Effective in all measurement modes	Language setting	: Japanese or English
		System reset	: Clears all data including measurement conditions and measurement data. Clears all saved measurement data. Clears all saved condition setting data including panels.

SG032(3) 様式4

計測装置社外校正記録書

管理番号 Control Number	F-56	管理部門 Division	リペアサービス Repair Service
計測装置名称 Product	LEAK CURRENT HITESTER S/N 051238345		
依頼先供給者名 Supplier	日本電計株式会社 Denki		
校正実施日/期間 Date of Calibration	2018年11月5日 11/05/2018		
校正証明書の文書番号 No.	校正証明書：発行番号 UI-4518429 試験成績書：HIK-0021-03 / 170217R01 トレーサビリティ体系図：UTM-006-1-017-20 / UTM-006-1-119-18		
※受取った校正証明書の全ての文書番号を転記			
校正結果の確認 Verification Result	①校正担当者は校正結果の校正項目、公称値、測定値をレビューのうえ、許容範囲(合格範囲)内であることを確認し、校正結果にレ点チェックにより記録する。 (レ点チェックは、赤、青、緑などの色つきのペンを使用すること。) ②校正担当者は①の確認の結果、問題がなければ下記にチェックを記載する。 <input checked="" type="checkbox"/> 校正項目、公称値、測定値をレビューし、許容範囲(合格範囲)内であることを確認しました。		
合否判定 Result	pass	合格	fail
次回校正実施期限*1 Due Date	2019年11月 11/2019		

*1 次回校正実施期限は、前回校正実施年月に校正周期を加えた年月の月末までとする。

【作成】校正担当者 署名: Yuki Shimura 11/8/2018

【レビュー】計測装置管理責任者 署名: Yuki Shimura 11/8/2018

【承認】QO 署名: E.Ota 11/08/2018

Date of Issue 11/15/2018
 発行日: 2018年11月5日
 Control Number UI-4518429
 発行番号: UI-4518429

校正証明書

Calibration
Certification

顧客名: ジョンソン・エンド・ジョンソン株式会社 御中
 Customer Johnson and Johnson k.k.

事業所名:
 Plant

部門名:
 Division

製品名: リークカレントハイテスター
 Product Leak Current + Hi Tester

型番: 3156
 Product Code

製造者名: 日置電機
 Manufacturer HIOKI

製造番号: 051238345
 Manufacturing Number

管理番号: F-56
 Control Number

校正日: 2018年11月5日
 Date of Calibration 11/15/2018

上記の製品は、当社の管理規定にもとづき校正されています。
 使用した基準器は国際度量衡委員会(CIPM)/国際度量衡局(BIPM)に加盟する、国立研究開発法人産業技術総合研究所(AIST)等の国家、国際標準にトレーサブルである事を証明します。

ISO / IEC 17025 (RCL00030) 認定校正機関
 (直流電圧・電流/交流電圧・電流/直流抵抗)
 ISO9001 (JQA-QMA15393), ISO14001 (JQA-EM7241) 認証取得



ユウアイ電子株式会社

埼玉県川越市かし野台2-22-4

TEL:049-243-8611 FAX:049-242-1928

E-mail:calroom@yuai.co.jp



試験成績書

(1/1)

顧客名 ジョンソン・エンド・ジョンソン株式会社 御中
 Customer Johnson and Johnson K.K.

校正日 2018年11月5日
 Date of Calibration 11/05/2018
 温度 23 °C 湿度 47 %
 Temperature 23 °C Humidity 47 %

製品名 リークカレントハイテスター
 Product Leak Current Hi Tester
 製造者名 日置電機
 Manufacturer HIOKI
 管理番号 F-56
 Control Number

型番 3156
 Product Code
 製造番号 051238345
 Manufacturing Number

校正者 金子正志
 Tester

承認者 松本
 Approved

判定 合格 Pass
 Result Pass
 埼玉県川越市かし野台2-22-20
 ユウアイ電子株式会社

上記の製品は、当社の管理規定にもとづき校正されています。使用した基準器は国際度量衡委員会(CIPM)/国際度量衡局(BIPM)に加盟する、国立研究開発法人産業技術総合研究所(AIST)等の国家、国際標準にトレーサブルです。

試験項目 Test Contents

直流電流測定 DC Current

レンジ Range	分解能 Resolution	試験点 Test Point	確度 Spec	下限値 Lower	測定値 measured	上限値 Upper	結果 Result
50 μA	0.01	10 μA	1 % fs	9.50 μA	9.96 μA	10.50 μA	合格 V Pass
500 μA	0.1	100 μA	1 % fs	95.0 μA	99.9 μA	105.0 μA	合格 V Pass
5 mA	0.001	1 mA	0.2 + 3	0.995 mA	0.998 mA	1.005 mA	合格 V Pass
25 mA	0.01	10 mA	0.2 + 3	9.95 mA	9.99 mA	10.05 mA	合格 V Pass
確度 ±(% rdg+dgt)							

交流電流測定 AC Current f=50Hz

レンジ Range	分解能 Resolution	試験点 Test Point	確度 Spec	下限値 Lower	測定値 measured	上限値 Upper	結果 Result
50 μA	0.01	10 μA	2 % fs	9.00 μA	10.01 μA	11.00 μA	合格 V
500 μA	0.1	100 μA	2.0 + 6	97.4 μA	100.1 μA	102.6 μA	合格 V
5 mA	0.001	1 mA	2.0 + 6	0.974 mA	0.998 mA	1.026 mA	合格 V
25 mA	0.01	10 mA	2.0 + 6	9.74 mA	9.98 mA	10.26 mA	合格 V
確度 ±(% rdg+dgt)							

ネットワーク 患者測定モードT-1 T-2間で測定。

data end

使用基準器 Reference Equipment

No	基準器名 Name of Equipment	製造番号 Manufacturing Number	管理番号 Control Number	基準器の有効期限 Due Date
1	デジタルマルチメータ DMM 34410A	MY53003789	U-1-017	19年10月 V
2	マルチプロダクトキャリブレーター 5500A Multi Product Calibrator	7455002	U-1-119	19年08月 V

5.0 PHYSICAL REQUIREMENTS



Attention: Factbook # FB003289

PHYSICAL REQUIREMENTS

Date: October 29, 2019

From: Jason Stivers


Re: 1000

Refer to Attachment for the support documentation that the J&J K.K. Medical Products Service and Repair Depot, Sukagawa, Japan has the necessary physical requirements to provide adequate space for the service and repair of the 1000.

The service center meets the electrical power requirements and has separate work spaces for product receiving/shipping, decontamination, service area, and electrical safety.

E-Sig in EpiCenter

Jason Stivers
Service Engineer, EES – Service Staff Engineer



6.0 OPERATING AGREEMENT



Attention: Factbook # FB003289

OPERATING AGREEMENT

Date: October 29, 2019

From: Jason Stivers

Re: 1000

EES and Johnson & Johnson Medical KK, Japan have entered into an operating agreement. The Operating Agreement is located in International Contract database and can be located under e-ICD # 1340355. The most recent version of that agreement is also attached to factbook.

Additionally, Johnson and Johnson K.K. Medical Products, has entered into a quality agreement with Ethicon Endo-Surgery, LLC. The most recent version of that agreement is in Adaptive under #100592785.

E-Sig in EpiCenter

Jason Stivers
Service Engineer, EES – Service Staff Engineer

SERVICE OPERATING AGREEMENT

THIS SERVICE OPERATING AGREEMENT (“Agreement”), effective as of the (February 7, 2019), between **ETHICON ENDO-SURGERY, Inc.**, (“EES”) an Ohio corporation, having offices at 4545 Creek Road, Cincinnati, Ohio 45242 and **JOHNSON AND JOHNSON K.K. Medical Company** with offices at 5-2 Nishi-kanda 3 chome Chiyoda-ku, Tokyo 101-0065, Japan.

RECITALS

WHEREAS, EES is engaged in the manufacturing, distribution, and sale of medical device products and employs the necessary qualified personnel and/or has agreements with service providers, including without limitation, SERVICE CENTER, to provide the services described hereto; and

WHEREAS, SERVICE CENTER has the capability and expertise to perform the Services described herein; and has been engaged to provide for the service and repair of affiliate units including but not limited to the purchase of spare parts and technical services;

NOW, THEREFORE, in consideration of the foregoing and the mutual agreements, provisions and covenants contained herein, and for other good and valuable consideration, the receipt and legal sufficiency whereof are hereby acknowledged, the parties hereto further agree as follows:

ARTICLE I DEFINITIONS

Section 1.1 Definitions

- (a) **Affiliate**. “Affiliate” of a party shall mean any entity or person that directly or indirectly controls, is controlled by or is under common control with such party. For purposes of this definition, “control” shall mean the possession, directly or indirectly, of the power to direct or cause the direction of the management and policies of an entity, whether through the ownership of voting securities, by contract or otherwise.
- (b) **Agreement**. “Agreement” shall mean this agreement, including the appendices, exhibits and attachments hereto, and any Modifications issued hereunder.

- (c) Management Services. “Services” shall mean (i) all services, tasks, functions and other responsibilities and activities associated to the service and repair of the EES products to be performed by SERVICE CENTER in accordance, including without limitation, the provision of providing spare parts, technician training, service bulletins; and (ii) all Implied Services. Services are described in more detail in Attachment 2, attached hereto and incorporated herein.
- (d) Technician In Training. “Technician In Training” shall mean an employee of SERVICE CENTER with at least 2 years of formal electro-mechanical training or equivalent experience.
- (e) Term. “Term” shall mean the period of time during which this Agreement is in effect in accordance herewith.

ARTICLE II TERM

Section 2.1 Term.

- (a) This Agreement shall be effective for an initial period of twelve (12) months from and after the Effective Date. At the end of this period, the Agreement shall automatically renew for additional successive twelve (12) month Terms, unless one party shall give the other party thirty (30) days' notice of its intention not to renew the Agreement at the end of its then-current term.
- (b) This is the entire Agreement. Previous agreements before the agreement hereto are superseded. This agreement can only be modified in writing, in an amendment signed by both parties

Section 2.2 Termination

- (a) Upon termination of business by either party, for one product or all business, the responsibility for the transfer of assets will be as follows:
- a. SERVICE CENTER will transfer all tooling, equipment and materials owned by EES or its independent contractor to the location determined by EES or its independent contractor.

-
- b. SERVICE CENTER will provide all original documentation, or copies including device master records, device history records, drawings, material inspection documentation, process and inspection instruction sheets for all EES products to EES.

ARTICLE III RETENTION OF SERVICE CENTER

Section 3.1 Performance of Services

- (a) EES hereby engages and retains SERVICE CENTER to perform the Services described in Article IV hereof and SERVICE CENTER hereby accepts and agrees to provide such Services to EES upon the terms and conditions set forth herein.
- (b) SERVICE CENTER shall, in its sole discretion, determine the corporate facilities to be used in rendering the Services and the personnel who will render such Services.
- (c) Nothing herein shall be deemed to restrict SERVICE CENTER or its directors, officers or employees from engaging in any business, or from contracting with other parties, including, without limitation, other subsidiaries of SERVICE CENTER, to provide such services contracted through this Agreement. Nothing herein shall be deemed to restrict EES or its directors, officers or employees from engaging in any business, or from contracting with other parties, to obtain services such as those contracted through this Agreement.
 - a. SERVICE CENTER and EES hereby agree that all Services for which EES engages SERVICE CENTER hereunder including but not limited to the service and repair of EES capital device units, will be done in accordance with the procedures and processes requirements as established in: Product Qualification at Service Centers per EES procedure WE001534, "Establishment and Certification of Service Depot, Field Service and Parts Depot Centers."

Section 3.2 Disclaimer, Limited Liability

- (a) Neither party makes any express or implied representations, warranties or guarantees relating to the subject matter of this Agreement.
- (b) EES will use reasonable efforts to make the Services available with substantially the same degree of care as it employees in making the same or similar services available for its own operations; provided, however, that EES shall not be liable to SERVICE CENTERS or any other person for any loss, damage or expense which may result there from or from any change in the manner in which EES renders the Services, so long as EES deems such change necessary or desirable in the conduct of its own operations.

ARTICLE IV SERVICES

SERVICE CENTER shall, at the request of EES, provide Services in the following areas as required by EES, or such other services as the parties may agree from time to time.

Section 4.1 Good Manufacturing Practices

- (a) SERVICE CENTER agrees to comply with Good Manufacturing Practices (GMP) for service and repair of all EES products as defined by the FDA Code of Federal Regulations, 21 CFR – Part 820, the applicable ISO 9000 and /or 13485, Quality Systems standards, and EES requirements.
- (b) These requirements include but are not limited to:
 - a. Maintaining a certified quality system is preferred. When the SERVICE CENTER does not maintain a certified quality system, the establishment process will ensure the FDA Code of Federal Regulations, 21 CFR – Part 820 is adhered to.
 - b. Maintain viable quality assurance functions separate from and not reporting to manufacturing.
 - c. Conduct ISO 9000, 13485, or similar training as required by the service center's Quality Management System of all personnel on a regular basis and at least once every two years.

- d. Conduct complaint training at least once every year.
 - e. Comply with the listed procedures and process requirements as well as with spare parts serial number traceability.
- (c) Maintain and comply with the Documentation Retention Policy. Each of SERVICE CENTER and EES will maintain all records for which it is responsible that are required to support compliance to FDA Code of Federal Regulations, 21, CFR – Part 820 in compliance with Worldwide Records and Information Management; WWRIM Standard RIMS-12.
- (d) Provide, within five business days, any information and/or records requested by EES.

Section 4.2 Product and Service Quality

- (a) In compliance to WE001534, SERVICE CENTER shall maintain at least one technician who received certified training by EES for each product line. Said technician(s) will be given a training certificate upon completion of the EES-provided training and will then be qualified to train other technicians at SERVICE CENTER on the service and repair of EES capital device products at EES. All associates working on EES products will be trained by EES or its designated independent contractor trained associate. EES will provide all technical service documentation required to service the EES products.
- (b) SERVICE CENTER will perform the activities of service and repair on EES products only if (i) the individual performing such services have a minimum of 2 years formal electro-mechanical training or equivalent experience, as determined by the center; and/or (ii) the individual performing such services is considered to be a Technician In Training and SERVICE CENTER employs an individual with a minimum of 2 years formal electro-mechanical training or equivalent experience, as determined by the center, which overviews the work performed by such Technician In Training.
- (c) SERVICE CENTER shall not commence servicing any device until the requirements of WE001534 have been met and documented, also per WE001534, in writing from EES.

- (d) Replacement/spare parts are only those listed in the service manual or service bulletins. All spare parts shall be purchased through EES. Serial numbers of spare parts used to repair a device shall be recorded in the service record.
- (e) SERVICE CENTER will use the medical device sector's complaint and service data bases to capture service records. SERVICE CENTER that have no access to these databases will provide EES with weekly service reports containing the following information
 - a. Product Code
 - b. Serial Number
 - c. Problem reported by the customer
 - d. Service performed
 - e. Date equipment was received at the Service Center
 - f. Date of service
 - g. Person who performed the Service
 - h. Was the customer problem verified?
 - i. Which part/parts caused the customer's reported problem
 - j. Which parts were related to the reported issue
 - k. Which parts were unrelated
 - l. Description(s), part number(s), and quantity of parts used.
- (f) All documentation related to the service and repair of EES products must be in English.
- (g) SERVICE CENTER shall comply with and/or execute all EES provisions and/or modification instructions regarding the equipment when carrying out service and repair work according to Service Manuals and Service Bulletins. SERVICE CENTER likewise shall check and record the function of equipment as per current equipment specifications once it has been repaired.
- (h) EES shall inform SERVICE CENTER of all hardware and software changes of equipment and their effect on the functionality of the equipment. EES shall also make the relevant technical modification documentation available to SERVICE CENTER. This shall be accomplished per the EES service bulletin process.

- (i) SERVICE CENTERS will be notified of product/process/part changes through EES service bulletins issued by the equipment service team.
- (j) Annual Business/Quality Reviews will be conducted between EES and SERVICE CENTER as described in Attachment 3, attached hereto and incorporated herein.
- (k) SERVICE CENTER has the options (or a combination thereof) listed below for communicating complaint and service information. This information will be communicated in accordance with Johnson & Johnson global medical device reporting mandates. SERVICE CENTER shall work closely with its Johnson & Johnson complaint management organizations to ensure reportable linkages between the service records and EES product issue records.
 - a. SERVICE CENTER may use a Johnson & Johnson complaint management system; or,
 - b. SERVICE CENTER may provide the information described in paragraph 4.2 (e) above using EES's service database/complaint system.
- (l) EES shall respond to escalations that require technical support

Section 4.3 End of Life

- (a) SERVICE CENTER will continue to service each device for seven (7) years after the end of manufacturing or a period determined by EES.
- (b) End of product notification will be communicated to SERVICE CENTER via a service bulletin.

ARTICLE V **RECORD RETENTION**

Section 5.1 Record Retention

- (a) SERVICE CENTER shall be responsible for retention of documentation generated in the performance of the Services. SERVICE CENTER shall comply with the Documentation Retention Policy. All records required to support compliance to

FDA Code of Federal Regulations, 21, CFR – Part 820 will be stored for fifteen (15) years after the last of the Services has been discontinued.

- (b) To the extent the Services under this Agreement involve the creation of records that demonstrate regulatory, legal and/or other compliance, the party creating such records and/or its subcontractors will comply with the provisions hereto, and will include such provisions in its contracts with any subcontractors it engages to perform the Services.

ARTICLE VI **COMPENSATION AND ASSET MANAGEMENT**

Section 6.1 Compensation for Services

- (a) Pricing for service and repair of Affiliate units shall be negotiated between EES and SERVICE CENTER.
- (b) SERVICE CENTER will place purchase orders to EES or any of its independent contractors for any service and/or spare parts to be provided. Compensation for services rendered will be paid to the company to which the SERVICE CENTER has placed the purchase order with.

Section 6.2 Asset Management

- (a) SERVICE CENTER, as the owner of the unique capital assets, has financial responsibility for maintaining them.
- (b) SERVICE CENTER agrees to keep and maintain service life history for all EES assets used by the SERVICE CENTER in the service and repair of EES capital products. The service history will be maintained in a Johnson & Johnson service database.

- (c) SERVICE CENTER agrees to maintain all assets required to service and repair EES products to ensure continuity of product service. SERVICE CENTER agrees to fund all maintenance expenses to refurbish EES related custom test fixtures/tooling.

ARTICLE VII REGULATORY AND SAFETY

Section 7.1 Regulatory and Safety.

- (a) SERVICE CENTER will provide access and support for the audit plan for EES/ Johnson & Johnson/ISO/FDA quality system audits. Written Corrective Action plans will be submitted to the EES Compliance organization per applicable audit procedures.
- (b) With respect to all environmental, safety and industrial hygiene matters related to the SERVICE CENTER activities in supply of services to EES, SERVICE CENTER shall (a) comply with all applicable laws and regulations issued by national, state and local authorities, (b) Inform EES promptly of any significant adverse events (e.g. fires, explosions, accidental discharges), (c) Inform EES promptly of any allegations of findings of violations of applicable laws or regulations, (d) allow EES to inspect the SERVICE CENTER facilities, such inspection to be at reasonable times and upon reasonable notice, and (e) implement promptly and corrective action which may be reasonably requested by EES including (without limitation) adhering to reasonable and significant elements of the Environmental, Safety and Industrial Hygiene Program adhered to EES in its own operations.
- (c) SERVICE CENTER agrees to notify EES immediately when contacted by a regulatory agency that has jurisdiction over EES products or business.

ARTICLE VIII MISCELLANEOUS

Section 8.1 Indemnity. SERVICE CENTER assumes all liability for and agrees to defend, indemnify and hold EES, its employees, officers, directors, shareholders and agents harmless from

and against all demands, liability, damages, costs and expenses, including attorneys' and expert witness fees ("Loss"), incurred by EES arising from or in connection with the Services, other than any loss caused by the gross negligence or willful misconduct of EES.

Section 8.2 Notices. All notices and other communications hereunder shall be in writing and shall be delivered by hand or through electronic mail or mailed by registered or certified mail (return receipt requested) or transmitted by facsimile to the parties at the following addresses (or at such other addresses for a party as shall be specified by like notice) and shall be deemed given on the date on which such notice is received:

If to EES:

Ethicon Endo-Surgery, Inc.
4545 Creek Road
Cincinnati, Ohio 45242
ATTN: Service and Repair Manager

If to (SERVICE CENTER):

Johnson and Johnson K.K. Medical Company
1 Meotozaka Ohkanbara Sukagawa Fukushima 962-8501 Japan
ATTN: Technical Services Manager

Section 8.3 Independent Contractor. This Agreement shall not be construed to appoint either SERVICE CENTER or EES as agent of the other, and the parties shall at all times be considered as an independent entity in the performance of this Agreement. Neither party shall be responsible in any way for any obligations or liability incurred or assumed by the other party, except as expressly set forth in this Agreement.

Section 8.4 Force Majeure. Neither party shall be in default of this Agreement or liable to the other party for any delay or default in performance where occasioned by any cause of any kind or extent beyond its control, including but not limited to, armed conflict or economic dislocation resulting there from; embargoes; shortages of labor, raw materials, production facilities or transportation; labor difficulties; civil disorders of any kind; action of any civil or military authorities (including priorities and allocations); fires; floods; and accidents. The dates on which obligations of a party are to be fulfilled shall be extended for a period equal to the time lost by reason of any delay arising directly or indirectly from;

- (a) Any of the foregoing causes, or

- (b) Inability of that party, because of causes beyond its reasonable control, to obtain instruction or information from the other party in time to perform its obligations by such dates.

Section 8.5 Entire Agreement. This Agreement constitutes the entire understanding between the parties with respect to the subject matter hereof and all prior agreements or understandings shall be deemed merged herein. No representations, warranties and certifications, express or implied, shall exist as between the parties except as stated herein.

Section 8.6 Amendments. No amendments, waivers or modifications hereof shall be made or deemed to have been made unless in writing executed by the party to be bound hereby.

Section 8.7 Severability. If any provisions in this Agreement or the application of such provision to any person or circumstance shall be invalid, illegal or unenforceable, the remainder of this Agreement or the application of such provision to persons or circumstances other than those to which it is held invalid, illegal or unenforceable shall not be affected thereby.

Section 8.8 Counterparts. This Agreement may be executed in any number of counterparts, each of which when so executed shall be deemed to be an original and all of which when taken together shall constitute this Agreement.

Section 8.9 Successors and Assigns. This Agreement shall not be assignable, in whole or in part, directly or indirectly, by any party hereto without the prior written consent of the other party hereto, and any attempt to assign any rights or obligations arising under this agreement without such consent shall be void. This Agreement shall be binding upon and inure to the benefit of the parties hereto and their respective successors and permitted assigns.

Section 8.10 Governing Law. This Agreement shall be governed by and construed in accordance with the laws of the State of Ohio without reference to the conflict of law doctrine of such state.

Section 8.11 No Third-Party Beneficiaries. This Agreement is solely for the benefit of the parties hereto and should not be deemed to confer upon third parties any remedy, claim, liability, reimbursement, claim of action or other right in excess of those existing without reference to this agreement.

Section 8.12 Confidentiality. The SERVICE CENTER shall treat as confidential all information, including customer data, which SERVICE CENTER acquires in the context of the co-operation

with EES and/or its affiliates and shall not make information accessible to third parties for the period of validity of the Agreement and after its expiration, with the exception of those employees of the SERVICE CENTER to whom information has to be revealed for the purpose of executing this Agreement. The SERVICE CENTER shall have such employees undertake in writing to maintain confidentiality.

IN WITNESS THEREOF, the undersigned have caused this Agreement to be duly executed and operable the Effective Date.

ETHICON ENDO-SURGERY, INC.

By: MARK COSIMI

Name: Mark Cosimi

Title: Director, Project Management

Date: _____

JOHNSON & JOHNSON K.K. Medical

Company

KOHEI SEKI

By: _____

Name: Kohei Seki

Title: Repair Services Manager

Date: _____

Digitally signed by MARK COSIMI
DN: c=US, o=JNJ, ou=Subscribers, cn=MARK
COSIMI, 0.9.2342.19200300.100.1.1=32153
Reason: I am approving this document.
Date: 2019.09.13 11:15:20 -04'00'

Digitally signed by KOHEI SEKI
DN: c=US, o=JNJ, ou=Subscribers, cn=KOHEI
SEKI, 0.9.2342.19200300.100.1.1=152816107
Reason: I attest to the accuracy and integrity of
this document.
Date: 2019.09.13 07:45:20 +09'00'
Adobe Reader version: 11.0.10

Attachment 1**Key Contacts**

The primary contacts are:

(SERVICE CENTER):

Kohei Seki, Repair Services Manager
Kenta Tomikoshi, Quality Assurance Manager

Ethicon Endo Surgery, Inc.:

Eric Smith, Sr Manager, Service Center
Shannon Gillespie, Service Center Manager
Milton Garrett, Director, Customer Quality
Robert Peters, Team Lead, Customer Quality

Attachment 2
Description of Services

Article I. DEFINITIONS

The terms used in this Attachment 2 shall have the meaning as set forth in Article 1 of the Agreement.

Article II. PRODUCT DESCRIPTION

- 1000 – MegaPower 1000 A Generator
- 0800, 0800S, 0830, 0830S, 0835, 0835S, 0840, 0840S, 0845, 0845S, 0846, 0846S, 0847, 0847S, 0848, 0848S - MEGADYNE MEGA SOFT Reusable Patient Return Electrodes
- GEN04 – Generator GEN04
- GEN11 – Generator GEN11
- RF60 – Generator RF-60

Article III. SCOPE

This agreement covers all aspects of the service and repair process including decontamination, product, evaluations, troubleshooting, repair, spare parts, testing, quality assurance, packaging, daily reporting, and database administration.

Article IV. SERVICE REQUIREMENTS

SERVICE CENTER is responsible for:

- Generating, maintaining and controlling internal process work instructions for servicing EES designated repairable products.
 - EES Service Manuals and Service Bulletins will address the following areas:
Service and Repair, testing, upgrades, and product certification and decontamination (when applicable).
- Return and packaging of replaced devices/parts to the EES Service Center along with the corresponding Service Report (when applicable).
- Complete copy of the service report for the complaint handling system
- Product specific letters of service to customers
- Storage and inventory of EES capital devices and spare parts

Article V. DATABASE ADMINISTRATION

The main interface to EES complaint and service databases will be using Unity and OneMD database software. It is the joint responsibility of both EES and SERVICE CENTER to create, and maintain, compliant/service documentation. SERVICE CENTER technicians will complete the necessary documentation for service and repair. When necessary, the SERVICE CENTER will ensure timely return of documents to EES for entry into the EES database software. When

the SERVICE CENTER technicians can enter their documentation into a service software database, they will be responsible for entering their own entries in the database.

Article VI. PRODUCT UPGRADES

Only approved product upgrades will be made to the product and the DHR. EES service bulletins or ECO changes will define approved product upgrades. All product service upgrades will be documented through the service manual or service bulletins. The upgrades will be conducted at time of service.

Article VII. PARTS AND ACCESSORY INVENTORY

SERVICE CENTER will store and provide inventory control for the parts and accessories designated by EES. These parts and accessories are owned by SERVICE CENTER. Non-conforming parts and accessories will be segregated from the inventory and must be dispositioned.

Article VIII. SERVICE PROCESS INSTRUCTIONS

Servicing is conducted using the current revision of the service manual and released service bulletins. The applicable service manual may be used as a reference guide for troubleshooting and repair. Service Manuals and Service Bulletins will be supplied by EES to change or aid the service and repair process. SERVICE CENTER is responsible for providing relevant and valid process documents. EES is responsible to conduct training on the repair process for each product serviced by the EES. SERVICE CENTER is responsible for being trained and performing re-certification training annually.

Article IX. TROUBLESHOOTING AND REPAIR

Planned maintenance, troubleshooting, repair, and testing will be performed utilizing approved test fixtures and calibrated equipment. Fixtures will be documented and have current preventive maintenance/calibration schedules completed.

Planned maintenance, troubleshooting, repair, and testing will be conducted by trained technicians only and the results of their activities, findings, service and repair will be recorded and reviewed in accordance with the Service and Repair Procedure. Calibrated equipment used in the service and repair process will be documented in the service record containing the equipment calibration number and calibration due date.

Article X. PARTS PROCUREMENT AND INSPECTION

All spare parts will only be purchased from the EES Part Depot. The parts planning, procurement, inspection, and warehousing will be performed by SERVICE CENTER when received by the service center.

Article XI. EQUIPMENT

The equipment used for servicing and validation of product are detailed in the product specific Validation Master Plan. SERVICE CENTER will only use equipment that has been successfully validated for service and repair.

Article XII. SERVICE LEVEL

The goal is to service a customer's EES device to manufacturer's specification within five business days from receipt of the device. SERVICE CENTER can set the required service timeliness based on the customer needs within their country or region.

Attachment 3**Business/Quality Reviews**

Business/Quality Reviews between EES and SERVICE CENTER will be held annually. The minutes will be reviewed by the service and quality organizations from both EES and SERVICE CENTER. The topics to be reviewed at a minimum will include:

1. Follow Ups from previous review
2. Service and Quality Agreements Update
3. Review of Service Metrics
4. Review of Specifications
5. Review of trending data by product
6. Service documentation review
7. Data entry into Service Database/Service reports
8. Project Review and Optimization
9. Loaner, Surplus Pools Review
10. Cycle Time for Complaints
11. Service Bulletins
12. Changes (Personnel)
13. Audits
14. Miscellaneous

7.0 START-UP ACTIVITIES



START-UP ACTIVITIES

Date: October 29, 2019

From: Jason Stivers

Attention: Factbook # FB003289

Re: 1000

In accordance with WE001534 Rev F, the following start-up activities have taken place to prepare the EES Service and Repair Depot, Cincinnati, Ohio to begin service and repair on the Megadyne Mega Power 1000 Electrosurgical Generator.

- Spare parts needed to perform repairs have been communicated to the service center and will be ordered by the center from the Ethicon parts depot upon factbook approval.
- Johnson & Johnson Medical KK has been added to the authorized supplier list. (EES CASL)
- Loaner pool processes have been established per business needs.
- The linkage to the Complaint Management/Service System has been established. The service center has been using the process for the routing of all service data to the Complaint Management/Service System and will include 1000 reporting as well.
- Equipment required to perform repairs and testing has been purchased and properly installed.
- The service center has been made aware of the requirements for repair tracking and expediting through the repair center, and monthly reporting requirements.
- Products have been added to the monthly reporting process.
- Product-related service documentation has been delivered and implemented.
- Established access to the Megadyne SharePoint site has been completed.

E-Sig in EpiCenter

Jason Stivers
Service Engineer, EES – Service Staff Engineer

8.0 SUPPLIER APPROVAL



Supplier Approval

Date: October 29, 2019

From: Jason Stivers

Re: 1000

Attention: Factbook # FB003289

Reference DOC000657 in Epi Center for supporting documentation that established Johnson and Johnson K.K. Medical Products as an authorized supplier. The service center will be maintained on the approved supplier list according to WE0652, Monitoring/Communication with Service Centers and Parts Depots. A snapshot of the Corporate Approved Supplier List is below showing Johnson and Johnson K.K. Medical Products, Sukagawa, Japan as currently being approved as of October 29, 2019.

Screenshot of an Excel spreadsheet titled "Approved Supplier List.xls" showing the "Approved Supplier List". The spreadsheet contains data from rows 389 to 408, columns A through H. The data includes various company names, addresses, and contact information for different countries.

	A	B	C	D	E	F	G	H
389	137001	137001-01	Jarden Plastics Solutions	20 Setar Way	Reedsville	PA	17084	USA
390	255058001	255058001-01	Jesus Manuel Duarte Buncen	Tapachula 356	Juarez	CH	32674	Mexico
391	1174621	1174621-01	JF Manufacturas/Maquinados El	Amozoc 7508	Juarez	CH		Mexico
392	430069833	430069833-01	John Keith Rhoden BDA Pacific Magnefor	15113 Sierra Bonito Ladr	Chino	CA	91710	USA
393	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson Customer & Logistics	23 Orchard Road	Skillman	NJ	08558	USA
394	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson De Chile S.A.	San Fransisco 253,	Santiago, RM			Chile 8330267
395	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson De Chile S.A.	San Fransisco 253	Santiago, RM			Chile 8330267
396	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson do Brasil Indústria e C	Rod. Presidente Dutra	Sao Paulo			Brazil
397	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson Health Care Systems	3955 E. Holmes Road,	Memphis	TN	38118	USA
398	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson K.K. Medical Compar	Sukagawa Plant, 1 Meio	Fukushima	962-8501	Japan	
399	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson K.K., Medical Compar	Sukagawa Plant, 1 Meio	Sukagawa City	Fukushima	962-8501	Japan
400	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson Medical Australia	Mar 1-5 Khartoum Road	North Ryde NSW 2113			Australia
401	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson Medical Brasil (J&J P	Avenida Julia Gaioli No	Aqua Chata	Guarulhos		Brazil
402	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson Medical Brasil (J&J P	Av Honorio Alvares Per	Sao Paulo			Brazil
403	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson Medical Canada	200 Whitehall Drive	Markham	ON	L3ROT5	Canada
404	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson Medical Caribbean - f	Los Frailes Industrial Pa	Guaynabo	Puerto Rico	00969	USA
405	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson Medical Caribbean - f	Los Frailes Industrial Pa	Guaynabo	Puerto Rico	00969	USA
406	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson Medical China	14F,Chend Kai Internati	Shanghai		200030	China
407	N/A - J&J Affiliate	N/A - J&J Affiliate	Johnson & Johnson Medical Columbia	Calle 20 # 68-C 57 Bod	Santa Fe De Bog	DC		Columbia
ANR	N/A - IR. I Affiliate	N/A - IR. I Affiliate	Johnson & Johnson Medical Columbia	Calle 20 # 68-C 57 Bod	Santa Fe De Bog	DC		Columbia

E-Sig in EpiCenter

Jason Stivers
Service Engineer, EES – Service Staff Engineer