



The Electrosurgical Authority®

DOCUMENT NUMBER: ENG-RPT-535

DOCUMENT TITLE: ZIP ACE Modified, 6-Pack Ship Test Report, T=0

DOCUMENT NOTES:

Document Information

Revision: 001

Vault: MEG-Rel

Status: Release

Document Type: ENG-RPT

Date Information

Effective Date: 22 Jun 2018

Expiration Date:

Release Date: 22 Jun 2018

Next Review Date:

Control Information

Author: TSKINNER

Owner: TSKINNER

Previous Number:

Change Number: 2018-ENG-DCO-043

Signature Manifest**Document Number:** ENG-RPT-535**Revision:** 001**Title:** ZIP ACE Modified, 6-Pack Ship Test Report, T=0

All dates and times are in Mountain Standard Time.

Zip ACE Mod T=0 Ship testing**Change Request**

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		07 Feb 2018, 10:05:07 AM	Approved

Collaboration

Name/Signature	Title	Date	Meaning/Reason
Joni Stegeman (JSTEGEMAN)	Ethicon Quality	14 Jun 2018, 06:18:53 PM	Complete
Paul Borgmeier (PBORGMEIER)		18 Jun 2018, 09:50:12 AM	Complete
Darlene Hull (DHULL)	Regulatory	18 Jun 2018, 10:56:54 AM	Complete
Tyler Skinner (TSKINNER)	Project Engineer	19 Jun 2018, 01:36:40 PM	Complete

Document Review

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		19 Jun 2018, 01:46:38 PM	Complete

RA-Approval

Name/Signature	Title	Date	Meaning/Reason
Darlene Hull (DHULL)	Regulatory	19 Jun 2018, 02:52:08 PM	Approved

QA-Approval

Name/Signature	Title	Date	Meaning/Reason
Joni Stegeman (JSTEGEMAN)	Ethicon Quality	19 Jun 2018, 04:06:01 PM	Approved

ENG-Approval

Name/Signature	Title	Date	Meaning/Reason
Paul Borgmeier (PBORGMEIER)		21 Jun 2018, 03:30:51 PM	Approved

Training Review

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		22 Jun 2018, 09:16:38 AM	Approved

Final Release

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		22 Jun 2018, 09:16:50 AM	Approved

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 1 of 37

Authored By: Tyler Skinner

1. ABSTRACT

Test Protocol ENG-PRT-441 was completed successfully. Testing was limited to Bubble Leak testing, Dye testing, Burst testing, Minimum Seal Width testing, and a Product Damage Inspection in order to verify the proposed 6-Pack shipping configuration. All tests were completed successfully.

These results demonstrate that the proposed 6-Pack shipping configuration (ME725M1C and ME725M1E) does not damage the Tyvek pouch seal and protects the product from damage. The results also provide confidence that the product will withstand the anticipated shipping environment and meet DMR requirements in ENG-DMR-012 after EO Sterilization.

2. REFERENCES

ENG-DMR-012	DMR, Smoke Evacuation Pencil and Accessories
ENG-RMF-045	Risk Analysis, Smoke Evacuation Accessories
ENG-PRT-441	ZIP ACE Modified, 6-Pack Ship Test, T=0
ME725M1C	Ace Blade 700, 2.5" Zip Pen, "C" Connector, 10 ft. Tubing
ME725M1E	Ace Blade 700, 2.5" Zip Pen, EC Connector, 10 ft. Tubing

3. OBJECTIVE

This Test Report documents that using the proposed 6-Pack shipping configuration to distribute Zip Pens has no effect on the packaging integrity of the EO sterile product.

4. APPENDICES

Appendix I – 2X EO Exposure
Appendix II – Preconditioning and Ship Testing
Appendix III – Bubble Leak Testing
Appendix IV – Dye Testing
Appendix V – Burst Testing
Appendix VI – Minimum Seal Width Testing
Appendix VII – Product Damage Inspection

5. RESULTS

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 2 of 37

5.1. 2X EO Sterilization

36 Zip Pen samples (6 boxes) SKU: ME725M1C Lot # 170323 were EO sterilized twice per Sterigenics Cycle 115. See Appendix I for EO sterilization documentation.

5.2. Ship Conditioning

36 Zip Pen Samples (6 boxes) SKU: ME725M1C Lot # 170323 were then ship conditioned/tested as outlined in ENG-PRT-441. See Appendix II for a complete summary of testing.

Tests were performed under typical warehouse conditions, which are:

Temperature: 23°C ±5°C

Relative Humidity: 50% ±35%

5.2.1. Preconditioning

Preconditioning followed the schedule below (see Appendix II for Cycle Data):

CONDITIONS	DURATION
Transition from ambient to -40°C	Based on Chamber Capability
Hold -40°C no humidity control	4 hours
Transition from -40°C to 55°C	Set time to 0:00 and set the standard deviation to 1°C
Transition from 55°C to 55°C and 95%RH	Set time to 0:00 and set the standard deviation to 1°C and 2% RH
Hold 55°C and 95%RH	4 hours
Transition from 55°C and 95% RH to 55°C and 15% RH	Set time to 0:00 and set the standard deviation to 1°C and 2% RH
Hold 55°C and 15%RH	4 hours
Transition to 23°C and 50%RH	Set time to 0:00 and set the standard deviation to 1°C and 2% RH
Hold 23°C and 50%RH	72 hours

ing - Manual (Drop Test)

The Manual Handling (Drop Test) was performed using a drop height of 15 in as outlined in ENG-PRT-441.

5.2.3. Vehicle Stacking (Compression Test)

The Vehicle Stacking (Compression Test) was performed using a computed load (L) of 200 lb as outlined in ENG-PRT-441.

5.2.4. Vehicle Vibration and Loose Load Vibration Tests

The Vehicle Vibration test was performed for 10 min as outlined in ENG-PRT-441. Following the Vehicle Vibration test, the Loose Load Vibration Test was performed for 40 min as outlined in ENG-PRT-441.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 3 of 37

5.2.5. Concentrated Impact Test

The Concentrated Impact Test was performed as outlined in ENG-PRT-441.

5.2.6. Manual - Handling (2nd Drop Test)

The Manual Handling (2nd Drop Test) was performed using a drop height of 15 in with the final drop at a height of 30 in as outlined in ENG-PRT-441.

5.2.7. Each box remained intact and did not break open during the test. See Appendix II.

5.3. Bubble Leak Testing

The Bubble Leak test was performed on 35 samples as outlined in ENG-PRT-441. There will no tears, holes, or open seals in any pouch. See Appendix III.

5.4. Dye Testing

The Dye Test was performed on 35 samples as outlined in ENG-PRT-441. There were no breaches in the seal and no signs of separation or degradation. See Appendix IV.

5.5. Burst Testing

The Burst Test was performed on 36 samples as outlined in ENG-PRT-441. All samples passed with a minimum of 23.1 in. H₂O and a mean of 26.625 in. H₂O. See Appendix V.

5.6. Minimum Seal Width Testing

The Minimum Seal Width Test was performed on 35 samples as outlined in ENG-PRT-441. The minimum seal width of all edges exceeded the passing criteria of 0.20" with an average of 0.31" and a minimum of 0.21". See Appendix VI.

5.7. Product Damage Inspection

Product Damage Inspection was performed on 35 samples as outlined in ENG-PRT-441. No damage to the electrode, coating, or any other part of the Zip Pen was observed. See Appendix VII.

6. DISCUSSION

6.1. 2X EO Sterilization

6.1.1. Products were EO sterilized as outlined in the test protocol.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 4 of 37

6.2. Ship Test Conditioning

6.2.1. Shipping and storage conditions reached all extremes required by the test protocol for all test groups.

6.2.2. The acceptance criteria were satisfied for the 6-Pack shipping configuration.

6.3. Bubble Leak Testing

6.3.1. The acceptance criteria were satisfied for the 6-Pack shipping configuration.

6.4. Dye Testing

6.4.1. The acceptance criteria were satisfied for the 6-Pack shipping configuration.

6.5. Burst Testing

6.5.1. The acceptance criteria were satisfied for the 6-Pack shipping configuration.

6.5.2. For Burst Testing a total of 36 samples were tested. This deviates from the required 35, however, this is acceptable as passing results were obtained for all samples and the testing of an extra sample provides an additional opportunity for a packaging failure.

6.6. Minimum Seal Width Testing

6.6.1. The acceptance criteria were satisfied for the 6-Pack shipping configuration.

6.7. Product Damage Inspection

6.7.1. The acceptance criteria were satisfied for the 6-Pack shipping configuration.

7. CONCLUSIONS

7.1. 2X EO Sterilization

7.1.1. Product was successfully exposed to EO sterilization.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 5 of 37

7.2. Ship Test Conditioning

7.2.1. The 6-Pack box meets ASTM D4169.

7.3. Bubble Leak Testing

7.3.1. The 6-Pack shipping configuration does not add additional risk of leaks in the product packaging.

7.4. Dye Testing

7.4.1. The 6-Pack shipping configuration does not add additional risk of breaches in the seal of the product packaging.

7.5. Burst Testing

7.5.1. The 6-Pack shipping configuration does not add additional risk of burst product packaging.

7.6. Minimum Seal Width Testing

7.6.1. The 6-Pack shipping configuration does not reduce the seal width of the product packaging.

7.7. Product Damage Inspection

7.7.1. The 6-Pack shipping configuration does not add additional risk of damage to the product.

8. RECOMMENDATIONS

Based on these test results, it is recommended that EO sterilization does not adversely affect the packaging of Zip Pens in the 6-Pack shipping configuration (ME725M1C and ME725M1E).

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 6 of 37

APPENDIX I – 2X EO EXPOSURE

1st Exposure

The user must ensure that they are using the correct/current revision of this document.
Document: OPER-FRM-016 Rev: 001 Effective: 12 Apr 2013 12:00 AM

EO STERILIZATION

PART#	LOT#	STERILE QTY.	# OF BOXES	DATE OUT
0012M	175317	18480	70	
0012M	175347	18744	71	
ACE12AM	175346	1848	7	
0014M	175403	7920	60	
0014M	175402	10,032	76	
0012M	175465	18480	70	
0014M	175413	8316	63	
0037H	175303	2699	54	
Dura	~	~	14	
ENG TEST	~		1	
ENG Test ZIP PEN			6	
PN: ME725MIC Lot # 5170323				TS 12-20-17
ACE14M	175405	264	2	

Confidential - Controlled copy on date printed

1st EO Exposure

PREPARERS SIGNATURE: [Signature] DATE: 11/9/17
P/O#: Sterigenics: 29511 Nelson: 29512 Full 11/10/17
COMMENTS: BOX TOTAL = 492 494
PALLETS = 4

REVISION HISTORY

REVISION	DOCUMENT CHANGE ORDER NUMBER	DESCRIPTION OF CHANGE	EFFECTIVE DATE
A	00-134-01	Initial Release	2000-09-21

Printed on: 17 Oct 2017, 08:50:45 am; Printed by: EMCKENNA

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT		Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0		Revision: 001
			Page 7 of 37



Certificate of Processing

STERIGENICS 5725 Harold Gatty Drive Salt Lake City UT 84116
TEL 801 328-9901 FAX 801 328-9951 www.sterigenics.com

R55480101

11/20/2017 15:44:19 GMT
Page 1 of 1

Customer Name: Megadyne Processing Facility: Salt Lake City Work Order # 2084963
P.O.# 29511 Sales Order # 1821030

Customer kits Cycle 115 Received Date/Time: 11/10/2017 15:12:02 GMT

SO Line #	Qty	UOM	Description #1	Description #2	Pallet ID	Customer Load Number	Customer Lot No.
101.000	1	PL	Customer kits	Cycle 115	47667197	29511	NA
102.000	1	PL	Customer kits	Cycle 115	47667198	29511	NA
103.000	1	PL	Customer kits	Cycle 115	47667199	29511	NA
104.000	1	PL	Customer kits	Cycle 115	47667200	29511	NA
	4	PL	Total				

Processing Summary

Op#	Operation Name	Location	Date/Time In (GMT)	Date/Time Out (GMT)	Total Time (Hours)
100.00	RECEIVING & LOAD PREP	UNPROC	11/10/2017 15:12:02	11/10/2017 17:39:10	2.45
110.00	WAREHOUSE HOLD	HOLD	11/10/2017 17:39:10	11/17/2017 17:08:38	167.49
120.00	PRECONDITIONING	P3LANE2	11/17/2017 17:08:38	11/18/2017 20:48:38	27.67
130.00	TRANSFER TO CHAMBER	WORKAISL	11/18/2017 20:48:38	11/18/2017 21:04:18	.26
200.00	CHAMBER	CHAMBER3	11/18/2017 21:04:18	11/19/2017 05:56:00	8.86
230.00	TRANSFER TO AERATION	WORKAISL	11/19/2017 05:56:00	11/19/2017 06:36:17	.67
300.00	AERATION	AERATION6	11/19/2017 06:36:17	11/20/2017 07:40:24	25.07
350.00	TRANSFER TO WAREHOUSE	PROC	11/20/2017 07:40:24	11/20/2017 07:50:21	.17
399.00	REVIEW	PROC	11/20/2017 07:50:21	11/20/2017 07:51:58	.03

Total Usage (to nearest whole number): EO 56 LB

Quality Test Summary

Op#	Quality Test Description	Min Spec	Max Spec	Result	Pass/Fail	User	Date / Time
110.00	8 HR Hold (Nov. - Apr.)	YES	YES	YES	Pass	AESPINO	11/17/2017 17:05:29 GMT
110.00	Place Bf's Per Specification	19 EA	19 EA	19 EA	Pass	ARMANDO ESPINO	11/17/2017 17:05:37 GMT
120.00	Precon. Room Temp. Specs. Met	YES	YES	YES	Pass	AESPINO	11/18/2017 20:54:38 GMT
120.00	Prec. Room RH Specs. Met	YES	YES	YES	Pass	JMREYES	11/18/2017 20:54:53 GMT
300.00	Aeration Temperature Specs Met	YES	YES	YES	Pass	JOSE REYES	11/20/2017 07:41:33 GMT
350.00	Remove Bf's Post Chamber	19 EA	19 EA	19 EA	Pass	RFREEMAN	11/20/2017 07:50:45 GMT
						RONALD FREEMAN	

The above products were processed according to the process specification requirements. All parameters reviewed were found to be in compliance with specifications.

Electronically Signed By: TONY VADNAIS
Reason: Work Order Completions

Date: 11/20/2017 15:43:13 GMT

ISO 9001 and ISO 13485 Registered

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 8 of 37



Megadyne Medical Products, Inc.
Work Order 2084963
2017-11-19



Sterigenics.
The Global Leader in
Contract Sterilization Services

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 9 of 37




WORK ORDER: 2084963 CYCLE: 115
SALES ORDER: 1821030 PURCHASE ORDER: 29511
FACILITY: Salt Lake City - Sterigenics CUSTOMER: Megadyne Medical Products, Inc.
5725 Harold Gatty Drive 11506 South State Street
Salt Lake City, UT 84116 Salt Lake City, UT 84020
United States
ITEM: 129_115_CH VERIFICATION RECIPE ID: 16585

Qty	Description #1	Description #2	Pallet ID	Customer Load #	Customer Lot #
1 PL	Customer kits	Cycle 115	47667197	29511	NA
1 PL	Customer kits	Cycle 115	47667198	29511	NA
1 PL	Customer kits	Cycle 115	47667199	29511	NA
1 PL	Customer kits	Cycle 115	47667200	29511	NA
4 PL					

2017-11-10 15:12	RECEIVING & LOAD PREP	UNPROC	2h27m	
	parts list total quantity between 3 and 7	4		PASS
2017-11-10 17:39	WAREHOUSE HOLD	HOLD	6d23h29m	
	duration at least 8h	6d23h29m		PASS
2017-11-17 17:08	PRECONDITIONING	P3LANE2	1d3h40m	
	duration between 1d and 3d	1d3h40m		PASS
	combi verifier			INFO
	- from start of PRECONDITIONING			INFO
	- until end of PRECONDITIONING			INFO
	- Precon Temperatures between 100 and 120 °F			INFO
	- Precon RHs between 45 and 75 %RH			INFO
	- location	P3LANE2		INFO
	- Precon Temperatures (PTSC)	109 .. 110 °F		INFO
	- Precon Temperatures (PTSM)	109 .. 110 °F		INFO
	- Precon RHs (PRHSC)	61 .. 67 %RH		INFO
	- Precon RHs (PRHSM)	61 .. 67 %RH		INFO
	- minimum 1d in spec	1d3h40m		PASS
	- maximum 30m out of spec (consecutive)	0m		PASS
	Precon Temperatures (PTSC) at most 120 °F	109 .. 110 °F		PASS
	Precon Temperatures (PTSM) at most 120 °F	109 .. 110 °F		PASS
	Precon RHs (PRHSC) at most 75 %RH	61 .. 67 %RH		PASS
	Precon RHs (PRHSM) at most 75 %RH	61 .. 67 %RH		PASS
2017-11-18 20:48	TRANSFER TO CHAMBER	WORKA/SL	16m	
	duration at most 1h	16m		PASS
2017-11-18 21:04	CHAMBER	CHAMBER3	8h52m	
	cycle check value must be 6115	6115		PASS

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 10 of 37



combi verifier		INFO
- from start of CHAMBER/VACUUM A		INFO
- until end of CHAMBER/RELEASE		INFO
- Sterilizer Temperatures between 110 and 130 °F		INFO
- location	CHAMBER3	INFO
- Sterilizer Temperatures (STSC)	117 .. 121 °F	INFO
- Sterilizer Temperatures (STSM)	117 .. 122 °F	INFO
- total time out of spec must be 0	0m	PASS
CURRENT TANK: LOT # UTLX902338K17E001266, TARE 304.0 LB		INFO
2017-11-18 21:06 VACUUM A	CHAMBER3	27m
Sterilizer Ctrl Pressure (SCP) between 1.5 and 2.5 inHg at the end of this phase	2.0 inHg	PASS
2017-11-18 21:34 LEAK TEST	CHAMBER3	5m
Sterilizer Ctrl Pressure (SCP) must not increase faster than 0.2 inHg per 5m (between begin and end)	0.1 inHg/5m	PASS
2017-11-18 21:39 NITROGEN DILUTION	CHAMBER3	46m
2017-11-18 21:39 NITROGEN	CHAMBER3	21m
Sterilizer Ctrl Pressure (SCP) between 23.5 and 24.5 inHg at the end of this phase	24.0 inHg	PASS
2017-11-18 22:00 EVACUATION	CHAMBER3	25m
Sterilizer Ctrl Pressure (SCP) between 1.5 and 2.5 inHg at the end of this phase	2.0 inHg	PASS
2017-11-18 22:25 HUMIDIFICATION (PRESSURE)	CHAMBER3	1m
Sterilizer Ctrl Pressure (SCP) increase between 1.0 and 2.0 inHg	1.5 inHg	PASS
2017-11-18 22:26 HUMIDITY DWELL (PRESS)	CHAMBER3	1h
duration between 50m and 1h15m	1h	PASS
Sterilizer Ctrl Pressure (SCP) between 3.0 and 4.0 inHg	3.3 .. 3.7 inHg	PASS
2017-11-18 23:26 GAS A (EO)	CHAMBER3	12m
duration between 5m and 55m	12m	PASS
2017-11-18 23:26 STERILANT	CHAMBER3	12m
Sterilizer Ctrl Pressure (SCP) between 14.0 and 15.0 inHg at the end of this phase	14.5 inHg	PASS
2017-11-18 23:38 GAS DWELL (EO)	CHAMBER3	4h1m
duration between 4h and 4h30m	4h1m	PASS
Sterilizer Temperatures (STSC) between 115 and 125 °F	120 .. 120 °F	PASS
Sterilizer Temperatures (STSM) between 115 and 125 °F	121 .. 121 °F	PASS
Sterilizer Ctrl Pressure (SCP) between 14.0 and 15.0 inHg	14.2 .. 14.7 inHg	PASS
2017-11-19 03:39 AFTER VACUUM	CHAMBER3	20m
Sterilizer Ctrl Pressure (SCP) between 1.5 and 2.5 inHg at the end of this phase	2.0 inHg	PASS
duration between 13m and 1h3m	20m	PASS
2017-11-19 03:59 GAS WASH A	CHAMBER3	1h25m
duration between 1h9m and 1h43m	1h25m	PASS
2017-11-19 03:59 RELEASE	CHAMBER3	19m
Sterilizer Ctrl Pressure (SCP) between 21.5 and 22.5 inHg at the end of this phase	22.0 inHg	PASS
2017-11-19 04:19 EVACUATION	CHAMBER3	23m
Sterilizer Ctrl Pressure (SCP) between 1.5 and 2.5 inHg at the end of this phase	2.0 inHg	PASS

WO2084963

page 2 of 3

printed on 2017-11-20 07:53:41

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 11 of 37



2017-11-19 04:42	RELEASE	CHAMBER3	19m	
Sterilizer Ctrl Pressure (SCP) between 21.5 and 22.5 inHg at the end of this phase		22.0 inHg		PASS
2017-11-19 05:02	EVACUATION	CHAMBER3	23m	
Sterilizer Ctrl Pressure (SCP) between 1.5 and 2.5 inHg at the end of this phase		2.0 inHg		PASS
2017-11-19 05:25	RELEASE	CHAMBER3	23m	
Sterilizer Ctrl Pressure (SCP) between 24.7 and 25.8 inHg at the end of this phase		24.8 inHg		PASS
2017-11-19 05:48	CYCLE COMPLETE	CHAMBER3	8m	
Sterilant Usage Total (SUT) at least 0 lb 5m before the start of this phase		56 lb		PASS
2017-11-19 05:56	TRANSFER TO AERATION	WORKA/SL	40m	
2017-11-19 06:36	AERATION	AERATION6	1d1h4m	
duration between 1d and 3d		1d1h4m		PASS
combi verifier				INFO
- from start of AERATION				INFO
- until end of AERATION				INFO
- Aeration Temperatures between 100 and 120 °F				INFO
- Aeration Circulation between 90 and 110				INFO
- location		AERATION6		INFO
- Aeration Temperatures (ATSC)		104 .. 111 °F		INFO
- Aeration Temperatures (ATSM)		104 .. 111 °F		INFO
- Aeration Circulation (AC1)		99 .. 101		INFO
- minimum 1d in spec		1d1h4m		PASS
- maximum 30m out of spec (consecutive)		0m		PASS
Aeration Temperatures (ATSC) at most 120 °F		104 .. 111 °F		PASS
Aeration Temperatures (ATSM) at most 120 °F		104 .. 111 °F		PASS
2017-11-20 07:40	TRANSFER TO WAREHOUSE	PROC	10m	
2017-11-20 07:50	REVIEW	PROC	2m	
8 HR Hold (Nov. - Apr.) (YES)		YES		PASS
Prec. Room RH Specs. Met (YES)		YES		PASS
Precon. Room Temp. Specs. Met (YES)		YES		PASS
Aeration Temperature Specs Met (YES)		YES		PASS
Remove BI's Post Chamber (19 EA)		19		PASS
Place BI's Per Specification (19 EA)		19		PASS
RESULT OF VERIFICATION		PASS		

It is customers responsibility to ensure acceptance of any result other than PASS noted on this report prior to release of the load.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 12 of 37

Confidential - Controlled copy on date printed

CERTIFICATE OF CONFORMANCE

Process Challenge Device (PCD)

Nelson Laboratories, LLC (NL) hereby certifies that the PCDs prepared conform to the requirements and quality specifications outlined in NL standard operating procedure (SOP) SOP0180.

Scope of Certification:

Product: NLI-PCD-010
Batch Number: BO 1616
Preparation Date: 24 Jul 2017
BI Manufacturer's Information: STERIS
Part #NA005, Lot #0397
Exp. 08 Aug 2018
BI Type: Spore Strip
BI Manufacturer's Spore Population: 1.4×10^6 CFU/BI
NL Verified Population: 1.9×10^6 CFU/BI
Organism: *Bacillus atrophaeus*



Tori Dieffenbacher electronically approved

Released by

Tori Dieffenbacher

02 Aug 2017 09:52

Date and Time

A PCD is a surrogate for the actual product and the appropriateness of a given PCD should be determined through comparative resistance testing. This testing is performed to provide biological monitors for development of a sterilization cycle process. The validated D value for the PCD is available upon request. To confirm authenticity of this document, please contact NL. The most recent signature reflects the latest version of the certificate.

P.O. Box 671830 | Murray, UT 84157-1830 U.S.A. • 8280 South Redwood Road | Salt Lake City, UT 84123-8800 U.S.A.

Rev. 1.3.3

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 13 of 37



CRM42221 Rcvd: 10 Jul 2017
Ambient Use By: 08 Aug 2018

Bo 1616

STERIS



Opened: NA

STERIS Corporation Biological Operations
***Bacillus atrophaeus* BIOLOGICAL INDICATOR CERTIFICATE OF PERFORMANCE**
Exhibit K

CATALOG NUMBER: NA005

PRODUCT NAME: Spordex Strip

PRODUCT LOT NUMBER: C17013

SUBASSEMBLY LOT NUMBER: 0397

EXPIRATION DATE: August 8, 2018

Bacillus atrophaeus NRRL B4418:

Mean Population Recovery* (CFU): 1.4×10^6

Survival Time (Min.): 11.2

D₂₀ Value** (Min.***): 2.7

Kill Time (Min.): 27.3

D_{2H} Value** (Min.***): 1.3

Survival Time (Min.): 5.4

Kill Time (Min.): 13.1

* Colony forming units determined after a preliminary heat treatment.

** 600±30mg/L EO using 100% EO, 54±1°C, 60±10% RH or 160±1°C Dry Heat

*** Determined at time of manufacture by fraction negative procedure after graded exposure to sterilization conditions.

Incubate at 30° to 35°C for seven days

Store at 2° to 24°C and between 30% to 80% RH. Do not use after the indicated expiration date.

Dispose of as you would any other microbiological waste (121°C for a minimum of 30 minutes).

This document certifies that the biological indicator product listed above meets STERIS' quality assurance specifications and the performance criteria suggested by the current revision of the United States Pharmacopeia. This certifies that the product listed above has been tested in compliance with current versions of, ISO 11138-1, ISO 11138-2, and USP.

Quality Systems Representative: [Signature]

Date: 4/18/17

Reviewed By: [Signature]

Date: 4/18/17

LIMITATION OF LIABILITY AND INDEMNITY

Nothing in this Certificate of Performance shall, or is intended to, alter, expand, or diminish the terms and conditions of sale governing your purchase of the Biological Indicator Product from STERIS. In no event, whether as a result of breach of warranty, or tort (including negligence and strict liability) shall STERIS or its suppliers be liable as a result of any statement or information contained in this Certificate of Performance. In addition, STERIS shall not be liable for any consequential or incidental damages including, without limitation, loss of use or damage to your products or equipment, cost of substitute products, or down time costs, allegedly caused by the Biological Indicators Products. The responsibility of STERIS for damages due to injuries or death caused by the Biological Indicator Product shall be limited to that portion of such damages as might be attributable to the negligence or strict liability or other tortious conduct of STERIS.

PC600035 Revision Level: AM BCN# 71402
Effectivity Date: 03/31/2017

Page 1 of 1

STERIS Corporation • 5960 Heisley Road • Mentor, OH 44060-1834 USA • 440-354-2600

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 14 of 37



Sponsor:
Tim Kessinger
MegaDyne Medical Products, Inc.
11506 S. State St.
Draper UT 84020

Biological Indicator (BI) Sterility Test Final Report

Study Number: 1004458-S01
Test Article: Sterigenics Run #2084963
Purchase Order: 29512
Study Received Date: 20 Nov 2017
Testing Facility: Nelson Laboratories, LLC, a Business Unit of Sterigenics International
6280 S. Redwood Rd.
Salt Lake City, UT 84123 U.S.A.
Test Procedure(s): Standard Test Protocol (STP) Number: STP0079 Rev 16
Deviation(s): None

Summary: This BI sterility test was conducted to verify if any viable organisms remained on BIs used to monitor routine and/or validation cycles. All test method acceptance criteria were met.

Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.

BI Information:

Type	Manufacturer	Lot #	Expiration Date	Species Name	Minimum Incubation Time	Incubation Temperature
Spore Strip	STERIS	0397	08 Aug 2018	<i>Bacillus atrophaeus</i>	7 days	30-35°C

Results:

Type	Number Tested	Number Positive	Number Negative
NLI-PCD-010	19	0	19
Environmental Control	1	0	1
Media Negative Control	1	0	1
Positive Control	1	1	0



Dania G. Cortes electronically approved for

Study Director

Derek L. Miller

27 Nov 2017 11:54 AM

Study Completion Date and Time

P.O. Box 571820 | Murray, UT 84157-1830 U.S.A. | 6280 South Redwood Road | Salt Lake City, UT 84123-6020 U.S.A.
www.nelsonlabs.com - Telephone 801 290 7500 - Fax 801 290 7558 - sales@nelsonlabs.com

Page 1 of 1

These results relate only to the test article listed in this report. Reports may not be reproduced except in their entirety. Subject to NLS terms and conditions of www.nelsonlabs.com Rev: 2.6.3

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 15 of 37

2nd Exposure

The user must ensure that they are using the correct/current revision of this document.
Document: OPER-FRM-016 Rev: 001 Effective: 12 Apr 2013 12:00 AM

Confidential - Controlled copy on date printed

EO STERILIZATION				
PART #	LOT #	STERILE QTY.	# OF BOXES	DATE OUT
✓ 0012AM	175837	8976	34	11/30/17
✓ 0012M	175803	17160	65	
✓ 0014M	175799	2772	21	
✓ 0014MD	175648	408	2	
✓ 0012AM	175878	18744	71	
✓ 0014AM	175985	7920	60	
✓ 0014M	175960	9108	69	
✓ ACE36H	175989	300	6	
✓ ACE14M	175959	180	2	
✓ 0037H	175961	2700	54	
ENG TEST ~ N/A				1
ENG TEST ZIPPER (2)				6

2nd EO Exposure

→ PN: ME725 MIL lot # SPT0323 TS 122017

PREPARERS SIGNATURE: [Signature] DATE: 11/30/17

P/O: 29591/29592

COMMENTS: steris nelson TOTAL # 391 boxes / 4 pallets

w/o 2098746
004
001
001
11/19/17

REVISION HISTORY

REVISION	DOCUMENT CHANGE ORDER NUMBER	DESCRIPTION OF CHANGE	EFFECTIVE DATE
A	00-134-01	Initial Release	2000-09-21

Printed on: 17 Oct 2017, 06:50:45 am; Printed by: :EMCKENNA.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 16 of 37



Certificate of Processing

STERIGENICS 5725 Harold Gatty Drive Salt Lake City UT 84116
TEL 801 328-9901 FAX 801 328-9951 www.sterigenics.com

RS5480101

12/10/2017 21:33:17 GMT
Page 1 of 1

Customer Name: Megadyne Processing Facility: Salt Lake City Work Order # 2098746
P.O.# 29591 Sales Order # 1835318
Customer kits Cycle 115 Received Date/Time: 11/30/2017 18:03:59 GMT

SO Line #	Qty	UOM	Description #1	Description #2	Pallet ID	Customer Load Number	Customer Lot No.
101.000	1	PL	Customer kits	Cycle 115	47944534	29591	NA
102.000	1	PL	Customer kits	Cycle 115	47944535	29591	NA
103.000	1	PL	Customer kits	Cycle 115	47944536	29591	NA
104.000	1	PL	Customer kits	Cycle 115	47944537	29591	NA
	4	PL	Total				

Processing Summary

Op#	Operation Name	Location	Date/Time In (GMT)	Date/Time Out (GMT)	Total Time (Hours)
100.00	RECEIVING & LOAD PREP	UNPROC	11/30/2017 18:03:59	11/30/2017 18:45:23	.59
110.00	WAREHOUSE HOLD	HOLD	11/30/2017 18:45:23	12/06/2017 15:14:25	140.49
120.00	PRECONDITIONING	P3LANE12	12/06/2017 15:14:25	12/07/2017 17:30:07	26.25
130.00	TRANSFER TO CHAMBER	WORKA1SL	12/07/2017 17:30:07	12/07/2017 17:39:57	.16
200.00	CHAMBER	CHAMBER3	12/07/2017 17:39:57	12/08/2017 02:33:00	8.86
230.00	TRANSFER TO AERATION	WORKA1SL	12/08/2017 02:33:00	12/08/2017 02:42:30	.16
300.00	AERATION	AERATION9	12/08/2017 02:42:30	12/09/2017 20:12:05	41.49
350.00	TRANSFER TO WAREHOUSE	PROC	12/09/2017 20:12:05	12/09/2017 20:48:22	.60
399.00	REVIEW	PROC	12/09/2017 20:48:22	12/09/2017 20:48:51	.01
Total Usage (to nearest whole number):			EO 55 LB		

Quality Test Summary

Op#	Quality Test Description	Min Spec	Max Spec	Result	Pass/Fail	User	Date / Time
110.00	Place Bfs Per Specification	19 EA	19 EA	19 EA	Pass	VMCINTYRE	12/06/2017 15:15:15 GMT
110.00	& HR Hold (Nov. - Apr.)	YES	YES	YES	Pass	VALERIE MCINTYRE	12/06/2017 15:15:22 GMT
120.00	Precun. Room Temp. Specs. Met	YES	YES	YES	Pass	VALERIE MCINTYRE	12/07/2017 17:33:57 GMT
120.00	Prec. Room RH Specs. Met	YES	YES	YES	Pass	JMREYES	12/07/2017 17:34:09 GMT
300.00	Aeration Temperature Specs Met	YES	YES	YES	Pass	JOSE REYES	12/09/2017 20:13:14 GMT
350.00	Remove Bfs Post Chamber	19 EA	19 EA	19 EA	Pass	ARMANDO ESPINO	12/09/2017 20:49:41 GMT
						AESPINO	
						ARMANDO ESPINO	

The above products were processed according to the process specification requirements. All parameters reviewed were found to be in compliance with specifications.

Electronically Signed By: TONY VADNAIS
Reason: Work Order Completions

Date: 12/10/2017 21:22:54 GMT

ISO 9001 and ISO 13485 Registered

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 17 of 37



Megadyne Medical Products, Inc.
Work Order 2098746
2017-12-08



Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 18 of 37




WORK ORDER: 2098746 CYCLE: 115
SALES ORDER: 1835318 PURCHASE ORDER: 29591
FACILITY: Salt Lake City - Sterigenics CUSTOMER: Megadyne Medical Products, Inc.
5725 Harold Gatty Drive 11506 South State Street
Salt Lake City, UT 84116 Salt Lake City, UT 84020
United States
ITEM: 129_115_CH VERIFICATION RECIPE ID: 16585

Qty	Description #1	Description #2	Pallet ID	Customer Load #	Customer Lot #
1 PL	Customer kits	Cycle 115	47944534	29591	NA
1 PL	Customer kits	Cycle 115	47944535	29591	NA
1 PL	Customer kits	Cycle 115	47944536	29591	NA
1 PL	Customer kits	Cycle 115	47944537	29591	NA
4 PL					

2017-11-30 18:03	RECEIVING & LOAD PREP	UNPROC	41m	
parts list total quantity between 3 and 7		4		PASS
2017-11-30 18:45	WAREHOUSE HOLD	HOLD	5d20h29m	
duration at least 8h		5d20h29m		PASS
2017-12-06 15:14	PRECONDITIONING	P3LANE12	1d2h16m	
duration between 1d and 3d		1d2h16m		PASS
combi verifier				INFO
- from start of PRECONDITIONING				INFO
- until end of PRECONDITIONING				INFO
- Precon Temperatures between 100 and 120 °F				INFO
- Precon RHs between 45 and 75 %RH				INFO
- location		P3LANE12		INFO
- Precon Temperatures (PTSC)		109 .. 110 °F		INFO
- Precon Temperatures (PTSM)		109 .. 111 °F		INFO
- Precon RHs (PRHSC)		59 .. 67 %RH		INFO
- Precon RHs (PRHSM)		59 .. 66 %RH		INFO
- minimum 1d in spec		1d2h16m		PASS
- maximum 30m out of spec (consecutive)		0m		PASS
Precon Temperatures (PTSC) at most 120 °F		109 .. 110 °F		PASS
Precon Temperatures (PTSM) at most 120 °F		109 .. 111 °F		PASS
Precon RHs (PRHSC) at most 75 %RH		59 .. 67 %RH		PASS
Precon RHs (PRHSM) at most 75 %RH		59 .. 66 %RH		PASS
2017-12-07 17:30	TRANSFER TO CHAMBER	WORKA1SL	10m	
duration at most 1h		10m		PASS
2017-12-07 17:39	CHAMBER	CHAMBER3	8h53m	
cycle check value must be 6115		6115		PASS

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 19 of 37



combi verifier			INFO
- from start of CHAMBER/VACUUM A			INFO
- until end of CHAMBER/RELEASE			INFO
- Sterilizer Temperatures between 110 and 130 °F			INFO
- location	CHAMBER3		INFO
- Sterilizer Temperatures (STSC)	118 .. 121 °F		INFO
- Sterilizer Temperatures (STSM)	119 .. 122 °F		INFO
- total time out of spec must be 0	0m		PASS
CURRENT TANK: LOT # UTLX902336L17E001930, TARE 305.0 LB			
2017-12-07 17:43	VACUUM A	CHAMBER3	27m
Sterilizer Ctl Pressure (SCP) between 1.5 and 2.5 inHg at the end of this phase		2.0 inHg	PASS
2017-12-07 18:11	LEAK TEST	CHAMBER3	5m
Sterilizer Ctl Pressure (SCP) must not increase faster than 0.2 inHg per 5m (between begin and end)		0.1 inHg/5m	PASS
2017-12-07 18:16	NITROGEN DILUTION	CHAMBER3	46m
2017-12-07 18:16	NITROGEN	CHAMBER3	21m
Sterilizer Ctl Pressure (SCP) between 23.5 and 24.5 inHg at the end of this phase		24.0 inHg	PASS
2017-12-07 18:37	EVACUATION	CHAMBER3	25m
Sterilizer Ctl Pressure (SCP) between 1.5 and 2.5 inHg at the end of this phase		2.0 inHg	PASS
2017-12-07 19:02	HUMIDIFICATION (PRESSURE)	CHAMBER3	1m
Sterilizer Ctl Pressure (SCP) increase between 1.0 and 2.0 inHg		1.5 inHg	PASS
2017-12-07 19:03	HUMIDITY DWELL (PRESS)	CHAMBER3	1h
duration between 50m and 1h15m		1h	PASS
Sterilizer Ctl Pressure (SCP) between 3.0 and 4.0 inHg		3.2 .. 3.8 inHg	PASS
2017-12-07 20:03	GAS A (EO)	CHAMBER3	12m
duration between 5m and 55m		12m	PASS
2017-12-07 20:03	STERILANT	CHAMBER3	12m
Sterilizer Ctl Pressure (SCP) between 14.0 and 15.0 inHg at the end of this phase		14.5 inHg	PASS
2017-12-07 20:15	GAS DWELL (EO)	CHAMBER3	4h1m
duration between 4h and 4h30m		4h1m	PASS
Sterilizer Temperatures (STSC) between 115 and 125 °F		120 .. 121 °F	PASS
Sterilizer Temperatures (STSM) between 115 and 125 °F		121 .. 122 °F	PASS
Sterilizer Ctl Pressure (SCP) between 14.0 and 15.0 inHg		14.2 .. 14.7 inHg	PASS
2017-12-08 00:16	AFTER VACUUM	CHAMBER3	20m
Sterilizer Ctl Pressure (SCP) between 1.5 and 2.5 inHg at the end of this phase		2.0 inHg	PASS
duration between 13m and 1h3m		20m	PASS
2017-12-08 00:36	GAS WASH A	CHAMBER3	1h25m
duration between 1h8m and 1h43m		1h25m	PASS
2017-12-08 00:36	RELEASE	CHAMBER3	19m
Sterilizer Ctl Pressure (SCP) between 21.5 and 22.5 inHg at the end of this phase		22.0 inHg	PASS
2017-12-08 00:55	EVACUATION	CHAMBER3	23m
Sterilizer Ctl Pressure (SCP) between 1.5 and 2.5 inHg at the end of this phase		2.0 inHg	PASS

WO2098746

page 2 of 3

printed on 2017-12-09 20:50:36

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .


Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 20 of 37

VeriCycle™ <small>Automated Process Verification</small>		Sterigenics	
2017-12-08 01:19	RELEASE	CHAMBER3	19m
Sterilizer Ctl Pressure (SCP) between 21.5 and 22.5 inHg at the end of this phase		22.0 inHg	PASS
2017-12-08 01:38	EVACUATION	CHAMBER3	23m
Sterilizer Ctl Pressure (SCP) between 1.5 and 2.5 inHg at the end of this phase		2.0 inHg	PASS
2017-12-08 02:01	RELEASE	CHAMBER3	23m
Sterilizer Ctl Pressure (SCP) between 24.7 and 25.8 inHg at the end of this phase		24.8 inHg	PASS
2017-12-08 02:24	CYCLE COMPLETE	CHAMBER3	9m
Sterilant Usage Total (SUT) at least 0 lb 5m before the start of this phase		55 lb	PASS
2017-12-08 02:33	TRANSFER TO AERATION	WORKAUSL	10m
2017-12-08 02:42	AERATION	AERATION9	1d17h30m
duration between 1d and 3d		1d17h30m	PASS
combi verifier			INFO
- from start of AERATION			INFO
- until end of AERATION			INFO
- Aeration Temperatures between 100 and 120 °F			INFO
- Aeration Circulation between 85 and 115			INFO
- location		AERATION9	INFO
- Aeration Temperatures (ATSC)		98 .. 112 °F	INFO
- Aeration Temperatures (ATSM)		98 .. 112 °F	INFO
- Aeration Circulation (AC1)		96 .. 102	INFO
- minimum 1d in spec		1d17h27m	PASS
- maximum 30m out of spec (consecutive)		3m	PASS
Aeration Temperatures (ATSC) at most 120 °F		98 .. 112 °F	PASS
Aeration Temperatures (ATSM) at most 120 °F		98 .. 112 °F	PASS
2017-12-09 20:12	TRANSFER TO WAREHOUSE	PROC	36m
2017-12-09 20:48	REVIEW	PROC	29s
8 HR Hold (Nov. - Apr.) (YES)		YES	PASS
Prec. Room RH Specs. Met (YES)		YES	PASS
Precon. Room Temp. Specs. Met (YES)		YES	PASS
Aeration Temperature Specs Met (YES)		YES	PASS
Remove BI's Post Chamber (19 EA)		19	PASS
Place BI's Per Specification (19 EA)		19	PASS
RESULT OF VERIFICATION		PASS	

It is customers responsibility to ensure acceptance of any result other than PASS noted on this report prior to release of the load.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 21 of 37

Confidential - Controlled copy on date printed



Nelson Laboratories, LLC | Because every day matters


CERTIFICATE OF CONFORMANCE


Process Challenge Device (PCD)

Nelson Laboratories, LLC (NL) hereby certifies that the PCDs prepared conform to the requirements and quality specifications outlined in NL standard operating procedure (SOP) SOP0180.

Scope of Certification:

Product:	NLI-PCD-010
Batch Number:	BO 1815
Preparation Date:	13 Nov 2017
BI Manufacturer's Information:	STERIS
	Part #NA005, Lot #2177A
	Exp. 05 Feb 2019
BI Type:	Spore Strip
BI Manufacturer's Spore Population:	2.5×10^6 CFU/BI
NL Verified Population:	3.1×10^6 CFU/BI
Organism:	<i>Bacillus atrophaeus</i>





Sammy Diphibane electronically approved

Released by

Sammy Diphibane

01 Dec 2017 11:18

Date and Time

A PCD acts as a surrogate for the actual product and the appropriateness of a given PCD should be determined through comparative resistance testing. This testing is performed to provide biological monitors for development of a sterilization cycle process. The validated D value for the PCD is available upon request. To confirm authenticity of this certificate, please contact NL. The most recent signature reflects the latest version of the certificate.

P.O. Box 571830 | Murray, UT 84157-1830 U.S.A. • 6280 South Redwood Road | Salt Lake City, UT 84123-8600 U.S.A.

Rev. 1.3.3

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 22 of 37



CRM43650 Rcvd: 18 Oct 2017
Ambient Use By: 05 Feb 2019



Opened: _____

STERIS Corporation Biological Operations
Bacillus atrophaeus BIOLOGICAL INDICATOR CERTIFICATE OF PERFORMANCE
Exhibit K

CATALOG NUMBER: NA005

PRODUCT NAME: Spordex Strip

PRODUCT LOT NUMBER: C17020

SUBASSEMBLY LOT NUMBER: 2177A

EXPIRATION DATE: February 5, 2019

Bacillus atrophaeus NRRL B4418:

Mean Population Recovery* (CFU): 2.5×10^6

Survival Time (Min.): 11.9

D₁₀ Value** (Min. ***): 2.7

Kill Time (Min.): 28.0

D₉₉ Value** (Min. ***): 1.6

Survival Time (Min.): 7.1

Kill Time (Min.): 16.6

* Colony forming units determined after a preliminary heat treatment.

** 600±30mg/L EO using 100% EO, 54±1°C, 60±10% RH or 160±1°C Dry Heat

*** Determined at time of manufacture by fraction negative procedure after graded exposure to sterilization conditions.

Incubate at 30° to 35°C for seven days

Store at 2° to 24°C and between 30% to 80% RH. Do not use after the indicated expiration date.

Dispose of as you would any other microbiological waste (121°C for a minimum of 30 minutes).

This document certifies that the biological indicator product listed above meets STERIS' quality assurance specifications and the performance criteria suggested by the current revision of the United States Pharmacopeia. This certifies that the product listed above has been tested in compliance with current versions of, ISO 11138-1, ISO 11138-2, and USP.

Quality Systems Representative: _____

Date: 9/25/17

Reviewed By: _____

Date: 9/25/17

LIMITATION OF LIABILITY AND INDEMNITY

Nothing in this Certificate of Performance shall, or is intended to, alter, expand, or diminish the terms and conditions of sale governing your purchase of the Biological Indicator Product from STERIS. In no event, whether as a result of breach of warranty, or tort (including negligence and strict liability) shall STERIS or its suppliers be liable as a result of any statement or information contained in this Certificate of Performance. In addition, STERIS shall not be liable for any consequential or incidental damages including, without limitation, loss of use or damage to your products or equipment, cost of substitute products, or down time costs, allegedly caused by the Biological Indicators Products. The responsibility of STERIS for

PC60003S Revision Level: AN ECN# 76462
Effectivity Date: 09/22/2017

Page 1 of 1

STERIS Corporation • 5960 Heisley Road • Mentor, OH 44060-1834 USA • 440-354-2600

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 23 of 37



Sponsor:
Tim Kessinger
MegaDyne Medical Products, Inc.
11506 S. State St.
Draper UT 84020

Biological Indicator (BI) Sterility Test Final Report

Study Number: 1008255-S01
Test Article: Sterigenics Run #2098746
Purchase Order: 29592
Study Received Date: 10 Dec 2017
Testing Facility: Nelson Laboratories, LLC, a Business Unit of Sterigenics International
6280 S. Redwood Rd.
Salt Lake City, UT 84123 U.S.A.
Test Procedure(s): Standard Test Protocol (STP) Number: STP0079 Rev 16
Deviation(s): None

Summary: This BI sterility test was conducted to verify if any viable organisms remained on BIs used to monitor routine and/or validation cycles. All test method acceptance criteria were met.

Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.

BI Information:

Type	Manufacturer	Lot #	Expiration Date	Species Name	Minimum Incubation Time	Incubation Temperature
Spore Strip	STERIS	2177A	05 Feb 2019	<i>Bacillus atrophaeus</i>	7 days	30-35°C

Results:

Type	Number Tested	Number Positive	Number Negative
NLI-PCD-010	19	0	19
Environmental Control	1	0	1
Media Negative Control	1	0	1
Positive Control	1	1	0



Derek L. Miller electronically approved

Study Director

Derek L. Miller

17 Dec 2017 03:46 PM

Study Completion Date and Time

Page 1 of 1

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 24 of 37

9. APPENDIX II - PRECONDITIONING AND SHIP TESTING

Megadyne Medical Products, Inc.	TEST PROTOCOL	<u>Document Number</u> XENG-PRT-441
	ZIP ACE Modified Product 6-Pack Ship Test Protocol	Revision: A
		Page 11 of 17


10. APPENDIX I - THERMAL CYCLE DATA

Appendix I: Thermal Cycle Data

Maximum Temperature (°C):	55.7C
Minimum Temperature (°C):	-40.1C
Maximum Temperature (%RH):	47.6%
Minimum Temperature (%RH):	43.9%
Chamber conditions held @ -40°C and no humidity control for a duration of 4 hours:	6:45 to 10:45 12/28/17 yes on 12-28-2017
Chamber conditions held @ 55°C and 95%RH for a duration of 4 hours:	12:45 to 16:45 12/28/17 yes on 12-28-2017
Chamber conditions held @ 55°C and 15%RH for a duration of 4 hours:	17:45 to 21:45 on 12/28/17 yes on 12-28-2017
Chamber conditions held @ 23°C and 50%RH for a duration of 72 hours:	22:30 on 12/28/17 to 23:59 12/31/17 yes, from 12-28-2017 to 12-31-2017.
Paul Valpreda Test Technician Name	Paul Valpreda Signature
	1-12-2018 Date
Tyler Skinner Engineer Name	Tyler Skinner Signature
	2-1-2018 Date
42012 Thermotron SN	5-31-2018 Calibration Due Date

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 25 of 37

Graph



Raw Data

Date	Time	Air Temp	Humidity
12/28/2017	4:00:03	55	4.5
	4:30:03	55.1	4.5
	5:00:03	55	4.5
	5:30:03	55	4.5
	6:00:03	55	4.5
	6:30:03	-4.8	3.9
	7:00:03	-40	19.1
	7:30:03	-40	34.9
	8:00:03	-40	38.9
	8:30:03	-40	40.6
	9:00:03	-40	41.9
	9:30:03	-40	43
	10:00:03	-40.1	43.9
	10:30:03	-40	44.4
	11:00:03	23.2	26.6
	11:30:03	23	23.8
	12:00:03	55.7	47.6
	12:30:03	55.2	94.1
	13:00:03	55.1	94.4

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 26 of 37

	13:30:03	55	94.5
	14:00:03	55.1	94.7
	14:30:03	55.1	94.6
	15:00:03	55.1	94.8
	15:30:03	55.1	94.7
	16:00:03	55	94.8
	16:30:03	55	57.1
	17:00:03	55	23.6
	17:30:03	55.1	17.6
	18:00:03	55.1	15
	18:30:03	55.1	15.1
	19:00:03	55.1	14.4
	19:30:03	55	14.8
	20:00:03	55	15.1
	20:30:03	55	14.9
	21:00:03	55	15
	21:30:03	55	14.9
	22:00:03	22.3	56.5
	22:30:03	23	50
	23:00:03	23	50
	23:30:03	23	48.7
12/29/2017	0:00:03	23	49.4
	0:30:03	23	49.9
	1:00:03	23	50.2
	1:30:03	23	50
	2:00:03	23	49.9
	2:30:03	23	50.2
	3:00:03	22.9	50.3
	3:30:03	23.1	50.3
	4:00:03	23	50.1
	4:30:03	23	49.9
	5:00:03	23	49.9
	5:30:03	23	50.1
	6:00:03	23	50.1
	6:30:03	23	50
	7:00:03	23	50
	7:30:03	23	50
	8:00:03	23	50

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 27 of 37

	8:30:03	23.1	49.9
	9:00:03	23	49.9
	9:30:03	23	50.1
	10:00:03	23.1	50
	10:30:03	23	50
	11:00:03	23	49.8
	11:30:03	23	50.1
	12:00:03	23	50
	12:30:03	23	50
	13:00:03	23	50
	13:30:03	23	50.1
	14:00:03	23	50.1
	14:30:03	23	50
	15:00:03	22.9	50
	15:30:03	23	50
	16:00:03	23	49.8
	16:30:03	23	49.9
	17:00:03	23	50.1
	17:30:03	23	49.9
	18:00:03	23	49.8
	18:30:03	22.9	49.8
	19:00:03	23	49.7
	19:30:03	23	49.9
	20:00:03	23	50
	20:30:03	23	49.9
	21:00:03	23	50.1
	21:30:03	23	49.9
	22:00:03	22.9	50.3
	22:30:03	23	50.1
	23:00:03	23	50.1
	23:30:03	23	49.9
12/30/2017	0:00:03	23	49.9
	0:30:03	23	50
	1:00:03	23	49.9
	1:30:03	23	49.9
	2:00:03	23	50.1
	2:30:03	23	50
	3:00:03	23	50.2

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 28 of 37

	3:30:03	23	50
	4:00:03	23	50
	4:30:03	23	49.7
	5:00:03	23	49.9
	5:30:03	23	50.3
	6:00:03	23	50.1
	6:30:03	22.9	50.4
	7:00:03	23	49.9
	7:30:03	23.1	49.8
	8:00:03	23	49.9
	8:30:03	23	50
	9:00:03	23	50
	9:30:03	23	50.1
	10:00:03	23	50
	10:30:03	22.9	49.9
	11:00:03	23	49.9
	11:30:03	23	50.1
	12:00:03	23	50.4
	12:30:03	23	49.9
	13:00:03	23	49.9
	13:30:03	23	50.1
	14:00:03	23	49.8
	14:30:03	23	49.9
	15:00:03	23	50.1
	15:30:03	23	49.9
	16:00:03	23	50
	16:30:03	23	50.1
	17:00:03	23	50
	17:30:03	23	49.9
	18:00:03	23	49.8
	18:30:03	22.9	50.1
	19:00:03	23	50.1
	19:30:03	23	49.9
	20:00:03	23	49.8
	20:30:03	23	50.2
	21:00:03	23.1	50.1
	21:30:03	23	49.8
	22:00:03	23	49.9

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 29 of 37

	22:30:03	23	49.7
	23:00:03	23	50
	23:30:03	23	49.9
12/31/2017	0:00:03	23	50.3
	0:30:03	23	50
	1:00:03	23	50
	1:30:03	23	50
	2:00:03	23	49.8
	2:30:03	23	50.2
	3:00:03	23	49.9
	3:30:03	23	50.3
	4:00:03	23	50.3
	4:30:03	23	49.7
	5:00:03	23	50
	5:30:03	23	50.1
	6:00:03	23	49.9
	6:30:03	23.1	50
	7:00:03	23	49.7
	7:30:03	23	49.9
	8:00:03	23	50.1
	8:30:03	23	49.9
	9:00:03	23	49.8
	9:30:03	23	50.1
	10:00:03	23	49.9
	10:30:03	23	50.3
	11:00:03	23	49.9
	11:30:03	23	49.8
	12:00:03	23	50.2
	12:30:03	23	50.3
	13:00:03	23	50
	13:30:03	23	50
	14:00:03	23	49.6
	14:30:03	23	49.9
	15:00:03	23	50
	15:30:03	23	50.3
	16:00:03	23	49.9
	16:30:03	23	50
	17:00:03	23	49.8

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 30 of 37

	17:30:03	23	49.9
	18:00:03	23	50
	18:30:03	23	50.2
	19:00:03	23	50.3
	19:30:03	23	50.1
	20:00:03	23	50
	20:30:03	22.9	50.2
	21:00:03	22.9	49.9
	21:30:03	23	49.9
	22:00:03	23.1	49.7
	22:30:03	23.5	52.2
	23:00:03	23.4	54.3
	23:30:03	23.2	55.6
1/1/2018	0:00:03	23	56.8

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 31 of 37

Megadyne Medical Products, Inc.	TEST PROTOCOL	<u>Document Number</u> XENG-PRT-441
	ZIP ACE Modified Product 6-Pack Ship Test Protocol	Revision: A
		Page 12 of 17

11. APPENDIX II – SHIPPING TEST

Preconditioning:

Start Date: 12-28-2018 Chamber Number: 01268
Completion Date: 12-31-2018 Last Calibration: 5-23-2017
Signature/Date: Paul Valpreda Calibration due: 5-31-2018
1-9-2018

Drop Test:

Catalog Number: ME725MIC Weight: 3.5 lbs. Drop Height: 15"

Drop Sequence	Orientation	Specific face, edge or corner	Initials/Date
1	Top	Face 1	PV 1-9-18
2	Edge	Edge 5-3	PV 1-9-18
3	Edge	Edge 6-3	PV 1-9-18
4	Corner	Corner 2-3-5	PV 1-9-18
5	Corner	Corner 4-3-6	PV 1-9-18
6	Bottom	Face 3	PV 1-9-18

Comments:

Signature: Paul Valpreda Date: 1-9-2018

Compression Test:

Catalog Number: ME725MIC Pounds Force: 200

Comments: Passed.

Signature: Paul Valpreda Date: 1-9-2018

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2018, 10:41:20 am; Printed by: TSKINNER.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 32 of 37

Megadyne Medical Products, Inc.	TEST PROTOCOL	<u>Document Number</u> XENG-PRT-441
	ZIP ACE Modified Product 6-Pack Ship Test Protocol	Revision: A
		Page 13 of 17

Appendix II Continued
Shipping Test Log Sheet

Vibration:

Low Frequency, 40 minutes, Initials: PV

High frequency 10 minutes, Initials: PV

Completion Date: 1-9-2018

Signature: Paul Valpreda Date: 1-9-2018

Concentrated Impact Test:

Completion Date: 1-9-2018

Signature: Paul Valpreda Date: 1-9-2018

Second Drop Test:

Catalog Number: ME725M1C Weight: 3.5 lbs. Drop Height: 15" + 30"

Drop Sequence	Orientation	Specific face, edge or corner	Initials/Date
1	Edge	Edge 4-6	PV 1-9-18
2	Face	Face 4	PV 1-9-18
3	Face	Face 6	PV 1-9-18
4	Corner	Corner 2-1-5	PV 1-9-18
5	Edge	Edge 2-1	PV 1-9-18
6	Bottom	Face 3, Increase height to 30 inches.	PV 1-9-18

Comments:

Signature: Paul Valpreda Date: 1-9-2018

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2018, 10:41:20 am; Printed by: TSKINNER.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 33 of 37

APPENDIX III – BUBBLE LEAK TESTING

Bubble Leak Test Log Sheet For Megadyne ACE Blade 700			
Catalog - ME725M1C			
LOT - S170323			
Sample	Pass	Fail	Comments
1	X		
2	X		
3	X		
4	X		
5	X		
6	X		
7	X		
8	X		
9	X		
10	X		
11	X		
12	X		
13	X		
14	X		
15	X		
16	X		
17	X		
18	X		
19	X		
20	X		
21	X		
22	X		
23	X		
24	X		
25	X		
26	X		
27	X		
28	X		
29	X		
30	X		
31	X		
32	X		
33	X		
34	X		
35	X		

No failures were observed - PV

Paul Valpreda 1/10/2018
OPERATOR NAME DATE

Paul Valpreda 1-10-2018
OPERATOR SIGNATURE DATE

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 34 of 37

APPENDIX IV – DYE TESTING

Dye Penetration Evaluation			
Megadyne ACE Blade 700 ME725M1C LOT S170323			
Sample #	PASS	FAIL	Comments
1	X		
2	X		
3	X		
4	X		
5	X		
6	X		
7	X		
8	X		
9	X		
10	X		
11	X		
12	X		
13	X		
14	X		
15	X		
16	X		
17	X		
18	X		
19	X		
20	X		
21	X		
22	X		
23	X		
24	X		
25	X		
26	X		
27	X		
28	X		
29	X		
30	X		
31	X		
32	X		
33	X		
34	X		
35	X		

No failures were observed - PV

Paul Valpreda 1/11/2018
Operator Name Date

Paul Valpreda 1-11-2018
Operator Signature Date

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 36 of 37

APPENDIX VI – MINIMUM SEAL WIDTH TESTING

Package Seal Width Evaluation Data Collection Form

MINIMUM SEAL WIDTH EVALUATION						
Sample	Part/Lot	Cavity	Front 1	Back 3	Right 2	Left 4
1	ME725M1C LOT - S170323	2			0.32	
2	ME725M1C LOT - S170323	1			0.32	
3	ME725M1C LOT - S170323	1			0.32	
4	ME725M1C LOT - S170323	1			0.32	
5	ME725M1C LOT - S170323	1			0.31	
6	ME725M1C LOT - S170323	1			0.32	
7	ME725M1C LOT - S170323	1			0.32	
8	ME725M1C LOT - S170323	1			0.33	
9	ME725M1C LOT - S170323	2			0.30	
10	ME725M1C LOT - S170323	1			0.31	
11	ME725M1C LOT - S170323	2			0.32	
12	ME725M1C LOT - S170323	2			0.31	
13	ME725M1C LOT - S170323	2			0.32	
14	ME725M1C LOT - S170323	1			0.31	
15	ME725M1C LOT - S170323	1			0.30	
16	ME725M1C LOT - S170323	2			0.31	
17	ME725M1C LOT - S170323	2			0.32	
18	ME725M1C LOT - S170323	2			0.32	
19	ME725M1C LOT - S170323	1			0.32	
20	ME725M1C LOT - S170323	2			0.31	
21	ME725M1C LOT - S170323	1			0.31	
22	ME725M1C LOT - S170323	1			0.32	
23	ME725M1C LOT - S170323	1			0.31	
24	ME725M1C LOT - S170323	1			0.27	
25	ME725M1C LOT - S170323	2			0.21	
26	ME725M1C LOT - S170323	1			0.32	
27	ME725M1C LOT - S170323	2			0.32	
28	ME725M1C LOT - S170323	2			0.31	
29	ME725M1C LOT - S170323	1		0.31		
30	ME725M1C LOT - S170323	2	0.34			
31	ME725M1C LOT - S170323	1			0.32	
32	ME725M1C LOT - S170323	2			0.31	
33	ME725M1C LOT - S170323	2			0.31	
34	ME725M1C LOT - S170323	2		0.26		
35	ME725M1C LOT - S170323	1			0.31	

 Burst Side

CALIBRATION INFORMATION

Calipers
Starrett MMP-1003

Serial Number: 723
Megadyne Number: 01039
Calibration Date: 5/18/2016
Calibration Due: 5/31/2018

Test-A-Pack Burst Test Equipment
Megadyne #01397
SN 2098
Cal date - 8-3-17
Cal due - 8-3-18

Program:
Burst Test
Flow = 9
Sensit = 1
Prefill = Y
in H2O

Test Fixture - Large Needle
Part #F100-1320-2
SN 478

Paul Valpreda
Operator Name
1/22/2018
Date


Operator Signature
1-22-2018
Date

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-535
	ZIP ACE Modified, 6-Pack Ship Test Report, T=0	Revision: 001
		Page 37 of 37

APPENDIX VII – PRODUCT DAMAGE INSPECTION

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-441 Rev: A Effective: 20 Sep 2017

Megadyne Medical Products, Inc.	TEST PROTOCOL	<u>Document Number</u> XENG-PRT-441
	ZIP ACE Modified Product 6-Pack Ship Test Protocol	Revision: A
		Page 17 of 17

15. APPENDIX VI – PRODUCT DAMAGE INSPECTION

Inspect the product per the protocol and enter the number of units that pass or fail in the box below.

Catalog #	Pass	Fail
Damage	35	0

Comments:

No failures were observed. No coating failures were found.

Signature: Paul Valpreda

Date: 1-9-2018

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 10 Jan 2018, 09:28:54 am; Printed by: TSKINNER.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 11:12:08 pm; Printed by: .