



The Electrosurgical Authority®

DOCUMENT NUMBER: ENG-RPT-418

DOCUMENT TITLE: Zip Pen Product Specification Verification Test Report

DOCUMENT NOTES:

Document Information

Revision: 002

Vault: MEG-Rel

Status: Release

Document Type: ENG-RPT

Date Information

Effective Date: 21 Feb 2018

Expiration Date:

Release Date: 21 Feb 2018

Next Review Date:

Control Information

Author: MGLASSETT

Owner: ENG-APPR

Previous Number:

Change Number: 2018-ENG-DCO-009

Signature Manifest**Document Number:** ENG-RPT-418**Revision:** 002**Title:** Zip Pen Product Specification Verification Test Report

All dates and times are in Mountain Standard Time.

ENG-RPT-418 Zip Product Verif Rpt**Change Request**

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		10 Jan 2018, 11:41:39 AM	Approved

Collaboration

Name/Signature	Title	Date	Meaning/Reason
Mark Glassett (MGLASSETT)		25 Jan 2018, 03:36:23 PM	Complete
Dave Shimkus (DSHIMKUS)		05 Feb 2018, 03:40:24 PM	Complete
Joni Stegeman (JSTEGEMAN)	Ethicon Quality	14 Feb 2018, 04:45:02 PM	Complete
Darlene Hull (DHULL)	Regulatory	15 Feb 2018, 11:51:21 AM	Complete
Mallory Schroeder (MSCHROEDER)	Engineer	16 Feb 2018, 09:10:46 AM	Complete

Document Review

Name/Signature	Title	Date	Meaning/Reason
Mark Glassett (MGLASSETT)		19 Feb 2018, 03:30:01 PM	Complete
Lucy Richards (LRICHARDS)		20 Feb 2018, 10:37:42 AM	Complete

RA-Approval

Name/Signature	Title	Date	Meaning/Reason
Darlene Hull (DHULL)	Regulatory	20 Feb 2018, 12:26:50 PM	Approved

QA-Approval

Name/Signature	Title	Date	Meaning/Reason
Joni Stegeman (JSTEGEMAN)	Ethicon Quality	21 Feb 2018, 08:47:02 AM	Approved

ENG-Approval

Name/Signature	Title	Date	Meaning/Reason
Paul Borgmeier (PBORGMEIER)		20 Feb 2018, 11:08:18 AM	Approved

Training Review

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		21 Feb 2018, 10:41:43 AM	Approved

Final Release

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		21 Feb 2018, 10:41:56 AM	Approved

Quick Approval

Approve Now

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		21 Feb 2018, 10:46:05 AM	Approved

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 1 of 24

Authored By: Mark Glassett

Revised By: Mallory Schroeder (M.S.)

TABLE OF CONTENTS

Page

1.	ABSTRACT	1
2.	OBJECTIVE	1
3.	RESULTS	1
4.	CONCLUSIONS	7
5.	RECOMMENDATIONS	7

1. ABSTRACT

Zip Project Catalog Numbers 2525-10, 2525-15, 2540, 2560, 2211 and 2220 were evaluated for conformance to requirements of the product specification ENG-PS-007 that were not verified under other protocols. The evaluation was done per of ENG-PRT-290. The Zip Pen and other accessories for the Zip Project met the requirements of the protocol for 19 design attributes required by the Product Specification.

2. OBJECTIVE

The objective of this test report is to document conformance to requirements of the product specification ENG-PS-007 for attributes that are not verified in other protocols.

3. RESULTS

3.1. BUTTON SIZE

- 3.1.1. Zip Pen button size measures 0.0699 square inches. The disposable pencil button size is 0.0697 square inches. The surface area of the two buttons is approximately the same (within .0002 inches square) and meets the requirements of the protocol. See calculations in Appendix I.

3.2. BRANDING

- 3.2.1. The Zip Pen Handle is gray pantone 427 with green inly that is pantone 356. The Megadyne name is molded into the handle and colored green

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 2 of 24

pantone 356. This design meets the requirements of the protocol. See verification in Appendix II.

3.3. TIP EXPOSURE

- 3.3.1. The average tip exposure of the Zip Pen is 0.87". The average tip exposure of the UltraVac is 0.93". This comparison shows that the Zip Pen has a .06" shorter exposure than the Ultra Vac on average. Note that a shorter exposure is more desirable than a longer exposure. The comparison data is as follows:

Product	Average Tip Exposure, inches	Standard Deviation	Range
Zip Pen 2525-10	0.87	.007	.864 - .887
Ultra Vac 2110-10	0.93	.086	.820 – 1.012

The requirement in the product specification is that the Zip Paen tip exposure be comparable to the Ultra Vac. Note that the Zip Pen standard deviation is very small indicating that the tip exposure is very consistent from one pencil to the next. The Ultra Vac has a much larger standard deviation (over ten times larger) indicating that it is less consistent from one pencil to the other. The standard deviation of the Ultra Vac is larger than the difference between the two averages. Statistically, the averages are not equal for a confidence interval of 90%. However, the range of values for the Zip Pen is within the range of values for the Ultra Vac. Given this fact, it can be concluded that the tip exposure for the two devices is comparable and meets the protocol. In order to make the requirement clearer for the Zip pen Product Specification, the Product Specification will be revised to be more specific. See data in Appendix III.

3.4. NOZZLE CLARITY

- 3.4.1. The nozzle design is clear in the area of the electrode and meets the requirements of the protocol. See verification in Appendix IV

3.5. CORD CONTAINMENT

- 3.5.1. The cord is contained within the tubing for the first 64 inches and meets the requirements of the protocol. See verification in Appendix V.





Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 3 of 24

3.6. TUBING SWIVEL


3.6.1. The tubing has a swivel that is 8 inches from the connection to the Zip Pen and meets the requirements of the protocol. See verification in Appendix VI.

3.7. TUBING ERGONOMICS

3.7.1. The photo's below show that the tubing can be positioned in the manner that is shown in the IFU and meets the requirements of the protocol. See photos below and verification in Appendix VII.

Actual Use	IFU
	<p>Trigger Method</p> 
	<p>Tube Over Hand Method</p> 

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 4 of 24



3.8. BUTTON LOCATION

3.8.1. The yellow button is located nearest to the electrode end of the Zip Pen and the blue button is located farthest from the electrode end of the Zip Pen. This meets the requirements of the protocol. See verification in Appendix VIII.

3.9. ULPA FILTER EFFICIENCY

3.9.1. A review of the specification sheet and drawing for the ULPA filter media shows that the filter efficiency exceeds the required efficiency of 99.999% for 0.1 to 0.2 micron particles and meets the requirements of the protocol. See verification in Appendix IX. Note that there was a typographical error on the particle size for the efficiency in the protocol. The efficiency should be 99.999% at 0.1 – 0.2 micron particle size. The protocol listed 0.01 – 0.02 and will be corrected at DCO.

3.10. ULPA FILTER FLUID TRAP

3.10.1. A review of the drawing and of an actual product shows that the ULPA filter has a fluid trap and meets the requirements of the protocol. See verification in Appendix X.

3.11. ULPA FILTER CONNECTOR

3.11.1. A review of the drawing and of an actual product shows that the ULPA filter has a threaded connector that mates with the Zip Pen and meets the requirements of the protocol. See verification in Appendix XI.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 5 of 24

3.12. CARBON FILTER CHARCOAL

3.12.1. A review of the Carbon filter data sheet and the actual product shows that it contains Charcoal and meets the requirements of the protocol. See verification in Appendix XII.

3.13. ZIP ELECTRODE INSERTION

3.13.1. Installation trials of 11 electrodes into 11 Zip Pens demonstrated that the electrode cannot be installed incorrectly. The installed electrodes met the maximum resistance requirement of less than 50 ohms. The Zip Pen meets the requirements of the protocol. See verification in Appendix XIII.

3.14. ZIP PEN PACKAGING

3.14.1. The Zip Pen shipping box is an RSC design. The drawing shows that the shipping box will contain one IFU per box and have 20 units in each box. The Zip Pen packaging meets the requirements of the protocol. See verification in Appendix XIV.

3.15. ZIP EXTENSION NOZZLE PACKAGING

3.15.1. The Extension Nozzle shipping box is an RSC design. The drawing shows that the shipping box will contain one IFU per box and have 10 units in each box. The Extension Nozzle packaging meets the requirements of the protocol. See verification in Appendix XV

3.16. ZIP PEN AND EXTENSION NOZZLE UNIT LABELING

3.16.1. The Zip Pen and Extension Nozzle unit labels contain the information required by the protocol. The Zip Pen and Extension Nozzle unit labels meet the requirements of the protocol. See verification in Appendix XVI.

3.17. ZIP PEN AND EXTENSION NOZZLE SHIPPER LABELING

3.17.1. The Zip Pen and Extension Nozzle shipping labels contain the information required by the protocol. The Zip Pen and Extension Nozzle shipping labels meet the requirements of the protocol. See verification in Appendix XVII.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 6 of 24

3.18. ULPA AND CARBON FILTER UNIT LABELING

3.18.1. The ULPA and Carbon Filter unit labels contain the information required by the protocol. The ULPA and Carbon Filter unit labels meet the requirements of the protocol. See verification in Appendix XVIII.

3.19. ULPA AND CARBON FILTER SHIPPER LABELING

3.19.1. The ULPA and Carbon Filter shipping labels contain the information required by the protocol with one exception. The ULPA Filter shipping label does not have the product trade name on the label. The label has been approved through the DCO process by management and is what Marketing requested. The product specification and protocol will be revised at DCO. The labels meet the requirements of the protocol with exception of one item that will be revised. See verification in Appendix XIX.

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 7 of 24

4. CONCLUSIONS

The products in the Zip Project were verified to meet the requirements of the protocol which shows that the products meet the Product Specification with one minor exception. For that one minor exception the Product Specification will be revised.

5. RECOMMENDATIONS

The following documentation for the Zip Project will be updated as a result of this test report: The Input/output conformance test matrix ENG-IOM-012 will be updated to show verification per this test report. The Product Specification ENG-PS-007 will be updated to clarify and correct the requirements identified in this report. The Smoke Evacuation Accessories Risk Analysis ENG-RMF-045 will be updated to show this test report in the appropriate verification areas.

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 8 of 24

Appendix I Button Size

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 19 of 33

Appendix I: BUTTON SIZE

Zip Pen Button Size Calculation from ENG-DWG-588

Width .226 PV 2-6-2015

Length .358 PV 2-6-2015

Calculation $\pi \frac{(.226)^2}{4} + [(.358 - .226) \times .226] = .0699$ MB 2-6-15

Two Dimensional Area: .0699 in²

Disposable Pencil Button Size Calculation

Catalog Number: 0039 Lot Number: 10664

Button Diameter: .298 PV 2-6-2015

Calculation: $\pi \frac{(.298)^2}{4} = .0697$ in² MB 2-6-15

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 9 of 24

Appendix II Branding

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 20 of 33

Appendix II BRANDING

Reference ENG-DWG-650

Verify that the Pen Body has the following:

Body Color is Gray Pantone 427: ☒ YES ☐ NO Initials/Date PV 1-26-2015

TPR inlay is Green Pantone 356: ☒ YES ☐ NO Initials/Date PV 1-26-2015

Megadyne Name: ☒ YES ☐ NO Initials/Date PV 1-26-2015

Confidential - Pre-Production Controlled on date printed

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 10 of 24

Appendix III Tip Exposure

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 21 of 33

Appendix III

TIP EXPOSURE

Zip Pen Catalog Number: 2525-10 Lot Number: 5140303

Ultra Vac Catalog Number: 2110-10 Lot Number: 5408

Zip Pen Sample Number	Measurement	Ultra Vac Sample Number	Measurement
Z1	.867	U1	.996
Z2	.877	U2	.977
Z3	.871	U3	.978
Z4	.868	U4	1.012
Z5	.878	U5	.986
Z6	.864	U6	.992
Z7	.887	U7	.976
Z8	.870	U8	.821
Z9	.869	U9	.821
Z10	.875	U10	.820
Z11	.883	U11	.825

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 11 of 24

Appendix IV: Nozzle Clarity
Appendix V: Cord Containment

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 22 of 33

Appendix IV

NOZZLE CLARITY

Reference Drawing ENG-DWG-594

Verify that the Nozzle meets the following:

Material is clear and not opaque: ☒ YES ☐ NO Initials/Date PV 1-26-2015

Comments: _____

Appendix V

CORD CONTAINMENT

Reference Drawing ENG-DWG-716

Verify that the Assembly meets the following:

Cord is assembled inside the tube for at least the first 64 inches:

☒ YES ☐ NO Initials/Date PV 1-26-2015

Comments: _____

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 12 of 24

Appendix VI Tubing Swivel

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 23 of 33

Appendix VI

TUBING SWIVEL

Reference Drawing ENG-DWG-716

Verify that the Assembly meets the following:

The assembly has a swivel component approximately 8 inches from the handle connection:

☒ YES ☐ NO

Initials/Date PV 1-26-2015

Comments: _____

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 13 of 24

Appendix VII: Tubing Ergonomics

Appendix VIII: Button Location

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 24 of 33

Appendix VII

TUBING ERGONOMICS

Reference IFU MKT-LBL-531

Verify that the Assembly meets the following:

The Tube/Cord can be configured as shown for the three methods in the IFU. Photograph the three configurations for the test report:

☒ YES ☐ NO Initials/Date PV 2-4-2015

Comments: _____

Appendix VIII

BUTTON LOCATION

Reference Drawing ENG-DWG-716

Verify that the Assembly meets the following:

The Yellow CUT Button is nearest to the electrode and the Blue COAG Button is farthest from the electrode:

☒ YES ☐ NO Initials/Date PV 1-27-2015

Comments: _____

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 14 of 24

Appendix IX: Filter Efficiency
Appendix X: ULPA Filter Fluid Trap

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 25 of 33

Appendix IX

ULPA FILTER EFFICIENCY

Reference Assembly Drawing ENG-DWG-1017 and
Component Drawing ENG-DWG-1015

Verify that the Assembly meets the following:

The assembly has a minimum filtration efficiency of 99.999% for ^{0.1 0.2}~~0.1 0.2~~ micron particles:

☒ YES ☐ NO Initials/Date PV 1-29-2015

Comments: _____

Appendix X

ULPA FILTER FLUID TRAP

Reference Assembly Drawing ENG-DWG-1017

Verify that the Assembly meets the following:

The assembly has a Fluid Trap:

☒ YES ☐ NO Initials/Date PV 1-28-2015

Comments: _____

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 15 of 24

Appendix XI: ULPA Filter Connector
Appendix XII: Charcoal Carbon Filter

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 26 of 33

Appendix XI

ULPA FILTER CONNECTOR

Reference Component Drawing ENG-DWG-1017

Verify that the Assembly meets the following:

The assembly has a threaded connector that is compatible with the Zip Pen Proximal Adapter ENG-DWG-1160:

☒ YES/NO

Initials/Date PV 1-28-2015

Comments: _____

Appendix XII

CHARCOAL CARBON FILTER

Obtain a sample of Carbon Filter 2220

Verify that the Filter meets the following:

The Filter connects to the Mega Vac and/or Mega Vac Plus

The filter contains charcoal

☒ YES/NO

Initials/Date PV 1-29-2015

Comments: _____

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Confidential - Controlled copy on date printed

Confidential - Pre-Production Controlled on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 16 of 24

Appendix XIII Zip Electrode Insertion

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 27 of 33

Appendix XIII

ZIP ELECTRODE INSERTION

Zip Pen Catalog Number: 2525-10 Lot Number: S140303

Electrode Catalog Number: 0014A Lot Number: 150088

Zip Pen Sample Number	Insertion Correct Yes/No	Continuity (Ohms)	Pass/Fail
C1	yes	0.8	PASS
C2	yes	0.8	PASS
C3	yes	0.7	PASS
C4	yes	1.7	PASS
C5	yes	1.2	PASS
C6	yes	0.7	PASS
C7	yes	1.0	PASS
C8	yes	0.7	PASS
C9	yes	9.2	PASS
C10	yes	0.9	PASS
C11	yes	0.7	PASS

Equipment Information	
Equipment	Fluke Multimeter, True RMS
Cert Number	497572
Asset Number	01372
Cal Date	9/12/2014
Cal Due	9/30/2015

M.S. 1/10/18 – Adding missing equipment & calibration information. Calibration record attached in Appendix XIII

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Confidential - Controlled copy on date printed

Confidential - Pre-Production Controlled on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 17 of 24

WESTERN STATES
CALIBRATION
The Standard Measure in all things

Certificate Number:
497572

Asset Number:
01372

CERTIFICATE OF CALIBRATION

Western States Calibration certifies this instrument has been cleaned, calibrated and inspected in accordance with said Instrument Calibration Procedure. This calibration was performed in accordance with requirements of ISO/IEC 17025, ISO 9001:2008, and ANSI/NCCL Z540-1 with measuring standards traceable to the National Institute of Standards and Technology. Units of measurement are stated according to the International System of Units (SI). The results reported on this certificate apply only to the item calibrated. All data is reported as raw data and uncorrected for uncertainty or environmental effects. It is the end users responsibility to determine fitness for use. Measurement uncertainty is not taken into account when determining In/Out of tolerance conditions. This certificate shall not be reproduced except in full, without the written approval of Western States Calibration. Measurement Uncertainty (EMU) is reported per measurement, when available, @ $k = 2$.

Attention: Alan Holt
Megadyne Medical Products
11506 South State Street

Draper, UT 84020

PO Number: 23845

Asset Number: 01372
Manufacturer S/N: 93480388
Instrument: Multimeter, True RMS,
Manufacturer: Fluke Corp
Model Number: 179
Asset Location: N/A

Received: 09/11/2014
Date Done: 09/12/2014
Date Due: 09/30/2015
Calibration Interval: 12 Months
Calibrated at Customer's Site: No

As Found: In Tolerance
As Returned: In Tolerance
Physical Damage: No

Characteristic Tested	Nominal Value	Tolerance	As Found	% Error	As Returned	% Error	TAR	EMU
AC Volts								
600.0 mV Range	300.0 mV~@ 45 Hz	296.7 to 303.3 mV	299.7	-9 %	Same	-9 %	24:1	170.8 μ V
6.000 V Range	5.000 V~@ 500 Hz	4.947 to 5.053 V	4.990	-19 %	Same	-19 %	14:1	4.6 mV
	5.000 V~@ 1 kHz	4.897 to 5.103 V	4.951	-48 %	Same	-48 %	26:1	4.6 mV
60.00 V Range	50.00 V~@ 45 Hz	49.47 to 50.53 V	49.90	-19 %	Same	-19 %	20:1	30.9 mV
	50.00 V~@ 1 kHz	48.97 to 51.03 V	50.02	2 %	Same	2 %	39:1	30.9 mV
600.0 V Range	300.0 V~@ 45 Hz	296.7 to 303.3 V	299.7	-9 %	Same	-9 %	20:1	203.7 mV
	500.0 V~@ 500 Hz	494.7 to 505.3 V	500.1	2 %	Same	2 %	17:1	366.0 mV
	500.0 V~@ 1 kHz	489.7 to 510.3 V	500.1	1 %	Same	1 %	33:1	366.0 mV
1000 V Range	1000 V~@ 45 Hz	987 to 1013 V	995	-38 %	Same	-38 %	21:1	925.2 mV
DC Volts								
6.000 V Range	5.000 V	4.992 to 5.008 V	4.999	-13 %	Same	-13 %	10:1	1.0 mV
600.0 V Range	300.0 V	299.5 to 300.5 V	299.9	-20 %	Same	-20 %	20:1	64.0 mV
1000 V Range	1000 V	996 to 1004 V	1000	0 %	Same	0 %	44:1	584.7 mV
1000 V Range	-1000 V	-1002 to -998 V	-1000	0 %	Same	0 %	31:1	584.7 mV
AC Volts Frequency								
99.99 Hz Range	45.00 Hz~@ 1 V	44.94 to 45.06 Hz	45.00	0 %	Same	0 %	49:1	5.9 mHz
99.99 kHz Range	50.00 kHz~@ 5 V	49.94 to 50.06 kHz	50.00	0 %	Same	0 %	48:1	6.0 Hz
DC Volts Frequency								
99.99 Hz Range	45.00 Hz~@ 3 V	44.94 to 45.06 Hz	45.00	0 %	Same	0 %	49:1	5.9 mHz
99.99 kHz Range	50.00 kHz~@ 30 V	49.94 to 50.06 kHz	50.00	0 %	Same	0 %	48:1	6.0 Hz
DC Millivolts								
600.0 mV Range	30.0 mV	29.8 to 30.2 mV	30.0	0 %	Same	0 %	38:1	58.1 μ V
	-300.0 mV	-300.5 to -299.5 mV	-299.8	40 %	Same	40 %	21:1	63.2 μ V
600.0 mV Range	600.0 mV	599.3 to 600.7 mV	599.6	-57 %	Same	-57 %	9.5:1	106.6 μ V
Resistance								
600.0 Ohm Range	500.0 Ohm	495.3 to 504.7 Ohm	500.1	2 %	Same	2 %	30:1	188.1 mOhm
50.00 MOhm Range	19.00 MOhm	18.68 to 19.32 MOhm	18.99	-3 %	Same	-3 %	10:1	35.7 kOhm
Continuity								
Beeper On 25 Ohm	Pass/Fail	Pass/Fail	PASS		Same			
Beeper Off 250 Ohm	Pass/Fail	Pass/Fail	PASS		Same			

Page 1 of 2

westerncal.com phone 801.466.1700 fax 801.484.5109 105 west 2950 south salt lake city, utah 84115.3433

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 18 of 24

Certificate Number: 497572							
Characteristic Tested	Nominal Value	Tolerance	As Found	% Error	As Returned	% Error	TAR EMU
Diode							
2.400 V Range	2.000 V	1.978 to 2.022 V	2.002	9 %	Same	9 %	>100:1 606.8 μ V
Temperature							
	0.0 °C	-1.0 to 1.0 °C	0.0	0 %	Same	0 %	5.3:1 0.23 °C
	-40.0 °C	-41.4 to -38.6 °C	-40.0	0 %	Same	0 %	7.4:1 0.23 °C
	400.0 °C	395.0 to 405.0 °C	399.8	-4 %	Same	-4 %	22:1 0.27 °C
Capacitance							
1000 nF Range	900 nF	866 to 934 nF	900	0 %	Same	0 %	15:1 1.2 nF
AC Current							
6.000 A Range	4.000 A~@ 45 Hz	3.957 to 4.043 A	3.999	-2 %	Same	-2 %	5.7:1 12.7 mA
10.00 A Range	9.00 A~@ 1 kHz	8.84 to 9.16 A	8.99	-6 %	Same	-6 %	7.9:1 24.9 mA
60.00 mA Range	3.00 mA~@ 45 Hz	2.92 to 3.08 mA	3.02	25 %	Same	25 %	31:1 6.4 μ A
	50.00 mA~@ 1 kHz	49.22 to 50.78 mA	49.94	-8 %	Same	-8 %	11:1 83.3 μ A
400.0 mA Range	400.0 mA~@ 1 kHz	393.7 to 406.3 mA	399.0	-16 %	Same	-16 %	7.2:1 1.0 mA
DC Current							
60.00 mA Range	3.00 mA	2.94 to 3.06 mA	3.01	17 %	Same	17 %	>100:1 5.8 μ A
	50.00 mA	49.47 to 50.53 mA	49.99	-2 %	Same	-2 %	30:1 21.1 μ A
400.0 mA Range	-400.0 mA	-404.3 to -395.7 mA	-399.9	2 %	Same	2 %	12:1 417.4 μ A
6.000 A Range	4.000 A	3.957 to 4.043 A	4.001	2 %	Same	2 %	14:1 3.7 mA
10.00 A Range	-9.00 A	-9.12 to -8.88 A	-9.00	0 %	Same	0 %	20:1 8.9 mA

(* - "Out of Tolerance" condition) (0 % Error of Limit - May represent a condition with an Asymmetrical Tolerance that can not be calculated.)
(Unit of Measure for both As Found and As Returned are the same as the Nominal Value.)

Standards Utilized (Make, Model & Description)	Trace Number	Asset Number	Date Due
Wavetek Corp 9100 Calibrator, Multifunction,	474841	WSC601	03/06/2015
ARCO ELECTRONICS SS-32 Capacitor, Standard,	413166	WSC657	06/14/2015
Description of Material or Parts Used	Quantity	Material or Parts ID	
Fuse	1	DMM44/100	
Battery	1	9 V	

Comments (i.e. adjustments, repairs, modifications, limitations and/or deviations from procedure):

Battery and blown fuse were replaced.

Calibration Procedure: 10224 R1

Quality Manager: Lon Miles

Degrees Fahrenheit: 73.0
Percent Humidity: 32.0

Calibration Technician: *Kirk Scriber* Kirk Scriber
09/12/2014

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 19 of 24

Appendix XIV Zip Pen Packaging

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 28 of 33

Appendix XIV

ZIP PEN PACKAGING

Reference Assembly Drawings 2525-10 and 2525-15

Verify that the Assembly meets the following:

The assembly box is an RSC style box:

☒ YES ☐ NO Initials/Date PV 1-29-2015

The assembly requires an IFU in each box:

☒ YES ☐ NO Initials/Date PV 1-29-2015

The assembly requires 20 units per box:

☒ YES ☐ NO Initials/Date PV 1-29-2015

Comments: _____

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 20 of 24

Appendix XV Zip Pen Extension Nozzle Packaging

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 29 of 33

Appendix XV

ZIP EXTENSION NOZZLE PACKAGING

Reference Assembly Drawing 2540 and 2560

Verify that the Assembly meets the following:

The assembly box is an RSC style box:

☒ YES ☐ NO

Initials/Date PV 1-29-2015

The assembly requires an IFU in each box:

☒ YES ☐ NO

Initials/Date PV 1-29-2015

The assembly requires 10 units per box:

☒ YES ☐ NO

Initials/Date PV 1-29-2015

Comments: _____

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 21 of 24

Appendix XVI Zip Pen and Extension Nozzle Unit Labeling

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 30 of 33

Appendix XVI

ZIP PEN AND EXTENSION NOZZLE UNIT LABELING

Reference Unit Label Drawings ENG-DWG-272, ENG-DWG-277, MKT-LBL-514 and MKT-LBL-517. Verify that the Labels meet the following:

Requirement	2525-10	2525-15	2540	2560
Label Part Number	X	X	X	X
Label Revision and Date	X	X	X	X
Bar code	X	X	X	X
Catalog Number and Symbol	X	X	X	X
Product Trade Name	X	X	X	X
Manufacturer symbol and name/info	X	X	X	X
Graphic of the product	X	X	X	X
E-Z Clean Logo	X	X	N/A	N/A
Lot Number and Symbol	X	X	X	X
Expiration Date and Symbol	X	X	X	X
CE Mark	X	X	X	X
EC Rep and symbol	X	X	X	X
Sterilization Symbol	X	X	X	X
Do Not Reuse Symbol	X	X	X	X
Consult Instructions Symbol	X	X	X	X
Rx Only Symbol	X	X	X	X
Do Not Use if Package is Damaged Symbol	X	X	X	X
The words "OPEN HERE"	X	X	N/A	N/A
Refer to Instructions for Use for Patent Information	X	X	N/A	N/A
Temperature and Humidity Symbol	X	X	X	X
Rated Voltage Symbol	X	X	N/A	N/A

Acceptance YES/NO Initials/Date PV 1-29-2015

Comments: _____

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Confidential - Controlled copy on date printed

Confidential - Pre-Production Controlled on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 22 of 24

Appendix XVII

Zip Pen and Extension Nozzle Shipping Labels

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 31 of 33

Appendix XVII

ZIP PEN AND EXTENSION NOZZLE SHIPPING LABELS

Reference Shipper Label Drawings MKT-LBL-499, MKT-LBL-510, MKT-LBL-515 and MKT-LBL-516. Verify that the Labels meet the following:

Requirement	2525-10	2525-15	2540	2560
Label Part Number	X	X	X	X
Label Revision and Date	X	X	X	X
Bar code	X	X	X	X
Catalog Number and Symbol	X	X	X	X
Manufacturer symbol and name/info	X	X	X	X
Product Trade Name	X	X	X	X
Graphic of the product	X	X	X	X
E-Z Clean Logo	X	X	N/A	N/A
CE Mark	X	X	X	X
EC Rep and symbol	X	X	X	X
Lot Number and Symbol	X	X	X	X
Expiration Date and Symbol	X	X	X	X
Sterilization Symbol	X	X	X	X
Do Not Reuse Symbol	X	X	X	X
Consult Instructions Symbol	X	X	X	X
Rx Only Symbol	X	X	X	X
Do Not Use if Package is Damaged Symbol	X	X	X	X
Temperature and Humidity Symbol	X	X	X	X
Quantity per box symbol	X	X	X	X

Acceptance YES/NO Initials/Date PV 1-29-2015

Comments: _____

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Confidential - Controlled copy on date printed

Confidential - Pre-Production Controlled on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 23 of 24

Appendix XVIII ULPA and Carbon Filter Unit Labeling

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 32 of 33

Appendix XVIII

ULPA AND CARBON FILTER UNIT LABELING

Reference Unit Label Drawings MKT-LBL-518 and MKT-LBL-522. Verify that the Labels meet the following:

Requirement	2211	2220
Label Part Number	X	X
Label Revision and Date	X	X
Bar code	X	X
Catalog Number and Symbol	X	X
Product Trade Name	X	N/A
Manufacturer symbol and name/info	X	X
Graphic of the product	X	X
Lot Number and Symbol	X	X
CE Mark	X	X
EC Rep and symbol	X	X
Consult Instructions Symbol	X	X
Rx Only Symbol	X	X
Temperature and Humidity Symbol	X	X

Acceptance ☒ YES ☐ NO Initials/Date PV 1-29-2015

Comments: _____

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-418
	Zip Pen Product Specification Verification	Revision: 002
		Page 24 of 24

Appendix XIX ULPA and Carbon Filter Shipper Labels

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-290 Rev: A Effective: 23 Jan 2015 11:41 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-290
	Zip Pen Product Specification Verification Protocol	Revision: A
		Page 33 of 33

Appendix XIX

ULPA AND CARBON FILTER SHIPPING LABELS

Reference Shipping Label Drawings MKT-LBL-519 and MKT-LBL-523. Verify that the Labels meet the following:

Requirement	2211	2220
Label Part Number	X	X
Label Revision and Date	X	X
Bar code	X	X
Catalog Number and Symbol	X	X
Manufacturer symbol and name/info	X	X
Product Trade Name		N/A
Graphic of the product	X	X
CE Mark	X	X
EC Rep and symbol	X	X
Lot Number and Symbol	X	X
Consult Instructions Symbol	X	X
Rx Only Symbol	X	X
Temperature and Humidity Symbol	X	X
Quantity per box symbol	X	X

Acceptance ☒ YES/NO Initials/Date PV 1-10-2018

Comments: PV-Original conductor of test corrected missing initials.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 23 Jan 2015, 02:14:56 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:58:17 pm; Printed by: .

Confidential - Controlled copy on date printed

Confidential - Pre-Production Controlled on date printed