

DOCUMENT NUMBER: ENG-RPT-452

DOCUMENT TITLE: Test Report: ISOS2 EO Sterilized Product Design Verification (t=0)

DOCUMENT NOTES:

#### Document Information

Revision: 001

Vault: MEG-Rel

Status: Release

Document Type: ENG-RPT

#### Date Information

Effective Date: 18 Apr 2016

Expiration Date:

Release Date: 18 Apr 2016

Next Review Date:

#### Control Information

Author: CCROOK

Owner: CCROOK

Previous Number:

Change Number: 2016-ENG-DCO-030

## Signature Manifest

**Document Number:** ENG-RPT-452**Revision:** 001**Title:** Test Report: ISOS2 EO Sterilized Product Design Verification (t=0)

All dates and times are in Mountain Standard Time.

### ISOS2 EO Compatibility Testing

#### Change Request

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		08 Apr 2016, 01:55:56 PM	Approved

#### Collaboration

Name/Signature	Title	Date	Meaning/Reason
Ihsan Samara (ISAMARA)	Quality Manager	11 Apr 2016, 02:35:22 PM	Complete
Haven McCall (HMCCALL)	VP RA/QA	12 Apr 2016, 07:45:24 AM	Complete
Paul Borgmeier (PBORGMEIER)		12 Apr 2016, 04:33:57 PM	Complete
Mike Hintze (MHINTZE)		14 Apr 2016, 05:07:29 PM	Complete
Christian Crook (CCROOK)		15 Apr 2016, 09:02:53 AM	Complete

#### Document Review

Name/Signature	Title	Date	Meaning/Reason
Christian Crook (CCROOK)		15 Apr 2016, 11:49:59 AM	Complete
Lucy Richards (LRICHARDS)		15 Apr 2016, 04:03:59 PM	Complete

#### RA-Approval

Name/Signature	Title	Date	Meaning/Reason
Haven McCall (HMCCALL)	VP RA/QA	18 Apr 2016, 09:11:56 AM	Approved

#### QA-Approval

Name/Signature	Title	Date	Meaning/Reason
Ihsan Samara (ISAMARA)	Quality Manager	18 Apr 2016, 11:41:06 AM	Approved

#### ENG-Approval

Name/Signature	Title	Date	Meaning/Reason
Paul Borgmeier (PBORGMEIER)		18 Apr 2016, 09:47:46 AM	Approved

#### Training Review

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		18 Apr 2016, 11:47:48 AM	Approved

#### Final Release

Printed on: 21 Jan 2020, 11:00:10 pm; Printed by: .

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		18 Apr 2016, 11:47:59 AM	Approved

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b> <b>Page 1 of 79</b>

Authored By: Christian Crook

	<b>TABLE OF CONTENTS</b>	<b>Page</b>
1.	ABSTRACT .....	1
2.	OBJECTIVE .....	1
3.	RESULTS AND DISCUSSION .....	1
4.	DISCUSSION .....	7
5.	CONCLUSIONS .....	12
6.	RECOMMENDATIONS .....	15
7.	REVISION HISTORY .....	15
<b>1.</b>	<b>ABSTRACT</b>	

The single use electrosurgical pencil part numbers under ENG-DWG-656 (ND PN: PB252SM1 LOT# S140520) and ENG-DWG-657 (ND PN: PR252SM1 LOT# S150012) was evaluated for conformance to requirements of IEC 60601-1:2012, IEC 60601-2-2: 2009, and ENG-PS-010 per test protocol ENG-PRT-302 after exposure to the validated EO sterilization cycle specified in ENG-WI-001, Sterilization Chart. The electrosurgical pencil meets the requirements of these standards for all items tested in the protocol after exposure to the EO sterilization cycle, shipping conditions, and humidity and temperature extremes.

## **2. OBJECTIVE**

This test report documents test results obtained while evaluating conformance to requirements of IEC 60601-1: 2012, IEC 60601-2-2: 2009, and ENG-PS-010 per test protocol ENG-PRT-302.

## **3. REFERENCES**

XENG-PRT-306	ISOS2 EO Sterilization Product Design Verification
ENG-WI-003	Sterilization Chart
ENG-RPT-443	ISOS2 Limited Design Verification Aged (t=5year), Thermal Cycled, Ship and Vibration Cycle Product
ENG-PS-010	Product Specification Disposable Electrosurgical Pencil
ENG-DWG-1222	ESU Pencil, Cable
ENG-DWG-1227	PLUG, CLAMSHELL, TOP, ESU PENCIL
ENG-DWG-1228	PLUG, CLAMSHELL, BOTTOM, ESU PENCIL

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 2 of 79</b>

## 4. APPENDICES

APPENDIX I	HANDPiece HIGH FREQUENCY DIELECTRIC WITHSTAND TESTING
APPENDIX II	CABLE LEAKAGE TEST
APPENDIX III	PLUG MAINS DIELECTRIC
APPENDIX IV	FLUID INGRESS
APPENDIX V	IEC 60601-1 TEMPERATURE TEST
APPENDIX VI	HANDPiece STRAIN RELIEF (ROTATIONAL) TEST
APPENDIX VII	ACTIVATION FORCE TESTING
APPENDIX VIII	ACTIVATION OVER TIME
APPENDIX IX	ELECTRODE WOBBLE
APPENDIX X	DURABILITY OF PRINT TO CHEMICAL IN USE ENVIRONMENT (ROCKER PENCIL ONLY)
APPENDIX XI	PENCIL BODY PEEL STRENGTH

## 5. RESULTS

### 5.1. PENCIL BUILD DOCUMENTATION, STERILIZATION

- 5.1.1. Pencil components were built at NewDeantronics under LOTS S150083 and S150084 on a certificate of compliance. Pencils were complete to finished product at Megadyne under work orders 151361, 151362, and 151363.

An additional LOT of ACE37H pencils for fluid ingress testing was built under Megadyne work order 160868.

Build documentation is included as attachments to the Info Card in MasterControl.

- 5.1.2. Final product under LOTS 151361 and 151363 was sterilized twice at Sterigenics on 4/21/2015 and 5/25/2015 using Cycle 115. Sterilization of product under LOT 160868 was completed at Nelson Labs on March 17, 2016 using exposure conditions equivalent to Sterigenics Cycle 115.

Sterilization documentation is included as an attachment to the Info Card in MasterControl.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b> <b>Page 3 of 79</b>

## 5.2. SHIPPING AND STORAGE CYCLE PRECONDITIONING

- 5.2.1. Shipping and storage conditioning were completed as instructed in section 9.1 of XENG-PRT-315. See Appendix I.

## 5.3. SHIPPING TEST

- 5.3.1. All shipping test drops and compressions were completed for both 0036 and ACE37H according to Section 10 of XENG-PRT-302. See Appendix I.

## 5.4. HANDPIECE HF DIELECTRIC TEST

- 5.4.1. Continuity is checked at the beginning of testing. All samples of both button and rocker pencils had adequate continuity to perform testing.
- 5.4.2. For each sample of button and rocker pencils, the test voltage was maintained at 30 seconds without any observed breakdown. The lowest HF voltage applied was 6.60 kV. The highest activation resistance measured after the test was  $24.6\ \Omega$ , well below the  $50\ \Omega$  maximum. Raw data are included in Appendix II.

## 5.5. HANDPIECE MAINS DIELECTRIC TEST

- 5.5.1. No breakdown was observed for any sample of button and rocker pencils subjected to the  $4,600\ V_{RMS}$  condition for 30 seconds. Raw data is included in Appendix II.

## 5.6. CABLE LEAKAGE TEST

- 5.6.1. Continuity is checked at the beginning of testing. All samples of demonstrated adequate continuity to perform testing.
- 5.6.2. The allowable leakage was calculated to be 120.96 mA. The greatest leakage measured was 39.5, the lowest was 33.2 mA. See Appendix III for raw data.

## 5.7. CABLE HF DIELECTRIC TEST

- 5.7.1. For each sample of cable, the test voltage was maintained for 30 seconds. The lowest HF voltage applied was 6.60 kV. No sample

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 4 of 79</b>

showed signs of breakdown during the test. Raw data are included in Appendix III.

#### 5.8. CABLE MAINS DIELECTRIC TEST

- 5.8.1. A voltage of 4,596 V<sub>RMS</sub> was applied to the cable samples for 5 minutes. No damage was observed, no sparks were observed, and the HI Pot alarm did not sound. All resistances were within established limits after testing. See Appendix III.

#### 5.9. PLUG MAINS DIELECTRIC TEST

- 5.9.1. No visible signs of damage were observed and the HI Pot did not alarm when samples were subjected to 4,600 V<sub>RMS</sub> for 30 seconds. See Appendix IV.

#### 5.10. HOLSTER HIGH FREQUENCY DIELECTRIC WITHSTAND TESTING

- 5.10.1. No change to component. NOT TESTED.

#### 5.11. HOLSTER MAINS FREQUENCY DIELECTRIC WITHSTAND TESTING

- 5.11.1. No change to component. NOT TESTED.

#### 5.12. FLUID INGRESS

- 5.12.1. Button switch pencils: Switches on all test samples were successfully operated 10 times with current less than 2.5 mA for 0.5 seconds after released. All samples passed.

- 5.12.2. Rocker switch pencils: Switches on all test samples were successfully operated 10 times with current less than 2.5 mA for 0.5 seconds after released. All samples passed. One sample displayed self activation. Raw data is included in Appendix V.

#### 5.13. IEC 60601-1 TEMPERATURE TEST

- 5.13.1. Pencil temperature remained below 41 °C for button and rocker pencils during activation. See Appendix V for raw data and analysis.

Sample and Test Description	Temperature		
	Before (°C)	After (°C)	Change (°C)
ISOS2, Rocker, 300 W CUT	22.9	24.4	1.5
ISOS2, Rocker, 120 W COAG	22.7	22.9	0.7

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 5 of 79

ISOS2, Button, 300 W CUT	22.2	24.9	2.7
ISOS2, Button, 120 W COAG	22.6	22.1	0.1

#### 5.14. PLUG STRAIN RELIEF (ROTATIONAL) TEST

5.14.1. Testing was not performed.

#### 5.15. HANDPIECE STRAIN RELIEF (ROTATIONAL) TEST

5.15.1. No sample, button or rocker pencil, showed any evidence that the cord had worked loose of the housing or showed any damage.

5.15.2. In all cases, the cable maintained acceptable continuity (less than  $50\Omega$ ) between the CUT wire and the Active wire with the CUT button depressed. The greatest resistance measured was  $21.0\ \Omega$ . Raw data is included in Appendix VII.

5.15.3. In all cases, the cable maintained acceptable continuity (less than  $50\Omega$ ) between the COAG wire and the Active wire with the COAG button depressed. The greatest resistance measured was  $33.4\ \Omega$ . Raw data is included in Appendix VII.

#### 5.16. ACTIVATION FORCE TESTING

5.16.1. Rocker Pencil 0036: All button forces were well within the specifications. The values range from 527.75 grams force to 381.75 grams force. Cpk values for the button forces are well over the generally accepted minimum of 1.33.

Button Activation Force 0036 LOT 151363			
	SPEC	CUT	COAG
MAX	700	527.75	458.83
MIN	250	446.41	381.75
AVG	-	484.98	417.65
SD	-	19.57	20.67
Cpk	-	3.66	2.70

5.16.2. Button Pencil ACE37H: All button forces were well within the specifications. The values range from 434.72 grams force. to 343.39 grams force. Cpk values for the button forces are well over the generally accepted minimum of 1.33.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 6 of 79

Button Activation Force ACE37H LOT 151361			
	SPEC	CUT	COAG
MAX	700	434.72	386.74
MIN	250	381.87	343.39
AVG	-	409.21	366.31
SD	-	13.93	10.95
Cpk	-	3.81	3.54

5.16.3. See Appendix VIII for raw data.

## 5.17. ACTIVATION OVER TIME

5.17.1. All pencils were activated 500 times. See Raw Data in Appendix IX.

5.17.2. Rocker Pencil 0036: All button forces were well within the specifications after 500 activations. The values range from 489.52 grams force to 364.58 grams force. Cpk values for the button forces are well over the generally accepted minimum of 1.33.

Button Activation Force 0036 LOT 151363			
	SPEC	CUT	COAG
MAX	700	489.52	448.60
MIN	250	414.14	364.58
AVG	-	448.34	398.34
SD	-	16.68	18.92
Cpk	-	3.96	2.61

5.17.3. Button Pencil ACE37H: All button forces were well within the specifications after 500 activations. The values range from 403.06 grams force. to 329.02 grams force. Cpk values for the button forces are well over the generally accepted minimum of 1.33.

Button Activation Force ACE37H LOT 151361			
	SPEC	CUT	COAG
MAX	700	403.06	367.99
MIN	250	356.79	329.02
AVG	-	379.60	349.84
SD	-	13.04	10.00
Cpk	-	3.31	3.33

## 5.18. ELECTRODE WOBBLE

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 7 of 79</b>

5.18.1. Electrode wobble in ISOS2 pencils is slightly higher than ISOS1 pencils. See Appendix X.

<b>Electrode Wobble</b>			
	<b>0036</b>	<b>ACE37H</b>	<b>Control (0035)</b>
<b>MAX</b>	0.155	0.151	0.123
<b>MIN</b>	0.130	0.127	0.094
<b>MEAN</b>	0.141	0.139	0.116
<b>SD</b>	0.006	0.007	0.005

#### 5.18.2. CORD LENGTH

5.18.3. Cord length was not tested in this protocol.

#### 5.19. DURABILITY OF PRINT TO CHEMICAL IN USE ENVIRONMENT (ROCKER PENCIL ONLY)

5.19.1. 0036: The ink remained visible after applying both saline and lactated ringers on all samples. See Appendix XI.

#### 5.20. RESISTANCE TO ISOPROPYL ALCOHOL

5.20.1. Not tested.

#### 5.21. PENCIL BODY PEEL STRENGTH

5.21.1. All peel strength values were greater than 10 lbs. The lowest value observed is 22.07 lbs. and the greatest value observed is 49.96 lbs. See Appendix XII.

### 6. DISCUSSION

#### 6.1. PENCIL BUILD DOCUMENTATION, STERILIZATION

6.1.1. Pencils are traceable and conform to the design.

6.1.2. Pencils were sterilized at an exposure between 50 and 55 kGy.

#### 6.2. SHIPPING AND STORAGE CYCLE PRECONDITIONING

6.2.1. Shipping and storage conditions were completed on the equipment below.

<b>Thermotron #1</b>
Model number - SE-600-3-3
Megadyne number - 01095

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 8 of 79</b>

Serial number - 28105
Calibration Date - 5/19/2015
Due Date - 5/31/2016
Shipping Weight - 1,680 lbs.

Sample preconditioning was performed to ensure that product was tested after being exposed to temperatures and humidity it would encounter in the field prior to testing.

### 6.3. SHIPPING TEST

- 6.3.1. Both catalog numbers 0036 and ACE37H were subjected to shipping conditions and storage conditions experienced by the product in the field.
- 6.3.2. Product was subjected to ship testing to ensure that product experienced the shock, vibration, and compression likely to occur in the field prior to testing, not to validate the packaging configuration. Packaging is not validated as part of this protocol.

### 6.4. HANDPIECE HIGH FREQUENCY DIELECTRIC WITHSTAND TESTING

- 6.4.1. There were no alarms or visible signs of damage to any hand piece during exposure to 120% of the rated accessory voltage. All samples satisfied the acceptance criteria.

### 6.5. HANDPIECE MAINS FREQUENCY DIELECTRIC WITHSTAND TESTING

- 6.5.1. There were no alarms or visible signs of damage to any hand piece. All samples satisfied the acceptance criteria.

### 6.6. CABLE LEAKAGE TEST

- 6.6.1. Both pencils are manufactured with the same cable, ENG-DWG-1222.
- 6.6.2. Leakage current was much less than 120.96 mA for all test samples. All acceptance criteria for ENG-DWG-1222 were satisfied.

### 6.7. CABLE HF DIELECTRIC TEST

- 6.7.1. The cable samples were undamaged by 6.6 kV for 30 seconds. 6.6 kV corresponds to 120% of the Rated Accessory Voltage of 5.5kV. All acceptance criteria for ENG-DWG-1222 were satisfied.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b> <b>Page 9 of 79</b>

## 6.8. CABLE MAINS DIELECTRIC TEST

6.8.1. All acceptance criteria for cable ENG-DWG-1222 were satisfied.

## 6.9. PLUG MAINS TEST

6.9.1. The plug is constructed from three (3) components; ENG-DWG-1222, ENG-DWG-1227, and ENG-DWG-1228; for both 0036 and ACE37H pencils, and all ISOS2 pencils.

6.9.2. All acceptance criteria for the plug design were satisfied.

## 6.10. HOLSTER HIGH FREQUENCY DIELECTRIC WITHSTAND TESTING

6.10.1. The component was not available for testing. This component will not be qualified as part of this testing.

## 6.11. HOLSTER MAINS FREQUENCY DIELECTRIC WITHSTAND TESTING

6.11.1. The component was not available for testing. This component will not be qualified as part of this testing.

## 6.12. FLUID INGRESS

6.12.1. Button switch pencils: All pencil buttons (ACE37H, Lot# 160868) operated correctly after exposure to saline. All samples satisfied the acceptance criteria.

6.12.2. Rocker switch pencils: All pencil buttons (0036, Lot# 151363) operated correctly after exposure to saline. All samples satisfied the acceptance criteria.

6.12.3. It should be noted that fluid ingress testing on rocker switch pencils was performed on two lots of button pencils. The second lot as indicated in section 5.1.1 above, passed. However, in the initial sample group (ACE37H, Lot# 151361), one pencil (sample #9 of ACE37H) self-activated during operation of the CUT button when connected to a generator. This failure was confirmed with a multimeter. The failure started in the CUT button and migrated to the COAG button over time.

The pencil was inspected under magnification. No obvious cause was found. The test method was examined and a discrepancy was found. In the case of this pencil, fluid ingress testing was the last of all tests performed. All pencils had been subjected to button force

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b> <b>Page 10 of 79</b>

testing and 500 button activations prior to fluid ingress. The 500 activations were performed using a new method, an Instron measurement system and fixtures. In all testing previously, 500 button activations had been performed by hand. The crosshead travel distance on the Instron was set to a distance that ensured that each button was pressed. Because button travel distance was controlled, instead of button force, wear may have been caused by buttons being over-depressed or due to poor alignment between the button and the pencil housing. Either of these cases could cause wear that would induce a fluid ingress failure not related to the design of the pencil, but to an unrelated test method.

This design adequately resists fluid ingress in conditions of normal use.

#### 6.13. IEC 60601-1 TEMPERATURE TEST

6.13.1. The acceptance criteria were satisfied. The temperature on the pencil did not exceed 41°C on any pencil or in any image. In three cases, the temperature rose only slightly. All acceptance criteria were satisfied.

#### 6.14. PLUG STRAIN RELIEF (ROTATIONAL) TEST

6.14.1. Because the design for retaining the cable in the plug and hand piece are similar, only the pencil hand piece strain relief test was performed. The hand piece test is more rigorous and is representative of plug performance as well. See ENG-PRT-306 section 10.5.2 for a complete discussion.

#### 6.15. HANDPIECE STRAIN RELIEF (ROTATIONAL) TEST

6.15.1. As discussed in ENG-PRT-306, the pencil and plug use nearly identical means of holding the plug. The pencil strain relief test is more strenuous than the plug test because it applies more rotation cycles. Only the strain relief for the hand piece was tested, but results are also valid for the plug.

6.15.2. The strain relief adequately prevents the cable from wearing and exposing conductors or forming short circuits during use. The strain relief design meets the acceptance criteria.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b> <b>Page 11 of 79</b>

## 6.16. ACTIVATION FORCE TESTING

6.16.1. All samples were within specification. The design satisfies the acceptance criteria of the protocol.

## 6.17. ACTIVATION OVER TIME

6.17.1. In previous testing, pencils had been activated 500 times by hand, each button was activated independently. Test samples used in this protocol were activated using the Instron Measurement System and a fixture. The test fixture allowed CUT and COAG buttons of five (5) pencils to activated simultaneously reducing the time required to activate all pencils. During testing, periodic continuity checks were made to ensure that pencils were being electrically activated during the test.

6.17.2. Button resistance is less than  $50\ \Omega$  in the depressed condition and greater than  $100,000\ \Omega$  after 500 activations in every case. Activation force is also within specification after 500 activations in every case. The design meets the acceptance criteria of the protocol after 500 activation cycles.

## 6.18. ELECTRODE WOBBLE

6.18.1. The protocol requirement for the electrode wobble as listed in section 31.18 of ENG-PRT-302 requires the wobble for the ISOS2 pencil to be less than that of current ISOS pencils. However, upon completion of the testing reported under XENG-RPT-452, it is noted that the ISOS2 pencil electrodes wobble (deviate) more than the LOT of current ISOS button pencils tested and do not meet the requirements of the protocol. The range difference in these means is  $0.025''$ . This is deemed not significant.

### 6.18.2.

#### 6.18.1.

## 6.19. CORD LENGTH

6.19.1. Cord length was tested previously as documented in Test Report ENG-RPT-433 for 15 foot cord pencils and Test Report ENG-RPT-438 for 10 ft. pencil.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 12 of 79</b>

## 6.20. DURABILITY OF PRINT (ROCKER PENCIL ONLY)

- 6.20.1. Both saline and lactated ringers are chemicals likely used to clean the pencil in the use environment. The ink was resistant to both solutions and did not smear or degrade during the test. The pencil design satisfies the acceptance criteria.
- 6.20.2. This alternative test for compliance with IEC 60601-1:2012 Clause 7.1.3 was used according to IEC 6060101-1:2012 Clause 4.5.

The risk associated with the button color is that the user will activate the CUT functionality instead of COAG if the ink is removed during manufacture or use. Button color is a redundant safety measure. Section 201.12.2 of IEC60601-2-2 requires that the CUT button be located closest to the electrode. As part of HF Dielectric testing, CUT and COAG functionality (i.e. the CUT button activates cut functionality, the COAG button activates coagulation functionality) is checked prior to the resistance measurement. See Appendix VIII.

## 6.21. RESISTANCE TO ISOPROPYL ALCOHOL

- 6.21.1. Isopropyl alcohol (70%) is used to clean pencils as needed in the manufacturing environment, prior to EO sterilization. The pencil will not be exposed to isopropyl alcohol under normal use. Resistance to isopropyl alcohol was completed as recorded in ENG-RPT-443.

## 6.22. PENCIL BODY PEEL STRENGTH

- 6.22.1. All test samples satisfy the acceptance criteria.

# 7. CONCLUSIONS

## 7.1. PENCIL BUILD DOCUMENTATION, STERILIZATION

- 7.1.1. Pencils represent the final product.
- 7.1.2. Pencils were correctly sterilized at Sterigenics using Cycle 115 per ENG-WI-003.

## 7.2. SHIPPING AND STORAGE CYCLE PRECONDITIONING

- 7.2.1. Both catalog numbers were subjected to temperature and humidity conditions conforming to ASTM D4169 prior to functional testing.

## 7.3. SHIPPING TEST

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b> <b>Page 13 of 79</b>

7.3.1. Both catalog numbers were subjected to shock, vibration, and compression conditions conforming to ASTM D4169 prior to functional testing.

#### 7.4. HANPIPE HIGH FREQUENCY DIELECTRIC WITHSTAND TESTING

7.4.1. The ISOS2 pencil hand piece is compliant with IEC 60601-2-2: 2009 Clause 201.8.8.3.103 after two exposures to the validated EO sterilization cycle.

#### 7.5. HANPIPE MAINS FREQUENCY DIELECTRIC WITHSTAND TESTING

7.5.1. The ISOS2 pencil hand piece is compliant with IEC 60601-2-2: 2009 Clause 201.8.8.3.104 after two exposures to the validated EO sterilization cycle.

#### 7.6. CABLE LEAKAGE TEST

7.6.1. The ISOS2 cable ENG-DWG-1222 is compliant with IEC 60601-2-2 Clause 201.8.8.3.101 and 201.8.8.3.102 after two exposures to the validated EO sterilization cycle.

#### 7.7. CABLE HF DIELECTRIC WITHSTAND TEST

7.7.1. The ISOS2 cable is compliant with IEC 60601-2-2 Clause 201.8.8.3.101 and 201.8.8.3.103 after two exposures to the validated EO sterilization cycle.

#### 7.8. CABLE MAINS DIELECTRIC WITHSTAND TEST

7.8.1. The ISOS2 pencil cable is compliant with IEC 60601-2-2 Clause 201.8.8.3.101 and 201.8.8.3.104 after two exposures to the validated EO sterilization cycle.

#### 7.9. PLUG MAINS TEST

7.9.1. The plug is constructed from three (3) components; ENG-DWG-1222, ENG-DWG-1227, and ENG-DWG-1228; is compliant with IEC60601-2-2 Clause 201.8.8.3.101, and Clause 201.8.8.3.104.

#### 7.10. HOLSTER HIGH FREQUENCY DIELECTRIC WITHSTAND TESTING

7.10.1. No change to component. Not tested.

#### 7.11. HOLSTER MAINS FREQUENCY DIELECTRIC WITHSTAND TESTING

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 14 of 79</b>

7.11.1. No change to component. Not tested.

#### 7.12. FLUID INGRESS

7.12.1. CUT and COAG switches of button and rocker pencils conform to IEC 60601-2-2 Clause 201.11.6.5 b after two exposures to conditions of the validated EO sterilization cycle.

#### 7.13. IEC 60601-1 TEMPERATURE TEST

7.13.1. Button and rocker switch pencils conform to IEC60601-1 Clause 11.1.2.2 after two exposures to the validated EO sterilization cycle.

#### 7.14. PLUG STRAIN RELIEF (ROTATIONAL) TEST

7.14.1. The hand piece strain relief is functionally equivalent to the hand piece strain relief and meets the requirements of ICE60601-2-2 Clause 201.8.4.2 after two exposures to the validated EO sterilization cycle.

#### 7.15. CABLE ANCHORAGE (HANDPIECE) TEST

7.15.1. The hand piece and plug strain relief designs meet the requirements of ICE60601-2-2 Clause 201.8.4.2 after two exposures to the validated EO sterilization cycle.

#### 7.16. ACTIVATION FORCE TESTING

7.16.1. The button activation force is for CUT and COAG buttons on both 0036 and ACE37H is within the specification of ENG-PS-010 section 4.2.1 after two exposures to the validated EO sterilization cycle.

#### 7.17. ACTIVATION OVER TIME

7.17.1. Button resistance and activation force demonstrate that the switch withstands 500 activations per switch after 500 activations and two exposures to the validated EO sterilization cycle.

#### 7.18. ELECTRODE WOBBLE

7.18.1. The marketing department was consulted and a determination was made that the ISOS2 pencil wobble is acceptable though not meeting the intent of the protocol. It is determined that this wobble was not an issue in the field and the values are accepted as is. The pencil

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b> <b>Page 15 of 79</b>

wobble is acceptable after two exposures to the validated EO sterilization cycle.

#### 7.19. CORD LENGTH

7.19.1. Cord length was verified previously and was not evaluated in this testing.

#### 7.20. DURABILITY OF PRINT TO CHEMICAL IN USE ENVIRONMENT (ROCKER PENCIL ONLY)

7.20.1. The pencil design complies with IEC 60601-1:2012 Clause 7.1.3 after two exposures to the validated EO sterilization cycle. See also Clause 4.5.

#### 7.21. RESISTANCE TO ISOPROPYL ALCOHOL

7.21.1. Not evaluated in this report. See ENG-RPT-443.

#### 7.22. PENCIL BODY PEEL STRENGTH

7.22.1. ISOS2 pencils meet the requirement of product specification ENG-PS-010 requirement 4.2.9 after two exposures to the validated EO sterilization cycle.

### 8. RECOMMENDATIONS

8.1. No further testing is required.


### 9. REVISION HISTORY

Revision	Document Change Order Number	Description Of Change	Effective Date
See Master Control for Revision History			

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001 Page 16 of 79

## APPENDIX I – SHIPPING AND STORAGE CONDITIONING





CH 2/8/16

Paul Valpreda 2-8-2016

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001 Page 17 of 79

Confidential - Controlled copy on date printed



Paul Valpreda 2-8-2016

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 18 of 79

Conditioning Chamber Calibration Information	Chamber: 01268 Last Calibration: 5-19-2015	Cal Due: 5-31-2016	Pass/Fail
Maximum Temperature (°C):	70.2		
Minimum Temperature (°C):	-40.6		
Maximum Humidity (%RH):	95.1		
Minimum Humidity (%RH):	2.5		
Chamber conditions held at -40°C with no humidity control for a duration of 4 hours:	9:57 AM TO 2:11 PM on 3-17-2016 Duration of 4 hrs. 14 minutes		Pass
Chamber conditions held at 70°C and 95% RH for a duration of 4 hours:	14:05 TO 18:40 on 3/21/2016 Duration of 4 hrs. 35 minutes		Pass
Chamber conditions held at 70°C and 15% RH for a duration of 4 hours:	19:45 TO 23:55 on 3/21/2016 Duration of 4 hrs. 10 minutes		Pass
Chamber conditions held at 23°C and 50% RH for a duration of 72 hours:	3/22/16 12:50 AM TO 3/25/16 12:50 AM Duration of 72 hrs.		Pass

Paul Valpreda  
Test Technician Name

PAUL VALPREDA 3-25-2016  
Signature Date


Christian Crook  
Engineer Name

CC 3/25/2016  
Signature Date

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-452
ISOS2 Functional Verification After EO Exposure	Revision: 001	Page 19 of 79

### Low Temperature Conditioning.



<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>						<b>Document Number</b> <b>ENG-RPT-452</b>		
	<b>ISOS2 Functional Verification After EO Exposure</b>						<b>Revision: 001</b>		
							<b>Page 20 of 79</b>		


Low Temperature conditioning data.				Air				Air			
Date	Time	Temp	Humidity	3/17/2016	10:41:05	-40	51.8	3/17/2016	12:03:05	-40	50.3
3/17/2016	9:22:05	51.8	2.8	3/17/2016	10:42:05	-39.9	51.7	3/17/2016	12:04:05	-40	50.3
3/17/2016	9:23:05	45.6	3.5	3/17/2016	10:43:05	-40	51.7	3/17/2016	12:05:05	-40	50.3
3/17/2016	9:24:05	34.1	3.5	3/17/2016	10:44:05	-40	51.6	3/17/2016	12:06:05	-40.1	50.2
3/17/2016	9:25:05	24.3	3.1	3/17/2016	10:45:05	-40	51.6	3/17/2016	12:07:05	-40	50.2
3/17/2016	9:26:05	14.9	2.5	3/17/2016	10:46:05	-40	51.5	3/17/2016	12:08:05	-40	50.2
3/17/2016	9:27:05	8.1	2.5	3/17/2016	10:47:05	-40	51.5	3/17/2016	12:09:05	-40	50.2
3/17/2016	9:28:05	2	2.7	3/17/2016	10:48:05	-40.1	51.4	3/17/2016	12:10:05	-40	50.2
3/17/2016	9:29:05	-3.3	2.9	3/17/2016	10:49:05	-40.1	51.3	3/17/2016	12:11:05	-40	50.2
3/17/2016	9:30:05	-7.8	3.2	3/17/2016	10:50:05	-40	51.3	3/17/2016	12:12:05	-40	50.2
3/17/2016	9:31:05	-12	3.6	3/17/2016	10:51:05	-40	51.2	3/17/2016	12:13:05	-40	50.2
3/17/2016	9:32:05	-15.8	3.9	3/17/2016	10:52:05	-40	51.2	3/17/2016	12:14:05	-40	50.2
3/17/2016	9:33:05	-19.4	4.3	3/17/2016	10:53:05	-40.1	51.1	3/17/2016	12:15:05	-40.1	50.2
3/17/2016	9:34:05	-22.9	4.6	3/17/2016	10:54:05	-40	51.1	3/17/2016	12:16:05	-40	50.1
3/17/2016	9:35:05	-25.8	5	3/17/2016	10:55:05	-40.1	51	3/17/2016	12:17:05	-40	50.1
3/17/2016	9:36:05	-29	5.2	3/17/2016	10:56:05	-40.1	51	3/17/2016	12:18:05	-40	50.1
3/17/2016	9:37:05	-32	5.5	3/17/2016	10:57:05	-40.1	51	3/17/2016	12:19:05	-40	50.1
3/17/2016	9:38:05	-34.8	5.8	3/17/2016	11:00:05	-40	50.9	3/17/2016	12:22:05	-40	50.1
3/17/2016	9:39:05	-36.4	6.1	3/17/2016	11:01:05	-40	50.9	3/17/2016	12:23:05	-40.1	50
3/17/2016	9:40:05	-38.2	6.4	3/17/2016	11:02:05	-40.1	50.9	3/17/2016	12:24:05	-40	50
3/17/2016	9:41:05	-39.1	6.7	3/17/2016	11:03:05	-40	50.8	3/17/2016	12:25:05	-39.9	50
3/17/2016	9:42:05	-39.6	6.9	3/17/2016	11:04:05	-39.9	50.8	3/17/2016	12:26:05	-40	50
3/17/2016	9:43:05	-34	7.3	3/17/2016	11:05:05	-40.1	50.8	3/17/2016	12:27:05	-40	50
3/17/2016	9:44:05	-30.3	8.9	3/17/2016	11:06:05	-40	50.8	3/17/2016	12:28:05	-40	50
3/17/2016	9:45:05	-29.9	13.9	3/17/2016	11:07:05	-40	50.8	3/17/2016	12:29:05	-40	50
3/17/2016	9:46:05	-29.5	20.3	3/17/2016	11:08:05	-40	50.8	3/17/2016	12:30:05	-40	50
3/17/2016	9:47:05	-29.1	26.8	3/17/2016	11:09:05	-40	50.7	3/17/2016	12:31:05	-40	50
3/17/2016	9:48:05	-28.6	32.6	3/17/2016	11:10:05	-40	50.7	3/17/2016	12:32:05	-40	50
3/17/2016	9:49:05	-28	37.9	3/17/2016	11:11:05	-40.1	50.7	3/17/2016	12:33:05	-40	49.9
3/17/2016	9:50:05	-27.5	42.5	3/17/2016	11:12:05	-40	50.7	3/17/2016	12:34:05	-40	50
3/17/2016	9:51:05	-23.5	46.5	3/17/2016	11:13:05	-40	50.7	3/17/2016	12:35:05	-40	50
3/17/2016	9:52:05	-26.1	52.3	3/17/2016	11:14:05	-40	50.7	3/17/2016	12:36:05	-40.2	50
3/17/2016	9:53:05	-29.7	54.2	3/17/2016	11:15:05	-40	50.7	3/17/2016	12:37:05	-39.9	49.9
3/17/2016	9:54:05	-35	53.2	3/17/2016	11:16:05	-40	50.7	3/17/2016	12:38:05	-39.9	49.9
3/17/2016	9:55:05	-37.1	52.1	3/17/2016	11:17:05	-40	50.7	3/17/2016	12:39:05	-40	49.9
3/17/2016	9:56:05	-38.4	51.6	3/17/2016	11:18:05	-40.1	50.7	3/17/2016	12:40:05	-39.9	49.9
3/17/2016	9:57:05	-39.5	51.3	3/17/2016	11:19:05	-40	50.7	3/17/2016	12:41:05	-40	49.9
3/17/2016	9:58:05	-40	51.2	3/17/2016	11:20:05	-40	50.6	3/17/2016	12:42:05	-40	49.9
3/17/2016	9:59:05	-40.4	51.1	3/17/2016	11:21:05	-40.1	50.6	3/17/2016	12:43:05	-40.1	49.9
3/17/2016	10:00:05	-40.5	51	3/17/2016	11:22:05	-40.1	50.6	3/17/2016	12:44:05	-40.1	49.9
3/17/2016	10:01:05	-40.6	51	3/17/2016	11:23:05	-40	50.6	3/17/2016	12:45:05	-40.1	49.9
3/17/2016	10:02:05	-40.3	51.1	3/17/2016	11:24:05	-40.1	50.6	3/17/2016	12:46:05	-40	49.9
3/17/2016	10:03:05	-40.4	51.1	3/17/2016	11:25:05	-40	50.6	3/17/2016	12:47:05	-40.1	49.9
3/17/2016	10:04:05	-40.3	51.1	3/17/2016	11:26:05	-40.1	50.6	3/17/2016	12:48:05	-40.1	49.9
3/17/2016	10:05:05	-40.4	51.2	3/17/2016	11:27:05	-40.1	50.5	3/17/2016	12:49:05	-39.9	49.9
3/17/2016	10:06:05	-40.2	51.3	3/17/2016	11:28:05	-40	50.5	3/17/2016	12:50:05	-39.9	49.8
3/17/2016	10:07:05	-40	51.3	3/17/2016	11:29:05	-40	50.5	3/17/2016	12:51:05	-40	49.9
3/17/2016	10:08:05	-40.1	51.4	3/17/2016	11:30:05	-40.1	50.5	3/17/2016	12:52:05	-40	49.8
3/17/2016	10:09:05	-40	51.5	3/17/2016	11:31:05	-40	50.5	3/17/2016	12:53:05	-40.1	49.8
3/17/2016	10:10:05	-40	51.5	3/17/2016	11:32:05	-40	50.5	3/17/2016	12:54:05	-40	49.8
3/17/2016	10:11:05	-40.2	51.5	3/17/2016	11:33:05	-40	50.5	3/17/2016	12:55:05	-39.9	49.8
3/17/2016	10:12:05	-39.9	51.6	3/17/2016	11:34:05	-40	50.5	3/17/2016	12:56:05	-40.1	49.8
3/17/2016	10:13:05	-40.1	51.6	3/17/2016	11:35:05	-39.9	50.5	3/17/2016	12:57:05	-40	49.8
3/17/2016	10:14:05	-40	51.6	3/17/2016	11:36:05	-39.9	50.5	3/17/2016	12:58:05	-40	49.7
3/17/2016	10:15:05	-40.1	51.6	3/17/2016	11:37:05	-40	50.5	3/17/2016	12:59:05	-40	49.7
3/17/2016	10:16:05	-40	51.7	3/17/2016	11:38:05	-40	50.5	3/17/2016	13:00:05	-40	49.7
3/17/2016	10:17:05	-39.9	51.7	3/17/2016	11:39:05	-40	50.4	3/17/2016	13:01:05	-40	49.7
3/17/2016	10:18:05	-40.1	51.7	3/17/2016	11:40:05	-40	50.4	3/17/2016	13:02:05	-40	49.7
3/17/2016	10:19:05	-40.1	51.7	3/17/2016	11:41:05	-40	50.4	3/17/2016	13:03:05	-40	49.7
3/17/2016	10:20:05	-40	51.7	3/17/2016	11:42:05	-40	50.4	3/17/2016	13:04:05	-40	49.7
3/17/2016	10:21:05	-40	51.7	3/17/2016	11:43:05	-40	50.4	3/17/2016	13:05:05	-40.1	49.7
3/17/2016	10:22:05	-40	51.8	3/17/2016	11:44:05	-40.1	50.4	3/17/2016	13:06:05	-40	49.7
3/17/2016	10:23:05	-40	51.8	3/17/2016	11:45:05	-40	50.4	3/17/2016	13:07:05	-40	49.7
3/17/2016	10:24:05	-40	51.8	3/17/2016	11:46:05	-40	50.4	3/17/2016	13:08:05	-39.9	49.7
3/17/2016	10:25:05	-40.1	51.8	3/17/2016	11:47:05	-40	50.4	3/17/2016	13:09:05	-40	49.7
3/17/2016	10:26:05	-40	51.8	3/17/2016	11:48:05	-40	50.4	3/17/2016	13:10:05	-40	49.7
3/17/2016	10:27:05	-40.1	51.9	3/17/2016	11:49:05	-40.1	50.4	3/17/2016	13:11:05	-40	49.7
3/17/2016	10:28:05	-40	51.8	3/17/2016	11:50:05	-40.1	50.4	3/17/2016	13:12:05	-39.9	49.6
3/17/2016	10:29:05	-40	51.8	3/17/2016	11:51:05	-40	50.4	3/17/2016	13:13:05	-40	49.7
3/17/2016	10:30:05	-40	51.8	3/17/2016	11:52:05	-40	50.4	3/17/2016	13:14:05	-40	49.6
3/17/2016	10:31:05	-40.1	51.8	3/17/2016	11:53:05	-40	50.3	3/17/2016	13:15:05	-40.1	49.6
3/17/2016	10:32:05	-40	51.9	3/17/2016	11:54:05	-40.1	50.3	3/17/2016	13:16:05	-39.9	49.6
3/17/2016	10:33:05	-40.1	51.9	3/17/2016	11:55:05	-40	50.4	3/17/2016	13:17:05	-40	49.6
3/17/2016	10:34:05	-40	51.8	3/17/2016	11:56:05	-40	50.3	3/17/2016	13:18:05	-40	49.6
3/17/2016	10:35:05	-39.9	51.8	3/17/2016	11:57:05	-40	50.3	3/17/2016	13:19:05	-40	49.6
3/17/2016	10:36:05										

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>						<b>Document Number</b> <b>ENG-RPT-452</b>		
	<b>ISOS2 Functional Verification After EO Exposure</b>						<b>Revision: 001</b>		
							<b>Page 21 of 79</b>		

3/17/2016	13:25:05	-40	49.6		3/17/2016	13:27:05	-39.9	49.6	
3/17/2016	13:26:05	-40	49.5		3/17/2016	13:28:05	-39.9	49.5	
3/17/2016	13:30:05	-40	49.6						
3/17/2016	13:31:05	-39.9	49.6						
3/17/2016	13:32:05	-40	49.5						
3/17/2016	13:33:05	-40	49.5						
3/17/2016	13:34:05	-40	49.5						
3/17/2016	13:35:05	-40	49.5						
3/17/2016	13:36:05	-40.1	49.5						
3/17/2016	13:37:05	-40	49.5						
3/17/2016	13:38:05	-40	49.5						
3/17/2016	13:39:05	-40	49.5						
3/17/2016	13:40:05	-39.9	49.5						
3/17/2016	13:41:05	-40	49.5						
3/17/2016	13:42:05	-39.9	49.5						
3/17/2016	13:43:05	-40	49.5						
3/17/2016	13:44:05	-40	49.5						
3/17/2016	13:45:05	-40	49.4						
3/17/2016	13:46:05	-40	49.4						
3/17/2016	13:47:05	-40	49.4						
3/17/2016	13:48:05	-40	49.4						
3/17/2016	13:49:05	-40	49.4						
3/17/2016	13:50:05	-40.1	49.4						
3/17/2016	13:51:05	-40.1	49.3						
3/17/2016	13:52:05	-40	49.4						
3/17/2016	13:53:05	-39.9	49.3						
3/17/2016	13:54:05	-40	49.4						
3/17/2016	13:55:05	-40.1	49.3						
3/17/2016	13:56:05	-40	49.3						
3/17/2016	13:57:05	-40.1	49.3						
3/17/2016	13:58:05	-40	49.3						
3/17/2016	13:59:05	-40.1	49.3						
3/17/2016	14:00:05	-40	49.3						
3/17/2016	14:01:05	-40	49.3						
3/17/2016	14:02:05	-40	49.2						
3/17/2016	14:03:05	-40	49.2						
3/17/2016	14:04:05	-40.1	49.3						
3/17/2016	14:05:05	-39.9	49.2						
3/17/2016	14:06:05	-40	49.2						
3/17/2016	14:07:05	-40	49.2						
3/17/2016	14:08:05	-39.9	49.2						
3/17/2016	14:09:05	-40	49.2						
3/17/2016	14:10:05	-40	49.2						
3/17/2016	14:11:05	-40	49.2						
3/17/2016	14:12:05	-36.9	49.5						
3/17/2016	14:13:05	-24.3	51.3						
3/17/2016	14:14:05	-15.4	55.7						
3/17/2016	14:15:05	-8.1	56.9						
3/17/2016	14:16:05	-1.9	51.6						
3/17/2016	14:17:05	4	43.2						
3/17/2016	14:18:05	9.6	35.7						
3/17/2016	14:19:05	14.9	28.9						
3/17/2016	14:20:05	20.2	23.6						
3/17/2016	14:21:05	25.8	37.7						
3/17/2016	14:22:05	31.7	42.9						
3/17/2016	14:23:05	37.6	41						
3/17/2016	14:24:05	43.1	38.8						
3/17/2016	14:25:05	48.7	36.9						
3/17/2016	14:26:05	51.5	36.5						
3/17/2016	14:27:05	53.8	38						

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-452
ISOS2 Functional Verification After EO Exposure	Revision: 001	Page 22 of 79

### High Temperature Conditioning.




Megadyne Medical Products, Inc.	TEST REPORT				Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure				Revision: 001
					Page 23 of 79

Date	Time	Air Temp	Humidity				
3/21/2016	13:35:05	55	51.2	3/21/2016	19:05:05	70	45.4
3/21/2016	13:40:05	55	51	3/21/2016	19:10:05	70	37.4
3/21/2016	13:45:05	55.1	50.9	3/21/2016	19:15:05	70	28.9
3/21/2016	13:50:05	55	50.7	3/21/2016	19:20:05	70	22.1
3/21/2016	13:55:05	55.1	50.6	3/21/2016	19:25:05	70.1	19.3
3/21/2016	14:00:05	54.2	5.6	3/21/2016	19:30:05	70.1	20.1
3/21/2016	14:05:05	70.2	5.1	3/21/2016	19:35:05	69.9	20.1
3/21/2016	14:10:05	70.2	35.4	3/21/2016	19:40:05	69.9	18.7
3/21/2016	14:15:05	70.2	72.8	3/21/2016	19:45:05	69.9	15.8
3/21/2016	14:20:05	70	89.1	3/21/2016	19:50:05	70.2	14.8
3/21/2016	14:25:05	70	92.7	3/21/2016	19:55:05	70.1	14.9
3/21/2016	14:30:05	70.2	93.6	3/21/2016	20:00:05	70.2	14.9
3/21/2016	14:35:05	70.1	94.2	3/21/2016	20:05:05	70	14.9
3/21/2016	14:40:05	70.1	94.3	3/21/2016	20:10:05	70.1	14.7
3/21/2016	14:45:05	70.1	94.4	3/21/2016	20:15:05	70.1	14.4
3/21/2016	14:50:05	70.1	94.6	3/21/2016	20:20:05	70.1	14.6
3/21/2016	14:55:05	70.1	94.7	3/21/2016	20:25:05	70	15.1
3/21/2016	15:00:05	70.1	94.6	3/21/2016	20:30:05	69.9	14.8
3/21/2016	15:05:05	70.1	94.8	3/21/2016	20:35:05	70	15.1
3/21/2016	15:10:05	70.1	94.8	3/21/2016	20:40:05	70	14.9
3/21/2016	15:15:05	70.1	94.8	3/21/2016	20:45:05	70.1	15
3/21/2016	15:20:05	70.1	94.8	3/21/2016	20:50:05	69.9	14.9
3/21/2016	15:25:05	70	94.9	3/21/2016	20:55:05	70	14.9
3/21/2016	15:30:05	70	94.9	3/21/2016	21:00:05	70	15.2
3/21/2016	15:35:05	70	94.9	3/21/2016	21:05:05	70	15.2
3/21/2016	15:40:05	70.1	95	3/21/2016	21:10:05	70	14.9
3/21/2016	15:45:05	70	94.9	3/21/2016	21:15:05	70	14.9
3/21/2016	15:50:05	70	94.9	3/21/2016	21:20:05	70	15.1
3/21/2016	15:55:05	70	94.9	3/21/2016	21:25:05	70	15.2
3/21/2016	16:00:05	70	94.9	3/21/2016	21:30:05	70	15.2
3/21/2016	16:05:05	70	95	3/21/2016	21:35:05	70	15
3/21/2016	16:10:05	70	95	3/21/2016	21:40:05	70	14.9
3/21/2016	16:15:05	70	95	3/21/2016	21:45:05	69.9	14.9
3/21/2016	16:20:05	70.1	95	3/21/2016	21:50:05	70	14.9
3/21/2016	16:25:05	70	95	3/21/2016	21:55:05	70	15
3/21/2016	16:30:05	70	95	3/21/2016	22:00:05	70	15.1
3/21/2016	16:35:05	70	94.9	3/21/2016	22:05:05	70	15
3/21/2016	16:40:05	70	94.9	3/21/2016	22:10:05	70	14.8
3/21/2016	16:45:05	70	95	3/21/2016	22:15:05	70	14.9
3/21/2016	16:50:05	70	95	3/21/2016	22:20:05	70	15.1
3/21/2016	16:55:05	70	95	3/21/2016	22:25:05	70	15
3/21/2016	17:00:05	70	95.1	3/21/2016	22:30:05	70	15.2
3/21/2016	17:05:05	70	95.1	3/21/2016	22:35:05	70	15
3/21/2016	17:10:05	70	95.1	3/21/2016	22:40:05	70	15
3/21/2016	17:15:05	70	95.1	3/21/2016	22:45:05	70.1	15.2
3/21/2016	17:20:05	70	95.1	3/21/2016	22:50:05	70	14.9
3/21/2016	17:25:05	70	95	3/21/2016	22:55:05	70	15.1
3/21/2016	17:30:05	70	95.1	3/21/2016	23:00:05	70	15.1
3/21/2016	17:35:05	70	94.9	3/21/2016	23:05:05	70	14.9
3/21/2016	17:40:05	70	94.9	3/21/2016	23:10:05	70	14.9
3/21/2016	17:45:05	70	95	3/21/2016	23:15:05	70	15
3/21/2016	17:50:05	70	95	3/21/2016	23:20:05	70	15
3/21/2016	17:55:05	70	95.1	3/21/2016	23:25:05	70	14.9
3/21/2016	18:00:05	70	95.1	3/21/2016	23:30:05	70	15.1
3/21/2016	18:05:05	70	95	3/21/2016	23:35:05	70	15.1
3/21/2016	18:10:05	70	95.1	3/21/2016	23:40:05	70	15
3/21/2016	18:15:05	70	95	3/21/2016	23:45:05	70	15
3/21/2016	18:20:05	70	95.1	3/21/2016	23:50:05	70	15
3/21/2016	18:25:05	70	95.1	3/21/2016	23:55:05	70	15
3/21/2016	18:30:05	70	95.1	3/22/2016	0:00:05	61.8	13.3
3/21/2016	18:35:05	70	95.1	3/22/2016	0:05:05	49.1	16.1
3/21/2016	18:40:05	70	95.1	3/22/2016	0:10:05	40.5	17.9
3/21/2016	18:45:05	69.5	75.5	3/22/2016	0:15:05	33.4	21.1
3/21/2016	18:50:05	69.9	65.1	3/22/2016	0:20:05	27.5	28.7
3/21/2016	18:55:05	69.9	58.6	3/22/2016	0:25:05	23.3	44.4
3/21/2016	19:00:05	70	52.4	3/22/2016	0:30:05	22	55.7
				3/22/2016	0:35:05	22.6	56.1
				3/22/2016	0:40:05	22.9	49.3

Confidential - Controlled copy on date printed

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	<b>TEST REPORT</b>	Document Number <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 24 of 79</b>



<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>						<b>Document Number</b> <b>ENG-RPT-452</b>			
	<b>ISOS2 Functional Verification After EO Exposure</b>						<b>Revision: 001</b>			
							<b>Page 25 of 79</b>			

Date	Time	Air Temp	Humidity	3/22/2016	6:25:05	23	50.3	3/22/2016	13:15:05	22.9	49.7
3/22/2016	0:00:05	61.8	13.3	3/22/2016	6:30:05	23	50	3/22/2016	13:20:05	22.9	50.4
3/22/2016	0:05:05	49.1	16.1	3/22/2016	6:35:05	23	49.7	3/22/2016	13:25:05	23	49.6
3/22/2016	0:10:05	40.5	17.9	3/22/2016	6:40:05	22.9	49.9	3/22/2016	13:30:05	23	50.3
3/22/2016	0:15:05	33.4	21.1	3/22/2016	6:45:05	23	49.8	3/22/2016	13:35:05	23	50.1
3/22/2016	0:20:05	27.5	28.7	3/22/2016	7:05:05	23	49.9	3/22/2016	13:40:05	23	49.8
3/22/2016	0:25:05	23.3	44.4	3/22/2016	7:10:05	23	50.3	3/22/2016	13:45:05	23.1	50.1
3/22/2016	0:30:05	22	55.7	3/22/2016	7:20:05	23	50.2	3/22/2016	13:50:05	23	49.8
3/22/2016	0:35:05	22.6	56.1	3/22/2016	7:25:05	23	49.7	3/22/2016	14:10:05	23	49.7
3/22/2016	0:40:05	22.9	49.3	3/22/2016	7:30:05	23	50.2	3/22/2016	14:15:05	23	50.6
3/22/2016	0:45:05	23	47.5	3/22/2016	7:35:05	23	49.9	3/22/2016	14:20:05	23	49.5
3/22/2016	0:50:05	23	49.4	3/22/2016	7:40:05	23	49.9	3/22/2016	14:25:05	23	50.5
3/22/2016	0:55:05	23	49.5	3/22/2016	7:45:05	23	50.5	3/22/2016	14:30:05	23	49.6
3/22/2016	1:00:05	23	49.6	3/22/2016	7:50:05	23	49.6	3/22/2016	14:35:05	23	50.3
3/22/2016	1:05:05	23	50.1	3/22/2016	7:55:05	22.9	50.1	3/22/2016	14:40:05	23	49.8
3/22/2016	1:10:05	23	49.5	3/22/2016	8:00:05	23	50.1	3/22/2016	14:45:05	23	50.1
3/22/2016	1:15:05	22.9	50.2	3/22/2016	8:05:05	23	49.8	3/22/2016	14:50:05	23	50
3/22/2016	1:20:05	23	49.8	3/22/2016	8:10:05	23	50.4	3/22/2016	15:00:05	23.1	50.5
3/22/2016	1:25:05	23	49.9	3/22/2016	8:15:05	23	49.6	3/22/2016	15:05:05	23	49.7
3/22/2016	1:30:05	23	49.8	3/22/2016	8:20:05	23	50.1	3/22/2016	15:10:05	23	50.3
3/22/2016	1:35:05	23	50.1	3/22/2016	8:25:05	23	50	3/22/2016	15:15:05	23	49.7
3/22/2016	1:40:05	23	49.9	3/22/2016	8:30:05	23	49.9	3/22/2016	15:20:05	23.1	50.3
3/22/2016	1:45:05	23	49.9	3/22/2016	8:35:05	23	50.5	3/22/2016	15:25:05	23	49.7
3/22/2016	1:50:05	23	49.8	3/22/2016	8:40:05	23	49.6	3/22/2016	15:30:05	23	50.5
3/22/2016	1:55:05	23	50.1	3/22/2016	8:45:05	23	50.4	3/22/2016	15:35:05	23.1	49.6
3/22/2016	2:00:05	23	49.7	3/22/2016	8:50:05	23.1	49.6	3/22/2016	15:40:05	22.9	50
3/22/2016	2:05:05	23	50.3	3/22/2016	8:55:05	23	50.5	3/22/2016	15:45:05	23	50.4
3/22/2016	2:10:05	23	49.7	3/22/2016	9:00:05	23	49.5	3/22/2016	15:50:05	23	49.7
3/22/2016	2:15:05	23	50.2	3/22/2016	9:05:05	23	50.4	3/22/2016	15:55:05	23	50.4
3/22/2016	2:20:05	23	49.8	3/22/2016	9:10:05	23	49.5	3/22/2016	16:00:05	23	49.6
3/22/2016	2:25:05	23	49.8	3/22/2016	9:15:05	23	50.3	3/22/2016	16:05:05	23	50.6
3/22/2016	2:30:05	23	49.9	3/22/2016	9:20:05	23	49.5	3/22/2016	16:10:05	23	49.5
3/22/2016	2:35:05	23	49.9	3/22/2016	9:25:05	23	50.4	3/22/2016	16:15:05	23	50.2
3/22/2016	2:40:05	23	49.8	3/22/2016	9:30:05	23	49.5	3/22/2016	16:20:05	23	49.8
3/22/2016	2:45:05	23	50.6	3/22/2016	9:35:05	23	50.4	3/22/2016	16:25:05	23	49.9
3/22/2016	2:50:05	23	49.4	3/22/2016	9:40:05	23	49.5	3/22/2016	16:30:05	23	50.3
3/22/2016	2:55:05	23	50.4	3/22/2016	9:45:05	23	50.3	3/22/2016	16:35:05	23	49.6
3/22/2016	3:00:05	23	49.8	3/22/2016	9:50:05	23	49.6	3/22/2016	16:40:05	23.1	50.2
3/22/2016	3:05:05	23	49.7	3/22/2016	9:55:05	23	50.3	3/22/2016	16:45:05	23.1	49.9
3/22/2016	3:10:05	23	50.6	3/22/2016	10:00:05	23	49.7	3/22/2016	16:50:05	23	49.8
3/22/2016	3:15:05	23	49.5	3/22/2016	10:05:05	23	50.3	3/22/2016	16:55:05	23	50.2
3/22/2016	3:20:05	23	50.1	3/22/2016	10:10:05	23	49.5	3/22/2016	17:00:05	23	49.8
3/22/2016	3:25:05	23	50.2	3/22/2016	10:15:05	23	50.2	3/22/2016	17:05:05	23	50
3/22/2016	3:30:05	23	49.5	3/22/2016	10:20:05	23	49.7	3/22/2016	17:10:05	23	50.2
3/22/2016	3:35:05	23	50.3	3/22/2016	10:25:05	23	50.3	3/22/2016	17:15:05	22.9	49.9
3/22/2016	3:40:05	23	50.4	3/22/2016	10:30:05	23	49.9	3/22/2016	17:20:05	23	49.9
3/22/2016	3:45:05	23	49.3	3/22/2016	10:35:05	22.9	49.9	3/22/2016	17:25:05	23	50.3
3/22/2016	3:50:05	23	50.4	3/22/2016	10:40:05	23	50.2	3/22/2016	17:30:05	23	49.6
3/22/2016	3:55:05	22.9	50.1	3/22/2016	10:45:05	23	49.7	3/22/2016	17:35:05	23	50.6
3/22/2016	4:00:05	23	49.7	3/22/2016	10:50:05	22.9	50.7	3/22/2016	17:40:05	23	49.5
3/22/2016	4:05:05	23	50.5	3/22/2016	10:55:05	23	49.5	3/22/2016	17:45:05	23	50.4
3/22/2016	4:10:05	23	49.4	3/22/2016	11:00:05	23	50.7	3/22/2016	17:50:05	23	49.5
3/22/2016	4:15:05	23	50	3/22/2016	11:05:05	23	49.5	3/22/2016	17:55:05	23	50.3
3/22/2016	4:20:05	23	50.5	3/22/2016	11:10:05	23	50.2	3/22/2016	18:00:05	23	50
3/22/2016	4:25:05	23	49.4	3/22/2016	11:15:05	23	49.7	3/22/2016	18:05:05	23	49.9
3/22/2016	4:30:05	23	50.1	3/22/2016	11:20:05	23	50.3	3/22/2016	18:10:05	23	50.5
3/22/2016	4:35:05	23	50.3	3/22/2016	11:25:05	22.9	49.7	3/22/2016	18:15:05	23	49.5
3/22/2016	4:40:05	23	50	3/22/2016	11:30:05	23	50.2	3/22/2016	18:20:05	23	50.7
3/22/2016	4:45:05	23	49.5	3/22/2016	11:35:05	23	49.6	3/22/2016	18:25:05	23	49.3
3/22/2016	4:50:05	23	50	3/22/2016	11:40:05	23	50.1	3/22/2016	18:30:05	23	50.5
3/22/2016	4:55:05	23	50.4	3/22/2016	11:45:05	23	49.9	3/22/2016	18:35:05	23	49.6
3/22/2016	5:00:05	23	50.3	3/22/2016	11:50:05	23	49.8	3/22/2016	18:40:05	22.9	49.9
3/22/2016	5:05:05	23	49.5	3/22/2016	11:55:05	23	50.4	3/22/2016	18:45:05	23	50.6
3/22/2016	5:10:05	23	50	3/22/2016	12:00:05	23	49.7	3/22/2016	18:50:05	23	49.5
3/22/2016	5:15:05	23	50.5	3/22/2016	12:05:05	23	50.6	3/22/2016	18:55:05	23	50.8
3/22/2016	5:20:05	23	50.2	3/22/2016	12:10:05	23	49.6	3/22/2016	19:00:05	23	49.5
3/22/2016	5:25:05	23	49.6	3/22/2016	12:15:05	23	50.7	3/22/2016	19:05:05	23	50.5
3/22/2016	5:30:05	23	49.7	3/22/2016	12:20:05	23	49.6	3/22/2016	19:10:05	22.9	49.4
3/22/2016	5:35:05	23	49.9	3/22/2016	12:25:05	23	50.7	3/22/2016	19:15:05	23	50.2
3/22/2016	5:40:05	23	50.2	3/22/2016	12:30:05	23	49.4	3/22/2016	19:20:05	23	50.3
3/22/2016	5:45:05	23	50.5	3/22/2016	12:35:05	23	50.5	3/22/2016	19:25:05	23	49.6
3/22/2016	5:50:05	23	50.1	3/22/2016	12:40:05	23	49.5	3/22/2016	19:30:05	23	50.7
3/22/2016	5:55:05	23	49.5	3/22/2016	12:45:05	23	50.2	3/22/2016	19:35:05	23	49.3
3/22/2016	6:00:05	23	49.9	3/22/2016	12:50:05	23	49.7	3/22/2016	19:40:05	23	50.3
3/22/2016	6:05:05	23	50	3/22/2016	12:55:05	23	50.1	3/22/2016	19:45:05	23	49.7
3/22/2016	6:10:05	22.9	50.1	3/22/2016	13:00:05	23	49.9	3/22/2016	19:50:05	23	50
3/22/2016	6:15:05	23.1	50.4	3/22/2016	13:05:05	23	49.8	3/22/2016	19:55:05	22.9	50.5
3/22/2016	6:20:05	23	50.2	3/22/2016	13:10:05	23	50.2	3/22/2016	20:		

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>							<b>Document Number</b> <b>ENG-RPT-452</b>		
	<b>ISOS2 Functional Verification After EO Exposure</b>							<b>Revision: 001</b>		
								<b>Page 26 of 79</b>		

3/22/2016	20:05:05	23	50.6	3/23/2016	2:55:05	23	49.9	3/23/2016	9:45:05	23	50.1
3/22/2016	20:10:05	23	49.6	3/23/2016	3:00:05	23	49.7	3/23/2016	9:50:05	23	49.7
3/22/2016	20:15:05	23	49.9	3/23/2016	3:05:05	23	49.8	3/23/2016	9:55:05	23	49.8
3/22/2016	20:20:05	23	50.5	3/23/2016	3:10:05	23	49.7	3/23/2016	10:00:05	23	50.2
3/22/2016	20:25:05	23	49.4	3/23/2016	3:15:05	23	50	3/23/2016	10:05:05	23	50.3
3/22/2016	20:30:05	23	50.4	3/23/2016	3:20:05	23	50.3	3/23/2016	10:10:05	23	50.1
3/22/2016	20:35:05	23	49.7	3/23/2016	3:25:05	23	50.3	3/23/2016	10:15:05	23	50.1
3/22/2016	20:40:05	23	50	3/23/2016	3:30:05	23	49.9	3/23/2016	10:20:05	23	49.9
3/22/2016	20:45:05	23	50.3	3/23/2016	3:35:05	23	50	3/23/2016	10:25:05	23	50
3/22/2016	20:50:05	23	49.6	3/23/2016	3:40:05	23	49.9	3/23/2016	10:30:05	23	50
3/22/2016	20:55:05	23	50.5	3/23/2016	3:45:05	23	50.1	3/23/2016	10:35:05	23	50.1
3/22/2016	21:00:05	22.9	49.4	3/23/2016	3:50:05	23	50.1	3/23/2016	10:40:05	23	50.1
3/22/2016	21:05:05	23	50.5	3/23/2016	3:55:05	23	49.9	3/23/2016	10:45:05	23	49.8
3/22/2016	21:10:05	23	49.7	3/23/2016	4:00:05	22.9	50.1	3/23/2016	10:50:05	23	50
3/22/2016	21:15:05	23	50	3/23/2016	4:05:05	23	50.2	3/23/2016	10:55:05	23	50.2
3/22/2016	21:20:05	23	50.2	3/23/2016	4:10:05	23	50.3	3/23/2016	11:00:05	23	50
3/22/2016	21:25:05	23	49.7	3/23/2016	4:15:05	23	50.2	3/23/2016	11:05:05	23	50
3/22/2016	21:30:05	23	50.7	3/23/2016	4:20:05	23.1	50.2	3/23/2016	11:10:05	23	50.2
3/22/2016	21:35:05	23	49.4	3/23/2016	4:25:05	23	50.1	3/23/2016	11:15:05	23	49.9
3/22/2016	21:40:05	23	50.4	3/23/2016	4:30:05	23	50.1	3/23/2016	11:20:05	23	49.8
3/22/2016	21:45:05	23	49.6	3/23/2016	4:35:05	23.1	50.2	3/23/2016	11:25:05	22.9	49.9
3/22/2016	21:50:05	23	50	3/23/2016	4:40:05	23	50.2	3/23/2016	11:30:05	23	50
3/22/2016	21:55:05	23	50.5	3/23/2016	4:45:05	23	50.1	3/23/2016	11:35:05	23	49.9
3/22/2016	22:00:05	23	49.5	3/23/2016	4:50:05	23	50.2	3/23/2016	11:40:05	23	50.1
3/22/2016	22:05:05	23	50.1	3/23/2016	4:55:05	23	50.1	3/23/2016	11:45:05	22.9	50.1
3/22/2016	22:10:05	23	50.1	3/23/2016	5:00:05	23	49.7	3/23/2016	11:50:05	23.1	49.9
3/22/2016	22:15:05	23	49.6	3/23/2016	5:05:05	23	49.9	3/23/2016	11:55:05	22.9	49.8
3/22/2016	22:20:05	23	50.5	3/23/2016	5:10:05	23	50	3/23/2016	12:00:05	22.9	49.8
3/22/2016	22:25:05	23	49.4	3/23/2016	5:15:05	23	50.2	3/23/2016	12:05:05	23	49.9
3/22/2016	22:30:05	23	50.2	3/23/2016	5:20:05	23	50.1	3/23/2016	12:10:05	23	50
3/22/2016	22:35:05	22.9	50.3	3/23/2016	5:25:05	23	50.1	3/23/2016	12:15:05	22.9	49.9
3/22/2016	22:40:05	22.9	49.6	3/23/2016	5:30:05	23	50.1	3/23/2016	12:20:05	23	49.9
3/22/2016	22:45:05	23	50.6	3/23/2016	5:35:05	23	50.3	3/23/2016	12:25:05	23.1	49.8
3/22/2016	22:50:05	23	49.5	3/23/2016	5:40:05	23	50	3/23/2016	12:30:05	23	49.7
3/22/2016	22:55:05	23	49.9	3/23/2016	5:45:05	23	50	3/23/2016	12:35:05	23	50
3/22/2016	23:00:05	23	50.5	3/23/2016	5:50:05	23	50	3/23/2016	12:40:05	23	50.2
3/22/2016	23:05:05	23	49.5	3/23/2016	5:55:05	23	49.8	3/23/2016	12:45:05	23	50.3
3/22/2016	23:10:05	23	50.1	3/23/2016	6:00:05	23	50.1	3/23/2016	12:50:05	23	50.1
3/22/2016	23:15:05	23	50.5	3/23/2016	6:05:05	23	50	3/23/2016	12:55:05	23	50.1
3/22/2016	23:20:05	23	49.5	3/23/2016	6:10:05	23	50	3/23/2016	13:00:05	23	50
3/22/2016	23:25:05	23	50.3	3/23/2016	6:15:05	23	50.1	3/23/2016	13:05:05	23	50
3/22/2016	23:30:05	23	50	3/23/2016	6:20:05	23	50.3	3/23/2016	13:10:05	23	49.9
3/22/2016	23:35:05	23	49.6	3/23/2016	6:25:05	23	50	3/23/2016	13:15:05	23	49.9
3/22/2016	23:40:05	23	50.2	3/23/2016	6:30:05	23	49.9	3/23/2016	13:20:05	23	49.9
3/22/2016	23:45:05	23	50.3	3/23/2016	6:35:05	22.9	50.1	3/23/2016	13:25:05	23	49.9
3/22/2016	23:50:05	23	49.8	3/23/2016	6:40:05	23	50.1	3/23/2016	13:30:05	23	50
3/22/2016	23:55:05	23	49.7	3/23/2016	6:45:05	23	50	3/23/2016	13:35:05	23	50.1
3/23/2016	0:00:05	23	50	3/23/2016	6:50:05	23.1	50.1	3/23/2016	13:40:05	23	50
3/23/2016	0:05:05	23	50.4	3/23/2016	6:55:05	22.9	50.1	3/23/2016	13:45:05	23	49.8
3/23/2016	0:10:05	23	50.2	3/23/2016	7:00:05	23	49.7	3/23/2016	13:50:05	23	49.9
3/23/2016	0:15:05	23	49.5	3/23/2016	7:05:05	23	49.8	3/23/2016	13:55:05	23	49.9
3/23/2016	0:20:05	23	49.7	3/23/2016	7:10:05	22.9	50.1	3/23/2016	14:00:05	23	50
3/23/2016	0:25:05	23	50.5	3/23/2016	7:15:05	23	50.1	3/23/2016	14:05:05	23	49.8
3/23/2016	0:30:05	23	50	3/23/2016	7:20:05	23	50.2	3/23/2016	14:10:05	23	49.9
3/23/2016	0:35:05	23	49.7	3/23/2016	7:25:05	23	50.1	3/23/2016	14:15:05	23	50.1
3/23/2016	0:40:05	23	49.7	3/23/2016	7:30:05	22.9	50.1	3/23/2016	14:20:05	23	50.3
3/23/2016	0:45:05	23	50	3/23/2016	7:35:05	23	49.9	3/23/2016	14:25:05	23	50.2
3/23/2016	0:50:05	22.9	50.3	3/23/2016	7:40:05	23	50	3/23/2016	14:30:05	22.9	50.1
3/23/2016	0:55:05	23.1	50.5	3/23/2016	7:45:05	23	50.1	3/23/2016	14:35:05	23	49.8
3/23/2016	1:00:05	23	49.9	3/23/2016	7:50:05	23	50.1	3/23/2016	14:40:05	23	49.9
3/23/2016	1:05:05	23	49.6	3/23/2016	7:55:05	23	49.9	3/23/2016	14:45:05	23	49.8
3/23/2016	1:10:05	23	50	3/23/2016	8:00:05	23	50	3/23/2016	14:50:05	23	50.1
3/23/2016	1:15:05	23	50.4	3/23/2016	8:05:05	23	50	3/23/2016	14:55:05	23	50.2
3/23/2016	1:20:05	23	50.4	3/23/2016	8:10:05	23	49.9	3/23/2016	15:00:05	23	50.2
3/23/2016	1:25:05	23	49.9	3/23/2016	8:15:05	23	50.1	3/23/2016	15:05:05	23	50.1
3/23/2016	1:30:05	23	49.7	3/23/2016	8:20:05	23	50.2	3/23/2016	15:10:05	23	50
3/23/2016	1:35:05	23	49.8	3/23/2016	8:25:05	23.1	50.2	3/23/2016	15:15:05	23	50
3/23/2016	1:40:05	23	50.3	3/23/2016	8:30:05	23	50.1	3/23/2016	15:20:05	23	50.1
3/23/2016	1:45:05	23	50.5	3/23/2016	8:35:05	23	49.9	3/23/2016	15:25:05	23	50.1
3/23/2016	1:50:05	23	49.8	3/23/2016	8:40:05	23	49.7	3/23/2016	15:30:05	23	50.2
3/23/2016	1:55:05	23	49.7	3/23/2016	8:45:05	23	49.9	3/23/2016	15:35:05	23	50
3/23/2016	2:00:05	23	49.6	3/23/2016	8:50:05	23	50.1	3/23/2016	15:40:05	23	50.1
3/23/2016	2:05:05	23	50	3/23/2016	8:55:05	23	50.1	3/23/2016	15:45:05	23	50
3/23/2016	2:10:05	23	50.1	3/23/2016	9:00:05	23	50.3	3/23/2016	15:50:05	23	50
3/23/2016	2:15:05	23	50.4	3/23/2016	9:05:05	23	50.1	3/23/2016	15:55:05	23	49.9
3/23/2016	2:20:05	23	50.3	3/23/2016	9:10:05	23	49.9	3/23/2016	16:00:05	23	49.7
3/23/2016	2:25:05	23	50	3/23/2016	9:15:05	23	49.7	3/23/2016	16:05:05	23	50
3/23/2016	2:30:05	23	49.6	3/23/2016	9:20:05	23	50.1	3/23/2016	16:10:05	23	50.4
3/23/2016	2:3										

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>							<b>Document Number</b> <b>ENG-RPT-452</b>			
	<b>ISOS2 Functional Verification After EO Exposure</b>							<b>Revision: 001</b>			
								<b>Page 27 of 79</b>			

3/23/2016	16:35:05	23	49.9	3/23/2016	23:25:05	23	50.3	3/24/2016	6:15:05	22.9	49.7
3/23/2016	16:40:05	23	50	3/23/2016	23:30:05	23	50.2	3/24/2016	6:20:05	23	50.1
3/23/2016	16:45:05	23	49.9	3/23/2016	23:35:05	23	50.1	3/24/2016	6:25:05	23	50.2
3/23/2016	16:50:05	23	50.1	3/23/2016	23:40:05	23	50.3	3/24/2016	6:30:05	23	50.2
3/23/2016	16:55:05	23.1	50.1	3/23/2016	23:45:05	23	50	3/24/2016	6:35:05	23	50.2
3/23/2016	17:00:05	23	50.1	3/23/2016	23:50:05	23	49.8	3/24/2016	6:40:05	23	49.9
3/23/2016	17:05:05	23	50.2	3/23/2016	23:55:05	23	50.1	3/24/2016	6:45:05	23	50
3/23/2016	17:10:05	23	50.1	3/24/2016	0:00:05	23	50	3/24/2016	6:50:05	23	50.2
3/23/2016	17:15:05	23	49.7	3/24/2016	0:05:05	23	49.7	3/24/2016	6:55:05	23	50.1
3/23/2016	17:20:05	22.9	50	3/24/2016	0:10:05	23	49.9	3/24/2016	7:00:05	23	50
3/23/2016	17:25:05	22.9	50.1	3/24/2016	0:15:05	23	50	3/24/2016	7:05:05	23	49.9
3/23/2016	17:30:05	23	50.3	3/24/2016	0:20:05	23	50	3/24/2016	7:10:05	23	50
3/23/2016	17:35:05	23	50	3/24/2016	0:25:05	23	49.5	3/24/2016	7:15:05	23.1	50
3/23/2016	17:40:05	23	50.1	3/24/2016	0:30:05	23	50	3/24/2016	7:20:05	23	50.2
3/23/2016	17:45:05	23	50	3/24/2016	0:35:05	23	50.1	3/24/2016	7:25:05	23	50.1
3/23/2016	17:50:05	23	50.1	3/24/2016	0:40:05	23	50.2	3/24/2016	7:30:05	23	49.9
3/23/2016	17:55:05	23	50	3/24/2016	0:45:05	23	50.3	3/24/2016	7:35:05	23	49.8
3/23/2016	18:00:05	23	50.1	3/24/2016	0:50:05	23	50	3/24/2016	7:40:05	23	49.9
3/23/2016	18:05:05	23	50	3/24/2016	0:55:05	23	49.9	3/24/2016	7:45:05	23	49.9
3/23/2016	18:10:05	23	50.1	3/24/2016	1:00:05	23	50.4	3/24/2016	7:50:05	23	49.9
3/23/2016	18:15:05	23	50.1	3/24/2016	1:05:05	23	50.3	3/24/2016	7:55:05	23	49.9
3/23/2016	18:20:05	23	50.3	3/24/2016	1:10:05	23	49.9	3/24/2016	8:00:05	23	49.9
3/23/2016	18:25:05	23	50	3/24/2016	1:15:05	23	49.7	3/24/2016	8:05:05	23	50.1
3/23/2016	18:30:05	23	49.7	3/24/2016	1:20:05	23	49.8	3/24/2016	8:10:05	23	50
3/23/2016	18:35:05	23	50	3/24/2016	1:25:05	23	49.9	3/24/2016	8:15:05	23	50.2
3/23/2016	18:40:05	23	49.9	3/24/2016	1:30:05	23	49.7	3/24/2016	8:20:05	22.9	50.3
3/23/2016	18:45:05	23	49.9	3/24/2016	1:35:05	23	49.8	3/24/2016	8:25:05	23	50.1
3/23/2016	18:50:05	23	49.9	3/24/2016	1:40:05	23	49.9	3/24/2016	8:30:05	23	50
3/23/2016	18:55:05	23	49.7	3/24/2016	1:45:05	23	50.1	3/24/2016	8:35:05	23	50.2
3/23/2016	19:00:05	23	49.7	3/24/2016	1:50:05	23	49.9	3/24/2016	8:40:05	23	50.2
3/23/2016	19:05:05	22.9	50	3/24/2016	1:55:05	23	50	3/24/2016	8:45:05	23	50.1
3/23/2016	19:10:05	23	50	3/24/2016	2:00:05	22.9	50.2	3/24/2016	8:50:05	23	49.9
3/23/2016	19:15:05	23	50.4	3/24/2016	2:05:05	23	50.2	3/24/2016	8:55:05	23	49.7
3/23/2016	19:20:05	23	50.1	3/24/2016	2:10:05	23	50.2	3/24/2016	9:00:05	23.1	49.9
3/23/2016	19:25:05	23.1	50.2	3/24/2016	2:15:05	23	50.1	3/24/2016	9:05:05	23	50.2
3/23/2016	19:30:05	23	50.2	3/24/2016	2:20:05	23	50.1	3/24/2016	9:10:05	23	49.9
3/23/2016	19:35:05	23	50	3/24/2016	2:25:05	23	50.1	3/24/2016	9:15:05	23	49.8
3/23/2016	19:40:05	23	49.7	3/24/2016	2:30:05	23	50	3/24/2016	9:20:05	23	49.9
3/23/2016	19:45:05	23	49.8	3/24/2016	2:35:05	23	49.9	3/24/2016	9:25:05	23	50.1
3/23/2016	19:50:05	23	50.3	3/24/2016	2:40:05	23	49.8	3/24/2016	9:30:05	23	49.9
3/23/2016	19:55:05	23.1	50.2	3/24/2016	2:45:05	23	49.8	3/24/2016	9:35:05	23	50
3/23/2016	20:00:05	23	50.2	3/24/2016	2:50:05	23.1	50.1	3/24/2016	9:40:05	23	49.9
3/23/2016	20:05:05	23	50.2	3/24/2016	2:55:05	22.9	49.9	3/24/2016	9:45:05	23	50
3/23/2016	20:10:05	23	49.9	3/24/2016	3:00:05	23	49.8	3/24/2016	9:50:05	23	50.1
3/23/2016	20:15:05	23	50.2	3/24/2016	3:05:05	23	50	3/24/2016	9:55:05	23	50
3/23/2016	20:20:05	23	50.3	3/24/2016	3:10:05	23	49.9	3/24/2016	10:00:05	23	50
3/23/2016	20:25:05	23	50.1	3/24/2016	3:15:05	23	50	3/24/2016	10:05:05	23.1	50.1
3/23/2016	20:30:05	23	50.3	3/24/2016	3:20:05	23	49.9	3/24/2016	10:10:05	23	50
3/23/2016	20:35:05	23	50.2	3/24/2016	3:25:05	23	49.9	3/24/2016	10:15:05	23	50
3/23/2016	20:40:05	23	50.2	3/24/2016	3:30:05	23	49.8	3/24/2016	10:20:05	23	50.1
3/23/2016	20:45:05	23	50.3	3/24/2016	3:35:05	22.9	49.8	3/24/2016	10:25:05	23	50.2
3/23/2016	20:50:05	23	50.1	3/24/2016	3:40:05	23	49.9	3/24/2016	10:30:05	23	50
3/23/2016	20:55:05	23	49.9	3/24/2016	3:45:05	23	49.7	3/24/2016	10:35:05	23	50
3/23/2016	21:00:05	23.1	50.1	3/24/2016	3:50:05	23	50	3/24/2016	10:40:05	23	50.1
3/23/2016	21:05:05	23	50.2	3/24/2016	3:55:05	23	50	3/24/2016	10:45:05	23.1	49.9
3/23/2016	21:10:05	23	50.1	3/24/2016	4:00:05	23	50	3/24/2016	10:50:05	23	49.7
3/23/2016	21:15:05	23	49.6	3/24/2016	4:05:05	23	50.1	3/24/2016	10:55:05	23	49.9
3/23/2016	21:20:05	23	49.9	3/24/2016	4:10:05	23	50	3/24/2016	11:00:05	23	49.6
3/23/2016	21:25:05	23	49.8	3/24/2016	4:15:05	23	49.7	3/24/2016	11:05:05	23	49.8
3/23/2016	21:30:05	23	49.8	3/24/2016	4:20:05	23	49.7	3/24/2016	11:10:05	23	50.3
3/23/2016	21:35:05	23.1	49.9	3/24/2016	4:25:05	23	49.7	3/24/2016	11:15:05	23.1	50.3
3/23/2016	21:40:05	22.9	50.4	3/24/2016	4:30:05	23	49.8	3/24/2016	11:20:05	23	50.2
3/23/2016	21:45:05	23.1	50.2	3/24/2016	4:35:05	23	50.1	3/24/2016	11:25:05	23	49.7
3/23/2016	21:50:05	23	50.3	3/24/2016	4:40:05	23	50.4	3/24/2016	11:30:05	23	49.8
3/23/2016	21:55:05	23	50.1	3/24/2016	4:45:05	23	50.2	3/24/2016	11:35:05	23	49.8
3/23/2016	22:00:05	23	50.1	3/24/2016	4:50:05	23	50.2	3/24/2016	11:40:05	23	49.8
3/23/2016	22:05:05	23	50.4	3/24/2016	4:55:05	23	50.2	3/24/2016	11:45:05	23	49.9
3/23/2016	22:10:05	23	50.1	3/24/2016	5:00:05	23.1	50.2	3/24/2016	11:50:05	22.9	49.9
3/23/2016	22:15:05	23	49.7	3/24/2016	5:05:05	23	50	3/24/2016	11:55:05	23	50.1
3/23/2016	22:20:05	23	49.6	3/24/2016	5:10:05	23	49.7	3/24/2016	12:00:05	23	50.3
3/23/2016	22:25:05	23	49.9	3/24/2016	5:15:05	23	49.9	3/24/2016	12:05:05	23	50.2
3/23/2016	22:30:05	23	50.2	3/24/2016	5:20:05	22.9	49.7	3/24/2016	12:10:05	23	50
3/23/2016	22:35:05	23	50.3	3/24/2016	5:25:05	23	49.8	3/24/2016	12:15:05	22.9	49.9
3/23/2016	22:40:05	23	50.3	3/24/2016	5:30:05	23	50	3/24/2016	12:20:05	23	49.9
3/23/2016	22:45:05	23	50.1	3/24/2016	5:35:05	22.9	50	3/24/2016	12:25:05	23	49.7
3/23/2016	22:50:05	23	49.7	3/24/2016	5:40:05	23	49.9	3/24/2016	12:30:05	23.1	49.7
3/23/2016	22:55:05	23	50	3/24/2016	5:45:05	23	50	3/24/2016	12:35:05	23.1	49.9
3/23/2016	23:00:05	23	49.8	3/24/2016	5:50:05	23.1	50	3/24/2016	12:40:05	23</td	

Megadyne Medical Products, Inc.	TEST REPORT							Document Number <b>ENG-RPT-452</b>			
	ISOS2 Functional Verification After EO Exposure							Revision: 001			
								Page 28 of 79			

3/24/2016	13:05:05	22.9	50.3	3/24/2016	18:05:05	23	49.8	3/24/2016	23:05:05	23	50.2
3/24/2016	13:10:05	23	50.3	3/24/2016	18:10:05	23	49.7	3/24/2016	23:10:05	23.1	50.2
3/24/2016	13:15:05	23	50.2	3/24/2016	18:15:05	23	49.8	3/24/2016	23:15:05	23	50.1
3/24/2016	13:20:05	23	50.3	3/24/2016	18:20:05	23	49.8	3/24/2016	23:20:05	23	50.1
3/24/2016	13:25:05	23	49.8	3/24/2016	18:25:05	22.9	49.9	3/24/2016	23:25:05	23	50
3/24/2016	13:30:05	23	49.7	3/24/2016	18:30:05	23	50	3/24/2016	23:30:05	23	50
3/24/2016	13:35:05	23	49.8	3/24/2016	18:35:05	23	49.9	3/24/2016	23:35:05	23.1	50
3/24/2016	13:40:05	23	49.9	3/24/2016	18:40:05	23.1	50.1	3/24/2016	23:40:05	23.1	50.1
3/24/2016	13:45:05	23	49.9	3/24/2016	18:45:05	23	50.4	3/24/2016	23:45:05	23	50.2
3/24/2016	13:50:05	23	50.1	3/24/2016	18:50:05	23	50	3/24/2016	23:50:05	23	50.1
3/24/2016	13:55:05	23	50.1	3/24/2016	18:55:05	23	49.7	3/24/2016	23:55:05	23	50.2
3/24/2016	14:00:05	23	50.1	3/24/2016	19:00:05	23	49.9	3/25/2016	0:00:05	23	50.2
3/24/2016	14:05:05	23	50	3/24/2016	19:05:05	23	50.3	3/25/2016	0:05:05	23	50.1
3/24/2016	14:10:05	23	49.9	3/24/2016	19:10:05	23	50.3	3/25/2016	0:10:05	23	49.8
3/24/2016	14:15:05	23	49.9	3/24/2016	19:15:05	23	50.2	3/25/2016	0:15:05	23	49.8
3/24/2016	14:20:05	23	50	3/24/2016	19:20:05	23	50.2	3/25/2016	0:20:05	23	49.9
3/24/2016	14:25:05	22.9	50	3/24/2016	19:25:05	23.1	49.9	3/25/2016	0:25:05	23	49.7
3/24/2016	14:30:05	23	50.1	3/24/2016	19:30:05	23	49.7	3/25/2016	0:30:05	23	49.7
3/24/2016	14:35:05	23	50	3/24/2016	19:35:05	23	49.7	3/25/2016	0:35:05	23	49.9
3/24/2016	14:40:05	23	49.8	3/24/2016	19:40:05	23.1	50.2	3/25/2016	0:40:05	23	50
3/24/2016	14:45:05	23	50.2	3/24/2016	19:45:05	23	50.2	3/25/2016	0:45:05	23.1	50.1
3/24/2016	14:50:05	23	50.1	3/24/2016	19:50:05	23	50.3	3/25/2016	0:50:05	23.2	50.8
3/24/2016	14:55:05	23	50.4	3/24/2016	19:55:05	23.1	49.9	3/25/2016	0:55:05	23.3	51.5
3/24/2016	15:00:05	23	50	3/24/2016	20:00:05	23	49.7	3/25/2016	1:00:05	23.3	52.1
3/24/2016	15:05:05	23	49.7	3/24/2016	20:05:05	23	49.9	3/25/2016	1:05:05	23.3	52.8
3/24/2016	15:10:05	23	49.8	3/24/2016	20:10:05	23	50.1	3/25/2016	1:10:05	23.3	53.5
3/24/2016	15:15:05	23	50	3/24/2016	20:15:05	23	50.1	3/25/2016	1:15:05	23.2	54.4
3/24/2016	15:20:05	23	49.9	3/24/2016	20:20:05	23	50.1	3/25/2016	1:20:05	23.2	55.2
3/24/2016	15:25:05	23	49.9	3/24/2016	20:25:05	23	50.2	3/25/2016	1:25:05	23.1	55.9
3/24/2016	15:30:05	23	49.9	3/24/2016	20:30:05	23.1	50.3	3/25/2016	1:30:05	23.1	56.5
3/24/2016	15:35:05	23	50.3	3/24/2016	20:35:05	23	50.1	3/25/2016	1:35:05	23	56.9
3/24/2016	15:40:05	23	50.3	3/24/2016	20:40:05	23	50.1	3/25/2016	1:40:05	23	57.2
3/24/2016	15:45:05	23	50.1	3/24/2016	20:45:05	23	49.9	3/25/2016	1:45:05	22.9	57.6
3/24/2016	15:50:05	23	50.1	3/24/2016	20:50:05	23	50	3/25/2016	1:50:05	22.9	57.9
3/24/2016	15:55:05	23	49.8	3/24/2016	20:55:05	23	50.1	3/25/2016	1:55:05	22.8	58.3
3/24/2016	16:00:05	23	50.2	3/24/2016	21:00:05	23	50.3	3/25/2016	2:00:05	22.7	58.7
3/24/2016	16:05:05	23	49.8	3/24/2016	21:05:05	23	50.3	3/25/2016	2:05:05	22.6	59.2
3/24/2016	16:10:05	23	49.7	3/24/2016	21:10:05	23	50.2	3/25/2016	2:10:05	22.6	59.5
3/24/2016	16:15:05	23	49.8	3/24/2016	21:15:05	23	50.2	3/25/2016	2:15:05	22.5	59.8
3/24/2016	16:20:05	23	50	3/24/2016	21:20:05	23	50.4	3/25/2016	2:20:05	22.5	60.1
3/24/2016	16:25:05	23	50	3/24/2016	21:25:05	23	50	3/25/2016	2:25:05	22.4	60.5
3/24/2016	16:30:05	23	50.1	3/24/2016	21:30:05	23	49.7	3/25/2016	2:30:05	22.3	60.5
3/24/2016	16:35:05	23	50.3	3/24/2016	21:35:05	23	49.8	3/25/2016	2:35:05	22.3	61.1
3/24/2016	16:40:05	23	50.3	3/24/2016	21:40:05	23	49.9	3/25/2016	2:40:05	22.2	61.2
3/24/2016	16:45:05	23.1	49.7	3/24/2016	21:45:05	23	49.9	3/25/2016	2:45:05	22.1	61.5
3/24/2016	16:50:05	23	49.9	3/24/2016	21:50:05	23	49.9	3/25/2016	2:50:05	22.1	61.7
3/24/2016	16:55:05	23	49.8	3/24/2016	21:55:05	23	50	3/25/2016	2:55:05	22	61.7
3/24/2016	17:00:05	23	49.6	3/24/2016	22:00:05	23	50	3/25/2016	3:00:05	21.9	61.9
3/24/2016	17:05:05	22.9	50	3/24/2016	22:05:05	23.1	49.9	3/25/2016	3:05:05	21.9	62.3
3/24/2016	17:10:05	23	50	3/24/2016	22:10:05	23	50	3/25/2016	3:10:05	21.8	62.7
3/24/2016	17:15:05	23	50.1	3/24/2016	22:15:05	23	50	3/25/2016	3:15:05	21.7	62.8
3/24/2016	17:20:05	22.9	50.2	3/24/2016	22:20:05	23	50	3/25/2016	3:20:05	21.7	62.9
3/24/2016	17:25:05	23	50.3	3/24/2016	22:25:05	23	50	3/25/2016	3:25:05	21.6	63.2
3/24/2016	17:30:05	23	50	3/24/2016	22:30:05	22.9	50	3/25/2016	3:30:05	21.6	63.5
3/24/2016	17:35:05	23	50.1	3/24/2016	22:35:05	23	50	3/25/2016	3:35:05	21.5	63.5
3/24/2016	17:40:05	23	50.2	3/24/2016	22:40:05	22.9	49.9	3/25/2016	3:40:05	21.4	63.7
3/24/2016	17:45:05	23	49.9	3/24/2016	22:45:05	23	50.1	3/25/2016	3:45:05	21.4	63.9
3/24/2016	17:50:05	23	49.9	3/24/2016	22:50:05	23	50.4	3/25/2016	3:50:05	21.3	64.1
3/24/2016	17:55:05	23	49.6	3/24/2016	22:55:05	23	50.2				
3/24/2016	18:00:05	23	49.8	3/24/2016	23:00:05	23	50.1				

Confidential - Controlled copy on date printed

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001 Page 29 of 79

**Appendix I**  
**Shipping Test Log Sheet**

Preconditioning:

Start Date: 3-21-2016 Chamber Number: 01268  
Completion Date: 3-25-2016 Last Calibration: 5-19-2015  
Signature/Date: Paul Valprede 3-25-16 Calibration due: 5-31-2016

Drop Test:

Catalog ACE3714 Weight 12 lbs. Drop Height: 15 inches

Drop	Orientation	Specific face, edge or	Initials/Date
1	Top	Face 1	PV 3-25-2016
2	Edge	Edge 5-3	PV 3-25-2016
3	Edge	Edge 6-3	PV 3-25-2016
4	Corner	Corner 2-3-5	PV 3-25-2016
5	Corner	Corner 4-3-6	PV 3-25-2016
6	Bottom	Face 3	PV 3-25-2016

Comments: \_\_\_\_\_

Signature: Paul Valprede Date: 3-25-2016

Compression Test:

Catalog ACE3714 Pounds Force 220 lbs.

Comments: LOT # 160868

Signature: Paul Valprede Date: 3-25-2016



3-28-2016

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 30 of 79

**Appendix I Continued**  
**Shipping Test Log Sheet**

Vibration:

Low Frequency, 40 minutes, Initials PV High frequency 10 minutes, Initials PV

Completion Date: 3-25-2016

Signature: Paul Valpreda Date: 3-25-2016

Second Drop Test:

Catalog AZE3714 Weight 12 lbs. Drop Height: 15 + 30 inches.

Drop	Orientation	Specific face, edge or	Initials/Date
1	Edge	Edge 4-6	PV 3-25-2016
2	Face	Face 4	PV 3-25-2016
3	Face	Face 6	PV 3-25-2016
4	Corner	Corner 2-1-5	PV 3-25-2016
5	Edge	Edge 2-1	PV 3-25-2016
6	Bottom	Face 3, Increase height to 2 times specified height.	PV 3-25-2016

Comments: \_\_\_\_\_

Signature: Paul Valpreda Date: 3-25-2016.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 31 of 79</b>

## **APPENDIX II – HANDPIECE HIGH FREQUENCY DIELECTRIC WITHSTAND TESTING**

Sample	Configuration	High Frequency Max Value (kHz)	Mains P/F	Button Activation Ω P/F
1	0036 - 151363	6.88	PASS	PASS
2	0036 - 151363	6.60	PASS	PASS
3	0036 - 151363	6.72	PASS	PASS
4	0036 - 151363	6.76	PASS	8.6, 1.8 - PASS
5	0036 - 151363	6.60	PASS	8.6, 2.7 - PASS
6	0036 - 151363	6.60	PASS	1.7, 1.5 - PASS
7	0036 - 151363	6.60	PASS	2.3, 4.6 - PASS
8	0036 - 151363	6.88	PASS	1.9, 1.5 - PASS
9	0036 - 151363	6.72	PASS	1.7, 3.7 - PASS
10	0036 - 151363	6.72	PASS	PASS
11	0036 - 151363	6.60	PASS	4.5, 4.5 - PASS
12	0036 - 151363	6.68	PASS	10.5, 3.1 - PASS
13	0036 - 151363	6.60	PASS	PASS
14	0036 - 151363	6.60	PASS	2.0, 1.5 - PASS
15	0036 - 151363	6.60	PASS	1.6, 1.4 - PASS
16	0036 - 151363	6.88	PASS	4.5, 4.2 - PASS
17	0036 - 151363	6.72	PASS	1.6, 1.5 - PASS
18	0036 - 151363	6.60	PASS	1.5, 1.6 - PASS
19	0036 - 151363	6.68	PASS	3.2, 5.5 - PASS
20	0036 - 151363	6.64	PASS	1.7, 1.5 - PASS
21	0036 - 151363	6.60	PASS	1.7, 3.6 - PASS
22	0036 - 151363	6.72	PASS	1.6, 1.5 - PASS
23	0036 - 151363	6.60	PASS	1.6, 1.9 - PASS
24	0036 - 151363	6.64	PASS	13.7, 18.5 - PASS
25	0036 - 151363	6.64	PASS	14.0, 24.6 - PASS
26	0036 - 151363	6.60	PASS	1.6, 8.1 - PASS
27	0036 - 151363	6.64	PASS	1.5, 1.4 - PASS
28	0036 - 151363	6.76	PASS	2.0, 11.4 - PASS
29	0036 - 151363	6.72	PASS	5.6, 1.9 - PASS
30	0036 - 151363	6.80	PASS	5.8, 4.3 - PASS

CALIBRATION INFORMATION	
Fluke 79 Series II Multimeter	Multimeter
Serial Number: 57040475	Serial Number: 01010
Megadyne: N/A	Calibration Date: 4/10/2015
Calibration Due: 4/30/2016	Calibration Due: 4/30/2016
Oscilloscope	Generator
Tektronix TDS 3012B	Mega Power 1000
Serial Number: B010635	Serial Number: 10333001
Megadyne Number: 01142	Megadyne Number: N/A
Calibration Date: 10/16/2014	Calibration Due: 10/13/2015
Hipot Test Generator	Hipotronics Model HD 100 Series
Megadyne Number: 01037	Calibration Due: 8/26/2014
Calibration Due: 8/31/2015	Calibration Due: 8/31/2015
High Voltage Probe	Tektronix P6015A High Voltage Probe
Serial Number: B043063	Serial Number: 01138
Megadyne Number: 01138	Calibration Date: 8/6/2014
Calibration Date: 8/6/2014	Calibration Due: 8/31/2015
Calibration Due: 8/31/2015	Inductive Current Coil
Pearson Current Monitor, Model 2700	RMS Voltmeter
Serial Number: 010955	Fluke 89204 True RMS Voltmeter
Megadyne Number: 01288	Serial Number: 4220005
Calibration Date: 1/14/2015	Megadyne Number: 01254
Calibration Due: 1/14/2016	Calibration Date: 10/16/2014
Calibration Due: 10/16/2015	Calibration Due: 10/16/2015

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT			Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure			Revision: 001
				Page 32 of 79

**High Frequency Dielectric Withstand Testing**  
**ISOS2 ACE37H - LOT 151361 EOx2**  
**Data Collection Form**

CALIBRATION INFORMATION				
		Fluke 79 Series II Multimeter		
		Serial Number: 157040475		
		Megadyne Number: 01010		
		Calibration Date: 4/10/2015		
		Calibration Due: 4/30/2016		
		Generator		
		Mega Power 1000		
		Serial Number: 10353001		
		Megadyne Number: N/A		
		Calibration Date: N/A		
		Oscilloscope		
		Tektronix TDS 3012B		
		Serial Number: B010635		
		Megadyne Number: 01142		
		Calibration Date: 10/16/2014		
		Calibration Due: 10/31/2015		
		Hirox Test Generator		
		Hirox Model HD 100 Series		
		Serial Number: 010357		
		Megadyne Number: 01142		
		Calibration Date: 8/25/2015		
		Calibration Due: 8/31/2016		
		High Voltage Probe		
		Tektronix P6015A High Voltage Probe		
		Serial Number: B043063		
		Megadyne Number: 01138		
		Calibration Date: 8/6/2014		
		Calibration Due: 8/31/2015		
		Inductive Current Coil		
		Pearson Current Monitor, Model 2100		
		Serial Number: 106055		
		Megadyne Number: 01288		
		Calibration Date: 1/14/2015		
		Calibration Due: 1/14/2016		
		RMS Voltmeter		
		Fluke 8020A True RMS Voltmeter		
		Serial Number: 4220005		
		Megadyne Number: 01254		
		Calibration Date: 10/16/2014		
		Calibration Due: 10/16/2015		

Rated Accessory V <sub>peak</sub>	5,500	V <sub>peak</sub>	PV
Minimum Mains Test Value	4,596	V <sub>RM</sub>	PV
Actual Mains Test Value	4,600	V <sub>RM</sub>	PV
Mains Test Value Calculation:	(Rated Accessory V <sub>peak</sub> + 1000V <sub>peak</sub> ) / v2		

Paul Valpreda  
 Operator(s) Name

  
 Operator(s) Signature

9/4/2015  
 Date

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 33 of 79

## APPENDIX II – CABLE LEAKAGE TEST

### Cable Dielectric Withstand Testing

Cable with new inner jacket EOx2 on 0036 pencils LOT 151363

Data Collection Form

Sample	Configuration	d (mm)	L (cm)	f <sub>test</sub> (kHz)	U <sub>peak</sub> (V <sub>p-p</sub> ) <sup>2</sup>	LEAKAGE		
						Measured Leakage	Calculated Acceptable Leakage	P/F
1	New Inner Jacket Cable EOx2	2.8	30	400	400	35.1	120.96	PASS 14
2	New Inner Jacket Cable EOx2	2.8	30	400	400	35.3	120.96	PASS 14
3	New Inner Jacket Cable EOx2	2.8	30	400	400	37.1	120.96	PASS 14
4	New Inner Jacket Cable EOx2	2.8	30	400	400	38.6	120.96	PASS 15
5	New Inner Jacket Cable EOx2	2.8	30	400	400	35.5	120.96	PASS 14
6	New Inner Jacket Cable EOx2	2.8	30	400	400	34.7	120.96	PASS 14
7	New Inner Jacket Cable EOx2	2.8	30	400	400	37.7	120.96	PASS 15
8	New Inner Jacket Cable EOx2	2.8	30	400	400	33.8	120.96	PASS 14
9	New Inner Jacket Cable EOx2	2.8	30	400	400	33.5	120.96	PASS 14
10	New Inner Jacket Cable EOx2	2.8	30	400	400	39.5	120.96	PASS 16
11	New Inner Jacket Cable EOx2	2.8	30	400	400	33.8	120.96	PASS 14
12	New Inner Jacket Cable EOx2	2.8	30	400	400	36.5	120.96	PASS 14
13	New Inner Jacket Cable EOx2	2.8	30	400	400	33.3	120.96	PASS 14
14	New Inner Jacket Cable EOx2	2.8	30	400	400	35.8	120.96	PASS 14
15	New Inner Jacket Cable EOx2	2.8	30	400	400	34.6	120.96	PASS 15
16	New Inner Jacket Cable EOx2	2.8	30	400	400	35.9	120.96	PASS 15
17	New Inner Jacket Cable EOx2	2.8	30	400	400	34.5	120.96	PASS 14
18	New Inner Jacket Cable EOx2	2.8	30	400	400	34.6	120.96	PASS 14
19	New Inner Jacket Cable EOx2	2.8	30	400	400	38.8	120.96	PASS 15
20	New Inner Jacket Cable EOx2	2.8	30	400	400	34.6	120.96	PASS 14
21	New Inner Jacket Cable EOx2	2.8	30	400	400	33.4	120.96	PASS 14
22	New Inner Jacket Cable EOx2	2.8	30	400	400	33.2	120.96	PASS 14
23	New Inner Jacket Cable EOx2	2.8	30	400	400	35.6	120.96	PASS 15
24	New Inner Jacket Cable EOx2	2.8	30	400	400	34.3	120.96	PASS 14
25	New Inner Jacket Cable EOx2	2.8	30	400	400	37.2	120.96	PASS 15
26	New Inner Jacket Cable EOx2	2.8	30	400	400	37.7	120.96	PASS 15
27	New Inner Jacket Cable EOx2	2.8	30	400	400	37.8	120.96	PASS 15
28	New Inner Jacket Cable EOx2	2.8	30	400	400	33.7	120.96	PASS 14
29	New Inner Jacket Cable EOx2	2.8	30	400	400	37.9	120.96	PASS 15
30	New Inner Jacket Cable EOx2	2.8	30	400	400	39.2	120.96	PASS 16

7/22/2015

PV

d = smallest measured outer dimension of cable insulation

L = length of cable insulation up to a max. of 300 mm

f<sub>test</sub> = frequency of the pure cut signal

U<sub>peak</sub> = Oscilloscope test p-p voltage

Leakage Calculation -

$$I_{\text{leakage}} = 9.0 \times 10^{-6} \times d \times L \times f_{\text{test}} \times U_{\text{peak}} (\text{mA})$$

Leakage Calculation -

$$I_{\text{leakage}} = 9.0 \times 10^{-6} \times d \times L \times f_{\text{test}} \times U_{\text{peak}} (\text{mA})$$

Sample	Configuration	HIGH FREQUENCY		MAINS	
		Max Vpk (kV)	P/F	P/F	P/F
1	New Inner Jacket Cable EOx2	6.72	PASS	PASS	
2	New Inner Jacket Cable EOx2	6.68	PASS	PASS	
3	New Inner Jacket Cable EOx2	6.72	PASS	PASS	
4	New Inner Jacket Cable EOx2	6.76	PASS	PASS	
5	New Inner Jacket Cable EOx2	6.72	PASS	PASS	
6	New Inner Jacket Cable EOx2	6.72	PASS	PASS	
7	New Inner Jacket Cable EOx2	6.68	PASS	PASS	
8	New Inner Jacket Cable EOx2	6.76	PASS	PASS	
9	New Inner Jacket Cable EOx2	6.72	PASS	PASS	
10	New Inner Jacket Cable EOx2	6.68	PASS	PASS	
11	New Inner Jacket Cable EOx2	6.64	PASS	PASS	
12	New Inner Jacket Cable EOx2	6.60	PASS	PASS	
13	New Inner Jacket Cable EOx2	6.68	PASS	PASS	
14	New Inner Jacket Cable EOx2	6.68	PASS	PASS	
15	New Inner Jacket Cable EOx2	6.76	PASS	PASS	
16	New Inner Jacket Cable EOx2	6.68	PASS	PASS	
17	New Inner Jacket Cable EOx2	6.64	PASS	PASS	
18	New Inner Jacket Cable EOx2	6.68	PASS	PASS	
19	New Inner Jacket Cable EOx2	6.72	PASS	PASS	
20	New Inner Jacket Cable EOx2	6.72	PASS	PASS	
21	New Inner Jacket Cable EOx2	6.68	PASS	PASS	
22	New Inner Jacket Cable EOx2	6.68	PASS	PASS	
23	New Inner Jacket Cable EOx2	6.60	PASS	PASS	
24	New Inner Jacket Cable EOx2	6.72	PASS	PASS	
25	New Inner Jacket Cable EOx2	6.68	PASS	PASS	
26	New Inner Jacket Cable EOx2	6.64	PASS	PASS	
27	New Inner Jacket Cable EOx2	6.72	PASS	PASS	
28	New Inner Jacket Cable EOx2	6.76	PASS	PASS	
29	New Inner Jacket Cable EOx2	6.72	PASS	PASS	
30	New Inner Jacket Cable EOx2	6.64	PASS	PASS	

7/23/2015 9/3/2015

PV PV

Mains Test Value Calculation:	
$(\text{Rated Accessory } V_{\text{peak}} + 1000V_{\text{peak}}) / \sqrt{2}$	
Rated Accessory $V_{\text{peak}}$	5,900
Minimum Mains Test Value	4,596
Actual Mains Test Value	4,600

CALIBRATION INFORMATION			
Multimeter			
Fluke 79 Series II Multimeter			
Serial Number:	57040475	Megadyne Number:	01010
Calibration Date:	4/10/2015	Calibration Due:	4/30/2016
Generator			
Mega Power 1000			
Serial Number:	10353001	Megadyne Number:	N/A
Calibration Date:	N/A	Calibration Due:	N/A
Oscilloscope			
Tektronix TDS 3012B			
Serial Number:	B010635	Megadyne Number:	01142
Calibration Date:	10/16/2014	Calibration Due:	10/31/2015
Hitpot Test Generator			
Hoptronics Model HD 100 Series			
Megadyne Number:	01037	Calibration Date:	8/26/2014
Calibration Date:	8/31/2015	Calibration Due:	8/31/2015
High Voltage Probe			
Tektronix P6015A High Voltage Probe			
Serial Number:	B043063	Megadyne Number:	01138
Calibration Date:	8/6/2014	Calibration Due:	8/31/2015
Inductive Current Coil			
Pearson Current Monitor, Model 2100			
Serial Number:	109055	Megadyne Number:	01288
Calibration Date:	1/14/2015	Calibration Due:	1/14/2016
RMS Voltmeter			
Fluke 8920A True RMS Voltmeter			
Serial Number:	4222005	Megadyne Number:	01254
Calibration Date:	10/16/2014	Calibration Due:	10/16/2015

Paul Valpreda  
Operator Name

9/4/2015  
Date



W/16/15  
Date

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 34 of 79

### APPENDIX III – PLUG MAINS DIELECTRIC

#### Plug Dielectric Data collection Form

Sample	Configuration	MAINS
		P/ F
1	0036 - 151363 Rocker EOx2	PASS
2	0036 - 151363 Rocker EOx2	PASS
3	0036 - 151363 Rocker EOx2	PASS
4	0036 - 151363 Rocker EOx2	PASS
5	0036 - 151363 Rocker EOx2	PASS
6	0036 - 151363 Rocker EOx2	PASS
7	0036 - 151363 Rocker EOx2	PASS
8	0036 - 151363 Rocker EOx2	PASS
9	0036 - 151363 Rocker EOx2	PASS
10	0036 - 151363 Rocker EOx2	PASS
11	0036 - 151363 Rocker EOx2	PASS
12	0036 - 151363 Rocker EOx2	PASS
13	0036 - 151363 Rocker EOx2	PASS
14	0036 - 151363 Rocker EOx2	PASS
15	0036 - 151363 Rocker EOx2	PASS
16	0036 - 151363 Rocker EOx2	PASS
17	0036 - 151363 Rocker EOx2	PASS
18	0036 - 151363 Rocker EOx2	PASS
19	0036 - 151363 Rocker EOx2	PASS
20	0036 - 151363 Rocker EOx2	PASS
21	0036 - 151363 Rocker EOx2	PASS
22	0036 - 151363 Rocker EOx2	PASS
23	0036 - 151363 Rocker EOx2	PASS
24	0036 - 151363 Rocker EOx2	PASS
25	0036 - 151363 Rocker EOx2	PASS
26	0036 - 151363 Rocker EOx2	PASS
27	0036 - 151363 Rocker EOx2	PASS
28	0036 - 151363 Rocker EOx2	PASS
29	0036 - 151363 Rocker EOx2	PASS
30	0036 - 151363 Rocker EOx2	PASS

9/18/2015

PV

Mains Test Value Calculation:		
(Rated Accessory V <sub>peak</sub> +1000V <sub>peak</sub> ) / √2		
Rated Accessory V <sub>peak</sub>	5,500	V <sub>peak</sub>
Minimum Mains Test Value	4,596	V <sub>RMS</sub>
Actual Mains Test Value	4,600	V <sub>RMS</sub>

CALIBRATION INFORMATION
Multimeter
Fluke 79 Series II Multimeter
Serial Number: 57040475
Megadyne Number: 01010
Calibration Date: 4/10/2015
Calibration Due: 4/30/2016
Generator
Mega Power 1000
Serial Number: 10353001
Megadyne Number: N/A
Calibration Date: N/A
Calibration Due: N/A
Oscilloscope
Tektronix TDS 3012B
Serial Number: B010635
Megadyne Number: 01142
Calibration Date: 10/16/2014
Calibration Due: 10/31/2015
Hipot Test Generator
Hipotronics Model HD 100 Series
Megadyne Number: 01037
Calibration Date: 8/25/2015
Calibration Due: 8/31/2016
High Voltage Probe
Tektronix P6015A High Voltage Probe
Serial Number: B043063
Megadyne Number: 01138
Calibration Date: 8/6/2014
Calibration Due: 8/31/2015
Inductive Current Coil
Pearson Current Monitor, Model 2100
Serial Number: 109055
Megadyne Number: 01288
Calibration Date: 1/14/2015
Calibration Due: 1/14/2016
RMS Voltmeter
Fluke 8920A True RMS Voltmeter
Serial Number: 4220005
Megadyne Number: 01254
Calibration Date: 10/16/2014
Calibration Due: 10/16/2015

Paul Valpreda  
Operator Name

9/18/2015  
Date



Operator Signature



Date

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 35 of 79

## APPENDIX IV – FLUID INGRESS

Ingress of Fluids			
Data Collection Form			
Product Tested - 0036 LOT 151363 EOx2			
Sample	Configuration	INITIAL CONTINUITY CUT, COAG ( $\Omega$ )	FLUID INGRESS CUT & COAG - P/F
1	0036 - 151363	6.0, 2.3	PASS
2	0036 - 151363	11.4, 19.2	PASS
3	0036 - 151363	3.5, 9.7	PASS
4	0036 - 151363	1.5, 4.1	PASS
5	0036 - 151363	11.2, 8.4	PASS
6	0036 - 151363	7.3, 8.3	PASS
7	0036 - 151363	3.4, 8.7	PASS
8	0036 - 151363	4.3, 13.7	PASS
9	0036 - 151363	5.7, 8.1	PASS
10	0036 - 151363	17.1, 15.4	PASS
11	0036 - 151363	3.3, 9.6	PASS
12	0036 - 151363	3.0, 2.6	PASS
13	0036 - 151363	3.5, 2.8	PASS
14	0036 - 151363	4.9, 12.0	PASS
15	0036 - 151363	3.9, 5.3	PASS
16	0036 - 151363	7.7, 5.4	PASS
17	0036 - 151363	2.1, 2.2	PASS
18	0036 - 151363	3.7, 2.8	PASS
19	0036 - 151363	3.7, 13.9	PASS
20	0036 - 151363	7.1, 6.4	PASS
21	0036 - 151363	1.6, 29.8	PASS
22	0036 - 151363	2.5, 1.5	PASS
23	0036 - 151363	5.9, 2.6	PASS
24	0036 - 151363	5.8, 3.6	PASS
25	0036 - 151363	3.2, 16.9	PASS
26	0036 - 151363	5.4, 8.2	PASS
27	0036 - 151363	2.1, 10.9	PASS
28	0036 - 151363	2.5, 6.3	PASS
29	0036 - 151363	5.3, 24.3	PASS
30	0036 - 151363	2.5, 7.9	PASS

10/19/2015

PV

Paul Valpreda  
Operator Name(s) 10/19/2015  
Date Completed


  
Operator Signature 11/16/15  
Date Completed

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 36 of 79

Ingress of Fluids			
Data Collection Form			
Product Tested - ACE37H - 160868 EOx2			
Sample	Configuration	INITIAL CONTINUITY	FLUID INGRESS
		CUT, COAG ( $\Omega$ )	CUT & COAG - P/F
1	ACE37H - 160868	1.6, 1.6	PASS
2	ACE37H - 160868	1.6, 1.4	PASS
3	ACE37H - 160868	1.7, 1.7	PASS
4	ACE37H - 160868	1.6, 1.5	PASS
5	ACE37H - 160868	1.8, 1.6	PASS
6	ACE37H - 160868	1.8, 1.6	PASS
7	ACE37H - 160868	1.8, 2.8	PASS
8	ACE37H - 160868	1.5, 1.7	PASS
9	ACE37H - 160868	6.2, 2.2	PASS
10	ACE37H - 160868	2.0, 2.0	PASS
11	ACE37H - 160868	1.9, 2.1	PASS
12	ACE37H - 160868	3.2, 1.8	PASS
13	ACE37H - 160868	3.1, 10.0	PASS
14	ACE37H - 160868	4.0, 3.8	PASS
15	ACE37H - 160868	2.7, 6.3	PASS
16	ACE37H - 160868	1.8, 1.8	PASS
17	ACE37H - 160868	1.9, 2.3	PASS
18	ACE37H - 160868	1.6, 1.6	PASS
19	ACE37H - 160868	2.9, 1.6	PASS
20	ACE37H - 160868	1.8, 6.7	PASS
21	ACE37H - 160868	1.9, 7.4	PASS
22	ACE37H - 160868	1.7, 2.3	PASS
23	ACE37H - 160868	7.7, 3.0	PASS
24	ACE37H - 160868	2.4, 2.4	PASS
25	ACE37H - 160868	1.5, 1.6	PASS
26	ACE37H - 160868	1.4, 1.4	PASS
27	ACE37H - 160868	1.5, 2.6	PASS
28	ACE37H - 160868	3.2, 2.0	PASS
29	ACE37H - 160868	2.2, 1.6	PASS
30	ACE37H - 160868	2.0, 1.5	PASS

No failures were observed - PV

Paul Valpreda - PV  
Operator Name(s) 3/28/2016  
Date Completed  
  
Engineer Signature Date 3/28/2016

CALIBRATION INFORMATION	
Multimeter	Fluke 179 True RMS Multimeter
Serial Number:	33660168
Megadyne Number:	01507
Calibration Date:	1/29/2016
Calibration Due:	1/31/2017
Oscilloscope	Tektronix TDS 3012B
Serial Number:	B010635
Megadyne Number:	01142
Calibration Date:	10/12/2015
Calibration Due:	10/31/2016
Inductive Current Coil	Pearson Current Monitor, Model 2100
Serial Number:	109055
Megadyne Number:	01288
Calibration Date:	1/25/2016
Calibration Due:	1/31/2017
Waveform/Function Generator	BK Precision 4084AWG
Serial Number:	40849136407060000
Megadyne Number:	01458
Calibration Date:	10/20/2015
Calibration Due:	10/20/2016

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001 Page 37 of 79

Ingress of Fluids			
Data Collection Form			
Product Tested - ACE37H - 151361 EOx2			
Sample	Configuration	INITIAL CONTINUITY	FLUID INGRESS
		CUT, COAG (Ω)	CUT & COAG - P/F
1	ACE37H - 151361	3.8, 15.5	PASS
2	ACE37H - 151361	15.5, 6.5	PASS
3	ACE37H - 151361	2.7, 16.5	PASS
4	ACE37H - 151361	2.0, 19.5	PASS
5	ACE37H - 151361	5.5, 12.6	PASS
6	ACE37H - 151361	8.5, 7.5	PASS
7	ACE37H - 151361	5.5, 1.5	PASS
8	ACE37H - 151361	20.4, 30.9	PASS
9	ACE37H - 151361	2.5, 2.6	FAIL
10	ACE37H - 151361	8.1, 6.2	PASS
11	ACE37H - 151361	6.4, 1.3	PASS
12	ACE37H - 151361	10.1, 6.9	PASS
13	ACE37H - 151361	3.5, 6.5	PASS
14	ACE37H - 151361	12.5, 21.5	PASS
15	ACE37H - 151361	3.4, 12.9	PASS
16	ACE37H - 151361	9.4, 6.8	PASS
17	ACE37H - 151361	2.3, 4.2	PASS
18	ACE37H - 151361	3.3, 2.4	PASS
19	ACE37H - 151361	7.7, 2.4	PASS
20	ACE37H - 151361	20.4, 3.1	PASS
21	ACE37H - 151361	2.5, 20.6	PASS
22	ACE37H - 151361	2.3, 4.4	PASS
23	ACE37H - 151361	1.6, 1.5	PASS
24	ACE37H - 151361	2.5, 2.9	PASS
25	ACE37H - 151361	26.0, 2.7	PASS
26	ACE37H - 151361	4.7, 5.6	PASS
27	ACE37H - 151361	4.9, 6.1	PASS
28	ACE37H - 151361	1.2, 4.4	PASS
29	ACE37H - 151361	27.3, 21.4	PASS
30	ACE37H - 151361	4.9, 1.2	PASS

10/21/2015

PV

Paul Valpreda  
Operator Name(s) 10/21/2015  
  
Date Completed  
  
Date Completed

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

CALIBRATION INFORMATION	
Multimeter	
Fluke 79 Series II Multimeter	
Serial Number:	57040475
Megadyne Number:	01010
Calibration Date:	4/10/2015
Calibration Due:	4/30/2016
Oscilloscope	
Tektronix TDS 3012B	
Serial Number:	B010635
Megadyne Number:	01142
Calibration Date:	10/12/2015
Calibration Due:	10/31/2016
Inductive Current Coil	
Pearson Current Monitor, Model 2100	
Serial Number:	109055
Megadyne Number:	01288
Calibration Date:	1/14/2015
Calibration Due:	1/14/2016
RMS Voltmeter	
Fluke 8920A True RMS Voltmeter	
Serial Number:	4220005
Megadyne Number:	01254
Calibration Date:	10/20/2015
Calibration Due:	10/20/2016
Waveform/Function Generator	
BK Precision 4084AWG	
Serial Number:	40849136407060000
Megadyne Number:	01458
Calibration Date:	10/6/2015
Calibration Due:	10/6/2016

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	Revision: 001
		Page 38 of 79


## APPENDIX V – IEC 60601-1 TEMPERATURE TEST

*John S.M. CMC*  
2016-01-28

Button - Cut 300W - Pre

Analysis Label	Min	Max	Max - Min	Avg	Stdev
Image ARO1	20.7	25.5	4.8		
	22.2	22.8	0.7	22.5	0.2


Image Name Value  
Date 10/7/2015  
Time 12:44:02 PM  
File name 3a-b-pure cut-300 w.jpg



Button - Cut 300W - Post

Analysis Label	Min	Max	Max - Min	Avg	Stdev
Image ARO1	21.2	27.2	6		
	22.3	24.9	2.6	23.5	0.8


Image Name Value  
Date 10/7/2015  
Time 1:45:58 PM  
File name 3b-b-pure cut-300 w.jpg



Button - Coag 120W - Pre

Analysis Label	Min	Max	Max - Min	Avg	Stdev
Image ARO1	21.2	24.2	2.9		
	22.6	23.2	0.6	23	0.1

Image Name Value  
Date 10/8/2015  
Time 5:25:17 AM  
File name 4a-b-std coag-120 w.jpg



Button - Coag 120W - Post

Analysis Label	Min	Max	Max - Min	Avg	Stdev
Image ARO1	20.8	23.7	3		
	21.6	22.7	1.1	22.1	0.3





Image Name Value  
Date 10/8/2015  
Time 6:30:20 AM  
File name 4b-b-std coag-120 w.jpg



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 39 of 79</b>

*John M. Clark*  
2016-01-28



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 40 of 79

## APPENDIX VI – HANDPIECE STRAIN RELIEF (ROTATIONAL) TEST

PENCIL		
ROTATION ON 0036 - 151363, ISOS2 - 200 CYCLES		
Sample Number	Rotation CUT, COAG Ω P/F	Configuration
1	1.9, 4.6 - PASS	0036 - 151363
2	3.5, 2.5 - PASS	0036 - 151363
3	1.6, 5.9 - PASS	0036 - 151363
4	2.1, 5.9 - PASS	0036 - 151363
5	1.5, 8.7 - PASS	0036 - 151363
6	3.5, 3.9 - PASS	0036 - 151363
7	8.9, 2.6 - PASS	0036 - 151363
8	6.0, 1.8 - PASS	0036 - 151363
9	2.7, 21.0 - PASS	0036 - 151363
10	1.8, 2.0 - PASS	0036 - 151363
11	3.1, 18.6 - PASS	0036 - 151363
12	2.2, 2.6 - PASS	0036 - 151363
13	3.6, 6.1 - PASS	0036 - 151363
14	3.0, 3.5 - PASS	0036 - 151363
15	6.1, 4.3 - PASS	0036 - 151363
16	2.3, 1.9 - PASS	0036 - 151363
17	1.6, 1.9 - PASS	0036 - 151363
18	2.6, 9.1 - PASS	0036 - 151363
19	2.4, 1.8 - PASS	0036 - 151363
20	1.9, 9.0 - PASS	0036 - 151363
21	1.4, 7.1 - PASS	0036 - 151363
22	1.7, 1.5 - PASS	0036 - 151363
23	5.1, 2.2 - PASS	0036 - 151363
24	2.5, 13.6 - PASS	0036 - 151363
25	2.8, 6.1 - PASS	0036 - 151363
26	1.4, 4.0 - PASS	0036 - 151363
27	2.2, 2.7 - PASS	0036 - 151363
28	2.2, 2.2 - PASS	0036 - 151363
29	4.6, 6.9 - PASS	0036 - 151363
30	1.4, 1.8 - PASS	0036 - 151363
Mass: 80 g		Date: 10/6/2015
Resistance Pass Value (Ω): <50		

PENCIL		
ROTATION ON ACE37H - 151361, ISOS2 - 200 CYCLES		
Sample Number	Rotation CUT, COAG Ω P/F	Configuration
1	7.8, 3.3 - PASS	ACE37H - 151361
2	2.4, 1.5 - PASS	ACE37H - 151361
3	4.1, 1.2 - PASS	ACE37H - 151361
4	1.8, 5.5 - PASS	ACE37H - 151361
5	1.7, 1.5 - PASS	ACE37H - 151361
6	2.9, 1.4 - PASS	ACE37H - 151361
7	2.3, 1.5 - PASS	ACE37H - 151361
8	6.5, 16.7 - PASS	ACE37H - 151361
9	5.0, 2.3 - PASS	ACE37H - 151361
10	1.7, 4.9 - PASS	ACE37H - 151361
11	1.2, 1.1 - PASS	ACE37H - 151361
12	1.2, 9.2 - PASS	ACE37H - 151361
13	2.0, 1.9 - PASS	ACE37H - 151361
14	2.9, 6.4 - PASS	ACE37H - 151361
15	2.7, 3.2 - PASS	ACE37H - 151361
16	1.4, 1.4 - PASS	ACE37H - 151361
17	1.9, 6.3 - PASS	ACE37H - 151361
18	3.9, 2.2 - PASS	ACE37H - 151361
19	1.7, 1.0 - PASS	ACE37H - 151361
20	1.6, 4.2 - PASS	ACE37H - 151361
21	1.6, 1.7 - PASS	ACE37H - 151361
22	2.9, 1.9 - PASS	ACE37H - 151361
23	1.4, 4.5 - PASS	ACE37H - 151361
24	10.8, 1.6 - PASS	ACE37H - 151361
25	2.4, 6.9 - PASS	ACE37H - 151361
26	3.9, 1.5 - PASS	ACE37H - 151361
27	2.0, 2.1 - PASS	ACE37H - 151361
28	2.0, 2.7 - PASS	ACE37H - 151361
29	2.1, 2.0 - PASS	ACE37H - 151361
30	1.5, 1.2 - PASS	ACE37H - 151361
Mass: 80 g		Date: 10/2/2015
Resistance Pass Value (Ω): <50		

CALIBRATION INFORMATION		
Brass Weights		
N/A		
Serial Number	56206	
Megadyne Number	1085	
Calibration Date	1/30/2012	
Calibration Due	1/31/2022	

CALIBRATION INFORMATION		
Brass Weights		
N/A		
Serial Number	56206	
Megadyne Number	1085	
Calibration Date	1/30/2012	
Calibration Due	1/31/2022	

Paul Valpreda  
OPERATOR NAME  
DATE

  
11/6/15

Paul Valpreda  
OPERATOR NAME  
DATE

  
11/5/16

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.


Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 41 of 79

## APPENDIX VII – ACTIVATION FORCE TESTING


Wednesday, October 07, 2015

Button Force on ISOS2 0036 -  
151363.is\_comp

Specimen 1 to 24




Specimen 25 to 48



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 42 of 79



	Compressive load at Preset Point (Cursor ) (lbf)	Button Type	Note
1	1.00349	CUT	Sample 1
2	1.01155	COAG	Sample 1
3	1.07061	CUT	Sample 2
4	0.95786	COAG	Sample 2
5	1.00618	CUT	Sample 3
6	0.97719	COAG	Sample 3
7	1.02605	CUT	Sample 4
8	0.92108	COAG	Sample 4
9	0.99974	CUT	Sample 5
10	0.86980	COAG	Sample 5
11	1.06470	CUT	Sample 6
12	0.84162	COAG	Sample 6
13	1.16349	CUT	Sample 7
14	0.85370	COAG	Sample 7
15	1.05772	CUT	Sample 8
16	0.90282	COAG	Sample 8
17	1.10121	CUT	Sample 9
18	0.94175	COAG	Sample 9

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 43 of 79</b>

	Compressive load at Preset Point (Cursor ) (lbf)	Button Type	Note
19	1.07168	CUT	Sample 10
20	0.87329	COAG	Sample 10
21	0.98417	CUT	Sample 11
22	0.92162	COAG	Sample 11
23	1.07544	CUT	Sample 12
24	0.87947	COAG	Sample 12
25	1.11571	CUT	Sample 13
26	0.88860	COAG	Sample 13
27	1.07544	CUT	Sample 14
28	0.99813	COAG	Sample 14
29	1.08242	CUT	Sample 15
30	0.90121	COAG	Sample 15
31	1.08672	CUT	Sample 16
32	0.91652	COAG	Sample 16
33	1.15544	CUT	Sample 17
34	0.91759	COAG	Sample 17
35	1.06417	CUT	Sample 18
36	0.95303	COAG	Sample 18
37	1.06524	CUT	Sample 19
38	0.95813	COAG	Sample 19
39	1.07652	CUT	Sample 20
40	1.00833	COAG	Sample 20
41	1.04860	CUT	Sample 21
42	0.93209	COAG	Sample 21
43	1.11732	CUT	Sample 22
44	0.88457	COAG	Sample 22
45	1.12806	CUT	Sample 23
46	0.92698	COAG	Sample 23
47	1.05504	CUT	Sample 24
48	0.94631	COAG	Sample 24
49	1.04162	CUT	Sample 25
50	0.87437	COAG	Sample 25
51	1.07598	CUT	Sample 26
52	0.92027	COAG	Sample 26
53	1.06739	CUT	Sample 27
54	0.95410	COAG	Sample 27
55	1.04430	CUT	Sample 28

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 44 of 79</b>


	Compressive load at Preset Point (Cursor ) (lbf)	Button Type	Note
56	0.94927	COAG	Sample 28
57	1.11786	CUT	Sample 29
58	0.88860	COAG	Sample 29
59	1.03356	CUT	Sample 30
60	0.85316	COAG	Sample 30
Maximum	1.16349		
Minimum	0.84162		
Mean	0.99498		
Standard Deviation	0.09		

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	Revision: 001
		Page 45 of 79


Wednesday, October 14, 2015

Button Force on ISOS2 0036 - 151363  
post 500.is\_comp

Specimen 1 to 24




Specimen 25 to 48



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 46 of 79



	Compressive load at Preset Point (Cursor ) (lbf)	Button Type	Note
1	0.91303	CUT	Sample 1, post activation
2	0.89584	COAG	Sample 1, post activation
3	0.99490	CUT	Sample 2, post activation
4	0.98900	COAG	Sample 2, post activation
5	0.98631	CUT	Sample 3, post activation
6	0.93987	COAG	Sample 3, post activation
7	0.93396	CUT	Sample 4, post activation
8	0.86900	COAG	Sample 4, post activation
9	1.02443	CUT	Sample 5, post activation
10	0.89450	COAG	Sample 5, post activation
11	1.02014	CUT	Sample 6, post activation
12	0.84698	COAG	Sample 6, post activation
13	1.04860	CUT	Sample 7, post activation
14	0.81772	COAG	Sample 7, post activation
15	0.98148	CUT	Sample 8, post activation
16	0.86819	COAG	Sample 8, post activation
17	1.01692	CUT	Sample 9, post activation
18	0.84376	COAG	Sample 9, post activation

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 47 of 79

	Compressive load at Preset Point (Cursor ) (lbf)	Button Type	Note
19	0.98524	CUT	Sample 10, post activation
20	0.86766	COAG	Sample 10, post activation
21	0.98954	CUT	Sample 11, post activation
22	0.87141	COAG	Sample 11, post activation
23	0.98578	CUT	Sample 12, post activation
24	0.84430	COAG	Sample 12, post activation
25	1.00242	CUT	Sample 13, post activation
26	0.89450	COAG	Sample 13, post activation
27	0.93745	CUT	Sample 14, post activation
28	0.93558	COAG	Sample 14, post activation
29	0.96645	CUT	Sample 15, post activation
30	0.88725	COAG	Sample 15, post activation
31	1.01477	CUT	Sample 16, post activation
32	0.87517	COAG	Sample 16, post activation
33	1.03303	CUT	Sample 17, post activation
34	0.87652	COAG	Sample 17, post activation
35	1.00618	CUT	Sample 18, post activation
36	0.87571	COAG	Sample 18, post activation
37	0.99974	CUT	Sample 19, post activation
38	0.93101	COAG	Sample 19, post activation
39	0.95678	CUT	Sample 20, post activation
40	0.94229	COAG	Sample 20, post activation
41	1.07920	CUT	Sample 21, post activation
42	0.83007	COAG	Sample 21, post activation
43	0.96430	CUT	Sample 22, post activation
44	0.84967	COAG	Sample 22, post activation
45	0.98524	CUT	Sample 23, post activation
46	0.90443	COAG	Sample 23, post activation
47	0.99813	CUT	Sample 24, post activation
48	0.84752	COAG	Sample 24, post activation
49	0.93182	CUT	Sample 25, post activation
50	0.82202	COAG	Sample 25, post activation
51	0.99383	CUT	Sample 26, post activation
52	0.90309	COAG	Sample 26, post activation
53	1.00188	CUT	Sample 27, post activation
54	0.83947	COAG	Sample 27, post activation
55	0.98685	CUT	Sample 28, post activation

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 48 of 79


	Compressive load at Preset Point (Cursor ) (lbf)	Button Type	Note
56	0.88349	COAG	Sample 28, post activation
57	0.98900	CUT	Sample 29, post activation
58	0.89611	COAG	Sample 29, post activation
59	0.92537	CUT	Sample 30, post activation
60	0.80376	COAG	Sample 30, post activation
Maximum	1.07920		
Minimum	0.80376		
Mean	0.93331		
Standard Deviation	0.07		

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 49 of 79


Thursday, October 08, 2015

Button Force on ISOS2 ACE37H -  
151361.is\_comp

Specimen 1 to 24




Specimen 25 to 48



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 50 of 79

Specimen 49 to 60



	Compressive load at Preset Point (Cursor ) (lbf)	Button Type	Note
1	0.95383	CUT	Sample 1, pre-activation
2	0.83893	COAG	Sample 1, pre-activation
3	0.94900	CUT	Sample 2, pre-activation
4	0.82121	COAG	Sample 2, pre-activation
5	0.89692	CUT	Sample 3, pre-activation
6	0.83222	COAG	Sample 3, pre-activation
7	0.91947	CUT	Sample 4, pre-activation
8	0.82846	COAG	Sample 4, pre-activation
9	0.89719	CUT	Sample 5, pre-activation
10	0.78121	COAG	Sample 5, pre-activation
11	0.93692	CUT	Sample 6, pre-activation
12	0.78256	COAG	Sample 6, pre-activation
13	0.87974	CUT	Sample 7, pre-activation
14	0.79678	COAG	Sample 7, pre-activation
15	0.91517	CUT	Sample 8, pre-activation
16	0.81155	COAG	Sample 8, pre-activation
17	0.90605	CUT	Sample 9, pre-activation
18	0.78658	COAG	Sample 9, pre-activation

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 51 of 79

	Compressive load at Preset Point (Cursor ) (lbf)	Button Type	Note
19	0.95490	CUT	Sample 10, pre-activation
20	0.80215	COAG	Sample 10, pre-activation
21	0.91759	CUT	Sample 11, pre-activation
22	0.84027	COAG	Sample 11, pre-activation
23	0.92349	CUT	Sample 12, pre-activation
24	0.81799	COAG	Sample 12, pre-activation
25	0.90927	CUT	Sample 13, pre-activation
26	0.85262	COAG	Sample 13, pre-activation
27	0.95839	CUT	Sample 14, pre-activation
28	0.80215	COAG	Sample 14, pre-activation
29	0.89584	CUT	Sample 15, pre-activation
30	0.84564	COAG	Sample 15, pre-activation
31	0.85504	CUT	Sample 16, pre-activation
32	0.78631	COAG	Sample 16, pre-activation
33	0.89155	CUT	Sample 17, pre-activation
34	0.77074	COAG	Sample 17, pre-activation
35	0.88296	CUT	Sample 18, pre-activation
36	0.79490	COAG	Sample 18, pre-activation
37	0.86954	CUT	Sample 19, pre-activation
38	0.79088	COAG	Sample 19, pre-activation
39	0.90202	CUT	Sample 20, pre-activation
40	0.79303	COAG	Sample 20, pre-activation
41	0.84188	CUT	Sample 21, pre-activation
42	0.75705	COAG	Sample 21, pre-activation
43	0.90792	CUT	Sample 22, pre-activation
44	0.82417	COAG	Sample 22, pre-activation
45	0.85692	CUT	Sample 23, pre-activation
46	0.77235	COAG	Sample 23, pre-activation
47	0.85987	COAG	Sample 24, pre-activation
48	0.81799	COAG	Sample 24, pre-activation
49	0.92886	CUT	Sample 25, pre-activation
50	0.83303	COAG	Sample 25, pre-activation
51	0.88349	CUT	Sample 26, pre-activation
52	0.80537	COAG	Sample 26, pre-activation
53	0.88000	CUT	Sample 27, pre-activation
54	0.80886	COAG	Sample 27, pre-activation
55	0.90551	CUT	Sample 28, pre-activation

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.


<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 52 of 79</b>

	Compressive load at Preset Point (Cursor ) (lbf)	Button Type	Note
56	0.83329	COAG	Sample 28, pre-activation
57	0.87705	CUT	Sample 29, pre-activation
58	0.80054	COAG	Sample 29, pre-activation
59	0.90792	CUT	Sample 30, pre-activation
60	0.79839	COAG	Sample 30, pre-activation
Maximum	0.95839		
Minimum	0.75705		
Mean	0.85486		
Standard Deviation	0.05		


Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 53 of 79</b>

Wednesday, October 14, 2015      Button Force on ISOS2 ACE37H - 151361  
post 500.is\_comp

Specimen 1 to 24




Specimen 25 to 48



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 54 of 79

Specimen 49 to 60



	Compressive load at Preset Point (Cursor ) (lbf)	Button Type	Note
1	0.85235	CUT	Sample 1, post activation
2	0.76940	COAG	Sample 1, post activation
3	0.87705	CUT	Sample 2, post activation
4	0.75195	COAG	Sample 2, post activation
5	0.80242	CUT	Sample 3, post activation
6	0.76484	COAG	Sample 3, post activation
7	0.80081	CUT	Sample 4, post activation
8	0.76430	COAG	Sample 4, post activation
9	0.82658	CUT	Sample 5, post activation
10	0.72537	COAG	Sample 5, post activation
11	0.86470	CUT	Sample 6, post activation
12	0.75974	COAG	Sample 6, post activation
13	0.84457	CUT	Sample 7, post activation
14	0.75410	COAG	Sample 7, post activation
15	0.82524	CUT	Sample 8, post activation
16	0.77799	COAG	Sample 8, post activation
17	0.82658	CUT	Sample 9, post activation
18	0.77933	COAG	Sample 9, post activation

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 55 of 79

	Compressive load at Preset Point (Cursor ) (lbf)	Button Type	Note
19	0.82363	CUT	Sample 10, post activation
20	0.75652	COAG	Sample 10, post activation
21	0.80698	CUT	Sample 11, post activation
22	0.80913	COAG	Sample 11, post activation
23	0.85209	CUT	Sample 12, post activation
24	0.77209	COAG	Sample 12, post activation
25	0.83141	CUT	Sample 13, post activation
26	0.80806	COAG	Sample 13, post activation
27	0.88484	CUT	Sample 14, post activation
28	0.76162	COAG	Sample 14, post activation
29	0.83490	CUT	Sample 15, post activation
30	0.80081	COAG	Sample 15, post activation
31	0.78658	CUT	Sample 16, post activation
32	0.76913	COAG	Sample 16, post activation
33	0.80564	CUT	Sample 17, post activation
34	0.72886	COAG	Sample 17, post activation
35	0.81826	CUT	Sample 18, post activation
36	0.79437	COAG	Sample 18, post activation
37	0.79195	CUT	Sample 19, post activation
38	0.74954	COAG	Sample 19, post activation
39	0.88484	CUT	Sample 20, post activation
40	0.77584	COAG	Sample 20, post activation
41	0.82309	CUT	Sample 21, post activation
42	0.75544	COAG	Sample 21, post activation
43	0.85316	CUT	Sample 22, post activation
44	0.81128	COAG	Sample 22, post activation
45	0.82470	CUT	Sample 23, post activation
46	0.74148	COAG	Sample 23, post activation
47	0.81933	CUT	Sample 24, post activation
48	0.77504	COAG	Sample 24, post activation
49	0.88860	CUT	Sample 25, post activation
50	0.78792	COAG	Sample 25, post activation
51	0.84135	CUT	Sample 26, post activation
52	0.77182	COAG	Sample 26, post activation
53	0.85370	CUT	Sample 27, post activation
54	0.77209	COAG	Sample 27, post activation
55	0.82631	CUT	Sample 28, post activation

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 56 of 79</b>

	Compressive load at Preset Point (Cursor ) (lbf)	Button Type	Note
56	0.80108	COAG	Sample 28, post activation
57	0.84752	CUT	Sample 29, post activation
58	0.78578	COAG	Sample 29, post activation
59	0.88698	CUT	Sample 30, post activation
60	0.76323	COAG	Sample 30, post activation
Maximum	0.88860		
Minimum	0.72537		
Mean	0.80407		
Standard Deviation	0.04		

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 57 of 79

## APPENDIX VIII – ACTIVATION OVER TIME

### Pencil Dome Switch Activation Force Test Product Tested - 0036 - 151363 Data Collection Form

Sample Number	Initial Continuity ( $\Omega$ ) CUT, COAG	Initial Force CUT	Initial Force COAG	Final Force CUT	Final Force COAG	Final Continuity ( $\Omega$ ) CUT, COAG		
1	1.9, 4.6	455.17541	458.83136	414.14344	406.34619	6.0, 2.3 - PASS		
2	3.5, 2.5	485.62053	434.47799	451.27905	448.60285	11.4, 19.2 - PASS		
3	1.6, 5.9	456.39557	443.24593	447.38269	426.31786	3.5, 9.7 - PASS		
4	2.1, 5.9	465.40845	417.79486	423.63713	394.17177	1.5, 4.1 - PASS		
5	1.5, 8.7	453.47443	394.53464	464.67363	405.73837	11.2, 8.4 - PASS		
6	3.5, 3.9	482.93979	381.75241	462.72772	384.18366	7.3, 8.3 - PASS		
7	8.9, 2.6	527.75018	387.23180	475.63696	370.91155	3.4, 8.7 - PASS		
8	6.0, 1.8	479.77372	409.51226	445.19184	393.80436	4.3, 13.7 - PASS		
9	2.7, 21.0	499.50045	427.17061	461.26715	382.72310	5.7, 8.1 - PASS		
10	1.8, 2.0	486.10587	396.11768	446.89735	393.56395	17.1, 15.4 - PASS		
11	3.1, 18.6	446.41200	418.03980	448.84779	395.26493	3.3, 9.6 - PASS		
12	2.2, 2.6	487.81138	398.92088	447.14228	382.96804	3.0, 2.6 - PASS		
13	3.6, 6.1	506.07754	403.06218	454.69006	405.73837	3.5, 2.8 - PASS		
14	3.0, 3.5	487.81138	452.74415	425.22017	424.3/195	4.9, 12.0 - PASS		
15	6.1, 4.3	490.97745	408.78198	438.37434	402.44983	3.9, 5.3 - PASS		
16	2.3, 1.9	492.92790	415.72648	460.29193	396.97043	7.7, 5.4 - PASS		
17	1.6, 1.9	524.09877	416.21182	468.57452	397.58278	2.1, 2.2 - PASS		
18	2.6, 9.1	482.69939	432.28713	456.39557	397.21537	3.7, 2.8 - PASS		
19	2.4, 1.8	483.18473	434.60046	453.47443	422.29903	3.7, 13.9 - PASS		
20	1.9, 9.0	488.30126	457.37079	433.98811	427.41555	7.1, 6.4 - PASS		
21	1.4, 7.1	475.63696	422.78891	489.51688	376.51342	1.6, 29.8 - PASS		
22	1.7, 1.5	506.80783	401.23420	437.39912	385.40383	2.5, 1.5 - PASS		
23	5.1, 2.2	511.67941	420.47105	446.89735	410.24255	5.9, 2.6 - PASS		
24	2.5, 13.6	478.55809	429.23899	452.74415	384.42860	5.8, 3.6 - PASS		
25	2.8, 6.1	472.47088	396.60756	422.66644	372.86200	3.2, 16.9 - PASS		
26	1.4, 4.0	488.05632	417.42745	450.79370	409.63473	5.4, 8.2 - PASS		
27	2.2, 2.7	484.15996	432.77248	454.44512	380.77719	2.1, 10.9 - PASS		
28	2.2, 2.2	473.68651	430.58163	447.62763	400.74432	2.5, 6.3 - PASS		
29	4.6, 6.9	507.05276	403.06218	448.60285	406.46866	5.3, 24.3 - PASS		
30	1.4, 1.8	468.81493	386.98687	419.74077	364.57940	2.5, 7.9 - PASS		
	MAX	527.75018	458.83136	489.51688	448.60285			
	MIN	446.41200	381.75241	414.14344	364.57940			
	AVG	484.97899	417.65288	448.34234	398.34315			
	SD	19.57490327	20.67146154	16.68182837	18.91835671			
	Cpk	3.661508178	2.703451551	3.963241414	2.613742781			

#### Calibration Information:

Instron	
System ID Number	4464C2820
Megadyne Number	01028
Calibration Date	5/1/2015
Calibration Due	5/1/2016

Multimeter	
Fluke 79 Series II Multimeter	
Serial Number:	57040475
Megadyne Number:	01010
Calibration Date:	4/10/2015
Calibration Due:	4/30/2016

PAUL VALPREDA  
Operator Name

Paul Valpreda  
Operator Signature

1-26-2016

Date

CHRISTY CROOK  
Engineer Name

Christy Crook  
Operator Signature

1/25/2018

Date

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 58 of 79

Pencil Dome Switch Activation Force Test  
Product Tested - ACE37H - 151361  
Data Collection Form

Sample Number	Initial Continuity ( $\Omega$ ) CUT, COAG	Initial Force CUT	Initial Force COAG	Final Force CUT	Final Force COAG	Final Continuity ( $\Omega$ ) CUT, COAG	
						CUT	COAG
1	7.8, 3.3	432.65001	380.53225	386.61946	348.99397	3.8, 15.5 - PASS	
2	2.4, 1.5	430.45916	372.49459	397.82319	341.07878	15.5, 6.5 - PASS	
3	4.1, 1.2	406.83607	377.48864	363.97159	346.92559	2.7, 16.5 - PASS	
4	1.8, 5.5	417.06458	375.78313	363.24130	346.68065	2.0, 19.5 - PASS	
5	1.7, 1.5	406.95854	354.35089	374.93038	329.02230	5.5, 12.6 - PASS	
6	2.9, 1.4	424.97976	354.96324	392.22132	344.61227	8.5, 7.5 - PASS	
7	2.3, 1.5	399.04335	361.41333	383.09051	342.05401	5.5, 1.5 - PASS	
8	6.5, 16.7	415.11413	368.11289	374.32257	352.89033	20.4, 30.9 - PASS	
9	5.0, 2.3	410.97737	356.78669	374.93038	353.49814	2.5, 2.6 - PASS	
10	1.7, 4.9	433.13535	363.84912	373.59228	343.15170	8.1, 6.2 - PASS	
11	1.2, 1.1	416.21182	381.14006	366.03997	367.01519	6.4, 1.3 - PASS	
12	1.2, 9.2	418.88802	371.03402	386.50152	350.21413	10.1, 6.9 - PASS	
13	2.0, 1.9	412.43793	386.74193	377.12123	366.52985	3.5, 6.5 - PASS	
14	2.9, 6.4	434.71839	363.84912	401.35667	345.46502	12.5, 21.5 - PASS	
15	2.7, 3.2	406.34619	383.57585	378.70427	363.24130	3.4, 12.9 - PASS	
16	1.4, 1.4	387.83962	356.66422	356.78669	348.87150	9.4, 6.8 - PASS	
17	1.9, 6.3	404.40028	349.60178	365.43216	330.60533	2.3, 4.2 - PASS	
18	3.9, 2.2	400.50392	360.56057	371.15649	360.32017	3.3, 2.4 - PASS	
19	1.7, 1.0	394.41671	358.73713	359.22248	339.98562	7.7, 2.4 - PASS	
20	1.6, 4.2	409.14939	359.71236	401.35667	351.91510	20.4, 3.1 - PASS	
21	1.6, 1.7	381.87034	343.39210	373.34734	342.66182	2.5, 20.6 - PASS	
22	2.9, 1.9	411.82558	373.83722	386.98687	367.99042	2.3, 4.4 - PASS	
23	1.4, 4.5	388.69237	350.33207	374.07763	336.32967	1.6, 1.5 - PASS	
24	10.8, 1.6	390.03047	371.03402	371.64184	351.55223	2.5, 2.9 - PASS	
25	2.4, 6.9	421.32381	377.85605	403.06218	357.39450	26.0, 2.7 - PASS	
26	3.9, 1.5	400.74432	365.30969	381.62994	350.09166	4.7, 5.6 - PASS	
27	2.0, 2.1	399.16128	366.89272	387.23180	350.21413	4.9, 6.1 - PASS	
28	2.0, 2.7	410.73243	377.97398	374.80791	363.36377	1.2, 4.4 - PASS	
29	2.1, 2.0	397.82319	363.11883	384.42860	356.42381	27.3, 21.4 - PASS	
30	1.5, 1.2	411.82558	362.14361	402.32736	346.19530	4.9, 1.2 - PASS	
		MAX	434.71839	386.74193	403.06218	367.99042	
		MIN	381.87034	343.39210	356.78669	329.02230	
		AVG	409.20533	366.30940	379.59875	349.84294	
		SD	13.9254158	10.95244773	13.03610145	10.00429473	
		Cpk	3.810905506	3.539829825	3.313842283	3.32666935	

Calibration Information:

Instron	
System ID Number	4464C2820
Megadyne Number	01028
Calibration Date	5/1/2015
Calibration Due	5/1/2016

Multimeter	
Fluke 79 Series II Multimeter	
Serial Number:	57040475
Megadyne Number:	01010
Calibration Date:	4/10/2015
Calibration Due:	4/30/2016

PAUL VALPREDA Paul Valpreda 1-26-2016  
Operator Name Operator Signature Date

Christian Crank 1/25/2016  
Engineer Name Engineer Signature Date


Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 59 of 79


### 500 Cycles on 0036 Pencils

Text Inputs: Method Author	P. Valpreda
Text Inputs: Test Operator	P. Valpreda

Specimen 1



Specimen 2




Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.


Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 60 of 79

Confidential - Controlled copy on date printed

### Specimen 3





### Specimen 4



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 61 of 79



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 62 of 79

	Specimen note	Load at Machine Peak Load (kgf)	Notes
1	0036 - 151363	-16.58334	Samples 1, 2, 3, 4, 5 - 500 Cycles
2	0036 - 151363	-25.24242	Samples 6, 7, 8, 9, 10 - 500 Cycles
3	0036 - 151363	-24.91583	Samples 11, 12, 13, 14, 15 - 500 Cycles
4	0036 - 151363	-24.67089	Samples 16, 17, 18, 19, 20 - 500 Cycles
5	0036 - 151363	-25.11088	Samples 21, 22, 23, 24, 25 - 500 Cycles
6	0036 - 151363	-25.15624	Samples 26, 27, 28, 29, 30 - 500 Cycles
Maximum		-16.58334	
Minimum		-25.24242	
Mean		-23.61326	
Standard Deviation		3.45001	

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 63 of 79


## 500 Cycles on ACE37H Pencils

Monday, October 12, 2015


ACE37H - 151361.is\_tcyclic

Text Inputs: Method Author	P. Valpreda
Text Inputs: Test Operator	P. Valpreda

Specimen 1





Specimen 2



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 64 of 79




Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.


Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 65 of 79

Confidential - Controlled copy on date printed

### Specimen 5



### Specimen 6



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 66 of 79</b>

	Specimen note 1	Load at Machine Peak Load (kgf)	Specimen label
1	ACE37H - 151361	-24.65728	Samples 1, 2, 3, 4, 5 - 500 Cycles
2	ACE37H - 151361	-23.69567	Samples 6, 7, 8, 9, 10 - 500 Cycles
3	ACE37H - 151361	-22.33035	Samples 11, 12, 13, 14, 15 - 500 Cycles
4	ACE37H - 151361	-21.59553	Samples 16, 17, 18, 19, 20 - 500 Cycles
5	ACE37H - 151361	-20.90607	Samples 21, 22, 23, 24, 25 - 500 Cycles
6	ACE37H - 151361	-21.27348	Samples 26, 27, 28, 29, 30 - 500 Cycles
Maximum		-20.90607	
Minimum		-24.65728	
Mean		-22.40973	
Standard Deviation		1.47827	

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 67 of 79

## APPENDIX IX – ELECTRODE WOBBLE

ELECTRODE DEVIATION WOBBLE ON E0x2 SAMPLES			
ND PN: ACE37H	LOT# 151361	ND PN: 0036	LOT# 151363
Sample #	Wobble	Sample #	Wobble
1	0.143	1	0.145
2	0.128	2	0.132
3	0.149	3	0.146
4	0.127	4	0.140
5	0.131	5	0.140
6	0.128	6	0.147
7	0.145	7	0.149
8	0.142	8	0.135
9	0.132	9	0.136
10	0.133	10	0.140
11	0.131	11	0.132
12	0.133	12	0.143
13	0.145	13	0.143
14	0.143	14	0.147
15	0.147	15	0.146
16	0.144	16	0.135
17	0.136	17	0.137
18	0.147	18	0.155
19	0.143	19	0.149
20	0.145	20	0.131
21	0.143	21	0.137
22	0.130	22	0.144
23	0.145	23	0.132
24	0.144	24	0.130
25	0.135	25	0.138
26	0.143	26	0.143
27	0.151	27	0.141
28	0.131	28	0.146
29	0.143	29	0.138
30	0.142	30	0.144
Button		Rocker	

Confidential - Controlled copy on date printed

Paul Valpreda

9/30/2015

Operator Name

Date

Engineer Signature

Date

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 68 of 79

ELECTRODE DEVIATION		
ND PN: 0035	LOT# 1410115	
Sample #	Wobble	
1	0.118	
2	0.117	
3	0.111	
4	0.121	
5	0.121	
6	0.111	
7	0.114	
8	0.119	
9	0.117	
10	0.120	
11	0.119	
12	0.112	
13	0.115	
14	0.115	
15	0.112	
16	0.113	
17	0.121	
18	0.119	
19	0.118	
20	0.094	
21	0.122	
22	0.118	
23	0.123	
24	0.115	
25	0.117	
26	0.114	
27	0.118	
28	0.116	
29	0.120	
30	0.116	
Button		

10/1/2015

PV

Paul Valpreda

Operator Name

10/1/2015

Date



Operator Signature



Date

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 69 of 79


## APPENDIX X – DURABILITY OF PRINT TO CHEMICAL IN USE ENVIRONMENT (ROCKER PENCIL ONLY)

PRINT DURABILITY EVALUATION			
MEGADYNE PN: 0036			
ND PN:			LOT # 151363
Sample #	Saline	Lactated Ringer's	Notes
1	Visible	Visible	
2	Visible	Visible	
3	Visible	Visible	
4	Visible	Visible	
5	Visible	Visible	
6	Visible	Visible	
7	Visible	Visible	
8	Visible	Visible	
9	Visible	Visible	
10	Visible	Visible	
11	Visible	Visible	
12	Visible	Visible	
13	Visible	Visible	
14	Visible	Visible	
15	Visible	Visible	
16	Visible	Visible	
17	Visible	Visible	
18	Visible	Visible	
19	Visible	Visible	
20	Visible	Visible	
21	Visible	Visible	
22	Visible	Visible	
23	Visible	Visible	
24	Visible	Visible	
25	Visible	Visible	
26	Visible	Visible	
27	Visible	Visible	
28	Visible	Visible	
29	Visible	Visible	
30	Visible	Visible	
Saline: Hospira, Part # NDC 0409-7983-09, LOT # 36-049-JT, EXP 1 DEC 2015			
Lactated Ringer's: Hospira, Part # NDC 0409-7953-09, LOT # 52-089-JT, EXP 1 APR 2017			

Paul Valpreda  
Operator Name

10/29/2015  
Date

  
Engineer Signature

  
Date

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 70 of 79

## APPENDIX XI – PENCIL BODY PEEL STRENGTH


Friday, October 23, 2015

Weld Strength on ISOS2 0036 -  
151363.is\_tens


Specimen label

0036 LOT 151363

Specimen 1 to 15



Specimen 16 to 30



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 71 of 79

	Specimen label	Maximum Load (lbf)	Specimen note 1
1	0036 LOT 151363	34.42	Proximal End
2	0036 LOT 151363	34.07	Proximal End
3	0036 LOT 151363	33.13	Proximal End
4	0036 LOT 151363	43.27	Proximal End
5	0036 LOT 151363	30.01	Proximal End
6	0036 LOT 151363	37.80	Proximal End
7	0036 LOT 151363	33.72	Proximal End
8	0036 LOT 151363	42.09	Proximal End
9	0036 LOT 151363	49.96	Proximal End
10	0036 LOT 151363	46.01	Proximal End
11	0036 LOT 151363	38.79	Proximal End
12	0036 LOT 151363	46.93	Proximal End
13	0036 LOT 151363	45.64	Proximal End
14	0036 LOT 151363	45.82	Proximal End
15	0036 LOT 151363	41.56	Proximal End
16	0036 LOT 151363	47.06	Proximal End
17	0036 LOT 151363	45.53	Proximal End
18	0036 LOT 151363	47.76	Proximal End
19	0036 LOT 151363	40.32	Proximal End
20	0036 LOT 151363	38.58	Proximal End
21	0036 LOT 151363	38.12	Proximal End
22	0036 LOT 151363	45.50	Proximal End
23	0036 LOT 151363	39.65	Proximal End
24	0036 LOT 151363	47.68	Proximal End
25	0036 LOT 151363	45.21	Proximal End
26	0036 LOT 151363	48.27	Proximal End
27	0036 LOT 151363	43.46	Proximal End
28	0036 LOT 151363	37.77	Proximal End
29	0036 LOT 151363	39.62	Proximal End
30	0036 LOT 151363	38.50	Proximal End
Maximum		49.96	
Mean		41.54	
Minimum		30.01	
Standard Deviation		5.26815	

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.


Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 72 of 79</b>

Friday, October 23, 2015


Weld Strength on ISOS2 0036 -  
151363\_Distal.is\_tens

Specimen label	0036 LOT 151363
----------------	-----------------

Specimen 1 to 15



Specimen 16 to 30



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 73 of 79</b>

	Specimen label	Maximum Load (lbf)	Specimen note 1
1	0036 LOT 151363	28.11	Distal End
2	0036 LOT 151363	25.88	Distal End
3	0036 LOT 151363	27.17	Distal End
4	0036 LOT 151363	28.03	Distal End
5	0036 LOT 151363	26.25	Distal End
6	0036 LOT 151363	28.05	Distal End
7	0036 LOT 151363	27.09	Distal End
8	0036 LOT 151363	27.03	Distal End
9	0036 LOT 151363	27.87	Distal End
10	0036 LOT 151363	25.10	Distal End
11	0036 LOT 151363	26.74	Distal End
12	0036 LOT 151363	26.28	Distal End
13	0036 LOT 151363	28.78	Distal End
14	0036 LOT 151363	29.66	Distal End
15	0036 LOT 151363	29.10	Distal End
16	0036 LOT 151363	29.83	Distal End
17	0036 LOT 151363	28.48	Distal End
18	0036 LOT 151363	25.13	Distal End
19	0036 LOT 151363	26.87	Distal End
20	0036 LOT 151363	25.72	Distal End
21	0036 LOT 151363	27.01	Distal End
22	0036 LOT 151363	26.55	Distal End
23	0036 LOT 151363	25.85	Distal End
24	0036 LOT 151363	27.22	Distal End
25	0036 LOT 151363	28.27	Distal End
26	0036 LOT 151363	27.70	Distal End
27	0036 LOT 151363	25.15	Distal End
28	0036 LOT 151363	26.07	Distal End
29	0036 LOT 151363	25.18	Distal End
30	0036 LOT 151363	26.87	Distal End
Maximum		29.83	
Mean		27.10	
Minimum		25.10	
Standard Deviation		1.32832	

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.


Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 74 of 79

Monday, October 26, 2015


Weld Strength on ISOS2 ACE37H -  
151361.is\_tens

Specimen label	ACE37H LOT 151361
----------------	-------------------

### Specimen 1 to 15



### Specimen 16 to 29



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	ISOS2 Functional Verification After EO Exposure	Revision: 001
		Page 75 of 79

	Specimen label	Maximum Load (lbf)	Specimen note 1
1	ACE37H LOT 151361	43.81	Proximal End, Sample 1
2	ACE37H LOT 151361	44.99	Proximal End, Sample 2
3	ACE37H LOT 151361	27.60	Proximal End, Sample 3
4	ACE37H LOT 151361	45.45	Proximal End, Sample 4
5	ACE37H LOT 151361	41.45	Proximal End, Sample 5
6	ACE37H LOT 151361	46.15	Proximal End, Sample 6
7	ACE37H LOT 151361	27.81	Proximal End, Sample 7
8	ACE37H LOT 151361	40.48	Proximal End, Sample 8
9	ACE37H LOT 151361	47.09	Proximal End, Sample 10
10	ACE37H LOT 151361	44.89	Proximal End, Sample 11
11	ACE37H LOT 151361	39.97	Proximal End, Sample 12
12	ACE37H LOT 151361	28.54	Proximal End, Sample 13
13	ACE37H LOT 151361	22.07	Proximal End, Sample 14
14	ACE37H LOT 151361	35.84	Proximal End, Sample 15
15	ACE37H LOT 151361	27.25	Proximal End, Sample 16
16	ACE37H LOT 151361	40.70	Proximal End, Sample 17
17	ACE37H LOT 151361	37.88	Proximal End, Sample 18
18	ACE37H LOT 151361	27.62	Proximal End, Sample 19
19	ACE37H LOT 151361	32.83	Proximal End, Sample 20

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 76 of 79</b>

	Specimen label	Maximum Load (lbf)	Specimen note 1
20	ACE37H LOT 151361	37.83	Proximal End, Sample 21
21	ACE37H LOT 151361	38.79	Proximal End, Sample 22
22	ACE37H LOT 151361	26.23	Proximal End, Sample 23
23	ACE37H LOT 151361	37.31	Proximal End, Sample 24
24	ACE37H LOT 151361	43.87	Proximal End, Sample 25
25	ACE37H LOT 151361	27.41	Proximal End, Sample 26
26	ACE37H LOT 151361	37.96	Proximal End, Sample 27
27	ACE37H LOT 151361	45.15	Proximal End, Sample 28
28	ACE37H LOT 151361	27.38	Proximal End, Sample 29
29	ACE37H LOT 151361	38.31	Proximal End, Sample 30
Maximum		47.09	
Mean		36.64	
Minimum		22.07	
Standard Deviation		7.48984	

Confidential - Controlled copy on date printed

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	Revision: 001
		Page 77 of 79


Monday, October 26, 2015

Weld Strength on ISOS2 ACE37H -  
151361\_Distal.is\_tens


Specimen label

ACE37H LOT 151361

### Specimen 1 to 15



### Specimen 16 to 29



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 78 of 79</b>

	Specimen label	Maximum Load (lbf)	Specimen note 1
1	ACE37H LOT 151361	26.09	Distal End, Sample 1
2	ACE37H LOT 151361	26.17	Distal End, Sample 2
3	ACE37H LOT 151361	25.02	Distal End, Sample 3
4	ACE37H LOT 151361	26.34	Distal End, Sample 4
5	ACE37H LOT 151361	29.69	Distal End, Sample 5
6	ACE37H LOT 151361	25.64	Distal End, Sample 6
7	ACE37H LOT 151361	26.42	Distal End, Sample 7
8	ACE37H LOT 151361	27.11	Distal End, Sample 8
9	ACE37H LOT 151361	28.70	Distal End, Sample 10
10	ACE37H LOT 151361	25.40	Distal End, Sample 11
11	ACE37H LOT 151361	26.01	Distal End, Sample 12
12	ACE37H LOT 151361	24.81	Distal End, Sample 13
13	ACE37H LOT 151361	23.65	Distal End, Sample 14
14	ACE37H LOT 151361	27.19	Distal End, Sample 15
15	ACE37H LOT 151361	24.97	Distal End, Sample 16
16	ACE37H LOT 151361	25.88	Distal End, Sample 17
17	ACE37H LOT 151361	26.85	Distal End, Sample 18
18	ACE37H LOT 151361	24.94	Distal End, Sample 19
19	ACE37H LOT 151361	25.56	Distal End, Sample 20

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>ENG-RPT-452</b>
	<b>ISOS2 Functional Verification After EO Exposure</b>	<b>Revision: 001</b>
		<b>Page 79 of 79</b>

	Specimen label	Maximum Load (lbf)	Specimen note 1
20	ACE37H LOT 151361	26.79	Distal End, Sample 21
21	ACE37H LOT 151361	23.70	Distal End, Sample 22
22	ACE37H LOT 151361	27.76	Distal End, Sample 23
23	ACE37H LOT 151361	25.96	Distal End, Sample 24
24	ACE37H LOT 151361	29.74	Distal End, Sample 25
25	ACE37H LOT 151361	24.72	Distal End, Sample 26
26	ACE37H LOT 151361	27.78	Distal End, Sample 27
27	ACE37H LOT 151361	28.38	Distal End, Sample 28
28	ACE37H LOT 151361	26.44	Distal End, Sample 29
29	ACE37H LOT 151361	25.58	Distal End, Sample 30
Maximum		29.74	
Mean		26.32	
Minimum		23.65	
Standard Deviation		1.53848	