



DOCUMENT NUMBER: ENG-RPT-340

DOCUMENT TITLE: Test Report - ULPA Replacement Filter & Carbon/Muffler

DOCUMENT NOTES:

Initial release into MasterControl

**Document Information**

Revision: 001

Vault: MEG-Rel

Status: Release

Document Type: ENG-RPT

**Date Information**

Effective Date: 07 Apr 2014

Expiration Date:

Release Date: 07 Apr 2014

Next Review Date:

**Control Information**

Author: PBORGMEIER

Owner: ENG-APPR

Previous Number: 1150751-01

Change Number: 14-056-01

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b> <b>MASTER DOCUMENT</b>	<b>Document Number</b> <b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision: A</b> <b>Effective Date:</b> <b>2014 APR 07</b>

4/4/14

Engineering Verification:

4-7-14  
D.C. Verification:

Authored By: Shawn Horner

## TABLE OF CONTENTS

1.	ABSTRACT .....	1
2.	IEC SHIPPING/STORAGE CONDITIONING .....	2
3.	ULPA FILTER INSPECTION AND TESTING.....	2
4.	CARBON/MUFFLER TESTING .....	10
5.	CONCLUSION.....	10
6.	REVISION HISTORY .....	10

### 1. ABSTRACT

#### 1.1. ULPA Replacement Filter #2211 Testing

The ULPA Replacement Filter, Catalog # 2211, passed the verification testing outlined in 1150751-10. The 2211 Filter is equivalent to or better than the current 2210 following the testing outlined below:

Product tested:  
2211 Lot # 61964  
2210 Lot # 5260

Test #	Description	Pass/Fail
A	1 <sup>st</sup> Article Feature Measurements	Complete
B	Force to insert Filter	Passed
C	Flow Test	Passed
D	Weld Joint Seal Pressure Test	Passed
E	Shear Force Test	Passed
1150095-01	Thermal Conditioning	Complete
F	Ship Test ASTM D4169-05	Passed

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>  <b>MASTER DOCUMENT</b>	<b>Document Number</b> <b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision: A</b>
		<b>Effective Date:</b> <i>2014 APR 07</i>
	<b>Page 2 of 30</b>	

## 1.2. Carbon/Muffler (Vendor PN DF-0) Catalog #2220

The Carbon/Muffler (Vendor PN DF-0), Catalog #2220, passed the verification testing outlined in 1150751-10 (ship test). The DF-0 Filter is identical to the #2220 product.

Test #	Description	Pass/Fail
1150095-01	Thermal Conditioning	Complete
A	Ship Test ASTM D4169-05	Passed

## 2. IEC SHIPPING/STORAGE CONDITIONING

All devices were conditioned per 1150095-10 starting on March 20, 2014.

## 3. ULPA FILTER INSPECTION AND TESTING

### 3.1. 1<sup>st</sup> Article Feature Measurements (Test A); See Appendix A for raw data.

3.1.1. Using Calipers and the ROI a dimensional comparison on the exhaust ports of the two designs were done by QA. See the first article inspection points report for required measurements

3.1.2. Using a threaded/barb connector and 22 mm port connector, the intake interfaces achieve their respective desired seals.

3.1.3. Results:

3.1.3.1. The 2211 parts passed 1<sup>st</sup> article inspection with one adjustment as an input to the drawing revision. The stem of the exhaust port has been consistently shorter than the currently defined length of  $0.430 \pm .005$ . As a result of this first article and inspection by the molder this dimension will be shifted to  $0.425 \pm .005$ .

The 2210 did not pass all the drawing dimensions at various points of the exhaust part inspection. No product drawings were provided for this part. The drawing dimensions were developed following review of various parts that fit into the LCD 10006 connection. Though the dimensions didn't match, the 2210 has a long history of use therefore it passes.


3.1.3.2. Intake connections connect to filters as intended, as evaluated by Marketing.

### 3.2. Force to insert Filter (Test B); See Appendix B for raw data.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number 1150751-01
	MASTER DOCUMENT	
	ULPA REPLACEMENT FILTER & CARBON/MUFFLER	Revision: A
		Effective Date: 2014 APR 07
		Page 3 of 30

- 3.2.1. Using the Instron and a feed rate of 7.8 in/min to a stroke length of 0.32 inches, the #2211 and 2210 filters were inserted into the Colder PN LC10006 3/8 NPT non-valved coupling body (Filter Connector fitting for the MV & MVP) while measuring the force to insert.

3.2.2. Results:



Sample #	2210 (lbf)	2211 (lbf)
1	11.11	6.23
2	12.81	4.32
3	12.17	4.67
4	12.19	5.32
5	5.86	4.84
6	6.01	6.76
7	8.99	4.81
8	12.21	3.67
9	5.74	5.24
10	7.95	5.92
11	6.91	4.95
Average	9.268182	5.157273
St Dev	2.891438	0.880296

Force to insert filter

Megadyne Medical Products, Inc.	<b>TEST REPORT</b>	<b>Document Number</b> <b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision: A</b>
		<b>Effective Date:</b> <b>2014 APR 07</b>
		<b>Page 4 of 30</b>

3.2.2.1. T-test resulted in the following:

**Two-Sample T-Test and CI: Control 2210, Test 2211 for insertion force of filter**

Two-sample T for Control 2210 vs Test 2211

	N	Mean	StDev	SE Mean
Control 2210	11	9.27	2.89	0.87
Test 2211	11	5.157	0.880	0.27


Difference = mu (Control 2210) - mu (Test 2211)

Estimate for difference: 4.111

95% lower bound for difference: 2.539

T-Test of difference = 0 (vs >): T-Value = 4.51 P-Value = 0.000 DF = 20

Both use Pooled StDev = 2.1372



The P-Value of 0 means that the insertion force for the two filters are different. The Control (2210) has a higher force for insertion than the Test (2211). Therefore, the 2211 passes insertion force test.

3.3. Flow Test (Test C); See Appendix C for raw data.

Megadyne Medical Products, Inc.	TEST REPORT <b>MASTER DOCUMENT</b>	Document Number <b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision: A</b>
		<b>Effective Date:</b> <b>2014 APR 07</b>
		<b>Page 5 of 30</b>

3.3.1. Using the MVP in an open mode and max flow condition, the air flow was measured at the filter face (filter only), between the filter and 2110-10 smoke pencil/tube set and 100 lpm box flow media evaluation .

3.3.2. Results:

Sample #	Filter Only		W/ 2210-09		@ 100 lpm Box Flow	
	2210 (lpm)	2211 (lpm)	2210 (lpm)	2211 (lpm)	2210 (lpm)	2211 (lpm)
1	76.49	80.06	69.60	70.34	90.50	101.70
2	76.70	80.32	68.52	70.25	89.85	102.10
3	76.89	80.22	69.62	70.90	90.90	101.80
4	76.68	80.33	68.36	70.77	90.10	102.00
5	75.96	80.44	67.45	71.45	90.02	101.90
6	77.25	80.53	68.56	70.55	91.60	101.80
7	77.34	79.96	69.91	70.01	91.70	102.10
8	77.31	79.45	69.57	70.33	90.07	102.00
9	77.16	80.47	69.45	70.45	91.20	102.20
10	76.55	79.85	67.23	69.61	90.70	101.90
11	76.43	79.94	67.84	69.60	91.00	102.20
Average	76.80	80.14	68.74	70.39	90.69	101.97
St Dev	0.44	0.33	0.95	0.54	0.65	0.17

3.3.2.1. A T-test of the max flow rates resulted in the following:

#### Two-Sample T-Test and CI: Control 2210, Test 2211 Filter Only

Two-sample T for Control 2210 vs Test 2211

	N	Mean	StDev	SE Mean
Control 2210	11	76.796	0.438	0.13
Test 2211	11	80.143	0.326	0.098

Difference = mu (Control 2210) - mu (Test 2211)  
Estimate for difference: -3.346  
95% upper bound for difference: -3.062  
T-Test of difference = 0 (vs <): T-Value = -20.32 P-Value = 0.000 DF = 20  
Both use Pooled StDev = 0.3861

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b> <b>MASTER DOCUMENT</b>	<b>Document Number</b> <b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision: A</b>
		<b>Effective Date:</b> <b>2014 APR 07</b>
		<b>Page 6 of 30</b>

The P-Value of 0 means that the flow measure with filter only is different. The Control (2210) has lower flow than the Test (2211). Therefore, the 2211 passes the flow test.

### **Two-Sample T-Test and CI: Control 2210, Test 2211 with 2110-09 Pencil**

Two-sample T for Control 2210 vs Test 2211

	N	Mean	StDev	SE Mean
Control 2210	11	68.737	0.952	0.29
Test 2211	11	70.387	0.544	0.16

Difference = mu (Control 2210\_1) - mu (Test 2211\_1)  
Estimate for difference: -1.650  
95% upper bound for difference: -1.080  
T-Test of difference = 0 (vs <): T-Value = -4.99 P-Value = 0.000 DF = 20  
Both use Pooled StDev = 0.7751

The P-Value of 0 means that the flow measure between the filter and 2110-09 pencil is different. The Control (2210) has lower flow than the Test (2211). Therefore, the 2211 passes the flow test.

### **Two-Sample T-Test and CI: Control 2210, Test 2211 using a MVP @ 100 lpm box flow**

Two-sample T for Control 2210 vs Test 2211

	N	Mean	StDev	SE Mean
Control 2210	11	90.695	0.646	0.19
Test 2211	11	101.973	0.168	0.051

Difference = mu (Control 2210) - mu (Test 2211)  
Estimate for difference: -11.278  
95% upper bound for difference: -10.931  
T-Test of difference = 0 (vs <): T-Value = -56.05 P-Value = 0.000 DF = 20  
Both use Pooled StDev = 0.4719

The P-Value of 0 means that the flow measure using a modified MVP to reach 100 lpm is different. The Control (2210) has lower flow than the Test (2211). Therefore, the 2211 passes the flow test. The 100 lpm box flow test resulted in no visible separation under microscope of the filter media from the perimeter capture surface.

Therefore the 2211 filter passes the various aspects of the flow test defined.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b> <b>MASTER DOCUMENT</b>	<b>Document Number</b> <b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision: A</b>
		<b>Effective Date:</b> 2014 APR 07
		<b>Page 7 of 30</b>

3.4. Weld Joint Seal Pressure Test (Test D) ; See Appendix D for raw data.


3.4.1. The filter inlets were sealed with a threaded connector with the ID sealed. The exhaust port was connected to a tube set. The system was pressurized to 7 psi and placed unit under water to check for weld joint seal leaks. The units passed, there was no leak/bubbles visible around the weld joint of the filter for a period of 30 seconds.

3.4.2. Results:


Sample #	No Leak 2211	Time / Pressure
1	Pass	30s/7psi
2	Pass	30s/7psi
3	Pass	30s/7psi
4	Pass	30s/7psi
5	Pass	30s/7psi
6	Pass	30s/7psi
7	Pass	30s/7psi
8	Pass	30s/7psi
9	Pass	30s/7psi
10	Pass	30s/7psi
11	Pass	30s/7psi
12	Pass	30s/7psi
13	Pass	30s/7psi
14	Pass	30s/7psi
15	Pass	30s/7psi
16	Pass	30s/7psi

Sample #	No Leak 2211	Time / Pressure
17	Pass	30s/7psi
18	Pass	30s/7psi
19	Pass	30s/7psi
20	Pass	30s/7psi
21	Pass	30s/7psi
22	Pass	30s/7psi
23	Pass	30s/7psi
24	Pass	30s/7psi
25	Pass	30s/7psi
26	Pass	30s/7psi
27	Pass	30s/7psi
28	Pass	30s/7psi
29	Pass	30s/7psi
30	Pass	30s/7psi
31	Pass	30s/7psi
32	Pass	30s/7psi

Megadyne Medical Products, Inc.	<b>TEST REPORT</b> <b>MASTER DOCUMENT</b> <b>ULPA REPLACEMENT FILTER &amp;</b> <b>CARBON/MUFFLER</b>	<b>Document Number</b> <b>1150751-01</b>
		<b>Revision: A</b>
		<b>Effective Date:</b> 2014 APR 07
		<b>Page 8 of 30</b>



Picture of Leak Test setup



Pressure Gauge read-out

The 2211 filter showed no signs of weld joint leak. Therefore it passes the leak test.

### 3.5. Shear Force Test (Test E) ; See Appendix E for raw data.

- 3.5.1. Using the Instron with the 100 lb load cell, a feed rate of 1.6 in/sec and a stroke length of 1.5 inches, the #2211 and 2210 filters were attached to the Colder fitting #LC10006 in the secured vice on the base of the Instron.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number 1150751-01
	MASTER DOCUMENT	Revision: A
	ULPA REPLACEMENT FILTER & CARBON/MUFFLER	Effective Date: 2014 APR 07
		Page 9 of 30

### 3.5.2. Results:

Sample #	2210 (lbf)	2211 (lbf)
1	19.03	26.55
2	18.42	19.22
3	19.25	27.87
4	19.95	35.30
5	17.37	19.78
6	17.30	37.80
7	17.53	38.01
8	17.96	25.72
9	16.67	34.07
10	17.12	33.91
11	19.36	32.00
Average	18.18	30.02
St Dev	1.08	6.67

#### 3.5.2.1. A T-test resulted in the following:

##### Two-Sample T-Test and CI: Control 2210, Test 2211 for Shear Test

Two-sample T for Control 2210 vs Test 2211

	N	Mean	StDev	SE Mean
Control 2210	11	18.18	1.08	0.33
Test 2211	11	30.02	6.67	2.0

Difference = mu (Control 2210-3) - mu (Test 2211-3)  
Estimate for difference: -11.84  
95% upper bound for difference: -8.33  
T-Test of difference = 0 (vs <): T-Value = -5.82 P-Value = 0.000 DF = 20  
Both use Pooled StDev = 4.7755

The P-Value of 0 means that the shear forces of the two filters are different. The Control (2210) has a lower shear force than the Test (2211). Therefore, the 2211 passes the shear test.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b> <b>MASTER DOCUMENT</b>	<b>Document Number</b> <b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision: A</b>
		<b>Effective Date:</b> <b>2014 APR 07</b>
		<b>Page 10 of 30</b>

3.6. Ship Test (Test F) ; See Appendix F for raw data.

- 3.6.1. The worst case packaging, 10 units in a case box, passed the ASTM D4169-05 ship test. The box remained intact and the product showed no visual signs of damage.

#### 4. CARBON/MUFFLER TESTING

4.1. Ship Test (Test A) ; See Appendix G for raw data.

- 4.1.1. The worst case packaging, 10 units in a case box, passed the ASTM D4169-05 ship test. The box remained intact and the product showed no visual signs of damage.

#### 5. CONCLUSION

5.1. ULPA Filter

- 5.1.1. The 2211 Filter passed all requirements outlined per the protocol.

5.2. Carbon Muffler Testing

- 5.2.1. The 2220 Carbon/Muffler passed the ship testing outlined per the protocol.

#### 6. REVISION HISTORY

REVISION	DOCUMENT CHANGE ORDER NUMBER	DESCRIPTION OF CHANGE	EFFECTIVE DATE
A	14-056-01	Initial Release	2014 APR 07

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>1150751-01</b>
	<b>MASTER DOCUMENT</b>	<b>Revision: A</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Effective Date:</b> <b>2014 APR 07</b>
		<b>Page 11 of 30</b>

# APPENDIX A

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 09:23:59 pm; Printed by: .

## FIRST ARTICLE INSPECTION REPORT

PART NO.:	2211 X5800119-01	REV.:	07	P.O.#		2014 APR 07
PART NAME:	INSPECTION POINT X5800119-01 REV 07 DETAIL L				Sheet 1 OF ____	
INSPECTOR:	m. Fisher		DATE:	3/24/2014		
SUPPLIER:			REASON:			
APPROVAL:						
INSPECTOR	m. Fisher		DATE:	3/24/14	ACCEPT:	REJECT: X
QUALITY ENG.	<i>John S. Salyer</i>		DATE:	2014 Mar 24	ACCEPT: ✓	REJECT: _____
PROD. ENG.	<i>M. Fisher</i>		DATE:	3/24/14	ACCEPT: ✓	REJECT: _____
ITEM	SPECIFICATION	1	RESULTS	2	3	PASS FAIL
M40	0.360 ± .005	.359	.363	.359		✓
M32	0.063 ± .005	.063	.063	.061		✓
M37	0.430± .005	.4285	.4265	.427		✓
M36	0.328± .005	.331	.327	.331		✓
M35	0.270± .005	.270	.273	.273		✓
M34	0.077± .005	.074 <del>.074</del> mf 31	.075	.0765		✓
M33	0.430± .005	.426	.424	.427		✓
M15-M14	1.427 ± .005 -0.933 ± .005 = 0.494 ± .010	.493	.494	.497		✓
M15-M13	1.427 ± .005 -0.873 ± .005 = 0.554 ± .010	.554	.554	.554		✓
COMMENTS: Note: M33 is being adjusted to .425 ± .005						

# FIRST ARTICLE INSPECTION REPORT


## MASTER DOCUMENT

PART NO.:	2211	REV.:	P.O.#	2014 APR 07		
PART NAME: INSPECTION POINT X5800119-01 REV 07 DETAIL L				Sheet 1 OF		
INSPECTOR:	m Fisher	DATE:	3/25/2014			
SUPPLIER:		REASON:				
APPROVAL:						
INSPECTOR	m Fisher	DATE:	3/25/14	ACCEPT:	REJECT:	X
QUALITY ENG.	Jay Sae	DATE:	2014 APR 01	ACCEPT:	✓	REJECT:
PROD. ENG.	M H	DATE:	3/28/14	ACCEPT:	X	REJECT:
ITEM	SPECIFICATION	1	RESULTS	2	3	PASS FAIL
M40	0.360 ± .005	3.58	1.3601	.358	X	
M32	0.063 ± .005	.0100	.055	.057	X	
M37	0.430± .005	.427	.430	.428	X	
M36	0.328± .005	.330	.329	.331	X	
M35	0.270± .005	.288	.288	.288	X	
M34	0.077± .005	.081	.087	.081	X	
M33	0.430± .005	.439	.442	.440	X	
M15-M14	1.427 ±.005 -0.933 ± .005 = 0.494 ± .010	.499	.500	.494	X	
M15-M13	1.427 ±.005 -0.873 ± .005 = 0.554 ± .010	.554	.556	.556	X	
COMMENTS: We have purchased this part for many years. The dimension outlined did come from 2210 filter measurements. The purpose of these measurements are to assist in qualification of a new part per protocol #. The intent is not to qualify the existing part 2210, but show compatibility to the new unit.						



MASTER DOCUMENT

2014 APR 07



DETAIL  
SCALE 4: 1

**SOLIDWORKS DOCUMENT**  
 B X5800119-01  
 DO NOT SCALE PRINT    SCALE 1:1    SHEET 2 OF 2


Printed on: 21 Jan 2020, 09:23:59 pm; Printed by: .

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b> <b>MASTER DOCUMENT</b>	<b>Document Number</b> <b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision: A</b> <b>Effective Date:</b> <b>2014 APR 07</b>
		<b>Page 15 of 30</b>

# APPENDIX B

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 09:23:59 pm; Printed by: .

Specimen 1 to 11 *6.1# 2210*

	Load at Machine Peak Load (lbf)	Notes	Comments
1	-11.11000		
2	-12.81000		
3	-12.17000		
4	-12.19000		
5	-5.86600		
6	-6.01100		
7	-8.99300		
8	-12.21000		
9	-5.74200		
10	-7.95200		
11	-6.91800		
Maximum	-5.74200		
Minimum	-12.81000		
Mean	-9.27018		
Standard Deviation	2.89		

*7.8 in/min  
stroke .32 inches*


By Paul Valpredo on 3-21-2014

Instron 4464  
System id - 4464 C2820  
Calibrated on 5-3-2013  
Due on 5-3-2014

2014 APR 07

Cat# 2211

Specimen 1 to 11



	Load at Machine Peak Load (lbf)	Notes	Comments
1	-6.23100		
2	-4.31700		
3	-4.67400		
4	-5.32100		
5	-4.84000		
6	-6.76200		
7	-4.81100		
8	-3.67200		
9	-5.23800		
10	-5.91700		
11	-4.95300		
Maximum	-3.67200		
Minimum	-6.76200		
Mean	-5.15782		
Standard Deviation	0.88		

7.8 in/min  
stroke .32 inches

By Paul Valpreda 3-21-2014

Instron 4464

System id - 4464 C2820

Calibrated on 5-3-2014  
2013

Due on 5-3-2014

**MASTER DOCUMENT**

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision:</b> A
		<b>Effective Date:</b> 2014 APR 07
		<b>Page 18 of 30</b>

# APPENDIX C

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 09:23:59 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST PROTOCOL	<u>Document Number</u> <u>XI150751-10</u>
	<b>MASTER DOCUMENT</b>	
	ULPA REPLACEMENT FILTER & CARBON/MUFFLER	Revision: 01
		Effective Date: 2014 APR 07
		Page 13 of 17

**Appendix 2**  
Flow Testing Log Sheet

2110-09  
Lot 5103

Enter the flow value and if dislodgement occurs.

Sample	2211	Dislodgement @ 100 lpm	2210	Dislodgement @ 100 lpm
1	80.06	70.34	101.7	76.49 69.60
2	80.32	70.25	102.1	76.70 68.52
3	80.22	70.90	101.8	76.89 69.62
4	80.33	70.17	102.0	76.68 68.36
5	80.44	71.45	101.9	75.96 67.45
6	80.53	70.55	101.8	77.25 68.56
7	79.96	70.01	102.1	77.34 69.91
8	79.45	70.33	102.0	77.31 69.57
9	80.47	70.45	102.2	77.16 69.45
10	79.85	69.61	101.9	76.55 67.23
11	79.94	69.60	102.2	76.43 67.84

Comments: \_\_\_\_\_  
\_\_\_\_\_

Paul Valpreda  
Inspected by:

3-21-2014  
Date completed

Flowmeter - Model 4040F

Serial number - 4040 0409 003

Megadyne number - 01272

Calibration date - April 2013

Recal Date - April 2014

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 09:23:59 pm; Printed by: .

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b>
	<b>MASTER DOCUMENT</b>	<b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision: A</b>
		<b>Effective Date: 2014 APR 07</b>
		<b>Page 20 of 30</b>

# APPENDIX D

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 09:23:59 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST PROTOCOL	<u>Document Number</u> <u>X1150751-10</u>
	ULPA REPLACEMENT FILTER & <del>MASTER DOCUMENT</del> CARBON/MUFFLER	Revision: 01 Effective Date: 2014 APR 07 Page 14 of 17

## Appendix 3

### Bubble Leak Test @ 7psi for 30 seconds Log Sheet

Catalog #\_\_\_\_\_

Sample	Pass	Fail	Comment
1	✓		
2	✓		
3	✓		
4	✓		
5	✓		
6	✓		
7	✓		
8	✓		
9	✓		
10	✓		
11	✓		
12	✓		
13	✓		
14	✓		
15	✓		
16	✓		

Sample	Pass	Fail	Comment
17	✓		
18	✓		
19	✓		
20	✓		
21	✓		
22	✓		
23	✓		
24	✓		
25	✓		
26	✓		
27	✓		
28	✓		
29	✓		
30	✓		
31	✓		
32	✓		

Signature:

Paul Valpride

Date: 3-21-2014

Megadyne Pressure gauge: Calibration 01/31/19  
#01263 Due 01/31/15

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>  <b>MASTER DOCUMENT</b>	<b>Document Number</b> <b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision: A</b>
		<b>Effective Date:</b> <b>2014 APR 07</b>
		<b>Page 22 of 30</b>


# APPENDIX E

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 09:23:59 pm; Printed by: .

Specimen 1 to 11  
**MASTER DOCUMENT** *6/1/2010*

2014 APR 07



	Load at Machine Peak Load (lbf)	Notes	Comments
1	-19.03000		
2	-18.42000		
3	-19.25000		
4	-19.95000		
5	-17.37000		
6	-17.30000		
7	-17.53000		
8	-17.96000		
9	-16.67000		
10	-17.12000		
11	-19.36000		
Maximum	-16.67000		
Minimum	-19.95000		
Mean	-18.17818		
Standard Deviation	1.08		

By Paul Valpreda on 3-21-2014

Instron 4464

System id - 4464 C2820

Calibrated on 5-3-2013


Due on 5-3-2014

# MASTER DOCUMENT

Specimen 1 to 11

C#2211

2014 APR 07



	Load at Machine Peak Load (lbf)	Notes	Comments
1	-26.55000		New Filter
2	-19.22000		New Filter
3	-27.87000		New Filter
4	-35.30000		New Filter
5	-19.78000		New Filter
6	-37.80000		New Filter
7	-38.01000		New Filter
8	-25.72000		New Filter
9	-34.07000		New Filter
10	-33.91000		New Filter
11	-32.00000		New Filter
Maximum	-19.22000		
Minimum	-38.01000		
Mean	-30.02091		
Standard Deviation	6.67		

By Paul Valpreda on 3-21-2014

Instron 4464

System id - 4464CZ820

Calibrated on 5-3-2013

Due on 5-3-2014

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number</b> <b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision: A</b>
		<b>Effective Date:</b> <b>2014 APR 07</b>
		<b>Page 25 of 30</b>

# APPENDIX F

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 09:23:59 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST PROTOCOL	<u>Document Number</u> <u>XI150751-10</u>
	<b>MASTER DOCUMENT</b>	
	ULPA REPLACEMENT FILTER & CARBON/MUFFLER	Revision: 01
	BOX 1	Effective Date: 2014 APR 07
		Page 16 of 17

## Appendix 5

### Shipping Test Log Sheet

#### CARBON FILTER / MUFFLER

Preconditioning:

Start Date: 3-21-2014 Chamber Number: 01095  
Completion Date: 3-24-2014 Last Calibration: 3-28-2013  
Signature/Date: Paul Valpreda 3-24-2014 Calibration due: 3-31-2014

Drop Test:

Catalog \_\_\_\_\_ Weight 8.65 lbs. Drop Height: 15.5 inches

Drop	Orientation	Specific face, edge or	Initials/Date
1	Top	Face 1	PV - 3-25-2014
2	Edge	Edge 5-3	PV - 3-25-2014
3	Edge	Edge 6-3	PV - 3-25-2014
4	Corner	Corner 2-3-5	PV - 3-25-2014
5	Corner	Corner 4-3-6	PV - 3-25-2014
6	Bottom	Face 3	PV - 3-25-2014

Comments: \_\_\_\_\_

Signature: Paul Valpreda Date: 3-25-2014

Compression Test:

Catalog \_\_\_\_\_ Pounds Force 454 lbs.

Comments: \_\_\_\_\_

Signature: Paul Valpreda Date: 3-25-2014

Megadyne Medical Products, Inc.	TEST PROTOCOL	<u>Document Number</u> <u>XI150751-10</u>
	<b>MASTER DOCUMENT</b>	
	ULPA REPLACEMENT FILTER & CARBON/MUFFLER	Revision: 01 Effective Date: 2014 APR 07 Page 17 of 17

## Appendix 5 Continued

### Shipping Test Log Sheet

Vibration:

Low Frequency, 40 minutes, Initials PV High frequency 10 minutes, Initials PV

Completion Date: 3-25-2014

Signature: Paul Valpreda Date: 3-27-2014

Second Drop Test:

Catalog \_\_\_\_\_ Weight 8.65 lbs . Drop Height: 15.5 + 30 inches.

Drop	Orientation	Specific face, edge or	Initials/Date
1	Edge	Edge 4-6	<u>PV - 3-25-2014</u>
2	Face	Face 4	<u>PV - 3-25-2014</u>
3	Face	Face 6	<u>PV - 3-25-2014</u>
4	Corner	Corner 2-1-5	<u>PV - 3-25-2014</u>
5	Edge	Edge 2-1	<u>PV - 3-25-2014</u>
6	Bottom	Face 3, Increase height to 30 inches.	<u>PV - 3-25-2014</u>

Comments: \_\_\_\_\_

Signature: Paul Valpreda Date: 3-25-2014

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b> <b>MASTER DOCUMENT</b>	<b>Document Number</b> <b>1150751-01</b>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision: A</b> <b>Effective Date:</b> <b>2014 APR 07</b>
		<b>Page 28 of 30</b>

# APPENDIX G

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 09:23:59 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST PROTOCOL <b>MASTER DOCUMENT</b>	Document Number <u>XI150751-10</u>
	ULPA REPLACEMENT FILTER & CARBON/MUFFLER Box 2	Revision: 01 Effective Date: 2014 APR 07 Page 16 of 17

## Appendix 5

### Shipping Test Log Sheet

#### ULPA FILTERS

Preconditioning:

Start Date: 3-21-2014 Chamber Number: 01095

Completion Date: 3-24-2014 Last Calibration: 3-28-2013

Signature/Date: Paul Valpreda 3-24-14 Calibration due: 3-31-2014

Drop Test:

Catalog \_\_\_\_\_ Weight 31bs Drop Height: 15.5 inches

Drop	Orientation	Specific face, edge or	Initials/Date
1	Top	Face 1	PV - 3-25-2014
2	Edge	Edge 5-3	PV - 3-25-2014
3	Edge	Edge 6-3	PV - 3-25-2014
4	Corner	Corner 2-3-5	PV - 3-25-2014
5	Corner	Corner 4-3-6	PV - 3-25-2014
6	Bottom	Face 3	PV - 3-25-2014

Comments: \_\_\_\_\_

Signature: Paul Valpreda Date: 3-25-2014

Compression Test:

Catalog \_\_\_\_\_ Pounds Force 120

Comments: \_\_\_\_\_

Signature: Paul Valpreda Date: 3-25-2014

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	<b>TEST PROTOCOL</b> <b>MASTER DOCUMENT</b>	<b>Document Number</b> <u>X1150751-10</u>
	<b>ULPA REPLACEMENT FILTER &amp; CARBON/MUFFLER</b>	<b>Revision:</b> 01
		<b>Effective Date:</b> 2014 APR 07
		<b>Page 17 of 17</b>

## Appendix 5 Continued

### Shipping Test Log Sheet

Vibration:

Low Frequency, 40 minutes, Initials PV High frequency 10 minutes, Initials PV

Completion Date: 3-25-2014

Signature: Paul Valpreda Date: 3-25-2014

Second Drop Test:

Catalog \_\_\_\_\_ Weight 3 lbs Drop Height: 15.5 + 30 inches

Drop	Orientation	Specific face, edge or	Initials/Date
1	Edge	Edge 4-6	<u>PV</u> 3-25-2014
2	Face	Face 4	<u>PV</u> 3-25-2014
3	Face	Face 6	<u>PV</u> 3-25-2014
4	Corner	Corner 2-1-5	<u>PV</u> 3-25-2014
5	Edge	Edge 2-1	<u>PV</u> 3-25-2014
6	Bottom	Face 3, Increase height to 30 inches.	<u>PV</u> 3-25-2014

Comments: \_\_\_\_\_

Signature: Paul Valpreda Date: 3-25-2014