



The Electrosurgical Authority®

DOCUMENT NUMBER: ENG-RPT-403

DOCUMENT TITLE: Zip Pen Flow Evaluation Test Report

DOCUMENT NOTES:

#### Document Information

Revision: 002

Vault: MEG-Rel

Status: Release

Document Type: ENG-RPT

#### Date Information

Effective Date: 21 Feb 2018

Expiration Date:

Release Date: 21 Feb 2018

Next Review Date:

#### Control Information

Author: MGLASSETT

Owner: MGLASSETT

Previous Number:

Change Number: 2018-ENG-DCO-006

**Signature Manifest****Document Number:** ENG-RPT-403**Revision:** 002**Title:** Zip Pen Flow Evaluation Test Report

All dates and times are in Mountain Standard Time.

**ENG-RPT-403 Zip Flow Testing****Change Request**

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		09 Jan 2018, 02:59:35 PM	Approved

**Collaboration**

Name/Signature	Title	Date	Meaning/Reason
Joni Stegeman (JSTEGEMAN)	Ethicon Quality	07 Feb 2018, 09:43:27 AM	Complete
Dave Shimkus (DSHIMKUS)		07 Feb 2018, 03:30:44 PM	Complete
Paul Borgmeier (PBORGMEIER)		12 Feb 2018, 03:56:59 PM	Complete
Darlene Hull (DHULL)	Regulatory	20 Feb 2018, 06:02:45 AM	Complete
Mallory Schroeder (MSCHROEDER)	Engineer	20 Feb 2018, 10:55:13 AM	Complete

**Document Review**

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		20 Feb 2018, 11:01:45 AM	Complete

**RA-Approval**

Name/Signature	Title	Date	Meaning/Reason
Darlene Hull (DHULL)	Regulatory	20 Feb 2018, 12:26:33 PM	Approved

**QA-Approval**

Name/Signature	Title	Date	Meaning/Reason
Joni Stegeman (JSTEGEMAN)	Ethicon Quality	20 Feb 2018, 06:36:09 PM	Approved

**ENG-Approval**

Name/Signature	Title	Date	Meaning/Reason
Paul Borgmeier (PBORGMEIER)		20 Feb 2018, 11:07:52 AM	Approved

**Training Review**

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		21 Feb 2018, 07:55:16 AM	Approved

**Final Release**

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		21 Feb 2018, 08:13:10 AM	Approved

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number ENG-RPT-403</b>
	<b>Zip Pen Flow Evaluation Test Report</b>	<b>Revision: 002</b>
		<b>Page 1 of 15</b>

Authored By: Mark Glassett

Revised By: Mallory Schroeder (M.S.)

<b>TABLE OF CONTENTS</b>	<b>Page</b>
1. ABSTRACT .....	1
2. OBJECTIVE .....	1
3. RESULTS .....	1
4. DISCUSSION .....	7
5. CONCLUSIONS .....	8
6. RECOMMENDATIONS .....	8

## 1. ABSTRACT

Zip Pen Catalog # 2525-15 samples were tested per the requirements of ENG-PRT-280 to evaluate flow characteristics with different smoke evacuator devices. The performance of the Zip Pen was compared with the predicate device 2110-10 UltraVac (Megadyne branded version of the IC Medical Penevac). The Zip Pen met the requirement of the protocol of having equal or greater flow than the predicate. The reaction of the Zip Pen under occlusion conditions was also observed. There was no damage to the product and the tubing does not collapse radially under occlusion.

## 2. OBJECTIVE

The objective of this test report is to document equivalence of flow rate of the Zip Pen with the predicate UltraVac. The testing also determined an upper flow rate value where the Zip Pen can be safely used. The report documents the reaction of the device under occlusion to insure it is not damaged and remains functional after an occlusion event.

## 3. RESULTS

### 3.1. Zip Pen Flow Test


The Zip Pen was tested for flow in comparison with the UltraVac using the following smoke evacuator devices set at maximum flow.

Megadyne Mega Vac Plus  
Buffalo Filter Viro Vac

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-403
	Zip Pen Flow Evaluation Test Report	Revision: 002
		Page 2 of 15


Conmed AER Defense  
Lina Safe Air SFR-0200

3.1.1. **Mega Vac Plus:** The 't' test analysis below shows that the Zip Pen has higher flow than the Ultra Vac when attached to the Mega Vac Plus and therefore meets the requirements of the protocol.




<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number ENG-RPT-403</b>
	<b>Zip Pen Flow Evaluation Test Report</b>	<b>Revision: 002</b>
		<b>Page 3 of 15</b>

- 3.1.2 **Buffalo Filter Viro Vac:** The 't' test analysis below shows that the Zip Pen has higher flow than the Ultra Vac when attached to the Buffalo Filter Viro Vac and therefore meets the requirements of the protocol.




Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-403
	Zip Pen Flow Evaluation Test Report	Revision: 002
		Page 4 of 15

- 3.1.3 **Conmed AER Defense:** The 't' test analysis below shows that the Zip Pen has higher flow than the Ultra Vac when attached to the Conmed AER Defense and therefore meets the requirements of the protocol.



<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number ENG-RPT-403</b>
	<b>Zip Pen Flow Evaluation Test Report</b>	<b>Revision: 002</b>
		<b>Page 5 of 15</b>

3.1.4 **LINA Safe Air SFR-0200:** The 't' test analysis below shows that the Zip Pen has higher flow than the Ultra Vac when attached to the LINA Safe Air SFR-0200 and therefore meets the requirements of the protocol.



### 3.1.5 Data Summary

Smoke Evacuator	2525-15 Zip Pen Flow Avg.	Ultra Vac Flow Avg.
Mega Vac Plus	71.7 lpm	67.0 lpm
Buffalo Filter Viro Vac	91.0 lpm	71.7 lpm
Conmed AER Defense	60.7 lpm	46.6 lpm
Lina Safe Air	89.0 lpm	81.9 lpm



<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number ENG-RPT-403</b>
	<b>Zip Pen Flow Evaluation Test Report</b>	<b>Revision: 002</b>
		<b>Page 6 of 15</b>

### 3.2. Maximum Flow Rate Evaluation

Prior to Zip Pen testing the maximum flow of each smoke evacuator with filter was checked. This data is shown in comparison to the Zip Pen average flow rate for that smoke evacuator. The measured values are as follows:

<b>Smoke Evacuator</b>	<b>Max Open Flow</b>	<b>2525-15 Zip Pen Flow Avg.</b>
Mega Vac Plus	78.5 lpm	71.7 lpm
Buffalo Filter Viro Vac	358.2 lpm	91.0 lpm
Conmed AER Defense	258.7 lpm	60.7 lpm
Lina Safe Air	100.6 lpm	89.0 lpm

The above values illustrate a large range of flow rates for the different smoke evacuator devices. They all use different technologies for the vacuum source but when the Zip Pen (or competitive smoke evacuation pencil) is attached the restriction of the device lowers the flow down to a much lower value. Note that the higher flow rate of the smoke evacuator does not correlate linearly with the flow rate of the Zip Pen. This no linearity is due to the ability of the smoke evacuator to produce vacuum pressure. The Mega Vac Plus and the Lina Safe Air can produce more vacuum pressure and therefore the drop in flow from the Max Open Flow to the flow with the Zip Pen is small. Conversely, the Buffalo Filter Viro Vac and the Conmed AER Defense only produce a small amount of vacuum and therefore the drop in flow from the Max Open Flow to the flow with the Zip Pen is large.

### 3.3. Tubing Occlusion Evaluation

For each smoke evacuator, the effect on the Zip Pen during occlusion was observed. The following comments were recorded.

**Mega Vac Plus:** Tubing collapsed and accorded , no kinking or blockages

**Buffalo filter Viro Vac:** Tubing collapsed slightly and accorded , no kinks or blocks

**Conmed AER Defense:** Movement in tubing barely noticeable when occluded

**Lina Safe Air SFR-0200:** Tubing accorded quite strongly, but did not kink or impede flow in any way.

Reviewing the comments with the technician, the tubing did not collapse radially but rather along the axis in an accordion motion. The accordion condition relaxes back to the normal tubing length when the occlusion is discontinued. This is the

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number ENG-RPT-403</b>
	<b>Zip Pen Flow Evaluation Test Report</b>	<b>Revision: 002</b>
		<b>Page 7 of 15</b>

expected result for convoluted tubing. There was no damage to the Zip Pen as a result of occlusion, therefore the Zip Pen passes the protocol.

#### 4. DISCUSSION

Zip Pen catalog number 2525-15 with 15 ft tubing Lot S140045 were used for this testing. The Zip Pen with 15 foot tubing is the worst case representation for the Zip Pen in terms of flow because this is the longest available tubing. These also represent the 2525-10 product that has 10 foot tubing. The comparison device was Ultra Vac 2110-10 Lot 5408 with ten foot tubing.

Prior to testing, the Zip Pen samples were gamma sterilized to a minimum of 50 kGy. The samples were also artificially aged at 55°C for 111 days to simulate three years shelf life. The samples were exposed to thermal cycling from -40°C to 70 °C to simulate worst case shipping environment. The samples were preconditioned for ship testing and then subjected to the ship test. These conditions were performed to simulate worst case sterilization and handling prior to reaching the customer. Documentation of this treatment for this lot of samples is shown in ENG-RPT-377.

##### 4.1. Zip Pen Flow Test

The Zip Pen was tested per the protocol for flow in comparison with the UltraVac using common smoke evacuators found in the market. The flow data with Zip Pen shows that it has higher flow characteristics of the two devices. This is expected because the Zip Pen has a slightly larger minimum cross sectional area than the Ultra VAC.

##### 4.2. Maximum Flow Rate Evaluation.

This evaluation was done to show that with the large variation of flow values from the smoke evacuator devices, the flow through the Zip Pen is still equal to or greater than the Ultra Vac predicate. The flow with the smoke evacuation pencil attached meets the requirement of the protocol.

##### 4.3. Tubing Occlusion Evaluation

The tubing occlusion evaluation demonstrates that the tubing remains patent during and after occlusion. There was no damage to the Zip Pen as a result of occlusion.

<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number ENG-RPT-403</b>
	<b>Zip Pen Flow Evaluation Test Report</b>	<b>Revision: 002</b>
		<b>Page 8 of 15</b>

## 5. CONCLUSIONS

This testing demonstrates that the Zip Pen catalog numbers 2525-10 and 2525-15 meet the requirement that the flow is equal to or greater than the Ultra Vac 2110-10 when used with the Megadyne Mega Vac Plus smoke evacuator.

This testing demonstrates that the Zip Pen has equivalent or greater flow rate than the Ultra Vac when used with competitive smoke evacuator devices.

This testing demonstrates that the Zip Pen can be used with smoke evacuator devices that have flow rates as high as 358 lpm.

This testing demonstrates that the Zip Pen is not damaged when fully occluded.

## 6. RECOMMENDATIONS

This test report will be archived in the Zip Pen Design History Record.

## 7. REVISION HISTORY DESCRIPTION

The purpose of revision 002 of this document is to address a typo in the Appendix data sheets. The operator mistakenly recorded 2525-10 instead of 2525-15 in the catalog number field. Appendix V provides further evidence that lot number S140045 is catalog number 2525-15.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-403
	Zip Pen Flow Evaluation Test Report	Revision: 002
		Page 9 of 15

**Appendix I:**  
**ZIP PEN FLOW TEST LOG SHEET**  
**Mega Vac Plus**

The user must ensure that they are using the correct/current revision of this document.  
Document: XENG-PRT-280 Rev: A Effective: 04 Nov 2014 11:45 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-280
	Zip Pen Flow Evaluation Protocol	Revision: A
		Page 6 of 9

M.S. 1/8/18 Fix incorrectly  
recorded catalog number

**Appendix I:**  
**Zip Pen Flow Test Data Sheet**  
**Mega Vac Plus**

2525-15

Zip Pen Catalog# 2525-10 Lot# S140045

Ultra Vac Catalog# 2110-10 Lot# 5408

Mega Vac Plus S/N 14211

ULPA Filter Lot# 5058

Maximum Flow of the Mega Vac Plus 78.48

Zip Pen Flow Test		Control Samples Flow Test	
Sample #	Zip Pen 2525-15	Sample #	Ultra Vac 2110-10
	Flow		Flow
1.	72.05	1.	65.50
2.	71.57	2.	67.95
3.	71.63	3.	66.75
4.	71.82	4.	67.50
5.	71.82	5.	66.95
6.	71.52	6.	66.75
7.	71.69	7.	66.50
8.	72.02	8.	66.90
9.	71.53	9.	68.10
10.	71.70	10.	66.85
11.	71.40	11.	67.30

Comments on Zip Pen Occlusion: Tubing collapsed and accorded, no  
kinking or blockages.

Test Performed by: Paul Valpreda Date: 11-5-2014

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 04 Nov 2014, 11:47:02 am; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:55:47 pm; Printed by: .

Confidential - Controlled copy on date printed

Confidential - Pre-Production Controlled on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-403
	Zip Pen Flow Evaluation Test Report	Revision: 002
		Page 10 of 15

## Appendix II ZIP PEN FLOW TEST LOG SHEET Buffalo Filter Viro Vac

The user must ensure that they are using the correct/current revision of this document.  
Document: XENG-PRT-280 Rev: A Effective: 04 Nov 2014 11:45 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-280
	Zip Pen Flow Evaluation Protocol	Revision: A
		Page 7 of 9

M.S. 1/8/18 Fix incorrectly recorded catalog number

### Appendix II: Zip Pen Flow Test Data Sheet Buffalo Filter Viro Vac

2525-15  
Zip Pen Catalog# 2525-10 Lot# 5140045  
Ultra Vac Catalog# 2110-10 Lot# 5408  
Viro Vac S/N V00325  
Maximum Flow of the Buffalo Filter Viro Vac 358.20

Zip Pen Flow Test		Control Samples Flow Test	
Sample #	Zip Pen 2525-15	Sample #	Ultra Vac 2110-10
	Flow		Flow
1.	90.4	1.	68.5
2.	91.2	2.	71.6
3.	90.4	3.	71.9
4.	91.8	4.	72.6
5.	90.3	5.	73.3
6.	90.8	6.	70.8
7.	90.4	7.	71.6
8.	95.0	8.	72.5
9.	89.7	9.	73.8
10.	90.8	10.	70.6
11.	90.2	11.	71.7

Comments on Zip Pen Occlusion: Tubing collapsed slightly and accordioned.  
No kinks or blocks.

Test Performed by: Paul Valpreda Date: 11-5-2014

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 04 Nov 2014, 11:47:02 am; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:55:47 pm; Printed by: .

Confidential - Controlled copy on date printed

Confidential - Pre-Production Controlled on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-403
	Zip Pen Flow Evaluation Test Report	Revision: 002
		Page 11 of 15

Appendix III  
ZIP PEN FLOW TEST LOG SHEET  
Conmed AER Defense

The user must ensure that they are using the correct/current revision of this document.  
Document: XENG-PRT-280 Rev: A Effective: 04 Nov 2014 11:45 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-280
	Zip Pen Flow Evaluation Protocol	Revision: A
		Page 8 of 9

M.S. 1/8/18 Fix incorrectly  
recorded catalog number

Appendix III:  
Zip Pen Flow Test Data Sheet  
Conmed AER Defense

Zip Pen Catalog# 2525-15 Lot# 5140045  
Ultra Vac Catalog# 2110-10 Lot# 5408  
AER Defense S/N 12DCA001  
Maximum Flow of the Conmed AER Defense 258.7

Zip Pen Flow Test		Control Samples Flow Test	
Sample #	Zip Pen 2525-15	Sample #	Ultra Vac 2110-10
	Flow		Flow
1.	60.9	1.	44.5
2.	60.2	2.	46.3
3.	60.3	3.	47.2
4.	60.2	4.	47.1
5.	60.2	5.	47.8
6.	60.9	6.	46.2
7.	60.6	7.	46.9
8.	63.4	8.	47.1
9.	60.1	9.	47.8
10.	61.1	10.	45.1
11.	60.1	11.	46.2

Comments on Zip Pen Occlusion: Movement in tubing barely noticable  
when occluded.

Test Performed by: Paul Valpreda Date: 11-5-2014

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 04 Nov 2014, 11:47:02 am; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:55:47 pm; Printed by: .

Confidential - Controlled copy on date printed

Confidential - Pre-Production Controlled on date printed



Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-403
	Zip Pen Flow Evaluation Test Report	Revision: 002
		Page 12 of 15

### Appendix IV ZIP PEN FLOW TEST LOG SHEET Lina Safe Air SFR-0200

The user must ensure that they are using the correct/current revision of this document.  
Document: XENG-PRT-280 Rev: A Effective: 04 Nov 2014 11:45 AM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-280
	Zip Pen Flow Evaluation Protocol	Revision: A Page 9 of 9

M.S. 1/8/18 Fix incorrectly recorded catalog number

### Appendix VI: Zip Pen Flow Test Data Sheet Lina Safe Air SFR-0200

2525-15  
Zip Pen Catalog# 2525-10 Lot# 5140045  
Ultra Vac Catalog# 2110-10 Lot# 5408  
Safe Air S/N 13035-3  
Maximum Flow of the Lina Safe Air SFR-0200 100.6

Zip Pen Flow Test		Control Samples Flow Test	
Sample #	Zip Pen 2525-15	Sample #	Ultra Vac 2110-10
	Flow		Flow
1.	89.0	1.	80.8
2.	89.0	2.	81.7
3.	88.2	3.	81.5
4.	89.6	4.	82.2
5.	88.0	5.	82.6
6.	90.1	6.	81.3
7.	89.3	7.	81.8
8.	89.2	8.	82.6
9.	89.2	9.	82.6
10.	88.8	10.	81.5
11.	88.6	11.	82.2

Comments on Zip Pen Occlusion: Tubing accordioned quite strongly, but did not kink or impede flow in any way.

Test Performed by: Paul Valpreda Date: 11-6-2014

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 04 Nov 2014, 11:47:02 am; Printed by: MGLASSETT.

### Appendix V

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:55:47 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-403
	Zip Pen Flow Evaluation Test Report	Revision: 002
		Page 13 of 15

2525-15 Lot Number S140045



### Certificate of Compliance

<p><b>Taiwan</b></p> <p>天德企業股份有限公司 New Deantronics Taiwan, Ltd. 新北市民權路中角段4段11號12樓 12F, No 51, Sec 4, Chong Yang Rd., Tu Chung Dist., New Taipei City 23675, Taiwan R.O.C. Tel: +886 2 2250-1728 Fax: +886 2 2250-0000</p>	<p><b>Customer Name:</b> MEGADYNE MEDICAL PRODUCTS, INC. <b>Invoice Number :</b> MD-441/14 <b>P. O. Number:</b> 23954 <b>Customer P/N:</b> X2525-15 <b>Drawing Number:</b> X2525-15 Rev.02 <b>New Deantronics P/N:</b> PB752SM1 <b>Lot Number:</b> S140045 <b>Expiration Date:</b> N/A <b>Quantity:</b> 40 Pieces <b>Carton Number:</b> #1 ~ #2 (2 Cartons)</p>	<p><b>Signed:</b> <i>D. Y. Chen</i> <b>Printed Name:</b> Da-Yu Chen <b>Title:</b> Q.A. Director</p>	<p><b>Date:</b> 03-05-2014</p>
---	---	---	--------------------------------

**USA**

New Deantronics Ltd.  
1990 North California Blvd.  
Suite 1040  
Walnut Creek, CA, 94596  
Tel: +1 (925) 250-0300  
Fax: +1 (925) 250-1788

#### Note: Aging sample for sealing test

<u>Materials</u>		<u>ND</u> <u>Mat'l Lot#</u>	<u>Megadyne</u> <u>Lot#</u>
Cable	C302600 Conductor : 7/0.16*3BC(26 AWG) Bare copper Insulation: PPE, one red, one blue & one white Outer jacket: PVC, Gray	131014Z	
Plug material	F505100, TOP, ABS PA707 F505200, Bottom, ABS PA707	130917 130917	
Terminal	T101702 Brass	131030	
PCB Overmold	R900801, TPR	130801	
Swivel, male	F916600, HDPE	S130301	
Swivel, Female	F917000, HDPE	S130301	
Connector Proximal	F916800, LDPE	131231	
Collet Terminal	T202502, Phosphor Bronze contact plated nickel	130904	
Button	F302500, NYLON 66	131217	
Tape	A100401, PTFE	131002	
PCB	H201602, PCB FR-4	140211	
Dome	M100600, 6mm ROUND W/DIMPLE & FEET	130820	
Pen Body	F105400, ABS+TPR	131217	



<b>Megadyne Medical Products, Inc.</b>	<b>TEST REPORT</b>	<b>Document Number ENG-RPT-403</b>
	<b>Zip Pen Flow Evaluation Test Report</b>	<b>Revision: 002</b>
		<b>Page 14 of 15</b>



<b>Taiwan</b>				
大陆金丰机械有限公司	Nozzle	F916300, PC	131217	
New Deantronics Taiwan, Ltd.	Collet Holder	F918100, PC	131217	
新北市土城區中央路4段31號10樓	Carriage	F916500, ABS	131217	
12F., No.51, Sec.4,	Snap Swivel, Male	F916400, HDPE	S130302	
Chong Yang Rd., Tu Cheng Dist.,	Snap Swivel, Female	F918300, HDPE	S130302	
New Taipei City 23675,	Tubing, Connector	F916700, HDPE	S130318	
Taiwan R.O.C.	Tubing, Convolutd	P305600, EVA 8.5"	130828	
Tel: +886 2 2258-1726	Tubing, Convolutd	P305700, EVA 56"	130828	
Fax: +886 2 2258-3800	Holster	F916900, HDPE	131220	
-----	Blade	G102700, Coated Megadyne P/N: 0012BN5 (Provided by Megadyne)	140121 (40pcs)	134585
<b>USA</b>	Paper Band	A900300	140116	
New Deantronics Ltd.	Tyvek	A000400	140116	
1990 North California Blvd.	Nylon Film	A000500	131217	
Suite 1040	Glue	S400900, Loctite 4061	130409	
Walnut Creek, CA 94596	Ink	S102000, Green, PMS356C	121127	
Tel: +1 (925) 250-8388	IFU PN and REV	P/N3000185-01 Rev.A	130910	
Fax: +1 (925) 250-1788				

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-403
	Zip Pen Flow Evaluation Test Report	Revision: 002
		Page 15 of 15

03/05/2014/03 03:52 PM

FAX No.

P. 005



中國生化科技股份有限公司

CHINA BIOTECH CORPORATION

誠信正信 誠意服務 誠心公益

TEL: 886-4-23597515 FAX: 886-4-23597080

台中市工業區33路10號

10, 33<sup>rd</sup> Road, Taichung Industrial Park,  
Taichung, Taiwan R.O.C 407

DATE: 2014/3/5

照射證明書  
CERTIFICATE OF IRRADIATION

REVISED

行政院原能會核准設立照射廠執照證號 IRRADIATION PLANT NO: 物字第 1100223 號

客戶名稱 CUSTOMER NAME: 大瑞企業(股)公司  
NEW DEANTRONICS TAIWAN LTD.  
照射日期 IRRADIATION RUN DATE: 2014/02/27

照射批號 IRRADIATION RUN NUMBER: NEW14134-J1

客戶產品已照射 MATERIALS PROCESSED:

箱數 CASE	內容 DESCRIPTION	客戶產品批號 LOT NO
1 (40PCS)	X2525-15 PB752SM1	S140045

總數 TOTAL 1 箱數 CASE

中國生化科技股份有限公司證明上述產品經本公司劑量偵測系統判讀，吸收劑量如下：  
China Biotech Corporation certifies that the material listed above (has described by its manufacturer)  
received the following doses within the precision limits of the dosimetry system employed

最低劑量 MINIMUM DOSAGE 51.7 kGy; 最高劑量 MAXIMUM DOSAGE 54.2 kGy

使用放射性同位素 ISOTOPE UTILIZED: 鈾 60 COBALT-60

客戶劑量要求 DOSE REQUIREMENT: 最低劑量 MIN 50.0 kGy; 最高劑量 MAX 60.0 kGy

確認者:   
CERTIFIED BY 品保部主管  
QUALITY ASSURANCE

確保人類健康 · 珍惜自然環境