

| | | |
|--|---------------------------------------|--|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | Document Number ENG-RMF-045 |
| | Smoke Evacuation Risk Analysis | Revision: D |
| | | Page 1 of 52 |

1. REFERENCES

| | |
|-------------------|--|
| EN ISO 14971:2019 | Medical devices - Application of risk management to medical devices |
| QA-SOP-015* | Risk Management of Medical Devices |
| ENG-RMF-044 | Mega Soft Pad Family Hazard Assessment Summary |
| ENG-DMR-012 | Megadyne Smoke Evacuation Pencil and Accessories, Device Master Record |
| 100583575 | Franchise Global Complaint Trending and Signal Detection PMS Procedure |
| RA-RPT-005 | Megadyne Family of products Clinical Evaluation Report |
| RA-RPT-007 | Megadyne ACE Blade & ACE Blade 700™ Soft Tissue Dissector Clinical Evaluation Report (Product Codes: ME7251C, ME7251E, ME725M1C, ME725M1E) |

*Note: This document follows QA-SOP-015 Rev 005. Per QA-SOP-015 section 2.4.1, "Risk documentation for devices implemented under revisions prior to Rev 006 of QA-SOP-015 can be maintained according to the version of QA-SOP-015 in effect when the document was created."

2. SCOPE

This Risk Analysis applies to Smoke Evacuation Accessory products. The products covered by this document are established in the Device Master Record ENG-DMR-012 and Megadyne's Design History File.

This Risk Analysis pertains to the Zip Pen products, including cat numbers 2525-10 (252510), 2525-10EC (252510EC), 2525-10BN (252510BN), 2525-10ECBN (252510ECBN), 2525-15EC (252515EC), and 2525-15 (252515), ME7251C, ME7251E, ME725M1C, ME725M1E. This document also pertains to the Extension Nozzle cat numbers 2540 (2540J) and 2560 (2560J), Charcoal Filter 2220 (2220J), and ULPA Filter 2211 (2211J). The catalog numbers in parentheses are new catalog numbers referring to the same product, due to transition of product information from Megadyne Medical Products to Ethicon Endo Surgery. This document will refer to the original catalog numbers.

3. PURPOSE

The document summarizes the risks associated with these devices and establishes how they have been addressed.

4. COMPLAINT HISTORY

| | | |
|--|---------------------------------------|--|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | Document Number ENG-RMF-045 |
| | Smoke Evacuation Risk Analysis | Revision: D |
| | | Page 2 of 52 |

4.1. Summary

This complaint history includes updated reported complaints received from JULY 2018 through JUNE 2020 completed as a periodic review of the risk management file. Complaint numbers are presented, as are complaint frequencies. The frequency of complaints is expressed as the number of complaints versus the number of units distributed during the above time frame. All of the data for the periodic review (JULY 2018 – JUNE 2020) was obtained from the PQSS DRB MMP Rates and Signals spreadsheet. The complaint classification codes, and their obtained Rates are presented in the below table.

| Smoke Evacuation- 252510, 2525-10, 252510EC, 2522510-EC | | | | |
|---|-------------------|---------------|------|----------------|
| Product Experience Code (PEC) | No. of Complaints | Opportunities | CPMO | Complaint rate |
| ACTIVATION FAILURE | 3 | 145657 | 20.6 | 0.0021% |
| SWITCH BUTTON ISSUE | 3 | 145657 | 20.6 | 0.0021% |
| COAG ISSUE | 3 | 145657 | 20.6 | 0.0021% |
| TISSUE EFFECT ISSUE | 2 | 145657 | 13.7 | 0.0014% |
| ASPIRATION ISSUE | 2 | 145657 | 13.7 | 0.0014% |
| DAMAGED PRODUCT | 2 | 145657 | 13.7 | 0.0014% |
| SELF ACTIVATION | 1 | 145657 | 6.9 | 0.0007% |
| FLAKING ISSUE | 1 | 145657 | 6.9 | 0.0007% |
| CAUTERY ISSUE | 1 | 145657 | 6.9 | 0.0007% |
| FLAME FLASH FIRE | 1 | 145657 | 6.9 | 0.0007% |
| NOT SPECIFIED | 1 | 145657 | 6.9 | 0.0007% |
| CUTTING ISSUE | 0 | 145657 | 0 | 0% |

| | | |
|--|---------------------------------------|--|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | Document Number ENG-RMF-045 |
| | Smoke Evacuation Risk Analysis | Revision: D |
| | | Page 3 of 52 |

| Smoke Evacuation ACE 700 - ME7251C, ME7251E, ME725M1C, ME725M1E | | | | |
|---|-------------------|---------------|-------|-------------------|
| Product Experience Code (PEC) | No. of Complaints | Opportunities | CPMO | Complaint Rate |
| ASPIRATION ISSUE | 4 | 23969.5 | 166.9 | 0.017% |
| DAMAGED PRODUCT | 4 | 23969.5 | 166.9 | 0.017% |
| CAUTERY ISSUE | 3 | 23969.5 | 125.2 | 0.013% |
| ACTIVATION ISSUE | 3 | 23969.5 | 125.2 | 0.013% |
| TISSUE EFFECT ISSUE | 2 | 23969.5 | 83.4 | 0.008% |
| NOT SPECIFIED | 2 | 23969.5 | 83.4 | 0.008% |
| HEMOSTASIS CONTROLLABLE | 2 | 23969.5 | 83.4 | 0.008% |
| COAG ISSUE | 2 | 23969.5 | 83.4 | 0.008% |
| ASSEMBLY DISASSEMBLY ISSUE | 1 | 23969.5 | 41.7 | 0.004% |
| SWITCH BUTTON ISSUE | 1 | 23969.5 | 41.7 | 0.004% |
| CUTTING ISSUE | 1 | 23969.5 | 41.7 | 0.004% |
| UNINTENDED THERMAL INJURY LESS THAN 2ND DEGREE | 1 | 23969.5 | 41.7 | 0.004% |
| UNRETRIEVED DEVICE FRAGMENT | 1 | 23969.5 | 41.7 | 0.004% |
| CUSTOMER DISSATISFACTION | 1 | 23969.5 | 41.7 | 0.004% |
| IGNITION OF DEBRIS OR FLAMMABLE CONDITION | 1 | 23969.5 | 41.7 | 0.004% |
| FLAME FLASH FIRE | 1 | 23969.5 | 41.7 | 0.004% |
| FLAKING ISSUE | 0 | 23969.5 | 0 | 0% |

| | | |
|--|---------------------------------------|--|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | Document Number ENG-RMF-045 |
| | Smoke Evacuation Risk Analysis | Revision: D |
| | | Page 4 of 52 |

| Smoke Evacuation Accessory – 2220, 2220J, 2140J, 2255J | | | | |
|--|-------------------|---------------|--------|-------------------|
| Product Experience Code (PEC) | No. of Complaints | Opportunities | CPMO | Complaint rate |
| PACKAGING IDENTIFICATION | 101 | 75223 | 1342.7 | 0.134% |
| PACKAGING DEVICE ISSUE | 5 | 75223 | 66.5 | 0.007% |
| NOT SPECIFIED | 1 | 75223 | 13.3 | 0.001% |

5. CLINICAL EVALUATION

The review of scientific literature is documented in RA-RPT-005 & RA-RPT-007. Megadyne believes, and the literature supports, the Smoke Evacuation accessories are safe products, with no unacceptable risks.

6. RISK SUMMARY

Risk analysis study of this product family was conducted according to QA-SOP-015 Risk Management of Medical Devices. A copy of the detailed Failure Mode and Effect Analysis (FMEA) for the Smoke Evacuation accessories follows in Section 8. All identified risks have been deemed acceptable or as low as possible (see Notes / Justification column of FMEA).

None of the risks identified in the Risk Analyses are of a nature that they are confined to Megadyne's Smoke Evacuation family of devices. All smoke evacuator manufacturers face approximately the same set of risks and have addressed them in very similar ways.

Megadyne's Smoke Evacuation devices present no risks to patients or users that are unacceptable or unreasonable when weighed against the benefits to the patient.

7. RISK BENEFIT ANALYSIS

Megadyne's Smoke Evacuation family of accessories is similar in many ways to existing smoke evacuation accessories that have a well-established history of safety and reliable performance through use in millions of cases in thousands of operating rooms throughout the world. Many of these devices have been in general distribution for decades. Some of the benefits provided by the Megadyne devices when compared with other competitive products as follows:

- the Zip-Pen Smoke Evacuation Pencil has ergonomic design for ease of surgeon use
- the smoke nozzle is contoured for improved surgeon visibility
- the Zip-Pen Smoke Evacuation Pencil incorporates swivels for tubing attachment to minimize drag on the surgeon's hand.

| | | |
|--|---------------------------------------|--|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | Document Number ENG-RMF-045 |
| | Smoke Evacuation Risk Analysis | Revision: D |
| | | Page 5 of 52 |

These devices are intended for use by or on the order of a physician who is trained in the associated risks of electrosurgery. The anticipated clinical benefits of using these high-frequency surgical devices as compared to other surgical techniques outweigh the residual and overall risks identified in the FMEA. As with other electrosurgical smoke evacuation devices on the market, these products meet international standards for safety. This family of devices presents no risks to patients or users that are unacceptable or unreasonable when weighed against the benefits to the patient. Risks of electrosurgery are reviewed with patients by the physician prior to through informed consent.

8. FMEA

The Failure Mode and Effect Analysis (FMEA) for the smoke evacuation product family are detailed in the following table:

DRAFT

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 6 of 52 |

Zip-Pen aFMEA

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|---|---|---|---|--|-----------------|---------------------------|-----|-----------------------------------|--------------|---|--------------------|---|-----------------|--------------------|----------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | Risk reduction considered? Yes/No | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 1-a | Break seal to open outer Pencil package | To deliver aseptic product to sterile field | Package difficult to open, defective package or package materials | Inner bag gets contaminated (recognized) and contaminates sterile field | Causes a delay, must re-establish sterile field | 1 | 1 | 1 | NA | DFS | Package designed with same sealant as existing packaging and with chevron tab for intuitive opening. See footnote 1 | Not Practical | Verification Plan: Verify that the package drawing indicates “OPEN HERE” on tab. Validation Plan: ENG-PRT-290, Ensure IFU includes opening instructions. | See ENG-DMR-012 | Report ENG-RPT-418 | 1 | 1 | 1 | NA | NA | NA |
| 1.5-a | Break seal to open outer Pencil package | To deliver aseptic product to sterile field | Package difficult to open, defective package or package materials | Inner bag gets contaminated (recognized) and contaminates sterile field | Causes a delay, must re-establish sterile field | 1 | 1 | 1 | NA | IFS | Package says “OPEN HERE” on tab. See footnote 1 | Not Practical | Verification Plan: Verify that the package drawing indicates “OPEN HERE” on tab. Validation Plan: ENG-PRT-290, Ensure IFU includes opening instructions | See ENG-DMR-012 | Report ENG-RPT-418 | 1 | 1 | 1 | NA | NA | NA |
| 2-a | Break seal to open outer Pencil package | To deliver aseptic product to sterile field | Package difficult to open, defective package or package materials | Inner bag gets contaminated (not recognized) and contaminates sterile field | Non-sterile product used on patient resulting in infection | 10 | 1 | 10 | Yes | DFS | Package designed with same sealant as existing packaging and with chevron tab for intuitive opening. See footnote 1 | Not Practical | NA | See ENG-DMR-012 | NA | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 2.5-a | Break seal to open outer Pencil package | To deliver aseptic product to sterile field | Package difficult to open, defective package or package materials | Inner bag gets contaminated (not recognized) and contaminates sterile field | Non-sterile product used on patient resulting in infection | 10 | 1 | 10 | Yes | IFS | Package says “OPEN HERE” on tab See footnote 1 | Not Practical | Verification Plan: Verify that the package drawing indicates “OPEN HERE” on tab, Ensure See ENG-PRT-290, IFU includes opening instructions | See ENG-DMR-012 | Report ENG-RPT-418 | 10 | 1 | 10 | No | N/A | See footnote 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 7 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Yes/No | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|---|---|---|--|--|-----------------|---------------------------|-----|---------------------------|--------------|---|--------------------|--|------------------------------------|--------------------|----------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 3-a | Break seal to open outer Pencil package | To deliver aseptic product to sterile field | Package difficult to open, defective package or package materials | Can break seal on wrong end, inner bag gets contaminated (recognized) and contaminates sterile field | Causes a delay, must re-establish sterile field | 1 | 1 | 1 | Yes | DFS | Package designed with same sealant as existing packaging and with chevron tab for intuitive opening. See footnote 1 | Not Practical | Validation Plan: ENG-PRT-290 | See ENG-DMR-012 and label drawings | Report ENG-RPT-418 | 1 | 1 | 1 | No | N/A | Risk is as low as possible. |
| 3.5-a | Break seal to open outer Pencil package | To deliver aseptic product to sterile field | Package difficult to open, defective package or package materials | Can break seal on wrong end, inner bag gets contaminated (recognized) and contaminates sterile field | Causes a delay, must re-establish sterile field | 1 | 1 | 1 | Yes | IFS | Package says “OPEN HERE” on tab See footnote 1 | Not Practical | Verification Plan: Verify that the package drawing indicates “OPEN HERE” on tab. | See ENG-DMR-012 and label drawings | Report ENG-RPT-418 | 1 | 1 | 1 | No | N/A | Risk is as low as possible. |
| 4-a | Break seal to open outer Pencil package | To deliver aseptic product to sterile field | Package difficult to open, defective package or package materials | Can break seal on wrong end, inner bag gets contaminated (not recognized) and contaminates sterile field | Non-sterile product used on patient resulting in infection | 10 | 1 | 10 | Yes | DFS | Package designed with same sealant as existing packaging and with chevron tab for intuitive opening. See footnote 1 | Not Practical | NA | See ENG-DMR-012 and label drawings | NA | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 4.5-a | Break seal to open outer Pencil package | To deliver aseptic product to sterile field | Package difficult to open, defective package or package materials | Can break seal on wrong end, inner bag gets contaminated (not recognized) and contaminates sterile field | Non-sterile product used on patient resulting in infection | 10 | 1 | 10 | Yes | IFS | Package says “OPEN HERE” on tab See footnote 1 | Not Practical | Verification Plan: Verify that the package drawing indicates “OPEN HERE” on tab. | See ENG-DMR-012 and label drawings | Report ENG-RPT-418 | 10 | 1 | 10 | No | N/A | See footnote 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 8 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Yes/No | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|--|---|---|--|---|-----------------|---------------------------|-----|---------------------------|--------------|---|--------------------|---|---|--------------------------------------|----------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 5-a | Scrub removes (picks) inside pouch from outer packaging | To deliver aseptic product to sterile field | The opening on the inside pouch is not oriented correctly | Can drop the product on floor | Non sterile product | 1 | 1 | 1 | NA | DFS | Ensure the flap is on the pouch side and not the top Tyvek side so the flap stays in the package when peeled open | Not Practical | Verify that the packaging design indicates correct pouch placement inside the outer package Validation Plan: ENG-PRT-290 | Bag is long enough to fold over creating a flap, Drawing shows flap with correct orientation See ENG-DMR-012 | Report ENG-RPT-418 | 1 | 1 | 1 | NA | NA | NA |
| 6-a | Pull Pencil out of inside bag | To deliver aseptic product to sterile field | Tubing unfurls | Drop product | Non sterile product | 1 | 1 | 1 | NA | DFS | Product is delivered to the sterile field inside of a secondary pouch to contain the tubing until technician is ready to connect tubing | Not Practical | Verify that the packaging design indicates correct pouch placement inside the outer package Validation Plan: ENG-PRT-290 | See ENG-DMR-012 | Report ENG-RPT-418 | 1 | 1 | 1 | NA | NA | NA |
| 7-a | Hand connecting ends to non-sterile person for connection to ESU and smoke box | Use ESU and smoke box | Tangled cable and tubing Cable not secured | Product contamination (recognized) | Causes a delay, need to replace the product | 1 | 1 | 1 | NA | DFS | The cable is contained in paper cohesive tape so it will not get tangled with the tubing | Not practical | NA | NA | NA | 1 | 1 | 1 | NA | NA | NA |
| 8-a | Hand connecting ends to non-sterile person for connection to ESU and smoke box | Use ESU and smoke box | Tangled cable and tubing Cable not secured | Product contamination (not recognized) | Non-sterile product used on patient resulting in infection | 10 | 5 | 50 | Yes | DFS | The cable is contained in paper cohesive tape so it will not get tangled with the tubing | Not practical | Ensure product is packaged per drawing Validation Plan: ENG-PRT-290 | See drawings 6020190-01 and 6020303-01, First Article | Report ENG-RPT-418 | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 9-a | Pencil Use | Evacuate smoke URS1006 | Tangled cable and tubing | Tubing is kinked or wrapped around the product, smoke not evacuated. | Possible user exposure to potential carcinogens, and infectious by-products | 5 | 1 | 5 | Yes | DFS | See footnote 1 and Designed similar to standard electrosurgical pencil | Not practical | Validation Plan: ENG-PRT-290 | NA | Usability test report ENG-RPT-418 | 5 | 1 | 5 | NA | NA | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 9 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------|---|---|--|---|-----------------|---------------------------|-----|-----------------|--------------|---|--------------------|---|--|--|----------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 9.5-a | Pencil Use | Evacuate smoke URS1006 | Tangled cable and tubing | Tubing is kinked or wrapped around the product, smoke not evacuated. | Possible user exposure to potential carcinogens, and infectious by-products | 5 | 1 | 5 | Yes | IFS | NA | Not practical | Verify the warning is captured in IFU Validation Plan: ENG-PRT-290 | Warning listed in Zip IFU, Zip Ace IFU | NA | 5 | 1 | 5 | NA | NA | See footnote 4 |
| 10-a | Pencil Use | Activation of pencil | Use error | Buzzing hemostat | Shock to clinician | 1 | 5 | 5 | Yes | DFS | See footnote 1 and Designed similar to standard electrosurgical pencil | Not practical | Validation plan: Review history of clinical use of pencil | NA | See footnote 2 | 1 | 5 | 5 | No | N/A | Risk is as low as possible. |
| 10.5-a | Pencil Use | Activation of pencil | Use error | Buzzing hemostat | Shock to clinician | 1 | 5 | 5 | Yes | IFS | NA | Not practical | Verify the caution to “Activate pencil at target tissue” is in IFU | Caution listed in Zip IFU, Zip Ace IFU | See footnote 2 | 1 | 5 | 5 | N | N/A | Risk is as low as possible. |
| 11-a | Pencil Use | Activation of pencil | Use error such as inappropriate contact of insulated shaft with off target tissue | Capacitive coupling Corona effect | Patient burn | 10 | 1 | 10 | Yes | DFS | See footnote 1 and Designed similar to standard electrosurgical pencil | Not practical | NA | NA | See footnote 2 | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 11.5-a | Pencil Use | Activation of pencil | Use error such as inappropriate contact of insulated shaft with off target tissue | Capacitive coupling Corona effect | Patient burn | 10 | 1 | 10 | Yes | IFS | NA | Not practical | Verify the caution to “Activate pencil at target tissue” is in IFU. | Caution listed in Zip IFU, Zip Ace IFU | See footnote 2 | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 12-a | Pencil Use | Between activations URS1002 | Use error such as failure to store device in holster | Contact with hot tip or inadvertent activation when electrode not in holster | Patient or user burn | 10 | 5 | 50 | Yes | DFS | Holster provided with each device. See footnote 1 Label according to IEC 60601-2-2 clause 201.7.9.2.2.101 a)5) | Not practical | Label according to IEC 60601-2-2 clause 201.7.9.2.2.101 a)5) | NA | See footnote 2 UL Test Report ENG-RPT-059 | 10 | 5 | 50 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 10 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Yes/No | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------|---|---|---|--|-----------------|---------------------------|-----|---------------------------|--------------|--|--------------------|---|---|--|----------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 12.5-a | Pencil Use | Between activations URS1002 | Use error such as failure to store device in holster | Contact with hot tip or inadvertent activation when electrode not in holster | Patient or user burn | 10 | 5 | 50 | Yes | IFS | Label according to IEC 60601-2-2 clause 201.7.9.2.2.101 a)5) | Not practical | Ensure labeling restricts use to physicians Ensure IFU states that when not in use, store active electrode in an electrically insulated container. Label according to IEC 60601-2-2 clause 201.7.9.2.2.101 a)5) | Labeling appropriately restricts the use of the device. Zip IFU, Zip Ace IFU warning meets IEC 60601-2-2 clause 201.7.9.2.2.101 a)5) | See footnote 2 UL Test Report ENG-RPT-059 | 10 | 5 | 50 | No | N/A | See footnote 3 and 4 |
| 13-a | Pencil Use | Insert pencil into ESU | Plug into ESU upside down in non-Megadyne unit. | No COAG activation | User dissatisfaction Delay in surgery | 3 | 5 | 15 | Yes | DFS | Application of industry standards to product design | Not practical | Ensure device meets requirements of IEC 60601-1 for safety of electrosurgical accessories | The Handpiece is tested and certified for compatibility with Megadyne ESU. See UL report ENG-RPT-058 | UL Test Report ENG-RPT-059 | 3 | 1 | 3 | No | N/A | See footnote 3 and 4 |
| 14-a | Pencil Use | Activation of pencil | Used in oxygen rich environment, with flammable anesthetics, materials, or gases, or by coming in contact with drape while activated. | Electrode starts fire in surgical field | Patient / User injury (burn) | 10 | 1 | 10 | Yes | IFS | See footnote 1 Label according to IEC 60601-2-2 clause 201.7.9.2.2.101 a)9) | Not practical | Ensure labeling restricts use to physicians Ensure IFU states that activated or hot electrodes should not be placed near or in contact with flammable materials or substances. Label according to IEC 60601-2-2 clause 201.7.9.2.2.101 a)9) | Labeling appropriately restricts the use of the device. Zip IFU, Zip Ace IFU warning meets IEC 60601-2-2 clause 201.7.9.2.2.101 a)9) | See footnote 2 UL Test Report ENG-RPT-059 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 14.5-a | Pencil Use | Activation of pencil | Used in oxygen rich environment, with flammable anesthetics, materials, or gases, or by coming in contact with drape while activated. | Electrode starts fire in surgical field with no property damage or patient injury | Delay of surgery, User dissatisfaction | 3 | 1 | 3 | Yes | IFS | See footnote 1 Label according to IEC 60601-2-2 clause 201.7.9.2.2.101 a)9) | Not practical | Ensure labeling restricts use to physicians Ensure IFU states that activated or hot electrodes should not be placed near or in contact with flammable materials or substances. Label according to IEC 60601-2-2 clause 201.7.9.2.2.101 a)9) | Labeling appropriately restricts the use of the device. Zip IFU, Zip Ace IFU warning meets IEC 60601-2-2 clause 201.7.9.2.2.101 a)9) | See footnote 2 UL Test Report ENG-RPT-059 | 3 | 1 | 3 | No | N/A | NA |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 11 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk considered? Yes/No | Options | Risk Control | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------|--|---|--|--------------------------|-----------------|---------------------------|-----|---|---------|---|--------------------|--|---|---|----------------------------------|---------------------------|-----|-------------------------------------|---|--|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 15-a | Pencil Use | User performs electrosurgery with disposable smoke pencil | Contamination from improper sanitizing and sterilization | Product used for more than 1 procedure | Patient injury | 10 | 5 | 50 | Yes | DFS | Product designed to ISO 15223 and IEC 60601-2-2 See footnote 1 | Not practical | NA | NA | UL Test Report ENG-RPT-059 | 10 | 5 | 50 | No | N/A | See footnote 3 and 4 |
| 15.5-a | Pencil Use | User performs electrosurgery with disposable smoke pencil | Contamination from improper sanitizing and sterilization | Product used for more than 1 procedure | Patient injury | 10 | 5 | 50 | Yes | IFS | NA | Not practical | Ensure pencil is designated as single use in IFU and on labeling. | Labeling appropriately restricts the use of the device. | UL Test Report ENG-RPT-059 | 10 | 5 | 50 | No | N/A | See footnote 3 and 4 |
| 16-a | Pencil Use | User performs surgery with small cross-sectional area with monopolar output. | Use error of monopolar energy | Unwanted tissue damage | Patient injury | 10 | 1 | 10 | Yes | DFS | See footnote 1 | Not practical | NA | NA | UL Test Report ENG-RPT-059 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 16.5-a | Pencil Use | User performs surgery with small cross-sectional area with monopolar output. | Use error of monopolar energy | Unwanted tissue damage | Patient injury | 10 | 1 | 10 | Yes | IFS | Label according to IEC 60601-2-2 clause 201.7.9.2.2.101 a)6) See footnote 1 | Not practical | Ensure labeling restricts use to physicians Ensure IFU states electrosurgery should not be used on body parts with relatively small cross-sectional area. Label according to IEC 60601-2-2 clause 201.7.9.2.2.101 a)6) | Labeling appropriately restricts the use of the device. Zip IFU, Zip Ace IFU warning meets IEC 60601-2-2 clause 201.7.9.2.2.101 a)6) | UL Test Report ENG-RPT-059 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 17-a | Pencil Use | Activation of electrosurgical pencil. | Failure to use the lowest power setting to achieve the desired effect | Excessive tissue damage | Patient injury | 10 | 1 | 10 | Yes | DFS | Megadyne ESU defaults to null power setting on start up, See footnote 1 | Not practical | NA | NA | UL Test Report ENG-RPT-059, See complaint analysis. | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 Occurrence based upon complaint analysis |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 12 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------|--|---|---------------------------------|----------------------------|-----------------|---------------------------|-----|-----------------|--------------|---|--------------------|--|--|---|----------------------------------|---------------------------|-----|-------------------------------------|---|---|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 17.5-a | Pencil Use | Activation of electrosurgical pencil, PRS 1701 – 1720, PRS 1901 – 1920 | Failure to use the lowest power setting to achieve the desired effect | Excessive tissue damage | Patient injury | 10 | 1 | 10 | Yes | IFS | Megadyne ESU defaults to null power setting on start up . Ensure IFU states that the lowest power setting to achieve the desired effect should be used. Label according to IEC 60601-2-2 See footnote 1 | Not practical | Ensure labeling restricts use to physicians Ensure IFU states that the lowest power setting to achieve the desired effect should be used. Label according to IEC 60601-2-2 clause 201.7.9.2.2.101 a)6) and ISO 15223. Validation Plan: Monitor customer complaints. | Labeling appropriately restricts the use of the device. Zip IFU, Zip Ace IFU warning meets IEC 60601-2-2 clause 201.7.9.2.2.101 a)7) | UL Test Report ENG-RPT-059, See complaint analysis. | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 Occurrence based upon complaint analysis |
| 18-a | Pencil Use | Activate electrode for tissue dissection. | User error. Cleaned PTFE with an abrasive cleaner. | Excessive eschar on tip. | Reduced performance of tip | 1 | 1 | 1 | NA | IFS | NA | NA | Ensure labeling restricts use to physicians Ensure IFU states not to use a scratch pad or abrasive material to clean the tip | Labeling appropriately restricts the use of the device and includes the warning | NA | 1 | 1 | 1 | NA | NA | NA |
| 19-a | Pencil Use | Activate electrode | Failure to inspect the tip for damage. | Damaged Electrode | Patient / User Injury | 10 | 1 | 10 | Yes | IFS | See footnote 1 Include caution statement in IFU | Not practical | Ensure labeling includes caution to inspect electrode Ensure labeling restricts use to physicians | Labeling appropriately restricts the use of the device and includes caution | See footnote 2 | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 20-a | Pencil Use | Activate bent electrode | Use error. Bent electrode past 60 degrees | Coating/electrode is damaged. | Patient / User injury | 10 | 1 | 10 | Yes | IFS | See footnote 1 Include caution statement in IFU | Not practical. | Ensure labeling restricts use to physicians Ensure IFU states not to bend electrode past 60 degrees. | Verified instructions are listed in the Zip IFU, Zip Ace IFU | See footnote 2 | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 21-a | Pencil Use | Activation away from intended target tissue | Use error. Activated electrode when not in contact with tissue | Unintended injury | Patient / User injury | 10 | 1 | 10 | Yes | IFS | See footnote 1 Include warning statement in IFU | Not practical. | Ensure labeling restricts use to physicians Ensure IFU states appropriate caution | Verified instructions are listed in the Zip IFU, Zip Ace IFU | See footnote 2 | 10 | 1 | 10 | No | N/A | See footnote 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 13 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------|--|---|--|---|-----------------|---------------------------|-----|----------------------------|--------------|---|--------------------|--|--|----------------------------|----------------------------------|---------------------------|-----|-------------------------------------|---|---|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 22-a | Pencil Use | Placement of pencil cable on surgical field | Use error. Failure to position cable correctly. Smoke pencil cables are coiled or fastened with metal clips | Capacitive coupling injury due to cable contact with the patient or other leads. | Patient / User injury | 10 | 1 | 10 | Yes | IFS | Design to IEC 60601-2-2 Include warning statement in IFU See footnote 1 | Not practical | Ensure labeling restricts use to physicians Ensure IFU cautions cable placement. | Verified instructions are listed in the Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 58-a | Pencil Use | Protect electrode from damage when pencil not in use | Use Error | Tip protector not removed | Delay of Surgery | 3 | 1 | 3 | Yes | IFS | Provide information to user See footnote #1 | Not Practical | Ensure IFU instructs the user to remove tip protector | Zip IFU, Zip Ace IFU includes instructions | NA | 3 | 1 | 3 | No | N/A | See footnote 4 Occurrence determined by complaint analysis |
| 59-a | Pencil Use | Intuitive Hand Grips | Pencil is non-intuitive | User unable to grip pencil as intended | User dissatisfaction | 1 | 3 | 3 | Yes | IFS | User Validation, Instructions in IFU. See footnote 1 | NA | ENG-PRT-453 Validation, Ensure Zip-Pen IFU provides visual guide | Zip IFU, Zip Ace IFU | ENG-RPT-557 | 1 | 3 | 3 | No | NA | NA |
| 60-a | Pencil Use | Evacuation of smoke | Nozzle becomes occluded | Poor smoke capture | Potential user exposure to carcinogens and infectious by-products | 5 | 1 | 5 | Yes | DFS | Angled nozzle design, wider than predicate UltraVac | NA | Demonstrate outcome of full occlusion in ENG-PRT-280 | Demonstrate outcome of full occlusion in ENG-RPT-403 | NA | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 Occurrence based upon complaint analysis |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 14 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk Evaluation Yes(No) | Risk Control | | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------|---|---|---|---|-----------------|------------------------------|-----|--|--------------|--|--|--|---|-----------------------------|----------|-------------------------------------|-----|---|---|---|--|
| Item# | Application | Intended Use or Characteris tic (Function) | Initiating Events / Circumstanc es / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measure s | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes | |
| 61-a | Pencil Use | Holster holds pencil | User places pencil with extension nozzle in holster | Product falls on the floor | Delay of procedure | 3 | 1 | 3 | Yes | DFS | NA | Not Practical | Test per protocol ENG- PRT-344 | Test report ENG-RPT-401 | Test report ENG-RPT- 401 | 3 | 1 | 3 | No | N/A | See footnote 3 and 4 occurrence based upon complaint analysis | |
| 62-a | Pencil Use | Activate electrode on targeted tissue | Pencil with tubing has too much drag | User Fatigue | Delay of procedure | 3 | 3 | 9 | Yes | DFS/ IFS | Swivels are part of the design to reduce drag as well as allowing for different grips, see IFU. | Not Practical | Zip-Pen IFU provides visual guide Test Protocol ENG-PRT- 290 | Zip IFU, Zip Ace IFU Test Report ENG-RPT-418 | Test Report ENG-RPT- 418 | 3 | 3 | 9 | NA | NA | occurrence based upon complaint analysis | |
| 68-a | Pencil Use | Proper use of device per Intended Use | Use of device on main vessels in central circulatory system | Unintended use of device | Bleeding / Hemorrhage | 10 | 1 | 10 | Yes | IFS | See footnote #1 | Specify intended user of device | Ensure IFU contains warning about device not being evaluated in the central circulatory system. | IFU contains required statement. Labeling has symbol to consult IFU. CER (RA-RPT-007) | See Footnote #2 | 10 | 1 | 10 | No | NA | See Footnote #4 | |
| 69-a | Pencil Use | Proper use of device per Intended Use | Use of device on the central nervous system | Unintended use of device | Nerve injury- permanent and significant sensory or motor loss that will affect acts of daily living | 10 | 1 | 10 | Yes | IFS | See footnote #1 | Specify intended user of device | Ensure IFU contains warning about device not being evaluated in the central nervous system | IFU contains required statement. Labeling has symbol to consult IFU. CER (RA-RPT-007) | See Footnote #2 | 10 | 1 | 10 | No | NA | See Footnote #4 | |
| 70-a | Pencil Use | Proper use of device per Intended Use | User activates energy device near patient's implant (Inadequate filtering/ shielding) | Interference with patient implanted devices (e.g.: pacemaker) during activation | Impedes implant from functioning correctly | 10 | 1 | 10 | Yes | IFS | See footnote #1 | Specify intended user of device | Ensure IFU contains warning about device being used with implants. | IFU contains required statement. Labeling has symbol to consult IFU. CER (RA-RPT-007) | See Footnote #2 | 10 | 1 | 10 | No | NA | See Footnote #4 | |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 15 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|--------------------|--|---|--|--|-----------------|---------------------------|-----|-----------------|--------------|---|--------------------|--|---|--|----------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 23-a | Button | Pencil activation PRS 1105 | Button too small | Sore finger | User dissatisfaction | 3 | 1 | 3 | NA | DFS | Design for adequate button size | NA | Verification protocol ENG-PRT-290 Validation Plan: Protocol ENG-PRT-453 | Test report ENG-RPT-418 | Report ENG-RPT-557 | 3 | 1 | 3 | NA | NA | See footnote 3 and 4 |
| 24-a | Pencil Body | Handle of electrosurgical device PRS 1103 PRS 1104 | Branding not on product | User not aware of product brand | None | 1 | 1 | 1 | NA | DFS | Design with branding on the product | NA | Verification protocol ENG-PRT-290 | Test report ENG-RPT-418 | NA | 1 | 1 | 1 | NA | NA | NA |
| 25-a | Cable | Connects the ESU Generator to the Zip pen | Cable in close proximity to other equipment or other equipment cables | Causes electromagnetic interference (EMI) or radio frequency interference (RFI) in other equipment | Disrupt other equipment -customer dissatisfaction | 3 | 3 | 9 | Yes | DFS | See footnote 1 | Not Practical | Ensure device meets requirements of IEC 60601-1 | Similar cable tested for EMC, see report ENG-RPT-052 | Equivalent product tested to IEC 60601-1 Clause 17, see report ENG-RPT-052 | 3 | 2 | 6 | No | N/A | See footnote 4 |
| 25.5-a | Cable | Connects the Zip pen to the ESU Generator | Cable in close proximity to other equipment or other equipment cables | Causes electromagnetic interference (EMI) or radio frequency interference (RFI) in other equipment | Disrupt other equipment - Customer Dissatisfaction | 3 | 3 | 9 | NA | IFS | Label product as required by standards | Not Practical | Ensure IFU instructs the user on properly positioning cables | IFU instructs user on proper routing of cables. Labelling appropriately restricts the use of the device. Zip IFU | UL EMC Test report ENG-RPT-052 | 3 | 3 | 9 | NA | N/A | See footnote 4 |
| 26-a | Handpiece Assembly | Allows electrode to be fully inserted Insulates electrode shaft | Practitioner doesn't fully insert electrode into handpiece. | Electrode not fully inserted into handpiece (exposed uninsulated electrode shaft) | Patient or user burns | 10 | 1 | 10 | Yes | IFS | Product designed to IEC 60601-2-2 See footnote 1 | Not Practical | Ensure device meets requirements of IEC 60601-1. Ensure IFU instructs the user to fully insert the electrode into the handpiece | Labeling appropriately restricts the use of the device according to standard – IEC 60601-2-2. IFU instructs the user on how to insert the electrode and states to fully insert electrode into the handpiece. | UL Report ENG-RPT-059 | 10 | 1 | 10 | No | N/A | See footnote 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 16 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk Evaluation Yes/No | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------------------------|---|--|--------------------------------------|--------------------------|-----------------|---------------------------|-----|---|--------------|--|--------------------|---|---|---|-------------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 27-a | Multivac Packaging for Pencil | To maintain sterile barrier protect ES pencil during shipping and storage condition | Ineffective packaging for this application | Product damaged | Patient injury. | 10 | 1 | 10 | Yes | IFS | Product designed to ISO 15223 and IEC 60601-2-2 See footnote 1 | Not Practical | Ensure IFU includes instructions to discard if pouch is damaged. | Labeling appropriately restricts the use of the device according to standards ISO 15223 | Labeling appropriately restricts the use of the device according to standards ISO 15223 | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 27.5-a | Multivac Packaging for Pencil | To maintain sterile barrier protect ES pencil during shipping and storage condition | Ineffective packaging for this application | Product damaged | User dissatisfaction | 1 | 1 | 1 | Yes | IFS | Product designed to ISO 15223 and IEC 60601-2-2 See footnote 1 | Not Practical | Ensure IFU includes instructions to discard if pouch is damaged. | Labeling appropriately restricts the use of the device according to standards ISO 15223 | Labeling appropriately restricts the use of the device according to standards ISO 15223 | 1 | 1 | 1 | No | N/A | NA |
| 28-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Using on patient with electronic implant | Interaction with electronic implants | Patient injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 29-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Use device in presence of flammable anesthetics or oxidizing gases or in close proximity to volatile solvents. | Explosion / Fire | User / Patient injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 30-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Use with hybrid trocar | Capacitive coupling | Patient / User injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 31-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Failure to inspect product prior to use | Damaged insulation or cables | Patient / User injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis | 10 | 1 | 10 | No | N/A | See footnote 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 17 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Yes/No | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------|---|---|---------------------------------|--|-----------------|---------------------------|-----|---------------------------|--------------|--|--------------------|--|---|--|----------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 32-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Use device without aspirating fluids from the area. | Activating in conductive fluids | Patient injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 33-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Activate device when not on target tissue resulting in capacitive coupling. | Alternate site injury | Patient injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 34-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Lack of use of protective equipment / lack of smoke evacuation | Smoke inhalation | Possible user exposure to potential carcinogens and infectious by-products | 5 | 1 | 5 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis | 5 | 1 | 5 | No | N/A | See footnote 4 |
| 35-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Attach device to ESU when unit is on. | Burn / Fire / electrical shock | User / Patient injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 36-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | User does not use ESU with CQM and a compatible monitoring neutral electrode, No CQM alarm. | Neutral electrode site burn | Patient Injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | See complaint analysis | 10 | 1 | 10 | No | N/A | See footnote 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 18 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|--------------------------------------|---|---|---------------------------------|--------------------------|-----------------|---------------------------|-----|-----------------|--------------|--|--------------------|--|---|--|----------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 36.2-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | User does not use ESU with CQM and a compatible monitoring neutral electrode, No CQM alarm. | Reduced electrosurgical effect | User Dissatisfaction | 3 | 1 | 3 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | See complaint analysis | 3 | 1 | 3 | No | N/A | See footnote 4 |
| 36.5-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Insufficient contact of neutral electrode | Neutral electrode site burn | Patient Injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | CQM Alarm | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 36.7-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Insufficient contact of neutral electrode | Reduced electrosurgical effect | User Dissatisfaction | 3 | 1 | 3 | Yes | IFS | Label product as required by standards | CQM Alarm | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis | 3 | 1 | 3 | No | N/A | See footnote 4 |
| 37-a | REDUNDANT CONTENT – SEE 17-A, 17.5-A | | | | | | | | | | | | | | | | | | | | |
| 38-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Allow eschar or other material to build up on active electrode. | Reduced electrosurgical effect | Patient injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 39-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Activating device away from target tissue (such as during cleaning) | Accidental activation | Patient / User injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis | 10 | 1 | 10 | No | N/A | See footnote 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 19 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Yes/No | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|--------------------------------------|---|---|--|---|-----------------|---------------------------|-----|---------------------------|--------------|--|--------------------|--|--|---|----------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 40-a | Zip-Pen IFU | Provide information for the user PRS 1701 – 1720, PRS 1901 – 1920 | Not following recommended ESU duty cycle and/or RF voltage rating | Hot electrode | Patient / User injury | 10 | 5 | 50 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059 | 10 | 5 | 50 | No | N/A | See footnote 4 |
| 63-a | Zip-Pen IFU | Provide information for user | Use Error due to user not reviewing manuals of equipment necessary for use of Zip Pen | Devices used incorrectly | Patient / User Injury | 10 | 1 | 10 | Yes | IFS | See footnote 1 Include caution statement in IFU | Not practical | Ensure labeling includes caution to review ESU, smoke evacuator, neutral, and active electrode use instructions. | Zip IFU, Zip Ace IFU, Labeling appropriately restricts the use of the device and includes caution 2 | NA | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 41-a | Pencil and Nozzle Extension Labeling | Provide information for the user PRS 1701 – 1720, PRS 1901-1920 | Product not labeled per Standards requirement | Use wrong product Use expired product Re-use the product | Patient Injury | 10 | 5 | 50 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 | See label drawings | Drawings meet Standards - ISO 15223 and IEC 60601-2-2 | 10 | 5 | 50 | No | N/A | See footnote 4 |
| 42-a | Pencil and Nozzle Extension IFUs | Use Instructions PRS 1001 URS 1001 | Intended use not stated in IFU | Improper use of device | Patient Injury | 10 | 1 | 10 | Yes | IFS | Provide proper Instructions, Warnings, and Cautions in IFU according to ISO 15223 and IEC 60601-2-2. | Not Practical | Ensure IFU is compliant with ISO 15223 and IEC 60601-2-2 | Labeling appropriately restricts the use of the device. See IFUs: Zip IFU, Zip Ace IFU (pencil) and nozzle extension | Drawings meet Standards - ISO 15223 and IEC 60601-2-2 | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 43-a | Nozzle Extension (4.0" and 6.5") | Install Nozzle and electrode | Use error, not following IFU | Incorrect Electrode to Nozzle size (Electrode too short) | Cannot do surgery without removing nozzle, delay of surgery | 3 | 1 | 3 | NA | IFS | Put compatibility statement in IFU | NA | Ensure compatibility is in IFU | Zip IFU, Zip Ace IFU | IFU 3000313-01 | 3 | 1 | 3 | NA | NA | NA |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 20 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|----------------------------------|--|---|---|--|-----------------|---------------------------|-----|-----------------|--------------|---|--------------------|--|--|---|----------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 44-a | Nozzle Extension (4.0" and 6.5") | Install Nozzle and electrode | Use error, not following IFU | Incorrect Electrode to Nozzle size (Electrode too long), Poor smoke capture | Possible user exposure to potential carcinogens and infectious by-products | 5 | 1 | 5 | NA | IFS | Put compatibility statement in IFU | NA | Ensure compatibility is in IFU | Zip IFU, Zip Ace IFU | IFU 3000313-01 | 5 | 1 | 5 | NA | NA | See footnote 4 |
| 45-a | Nozzle Extension (4.0" and 6.5") | Install Nozzle | Nozzle not correctly installed or fully installed | Nozzle loose or falls off | Nozzle falls into patient | 5 | 1 | 5 | Yes | DFS | Visual and tactile confirmation of correct placement | Not practical | NA | Usability Test Report ENG-RPT-418 identified that nozzle was acceptable and was easy to attach | Usability Test Report ENG-RPT-418 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 45.5-a | Nozzle Extension (4.0" and 6.5") | Install Nozzle | Nozzle not correctly installed or fully installed | Nozzle loose or falls off | Nozzle falls into patient | 5 | 1 | 5 | Yes | IFS | NA | Not practical | Ensure installation instruction is in IFU Usability Test Protocol ENG-PRT-290 | Zip IFU, Zip Ace IFU | Usability Test Report ENG-RPT-418 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 46-a | Nozzle Extension (4.0" and 6.5") | Install Nozzle | Nozzle not correctly installed or fully installed | Installed upside down or incorrectly | Nozzle falls into patient | 5 | 1 | 5 | Yes | DFS | Visual and tactile confirmation of correct placement | Not practical | NA | Usability Test Report ENG-RPT-418 identified that nozzle was acceptable and was easy to attach | Usability Test Report ENG-RPT-418 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 46.5-a | Nozzle Extension (4.0" and 6.5") | Install Nozzle | Nozzle not correctly installed or fully installed | Installed upside down or incorrectly | Nozzle falls into patient | 5 | 1 | 5 | Yes | IFS | NA | Not practical | Ensure installation instruction is in IFU Usability Test Protocol ENG-PRT-290 | Zip IFU, Zip Ace IFU | Usability Test Report ENG-RPT-418 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 47-a | Nozzle Extension (4.0" and 6.5") | Install Nozzle and electrode | Use error, not following IFU | Install Nozzle first and cannot get electrode installed properly | Delay of surgery | 1 | 1 | 1 | NA | IFS | IFU gives instructions to install electrode prior to extension Nozzle | NA | Ensure installation instruction is in IFU Usability Test Protocol ENG-PRT-290 | Usability Test Report ENG-RPT-418 identified that nozzle was acceptable and was easy to attach | Usability Test Report ENG-RPT-418 | 1 | 1 | 1 | NA | NA | NA |
| 48-a | Use pencil with extension nozzle | Perform electrosurgery with extension nozzle | Use error | Nozzle used for more than 1 procedure | Patient infection | 10 | 5 | 50 | Yes | IFS | Product designed to ISO 15223 and IEC 60601-2-2 See footnote 1 | Not practical | Verify labeling designates nozzle as single use | Labeling appropriately restricts the use of the device according to standards ISO 15223 | Labeling appropriately restricts the use of the device according to standards ISO 15223 | 10 | 5 | 50 | No | N/A | See footnote 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 21 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Yes/No | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|------------------------------------|--|--|--|--|-----------------|---------------------------|-----|---------------------------|--------------|---|---|---|--|----------------|----------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 49-a | Use pencil with extension nozzle | Attach nozzle | Use error Nozzle attached with foreign material | Decreases visibility Product damage Reduced smoke evacuation | Possible user exposure to potential carcinogens and infectious by-products | 5 | 1 | 5 | NA | IFS | NA | NA | Verify labeling restricts the use to physicians and includes caution to not use foreign material to secure nozzle | Labeling cautions not to secure nozzle with foreign material | NA | 5 | 1 | 5 | NA | NA | See footnote 4 |
| 50-a | Use ULPA filter | Use filter to filter smoke PRS 1201 PRS 1202 | Use error. Improper use of filter | Filter clogged (with fluid and/or tissue) Filter operation affected by environmental conditions | User dissatisfaction | 1 | 1 | 1 | NA | IFS | NA | NA | Ensure IFU cautions user not to use filter as primary fluid / tissue collector Ensure IFU includes shipping, storage, and use conditions | IFU | IFU 3000190-01 | 1 | 1 | 1 | NA | NA | NA |
| 51-a | Use ULPA filter | Use filter to filter smoke | Use error. Did not put cap on. | Biohazard materials spill out of ULPA filter during disconnection or disposal. | Potential for biohazard exposure | 10 | 1 | 10 | Yes | DFS | Connector designed above fluid level. Cap provided for used filter. | Clear material design provides visual indication of fluid presence. | Verify design provides connector above fluid level, the material is clear and the cap is provided | ENG-RPT-418 Usability Report | NA | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 51.5-a | Use ULPA filter | Use filter to filter smoke | Use error. Did not put cap on. | Biohazard materials spill out of ULPA filter during disconnection or disposal. | Potential for biohazard exposure | 10 | 1 | 10 | Yes | IFS | NA | NA | Ensure IFU cautions user to install cap prior to disposal. | IFU | IFU | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 52-a | Use filters (ULPA/Charcoal filter) | Use to filter smoke PRS 1114 | Use of damaged filters | Filter does not attach. Does not work properly. | User dissatisfaction | 1 | 1 | 1 | NA | IFS | Include instructions in IFU to not use damaged filter | NA | Ensure IFU cautions against use of damaged filters | See IFU | NA | 1 | 1 | 1 | NA | NA | NA |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 22 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|------------------------------------|---|--|---|--|-----------------|---------------------------|-----|-----------------|--------------|--|--------------------|---|---|---|----------------------------------|---------------------------|-----|-------------------------------------|---|-------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 53-a | Use filters (ULPA/Charcoal filter) | Use to filter smoke | Cleaning or otherwise sterilizing filters | Damaged filters | User dissatisfaction | 1 | 1 | 1 | NA | IFS | Include instructions in IFU to not clean or otherwise sterilize filters | NA | Ensure IFU cautions against cleaning or otherwise sterilizing filters | See IFU | NA | 1 | 1 | 1 | NA | NA | NA |
| 54-a | Charcoal Filter | Filter odors and residual chemicals PRS 1318 | Overuse of filter | Charcoal absorption depleted | Odor/chemical removal inadequate, User dissatisfaction | 1 | 1 | 1 | NA | IFS | Source from approved supplier Identify requirements for replacement in IFU | NA | Ensure IFU restricts to use by physician and includes instruction on conditions that require filter change. | Required information contained in the IFU | NA | 1 | 1 | 1 | NA | NA | NA |
| 55-a | ULPA and Charcoal Filter Labeling | Provide information for the user PRS 1721 – 1733, PRS 1921 – 1935 | Product not labeled per Standards requirement | Use wrong product Re-use the product | User injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan: Monitor customer complaints | See label drawings | UL Test Report ENG-RPT-059, See complaint analysis. | 10 | 1 | 10 | No | N/A | See footnote 4 |
| 56-a | RF Sensor | Detect Pencil Activation | Defective / damaged sensor Not properly installed/ connected to unit or pencil | Not sensing activation | Customer dissatisfaction | 1 | 1 | 1 | NA | DFS | Design, material selection, and in process inspection See footnote 1 | NA | Verify certificate of conformance is required on receipt | Drawing requires Certificate of Conformance | NA | 1 | 1 | 1 | NA | NA | NA |
| 56.5-a | RF Sensor | Detect Pencil Activation | Defective / damaged sensor Not properly installed/ connected to unit or pencil | Not sensing activation | Customer dissatisfaction | 1 | 1 | 1 | NA | IFS | Provide information to user See footnote 1 | NA | Verify IFU includes installation instruction | IFU ICM-001-0866 includes instruction | NA | 1 | 1 | 1 | NA | NA | NA |
| 57-a | Handpiece Assembly | Provide information for user | Supplier defect Use Error | Cut/Coag buttons fall out, Fractures or pieces breakoff, etc. | Patient/User Injury | 10 | 1 | 10 | Yes | DFS | Design, material selection, and in process inspection See footnote 1 | Not Practical | Ensure device meets IEC 60601-2-2 Verify certificate of conformance is required on receipt | Drawing requires Certificate of Conformance | UL Test Report ENG-RPT-059 See footnote 2 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 23 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|--------------------|---|---|---|--------------------------|-----------------|---------------------------|-----|-----------------|--------------|--|--------------------|---|---|--|----------------------------------|---------------------------|-----|-------------------------------------|---|---|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 57.2-a | Handpiece Assembly | Provide information for user | Supplier defect Use Error | Device damage such as Cracked or damaged components | Customer dissatisfaction | 1 | 1 | 1 | Yes | DFS | Design, material selection, and in process inspection See footnote 1 | Not Practical | Ensure device meets IEC 60601-2-2 Verify certificate of conformance is required on receipt | Drawing requires Certificate of Conformance | UL Test Report ENG-RPT-059 See footnote 2 | 1 | 1 | 1 | No | N/A | See footnote 3 and 4 |
| 57.5-a | Handpiece Assembly | Provide information for user | Supplier defect Use Error | Cut/Coag buttons fall out, Fractures or pieces breakoff, etc. | Patient / User Injury | 10 | 1 | 10 | Yes | IFS | Provide information to user See footnote 1 | Not Practical | Verify IFU includes proper use instruction including to inspect for damage | IFU & ICM-001-0866 includes instructions | UL Test Report ENG-RPT-059 See footnote 2 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 57.7-a | Handpiece Assembly | Provide information for user | Supplier defect Use Error | Device damage such as Cracked or damaged components | Customer dissatisfaction | 1 | 1 | 1 | Yes | IFS | Provide information to user See footnote 1 | Not Practical | Verify IFU includes proper use instruction including to inspect for damage | IFU & ICM-001-0866 includes instructions | UL Test Report ENG-RPT-059 See footnote 2 | 1 | 1 | 1 | No | N/A | See footnote 3 and 4 |
| 64-a | Electrode | Electrode removal | Unable to grasp electrode to remove electrode for the purpose of swapping in a different electrode type | Electrode stuck in Zip Pen | Delay of surgery | 3 | 5 | 15 | Yes | DFS | Use of common tools found in surgical suite | NA | ENG-PRT-452 | NA | ENG-RPT-559 | 3 | 1 | 3 | No | NA | See footnote 3 and 4, Occurrence based on complaints analysis |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 24 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|------------------------|---|--|---|------------------------------|-----------------|---------------------------|-----|-----------------|--------------|---|---------------------------------|--|---|---|----------------------------------|---------------------------|-----|-------------------------------------|---|--|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? | Additional Risk Control (if applicable) | Risk Benefit Analysis / Notes |
| 65-a | Electrode | Electrode removal | Re-inserting or adjusting electrode after removal with a tool | Significant coating / insulation damage, resulting in significant eschar build up and reduced surgical effect | Patient Injury | 10 | 1 | 10 | Yes | IFS | Label product as required by standards | Not Practical | Ensure label drawings meet Standards - ISO 15223 and IEC 60601-2-2 Validation Plan; Monitor customer complaints | See label drawings Zip IFU, Zip Ace IFU | UL Test Report ENG-RPT-059, See complaint analysis. | 10 | 1 | 10 | NA | NA | See footnote 4 Occurrence based on complaints analysis |
| 66-a | Electrode | Electrode removal | Re-inserting or adjusting electrode after removal with tool | Electrode bent past 60 degrees | Patient / User injury | 10 | 1 | 10 | Yes | IFS | See footnote 1 Include caution statement in IFU | Not practical. | Zip IFU 3000312-01, Zip Ace IFU 3000317-01 | Zip IFU, Zip Ace IFU | NA | 10 | 1 | 10 | No | N/A | See footnote 4 Occurrence based on complaints analysis |
| 67-a | Electrode | Electrode removal | Removing electrode with metal instrument | Activating while removing with metal instrument | User injury | 10 | 1 | 10 | Yes | IFS | Caution statement in IFU | Not practical. | Zip IFU 3000312-01, Zip Ace IFU 3000317-01 Caution 9 | Zip IFU, Zip Ace IFU | NA | 10 | 1 | 10 | No | N/A | See footnote 4 Occurrence based on complaints analysis |
| 71-a | Perform Electrosurgery | Proper use of device per Intended Use | Use of Electrosurgery - Do not activate in contact with Metal instrument or device | Unintended monopolar energy is applied to critical structure | Patient / User injury (burn) | 10 | 1 | 10 | Yes | IFS | See footnote 1 | Specify intended user of device | Ensure labeling restricts use to physicians familiar with electrosurgery Verify potential distraction for user relative to printed image on electrode insulation & Ensure IFU specifies electrode to be in contact or in close proximity to target tissue prior to activation | IFU contains required statement. Labeling has symbol to consult IFU. CER (SCN070741) | See Footnote 2 | 10 | 1 | 10 | No | NA | See footnote 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 25 of 52 |

Zip-Pen DESIGN FMEA

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|---------------|---|--|--|---|-----------------|------------------------------|-----|--------------------------------|--------------|---|-----------------------|--|--|--|--|------------------------------|-----|---|---|---------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 1-d | Circuit Board | Cut mode activates when distal switch depressed PRS 1101 PRS 1402 PRS 1403 PRS 1409 PRS 2001 | Incorrect tracings on circuit board. Incorrect cable lead position on circuit board | Coag mode activated when distal (cut) switch depressed (will not cut properly) | Practitioner dissatisfaction, wrong effect, patient tissue damage | 5 | 5 | 25 | Yes | DFS | Material and supplier selection and product validation | Not Practical | Check cut and coag continuity, prior to and after cable flex test and fluid ingress test during validation | Mechanical Test Report ENG-RPT- 329 Fluid Ingress Test Report ENG-RPT- 555 Fluid Ingress Test Report 500433702 | Tested to IEC 60601-2-2 clause 201.11.6.5 b | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 2-d | Circuit Board | Coag mode activates when proximal switch depressed PRS 1101 PRS 1402 PRS 1403 PRS 1409 PRS 2001 | Incorrect tracings on circuit board. Incorrect cable lead position on circuit board. | Cut mode activated when proximal (coag) switch depressed (will not coag properly) | Practitioner dissatisfaction, wrong effect | 5 | 5 | 25 | Yes | DFS | Material and supplier selection and product validation | Not Practical | Check cut and coag continuity, prior to and after cable flex test and fluid ingress test during validation | Mechanical Test Report ENG-RPT- 329 Fluid Ingress Test Report ENG-RPT- 555 Fluid Ingress Test Report 500433702 | Tested to IEC 60601-2-2 clause 201.11.6.5 b | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 3-d | Circuit Board | Coag mode activates when proximal switch depressed PRS 1101 PRS 1402 PRS 1403 PRS 1409 PRS 2001 | Cable leads not connected to circuit board properly | Both “cut” and “coag” modes activate simultaneousl y or unintentional/ self-activation | Practitioner dissatisfaction, wrong effect | 5 | 5 | 25 | Yes | DFS | Material and supplier selection and product validation | Not Practical | Check cut and coag continuity, prior to and after cable flex test and fluid ingress test during validation | Mechanical Test Report ENG-RPT- 329 Fluid Ingress Test Report ENG-RPT- 555 Fluid Ingress Test Report 500433702 | Tested to IEC 60601-2-2 clause 201.11.6.5 b | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 26 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|---------------|---|--|---|---|-----------------|---------------------------|-----|---|--------------|--|--------------------|--|--|--|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 3.5-d | Circuit Board | Coag mode activates when proximal switch depressed PRS 1101 PRS 1402 PRS 1403 PRS 1409 PRS 2001 | Cable leads not connected to circuit board properly | Self-activation | Patient injury | 10 | 5 | 50 | Yes | DFS | Material and supplier selection and product validation | Not Practical | Check cut and coag continuity, prior to and after cable flex test and fluid ingress test during validation | Mechanical Test Report ENG-RPT-329 Fluid Ingress Test Report ENG-RPT-555 Fluid Ingress Test Report 500433702 | Tested to IEC 60601-2-2 clause 201.11.6.5 b | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 4-d | Circuit Board | Deliver current from ESU to electrode. PRS 1101 PRS 1402 PRS 1403 PRS 1409 PRS 2001 | Incorrect tracings on circuit board. | No ES current at electrode (inoperable unit) | Practitioner dissatisfaction (inconvenience), Delay of surgery | 3 | 1 | 3 | NA | DFS | Material and supplier selection and product validation | NA | Check cut and coag continuity, prior to and after cable flex test and fluid ingress test during validation | Mechanical Test Report ENG-RPT-329 Fluid Ingress Test Report ENG-RPT-555 Fluid Ingress Test Report 500433702 | Tested to IEC 60601-2-2 clause 201.11.6.5 b | 3 | 1 | 3 | NA | NA | NA |
| 5-d | Dome Switches | Activate ESU to provide signal to electrode Provide tactile response to user PRS 1101 PRS 1301 PRS 2001 | Incorrect geometry (or material type) of dome switch | Switches too hard to depress | User dissatisfaction | 1 | 1 | 1 | NA | DFS | Material and supplier selection and product validation | NA | Verify button force is validated, include button force test in Mechanical Protocol | Mechanical Test Report ENG-RPT-329 Dome Switch Force Test Report 500433702 | Mechanical Test Report ENG-RPT-329 Use Test Report ENG-RPT-779 | 1 | 1 | 1 | NA | NA | NA |
| 6-d | Dome Switches | Activate ESU to provide signal to electrode Provide tactile response to user PRS 1101 PRS 1301 PRS 2001 | Incorrect geometry (or material type) of dome switch | Switches too soft (no tactile response, but operable) | User dissatisfaction | 1 | 1 | 1 | NA | DFS | Material and supplier selection and product validation | NA | Verify button force is validated, include button force test in Mechanical Protocol | Mechanical Test Report ENG-RPT-329 Dome Switch Force Test Report 500433702 | Mechanical Test Report ENG-RPT-329 Same design as disposable pencils, see footnote 2 Use Test Report ENG-RPT-779 | 1 | 1 | 1 | NA | NA | NA |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 27 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------------------------|---|--|--|--|-----------------|---------------------------|-----|--|--------------|--|------------------------------|--|---|--|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 7-d | Dome Switches | Activate ESU to provide signal to electrode Provide tactile response to user PRS 1101 PRS 1301 PRS 2001 | Incorrect geometry (or material type) of dome switch | Switches too soft (always activated) | Potential user or patient injury | 10 | 5 | 50 | Yes | DFS | Material and supplier selection and product validation | Not Practical | Verify button force is validated, include button force test in Mechanical Protocol | Mechanical Test Report ENG-RPT-329 Dome Switch Force Test Report 500433702 | Mechanical Test Report ENG-RPT-329 Same design as disposable pencils, see footnote 2 Use Test Report ENG-RPT-779 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 8-d | Dome Switches | Activate ESU to provide signal to electrode Provide tactile response to user PRS 1101 PRS 1301 PRS 2001 | Incorrect geometry (or material type) of dome switch | Switches fail during surgery | Patient Injury | 10 | 5 | 50 | Yes | DFS | Material and supplier selection and product validation | Not Practical | Verify button force is validated, include button force test in Mechanical Protocol | Mechanical Test Report ENG-RPT-329 Dome Switch Force Test Report 500433702 | Mechanical Test Report ENG-RPT-329 Same design as disposable pencils, see footnote 2 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 9-d | Electrical Connection/ Collet | Provide secure mechanical and electrical connection for electrode PRS 1101 PRS 1302 PRS 1306 PRS 2001 URS 1002 | Collet diameter out of tolerance | Electrode falls out (possibly during procedure) | Foreign object in patient (potential injury) | 10 | 5 | 50 | Yes | DFS | Material and supplier selection and product validation | Detect electrode under X-Ray | Validate collet retention force per standards. | Mechanical Test Report ENG-RPT-329 | Tested to IEC 60601-2-2 clause 201.15.4.1.012 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 10-d | Electrical Connection/ Collet | Provide secure mechanical and electrical connection for electrode PRS 1101 PRS 1107 PRS 1302 PRS 1303 PRS 1306 PRS 2001 URS 1002 | Loose Collet, electrode not installed properly | Electrode falls out (possibly during procedure), electrode falls on floor (becomes non-sterile), no electrical contact | Delay of surgery | 1 | 1 | 1 | NA | DFS | Design, material and supplier selection and product validation | NA | Validate collet retention force per standards | Mechanical Test Report ENG-RPT-329 | Tested to IEC 60601-2-2 clause 201.15.4.1.012 | 1 | 1 | 1 | NA | NA | NA |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 28 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------------------------|---|---|----------------------------------|----------------------------------|-----------------|---------------------------|-----|-----------------|--------------|---|--------------------|--|---|--|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 11-d | Electrical Connection/ Collet | Provide secure mechanical and electrical connection for electrode PRS 1101 PRS 1302 PRS 1303 PRS 1306 PRS 1319 PRS 1409 PRS 2001 | Insufficient contact area between electrode and electrical connector. | Pencil overheats and melts | Potential user or patient injury | 10 | 5 | 50 | Yes | DFS | Material and supplier selection and design, product validation | Not Practical | Check cut and coag continuity, prior to and after cable flex test and fluid ingress test during validation | Mechanical Test Report ENG-RPT-329 Fluid Ingress Test Report ENG-RPT-555 Thermal Test Report ENG-RPT-395 Fluid Ingress Test Report 500433702 | Mechanical Test Report ENG-RPT-329 Fluid Ingress Test Report ENG-RPT-555 Tested to IEC 60601-2-2 clause 201.11.6.5 b and IEC 60601-1 Clause 11.1.2.2 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 12-d | Electrical Connection/ Collet | Provide secure mechanical and electrical connection for electrode PRS 1101 PRS 1302 PRS 1303 PRS 1306 PRS 1409 PRS 2001 | Detached or damaged collet | No continuity | Practitioner dissatisfaction | 1 | 1 | 1 | NA | DFS | Material and supplier selection and design, product validation | NA | Check cut and coag continuity, prior to and after cable flex test and fluid ingress test during validation | Mechanical Test Report ENG-RPT-329 Fluid Ingress Test Report ENG-RPT-555 Fluid Ingress Test Report 500433702 | Mechanical Test Report ENG-RPT-329 Fluid Ingress Test Report ENG-RPT-555 Tested to IEC 60601-2-2 clause 201.11.6.5 b | 1 | 1 | 1 | NA | NA | NA |
| 13-d | Buttons / Rocker | User interface for ESU activation on the pencil PRS 1106 | Incorrect production or assembly | Both switches are the same color | User dissatisfaction | 1 | 1 | 1 | NA | DFS | Part designed for two shot color. Designed for oriented assembly. With yellow button nearest the active electrode | NA | Design does not allow buttons to be assembled incorrectly | Drawing 5800097-01, 6020190-01 and 6020191-01 | NA | 1 | 1 | 1 | NA | NA | NA |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 29 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|---------------------|--|---|---|---|-----------------|------------------------------|-----|---|--------------|---|-----------------------|---|--|---|--|------------------------------|-----|---|---|---------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 14-d | Buttons / Rocker | User interface for ESU activation on the pencil PRS 1101 PRS 2001 | Switch components binding | ES pencil remains activated after switch is released. | Potential user or patient injury | 10 | 5 | 50 | Yes | DFS | Design for proper clearances Material and supplier selection and product validation | Not Practical | Verify button force is validated | Mechanical test report ENG-RPT- 329 | Mechanical test report ENG-RPT-329 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 15-d | Pencil Nozzle | Provides port for smoke capture PRS 1101 PRS 1108 PRS 1304 PRS 1501 PRS 1502 PRS 1503 PRS 2001 URS 1007 | Incorrect design or material Not a gamma stable material Opaque nozzle Inadequate design for flow | Fractures leaving material in patient Falls off into patient | Foreign body reaction or additional surgical intervention, patient injury | 10 | 5 | 50 | Yes | DFS | Design, material selection, and product validation | Not Practical | Test for mechanical strength Test for flow volume | Mechanical test report ENG-RPT- 329 | Tested per ISO 10993-1:2009 clauses 6.2.2.2, 6.2.2.3, & 6.2.2.4 See test report ENG-RPT-327 See test report ENG-RPT-329 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 15.5-d | Pencil Nozzle | Provides port for smoke capture PRS 1101 PRS 1108 PRS 1304 PRS 1501 PRS 1502 PRS 1503 PRS 2001 URS 1007 | Incorrect design or material Not a gamma stable material Opaque nozzle Inadequate design for flow | Transparency obstructs surgeons view Air flow not adequate, poor smoke evacuation | Possible user exposure to potential carcinogens and infectious by-products | 5 | 1 | 5 | Yes | DFS | Design, material selection, and product validation | Not Practical | Test for mechanical strength Test for flow volume | Mechanical test report ENG-RPT- 329 | Tested per ISO 10993-1:2009 clauses 6.2.2.2, 6.2.2.3, & 6.2.2.4 See test report ENG-RPT-327 See test report ENG-RPT-329 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 16-d | Housing / Handle | “Splash proof” – protect internal components from splashed liquids PRS 1101 PRS 1409 PRS 2001 | Fluid enters housing, causing short circuits to internal components | Self- activation (without switches pressed) | Potential user or patient injury | 10 | 5 | 50 | Yes | DFS | Design specifies material and assembly to prevent fluid ingress | Not Practical | Ensure that material requirements and assembly are defined. The test for fluid ingress will be performed as part of product validation. ENG- PRT-439 | Material and assembly requirements are defined on drawing. Fluid Ingress Test Report ENG-RPT- 555 Fluid Ingress Test Report 500433702 | Tested to IEC 60601-2-2 clause 201.11.6.5 b | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 30 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|------------------|---|---|---------------------------------|----------------------------------|-----------------|---------------------------|-----|--|--------------|--|--------------------|---|---|---|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 63-d | Housing / Handle | “Splash proof” – protect internal components from splashed liquids PRS 1101 PRS 1409 PRS 2001 | Fluid enters housing | Device circuit non-functional | Delay of surgery | 3 | 1 | 3 | Yes | DFS | Design specifies material and assembly to prevent fluid ingress | Not practical | Ensure that material requirements and assembly are defined. The test for fluid ingress will be performed as part of product validation. ENG-PRT-439 | Material and assembly requirements are defined on drawing. Fluid Ingress Test Report ENG-RPT-555 Fluid Ingress Test Report 500433702 | Tested to IEC 60601-2-2 clause 201.11.6.5 b | 3 | 1 | 3 | No | NA | See footnote 3 and 4 |
| 17-d | Housing / Handle | Electrical connection from cable to electrode Electrically insulates external surfaces from internal conductors. PRS 1101 PRS 1407 PRS 1408 PRS 2001 | Insufficient amount and/or dielectric strength of housing material, or inner conductors too close to housing walls. | Burn holes through housing | Potential user or patient injury | 10 | 5 | 50 | Yes | DFS | Handpiece design, material and supplier selection and product validation | Not Practical | Ensure that material requirements and assembly are defined. The test for dielectric withstand will be performed as part of product validation. | Material and assembly requirements are defined on drawing. Dielectric Test Reports ENG-RPT-328 and ENG-RPT-412. Verified post EO sterilization in test report ENG-RPT-476 Dielectric Test Report 500433702 | IEC 60601-2-2 Clauses 201.8.8.3.103 and .104 UL Test Report ENG-RPT-059 See test reports ENG-RPT-328, ENG-RPT-412, and ENG-RPT-476. | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 18-d | Electrode | Connects to the distal end of the pencil body and provides electrosurgical activation at the tissue interface PRS 1101 PRS 1102 PRS 1601 PRS 1603 PRS 1604 PRS 1803 PRS 2001 | Electrode punctures the cap which then punctures package | Electrode punctures package | Patient infection | 10 | 5 | 50 | Yes | DFS | Application of industry standards to package design | Not Practical | Test pouches after sterilization and shipping | Shipping test reports ENG-RPT-330, ENG-RPT-413, ENG-RPT-503, ENG-RPT-535, and ENG-RPT-585 | Tested per ISO 11607-1 See test reports ENG-RPT-330, ENG-RPT-413, ENG-RPT-503, ENG-RPT-535, and ENG-RPT-585 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |

| | | | | | | | | | | |
|---------------------------------|--|--|--|--|--|--------------------------------|--|--|--|--|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | | | | Document Number ENG-RMF-045 | | | | |
| | Smoke Evacuation Accessories Risk Analysis | | | | | Revision: D | | | | |
| | | | | | | Page 31 of 52 | | | | |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------|---|---|---------------------------------|--------------------------|-----------------|---------------------------|-----|-----------------|------------------------------------|--|---------------------|--|--|---|----------------------------------|----------|---------------------------|----------------------------------|--|---|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Risk reduction considered (Yes/No) | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) |
| 19-d | Electrode | Connects to the distal end of the pencil body and provides electrosurgical activation at the tissue interface | SEE ANY AND ALL ELECTRODE RISK ANALYSIS DOCUMENTS, SUCH AS: E-Z Clean Electrodes Risk Analysis (ENG-RMF-008) ACE Blade User FMEA (ENG-RMF-009) ACE Process FMEA (PROP-RMF-003) ACE Blade Design FMEA (PROP-RMF-004) Resposable and LERIS Risk Analysis (ENG-RMF-010) LLETZ Loops FMEA (ENG-RMF-012) SS Electrode Risk Analysis (ENG-RMF-016) | | | | | | | | | | | | | | | | | | |
| 20-d | Cable | Connects the ESU plug connectors to the PCB of the Pencil Body Electrical and mechanical wire protection PRS 1101 PRS 1307 PRS 1308 PRS 1309 PRS 1404 PRS 1405 PRS 1406 PRS 2001 | Insufficient wire insulation Insufficient cable strain relief in pencil. | Shock or burn to user | Patient or user injury | 5 | 5 | 25 | Yes | DFS | Design, material selection, and product validation Design and test Cable Strain Relief to meet standard | Not Practical | Test dielectric properties of cable during validation Check strain relief properties of cable during validation | Dielectric Withstand Test Report ENG-RPT-328 and ENG-RPT-412 Mechanical Test Report ENG-RPT-329. Dielectric and Mechanical requirements met post EO sterilization. See test report ENG-RPT-452 | IEC 60601-2-2 Clauses 201.8.8.3.102, .103, .104 & 201.8.10.4.2 See test reports ENG-RPT-328, ENG-RPT-412 and ENG-RPT-329. Dielectric and Mechanical requirements met post EO sterilization. See test report ENG-RPT-452 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 32 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------|---|---|--|------------------------------------|-----------------|---------------------------|-----|--|--------------|---|--------------------|---|---|--|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 21-d | Cable | Connects the ESU plug connectors to the PCB of the Pencil Body Electrical and mechanical wire protection PRS 1101 PRS 1110 PRS 1307 PRS 1308 PRS 1309 PRS 2001 | Incorrect wire insulation Insufficient cable strain relief in pencil | Cable pulled from handpiece, possibly exposing bare conductors | Practitioner and/or patient injury | 5 | 5 | 25 | Yes | DFS | Design cable inside of tubing for first 64", material selection, and product validation Design and test Cable Strain Relief to meet standard | Not Practical | Check strain relief properties of cable during validation | Mechanical Test Report ENG-RPT-329, and ENG-RPT-418. Mechanical requirements met post EO sterilization. See test report ENG-RPT-452 | IEC 60601-2-2 Clause 201.8.10.4.2 See test reports ENG-RPT-329, and ENG-RPT-418. Mechanical requirements met post EO sterilization. See test report ENG-RPT-452 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 22-d | Cable | Connects the ESU plug connectors to the PCB of the Pencil Body Electrical and mechanical wire protection PRS 1101 PRS 1307 PRS 1308 PRS 1309 PRS 2001 | Incorrect geometry of cable / pencil housing junction and/or incorrect cable material resulting in damaged cord | Open circuit and / or intermittent connection | Delay of surgery | 5 | 5 | 25 | Yes | DFS | Design, material selection, and product validation | Not Practical | Check strain relief properties of cable during validation | Mechanical Test Report ENG-RPT-329. Mechanical requirements met post EO sterilization. See test report ENG-RPT-452 | IEC 60601-2-2 Clause 201.8.10.4.2 See test report ENG-RPT-329. Mechanical requirements met post EO sterilization. See test report ENG-RPT-452 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 23-d | Cable | Connects the ESU plug connectors to the PCB of the Pencil Body Electrical and mechanical wire protection | Cable in close proximity to other equipment or other equipment cables | Causes electromagnetic interference (EMI) or radio frequency interference (RFI) in other equipment | Disrupt other equipment | 5 | 5 | 25 | Yes | DFS | Cables are designed to meet requirements for HF applications | Not Practical | Ensure device meets requirements of IEC 60601-1. | Similar cable tested for EMC, see report ENG-RPT-052 | Tested to IEC 60601-1 Clause 17, see report ENG-RPT-052 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 33 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------|--|---|--|--|-----------------|---------------------------|-----|--|--------------|--|--------------------|--|---|---|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 24-d | Plug | Electrical connection between ESU and pencil cable PRS 1101 PRS 1401 PRS 1409 PRS 2001 | Cable leads not connected to plug pins | No ES current at electrode (inoperable unit) | Practitioner dissatisfaction (inconvenience) | 1 | 1 | 1 | NA | DFS | Design, material selection, and product validation | NA | Check cut and coag continuity, prior to and after plug strain relief test and fluid ingress test during validation | Mechanical Test Report ENG-RPT-329 Fluid Ingress Test Report ENG-RPT-555 | IEC 60601-2-2 Clauses 201.8.8.3.104 and 201.11.6.5 b See test reports ENG-RPT-329 and ENG-RPT-555 | 1 | 1 | 1 | NA | NA | NA |
| 25-d | Plug | Electrical connection between ESU and pencil cable PRS 1101 PRS 1310 PRS 1401 PRS 2001 | Plug pin(s) broke off in generator socket | Unable to plug in connector | Practitioner dissatisfaction (inconvenience) | 1 | 1 | 1 | NA | DFS | Design and material selection | NA | Ensure that material requirements and assembly are defined | Mechanical Test Report ENG-RPT-329 | IEC 60601-2-2 Clause 201.8.8.3.104 See test report ENG-RPT-329 | 1 | 1 | 1 | NA | NA | NA |
| 26-d | Plug | Electrical connection between ESU and pencil cable PRS 1101 PRS 1401 PRS 1409 PRS 2001 | Incorrect wire connections to plug pins Pin or cable resistance too high (does not activate ESU) | No ES current at electrode (inoperable unit) | Patient injury Practitioner dissatisfaction | 5 | 5 | 25 | Yes | DFS | Design, material selection, and product validation | Not Practical | Check cut and coag continuity, prior to and after cable flex test and fluid ingress test during validation | Mechanical Test Report ENG-RPT-329 Fluid Ingress Test Report ENG-RPT-555 | IEC 60601-2-2 Clauses 201.8.8.3.104 & 201.11.6.5 b See test reports ENG-RPT-329 and ENG-RPT-555 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 27-d | Plug | Electrical connection between ESU and pencil cable PRS 1101 PRS 1310 PRS 1401 PRS 1409 PRS 2001 | Cable leads not connected to correct plug pins Plug pins broke off in generator | cut switch activates coag or coag switch activates cut | Practitioner dissatisfaction, wrong effect | 5 | 5 | 25 | Yes | DFS | Design, material selection, and product validation | Not Practical | Check cut and coag continuity, prior to and after cable flex test and fluid ingress test during validation | Mechanical Test Report ENG-RPT-329 Fluid Ingress Test Report ENG-RPT-555 | IEC 60601-2-2 Clauses 201.8.8.3.104 & 201.11.6.5 b See test reports ENG-RPT-329 and ENG-RPT-555 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 34 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------------|---|--|--|--|-----------------|---------------------------|-----|--|--------------|---|--------------------|--|---|---|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 28-d | Plug | Electrical connection between ESU and pencil cable | Exposed pins of plug not fully inserted into the ES generator. | Exposed plug pins at ES generator. | User burn | 5 | 5 | 25 | Yes | DFS | Generators that comply with IEC 60601-1 are required to prevent electrical contact if pins are accessible | Not Practical | Ensure material, assembly, and inspection requirements are defined. | UL Testing to IEC 60601-1 and IEC 60601-2-2 | CB Test Report ENG-RPT-059 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 28.5-d | Plug | Electrical connection between ESU and pencil cable | Exposed pins of plug not fully inserted into the ES generator. | Exposed plug pins at ES generator. | User burn | 5 | 5 | 25 | Yes | IFS | Generators that comply with IEC 60601-1 are required to prevent electrical contact if pins are accessible | Not Practical | Ensure IFU instructs the user to fully insert the plug into the ESU | Zip IFU 3000312 01, Zip Ace IFU 3000317 01 | NA | 5 | 5 | 25 | No | N/A | See footnote 3 and 4 |
| 29-d | Plug | Electrical connection between ESU and pencil cable PRS 1101 PRS 1310 PRS 2001 URS 1004 | Wrong pin design or wrong ESU socket | Difficult to insert in ES generator Too loose – won’t stay plugged in | Delay of procedure | 1 | 1 | 1 | NA | DFS | Design plug to be compatible with Megadyne and other common ESU’s | NA | Test insertion / extraction force in Mechanical Protocol ENG-PRT-228 | Mechanical Test Report ENG-RPT-329 | See test report ENG-RPT-329 | 1 | 1 | 1 | NA | NA | NA |
| 30-d | Convoluted Tubing | Lightweight / strong tubing to convey smoke from handpiece to smoke evac unit PRS 1101 PRS 2001 URS 1008 URS 2001 URS 1005 | Improper tubing material specified | Disconnects from fitting Tubing web tears, No smoke collection | Possible user exposure to potential carcinogens and infectious by-products | 5 | 1 | 5 | Yes | DFS | Design, material selection, and product validation | NA | Test mechanical pull strength of tubing from connectors | Mechanical Test Report ENG-RPT-329 Tubing verified to maintain integrity post EO sterilization and after 3 yr. accelerated aging per ENG-PRT-466 | See test report ENG-RPT-329 and ENG-RPT-585 | 5 | 1 | 5 | NA | NA | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 35 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|---------------------------|--|--|---|--|-----------------|---------------------------|-----|--|--------------|---|--------------------|--|---|--|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 31-d | Swivel Connectors | Swivel to increase mobility of pencil body for ease of use PRS 1101 PRS 2001 URS 2001 URS 1005 | Incorrect tubing design or material | Disconnects from pencil or tubing, poor smoke collection | Possible user exposure to potential carcinogens and infectious by-products | 5 | 1 | 5 | Yes | DFS | Design and material selection | NA | Test mechanical pull strength of tubing from swivel | Mechanical Test Report ENG-RPT-329 Tubing verified to maintain integrity post EO sterilization and after 3 yr. accelerated aging per ENG-PRT-466 | See test report ENG-RPT-329 and ENG-RPT-585 | 5 | 1 | 5 | NA | NA | See footnote 3 and 4 |
| 32-d | Proximal Filter Connector | Designed for connection to a variety of smoke boxes PRS 1101 PRS 2001 URS 2001 | Customer does not use Megadyne smoke evacuator | Connector does not fit or disconnects easily, poor smoke collection | Possible user exposure to potential carcinogens and infectious by-products | 5 | 5 | 25 | Yes | DFS | Design universal Connector or provide adapter | Not Practical | Test with protocol | Mechanical Test Report ENG-RPT-329 | See test report ENG-RPT-329 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 33-d | Tubing Connector | Designed for connection to a variety of smoke boxes PRS 1101 PRS 2001 URS 2001 | Improper set-up | Connector compatible with other equipment, No smoke collection | Possible user exposure to potential carcinogens and infectious by-products, Customer dissatisfaction | 5 | 1 | 5 | NA | DFS | Design Connector to industry standards accepted | NA | Compare to predicate. | Mechanical Test Report ENG-RPT-329 | See test report ENG-RPT-329 | 5 | 1 | 5 | NA | NA | See footnote 3 and 4 |
| 34-d | Pencil Assembly | Sterile electrosurgical instrument PRS 1101 PRS 1501 PRS 1502 PRS 1503 PRS 1504 PRS 2001 | Materials not biocompatible | Patient reaction | Patient injury/ Infection/ reaction | 10 | 5 | 50 | Yes | DFS | Material selection | Not Practical | Ensure material requirements are defined. Biocompatibility testing | Biocompatibility Test Report ENG-RPT-337 | Tested per ISO 10993-1:2009 clauses 6.2.2.2, 6.2.2.3, & 6.2.2.4 See test report ENG-RPT-337 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |

| | | |
|---------------------------------|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | Revision: D |
| | | Page 36 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-----------------|--|---|--|-----------------------------|-----------------|---------------------------|-----|--|--------------|--|--------------------|---|--|--|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 35-d | Pencil Assembly | Sterile electrosurgical instrument PRS 1101 PRS 1402 PRS 1403 PRS 1404 PRS 1405 PRS 1407 PRS 1408 PRS 2001 | Materials have insufficient insulation properties after sterilization | Adversely affected by Gamma and Pencil housing and/or cable material breaks down | Patient or user gets burned | 10 | 5 | 50 | Yes | DFS | Design and material selection to withstand gamma | Not Practical | Test product for dielectric withstand after sterilization | Dielectric Test report ENG-RPT-328 and ENG-RPT-412 Dielectric Test Report 500433702 | IEC 60601-2-2 Clauses 201.8.8.3.102, .103, .104, & 201.8.10.4.2 See test report ENG-RPT-338 and ENG-RPT-412 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 35.1-d | Pencil Assembly | Sterile electrosurgical instrument PRS 1101 PRS 1402 PRS 1403 PRS 1404 PRS 1405 PRS 1407 PRS 1408 PRS 2001 | Materials have insufficient insulation properties after sterilization | Adversely affected by EtO and Pencil housing and/or cable material breaks down | Patient or user gets burned | 10 | 5 | 50 | Yes | DFS | Design and material selection to withstand gamma | Not Practical | Test product for dielectric withstand after sterilization | Dielectric Test report ENG-RPT-452 (Cable) and ENG-RPT-476 (Housing and Holster) Dielectric Test Report 500433702 | IEC 60601-2-2 Clauses 201.8.8.3.102, .103, .104, & 201.8.10.4.2 See test report ENG-RPT-452 (Cable) and ENG-RPT-476 (Housing and Holster) | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 36-d | Pencil Assembly | Sterile electrosurgical instrument PRS 1101 PRS 2001 | Sterilization Validation not performed on worst case configuration | Bioburden | Patient reaction | 10 | 5 | 50 | Yes | DFS | Design, material selection, and product validation | Not Practical | Test 15 ft. assembly as worst case for bioburden. | See Sterilization Validation test reports in DHF See DHF: Gamma Sterilization: New Dean Protocol P170829-E. EtO Sterilization: ENG-PRT-435, ENG-PRT-493, and ENG-PRT-489 | ISO 11137-2 See test reports in DHF See DHF: Gamma Sterilization: New Dean Report: R171106-E EtO Sterilization: ENG-RPT-545, ENG-RPT-593, and ENG-RPT-594 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 37 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|--------------------------|---|---|--|--|-----------------|---------------------------|-----|--|--------------|---|--------------------|--|---|---|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 37-d | Finished Device Assembly | Designed for smoke removal from the surgical site. PRS 1101 PRS 1311 PRS 2001 | Improper flow path design | Does not remove smoke adequately | Possible user exposure to potential carcinogens and infectious by-products | 5 | 5 | 25 | Yes | DFS | Design and product validation | Not Practical | Test with protocol | Mechanical Test Report ENG-RPT-329 | See test report ENG-RPT-329 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 38-d | Holster | Attaches securely to drape. Holds ES pencil when not in use Protects patient from burns when ES pencil is accidentally activated. | Inadequate design of drape attachment feature | Holster disconnects from drape and pencil falls on floor | Customer dissatisfaction Delay of surgery | 3 | 5 | 15 | Yes | DFS | Design to specification | Not Practical | Ensure design requirements are defined. | Holster design requirements are defined in drawing 5800100-01 | N/A | 3 | 1 | 3 | No | N/A | See footnote 3 and 4 |
| 39-d | Holster | Attaches securely to drape. Holds ES pencil when not in use Protects patient from burns when ES pencil is accidentally activated. PRS 1101 PRS 2001 | Inadequate dielectric strength of holster | Holster does not protect patient from accidental activations (when ES pencil stored in holster). | Patient burns | 10 | 5 | 50 | Yes | DFS | Design, material selection, and product validation | Not Practical | Test dielectric strength of holster per protocol | Dielectric Test report ENG-RPT-328, ENG-RPT-412, and ENG-RPT-476, Drawing: 5800100-01 | See test report ENG-RPT-328 and ENG-RPT-412 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 40-d | NOT USED | | | | | | | | | | | | | | | | | | | | |
| 41-d | Holster | Holster holds pencil PRS 1109 | Incorrect design | Product falls on the floor | Delay of procedure | 3 | 5 | 15 | Yes | DFS | Design for use and ensure color is Pantone 427, and Validation test | Not Practical | Test per protocol ENG-PRT-344 | Test report ENG-RPT-401, drawing 5800100-01 | Test report ENG-RPT-401 | 3 | 1 | 3 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 38 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk Evaluation (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|---------------------------|---|--|--|--------------------------|-----------------|---------------------------|-----|---|--------------|---------------------|--------------------|---|---|---|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 42-d | Pencil Multivac Packaging | To maintain sterile barrier protect ES pencil during shipping and storage condition PRS 1101 PRS 1102 PRS 1601 PRS 1603 PRS 1604 PRS 1803 PRS 2001 | Ineffective packaging for this application | Product unsterile (sterile barrier broken) Holes in packaging | Patient infection | 10 | 5 | 50 | Yes | DFS | Material Selection | Not Practical | Ensure material requirements are defined. Test according to protocol ENG-PRT-229, ENG-PRT-425, and ENG-PRT-441 | Material requirements are defined in drawing 2525-10, 2525-15, 2525-10EC, 2525-15EC, ME7251E, ME7251C, ME725M1C, and ME725M1E. See ENG-DMR-012 See test report ENG-RPT-330, ENG-RPT-413, ENG-RPT-503, ENG-RPT-546, and ENG-RPT-535 (EtO Sterilization) | Tested per ISO 11607-1, see test report ENG-RPT-330, ENG-RPT-413, ENG-RPT-546, and ENG-RPT-535 (EtO Sterilization). | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 39 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|---------------------------|---|--|---------------------------------|--------------------------|-----------------|---------------------------|-----|-----------------|--------------|---------------------|--------------------|---|---|---|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 42.5-d | Pencil Multivac Packaging | To maintain sterile barrier protect ES pencil during shipping and storage condition PRS 1101 PRS 1102 PRS 1601 PRS 1603 PRS 1604 PRS 1803 PRS 2001 | Ineffective packaging for this application | Product damaged | User dissatisfaction | 1 | 5 | 5 | Yes | DFS | Material Selection | Not Practical | Ensure material requirements are defined. Test according to protocol ENG-PRT-229, ENG-PRT-425, and ENG-PRT-441 | Material requirements are defined in drawing See ENG-DMR-0122525-10, 2525-15, 2525-10EC, 2525-15EC, ME7251E, ME7251C, ME725M1C, and ME725M1E. See test report ENG-RPT-330, ENG-RPT-413, ENG-RPT-503, ENG-RPT-546, and ENG-RPT-535 (EtO Sterilization) | Tested per ISO 11607-1, see test report ENG-RPT-330, ENG-RPT-413, ENG-RPT-546, and ENG-RPT-535 (EtO Sterilization). | 1 | 1 | 1 | No | N/A | NA |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 40 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-----------------|--|--|---------------------------------|-----------------------------------|-----------------|---------------------------|-----|--|--------------|---|--------------------|---|--|---|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 43-d | Inner bag/pouch | To contain and protect product | Inner bag/pouch opens unexpectedly | Product is contaminated | Unusable product delay of surgery | 3 | 5 | 15 | Yes | DFS | Design bag flap to be on the pouch side so flap stays in the package when peeled open | Not Practical | Ensure design requirements are defined. See ENG-PRT-425 and ENG-PRT-441 | Material requirements are defined in drawings See ENG-DMR-012 2525-10, 2525-15, 2525-10EC, 2525-15EC, 2525-10BN, and 2525-10ECBN. For ME7251E, and ME7251C See ENG-RPT-546. For ME725M1E and ME725M1C See ENG-RPT-535. | Test report ENG-RPT-546 and ENG-RPT-535 | 3 | 1 | 3 | No | N/A | See footnote 3 and 4 |
| 44-d | Inner bag/pouch | To contain and protect product PRS 1101 PRS 1501 PRS 1502 PRS 1503 PRS 2001 | Improper material selection | Bag material is cytotoxic | User injury | 10 | 5 | 50 | Yes | DFS | Vendor selection and part approval with testing for cytotoxicity | Not Practical | Test product for biocompatibility, see test protocol ENG-PRT-235 | Biocompatibility Report ENG-RPT-337 | Tested per ISO 10993-1:2009 clause 6.2.2.2, 6.2.2.3, & 6.2.2.4 See test report ENG-RPT-337 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 41 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------------------|--|---|---|------------------------------------|-----------------|---------------------------|-----|--|--------------|--|--------------------|---|--|---|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 45-d | Pencil Shipping Box | Protect product during shipment PRS 1101 PRS 1102 PRS 1601 PRS 1603 PRS 1604 PRS 1801 PRS 1802 PRS 1803 PRS 2001 | Box sized incorrectly – does not protect product adequately | Damaged pencil fails in use. | Practitioner and/or patient injury | 10 | 5 | 50 | Yes | DFS | Design, material selection, and product validation | Not Practical | Test according to protocol ENG-PRT-229 for 2525-10, ENG-PRT-327 for 2525-15, ENG-PRT-425 for ME7251C and ME7251E, and ENG-PRT-441 for ME725M1C and ME725M1E. Drawings 2525-10, 2525-15, 2525-10EC, 2525-15EC, 2525-10BN, 2525-10ECBN, ME7251E, ME7251C, ME725M1E, and ME725M1C. | Shipping test report ENG-RPT-330, ENG-RPT-413, ENG-RPT-503, ENG-RPT-546, and ENG-RPT-535. Drawings 2525-10, 2525-15, 2525-10EC, 2525-15EC, 2525-10BN, 2525-10ECBN, ME7251E, ME7251C, ME725M1E, and ME725M1C. | Tested per ISO 11607-1, see test reports ENG-RPT-330, ENG-RPT-413, ENG-RPT-503, ENG-RPT-546, and ENG-RPT-535. | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 46-d | Packaged Pencil Product | For electrosurgical use in the sterile field PRS 1101 PRS 1102 PRS 1201 PRS 1202 PRS 1601 PRS 1603 PRS 1604 PRS 1803 PRS 2001 | Transport and storage out of accepted limits | Pencil material breakdown chemically / physically from exposure to extreme shipping / storage conditions (Temp and/or Humidity) | Patient or user injury | 10 | 5 | 50 | Yes | DFS | Design, material selection, and product validation | Not Practical | Ensure material requirements are defined. Test according to protocol ENG-PRT-229 for 2525-10, ENG-PRT-327 for 2525-15, and ENG-PRT-425 for ME7251C. | Material requirements are defined in drawing 2525-10, 2525-15, 2525-10EC, 2525-15EC, 2525-10BN, 2525-10ECBN, ME7251E, and ME7251C. Test Reports ENG-RPT-330, ENG-RPT-413 and ENG-RPT-503 | Tested per ISO 11607-1, see test reports ENG-RPT-330, ENG-RPT-413 and ENG-RPT-503 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 42 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------------------|---|---|---|--|-----------------|---------------------------|-----|--|--------------|--|--------------------|--|--|--|----------------------------------|---------------------------|-----|--|---|---|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 46.2-d | Packaged Pencil Product | PRS 1720 PRS 1920 | Label print material not compatible with sterilization and shipping conditions | Product label information incomplete or illegible | Product not identifiable | 1 | 1 | 1 | Yes | IFS | Design material selection and product verification | Not Practical | Verify clarity and legibility post sterilization | After Gamma Sterilization see ENG-PRT-229. After EtO Sterilization see ENG-PRT-465. | After Gamma Sterilization see ENG-RPT-330. After EtO Sterilization see ENG-RPT-581. | 1 | 1 | 1 | No | N/A | See footnote 3 and 4 This is a potential regulatory compliance issue |
| 47-d | Extension Nozzles | Extend smoke capture when using extended electrodes in pockets or deep tissue areas PRS 1101 PRS 1111 PRS 1113 PRS 1201 PRS 1202 PRS 1305 PRS 1313 PRS 1501 PRS 1502 PRS 1503 PRS 2001 URS 1011 | Incorrect design or material Not a gamma stable material Opaque nozzle Wrong electrode exposure length | Fractures leaving material in patient Falls off into patient | Foreign body reaction or additional surgical intervention | 10 | 5 | 50 | Yes | DFS | Design, material selection, and product validation | Detect under X-ray | Test for mechanical strength | Mechanical test report ENG-RPT-329 | Tested per ISO 10993-1:2009 clause 6.2.2.2, 6.2.2.3, & 6.2.2.4 See test report ENG-RPT-337 See test report ENG-RPT-329 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 47.5-d | Extension Nozzles | Extend smoke capture when using extended electrodes in pockets or deep tissue areas PRS 1101 PRS 1111 PRS 1113 PRS 1201 PRS 1202 PRS 1305 PRS 1313 PRS 1501 PRS 1502 PRS 1503 PRS 2001 URS 1011 | Incorrect design or material Not a gamma stable material Opaque nozzle Wrong electrode exposure length | Transparency obstructs surgeon's view. Poor smoke capture | Possible user exposure to potential carcinogens and infectious by-products | 5 | 5 | 25 | Yes | DFS | Design, material selection, and product validation | Not Practical | Test for mechanical strength | Mechanical test report ENG-RPT-329 | Tested per ISO 10993-1:2009 clause 6.2.2.2, 6.2.2.3, & 6.2.2.4 See test report ENG-RPT-337 See test report ENG-RPT-329 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 43 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|----------------------------|--|--|--|--|-----------------|---------------------------|-----|--|--------------|-------------------------------|--------------------|--|--|---|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 48-d | Extension Nozzles | Extend smoke capture when using extended electrodes in pockets or deep tissue areas PRS 1101 PRS 1312 PRS 2001 | Improper flow path design | Does not remove smoke adequately | Possible user exposure to potential carcinogens and infectious by-products | 5 | 5 | 25 | Yes | DFS | Design and product validation | Not Practical | Flow test per protocol ENG-PRT-280 | Flow Test Report ENG-RPT-403 | See test report ENG-RPT-403 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 49-d | Extension Nozzle Packaging | To maintain sterile barrier protect pencil extension during shipping and storage condition PRS 1101 PRS 1102 PRS 1112 PRS 1602 PRS 1603 PRS 1604 PRS 1803 PRS 2001 | Ineffective packaging for this application | Product unsterile (sterile barrier broken) Holes in packaging | Patient infection | 10 | 5 | 50 | Yes | DFS | Material Selection | Not Practical | Ensure material requirements are defined. Test according to protocol ENG-PRT-229 | Material requirements are defined in drawing 2540 and 2560 Test with Protocol ENG-RPT-329 | Tested per ISO 11607-1, see test report ENG-RPT-329 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 49.5-d | Extension Nozzle Packaging | To maintain sterile barrier, protect pencil extension during shipping and storage condition PRS1101 PRS1102 PRS1112 PRS1602 PRS1603 PRS1604 PRS1803 PRS2001 | Ineffective packaging for this application | Product Damaged (noticed before surgery) | User dissatisfaction | 1 | 5 | 5 | Yes | DFS | Material Selection | Not Practical | Ensure material requirements are defined. Test according to protocol ENG-PRT-229 | Material requirements are defined in drawing 2540 and 2560 Test with Protocol ENG-RPT-329 | Tested per ISO 11607-1, see test report ENG-RPT-329 | 1 | 1 | 1 | No | N/A | NA |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 44 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|------------------------------------|---|---|--|------------------------------------|-----------------|---------------------------|-----|--|--------------|--|--------------------|--|--|---|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 49.7-d | Extension Nozzle Packaging | To maintain sterile barrier, protect pencil extension during shipping and storage condition PRS1101 PRS1102 PRS1112 PRS1602 PRS1603 PRS1604 PRS1803 PRS2001 | Ineffective packaging for this application | Product Damaged (noticed during surgery) | Delay of surgery | 3 | 5 | 15 | Yes | DFS | Material Selection | Not Practical | Ensure material requirements are defined. Test according to protocol ENG-PRT-229 | Material requirements are defined in drawing 2540 and 2560 Test with ReprotENG-RPT-329 | Tested per ISO 11607-1, see test report ENG-RPT-329 | 3 | 1 | 3 | No | N/A | NA |
| 50-d | Extension Nozzle Shipping/unit Box | Protect product during shipment PRS 1101 PRS 1102 PRS 1602 PRS 1603 PRS 1604 PRS 1801 PRS 1802 PRS 1803 PRS 2001 | Box sized incorrectly – does not protect product adequately | Product unsterile (sterile barrier broken) due to holes in packaging | Practitioner and/or patient injury | 10 | 5 | 50 | Yes | DFS | Design, material selection, and product validation | Not Practical | Test according to protocol ENG-PRT-229 Drawings 2540, 2560 | Shipping test report ENG-RPT-329 Drawings 2540, 2560 | Tested per ISO 11607-1, see test report ENG-RPT-329 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 50.5-d | Extension Nozzle Shipping/unit Box | Protect product during shipment PRS 1101 PRS 1102 PRS 1602 PRS 1603 PRS 1604 PRS 1801 PRS 1802 PRS 1803 PRS 2001 | Box sized incorrectly – does not protect product adequately | Product Damaged (noticed before surgery) | User dissatisfaction | 1 | 5 | 5 | Yes | DFS | Design, material selection, and product validation | Not Practical | Test according to protocol ENG-PRT-229 Drawings 2540, 2560 | Shipping test report ENG-RPT-329 Drawings 2540, 2560 | Tested per ISO 11607-1, see test report ENG-RPT-329 | 1 | 1 | 1 | No | N/A | NA |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 45 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|------------------------------------|---|---|--|--------------------------|-----------------|---------------------------|-----|--|--------------|--|--------------------|--|--|--|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 50.7-d | Extension Nozzle Shipping/unit Box | Protect product during shipment PRS 1101 PRS 1102 PRS 1602 PRS 1603 PRS 1604 PRS 1801 PRS 1802 PRS 1803 PRS 2001 | Box sized incorrectly – does not protect product adequately | Product Damaged (noticed during surgery) | Delay of surgery | 3 | 5 | 15 | Yes | DFS | Design, material selection, and product validation | Not Practical | Test according to protocol ENG-PRT-229 Drawings 2540, 2560 | Shipping test report ENG-RPT-329 Drawings 2540, 2560 | Tested per ISO 11607-1, see test report ENG-RPT-329 | 3 | 1 | 3 | No | N/A | NA |
| 51-d | Expiration Life | 3-year expiration date PRS 2002 | Unstable product | Does not meet DMR after 3 years | Patient injury | 10 | 5 | 50 | Yes | DFS | Design for 3-year real time aging | Not Practical | Test per protocol ENG-PRT-227 ENG-PRT-228 ENG-PRT-229 ENG-PRT-439 ENG-PRT-239 ENG-PRT-327 ENG-PRT-466 | Test report ENG-RPT-328 ENG-RPT-329 ENG-RPT-330 ENG-RPT-413 ENG-RPT-555 ENG-RPT-344 ENG-RPT-503 ENG-RPT-585 (after EtO sterilization) | Test report ENG-RPT-328 ENG-RPT-329 ENG-RPT-330 ENG-RPT-413 ENG-RPT-555 ENG-RPT-344 ENG-RPT-503 ENG-RPT-585 (after EtO sterilization) | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 52-d | ULPA Filter | Primary smoke filter and fluid trap PRS 1316 | Incorrect sizing of locking feature | Filter will not lock into box port | User dissatisfaction | 3 | 5 | 15 | Yes | DFS | Design port access feature to lock into smoke box port connection 1 st article measurements on Mold part feature | Not practical | Test according to Protocol ENG-PRT-238 force to insert and shear force tests and for 1 st article measurements. | See Report ENG-RPT-340 | See footnote 2 | 3 | 1 | 3 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 46 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|-------------|---|--|--|--|-----------------|---------------------------|-----|--|--------------|---|--------------------|---|------------------------|----------------|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 53-d | ULPA Filter | Primary smoke filter and fluid trap PRS 1316 PRS 1317 | Incorrect sizing of connecting features | Filter will not insert into box port Pencil will not connect to filter | User dissatisfaction | 1 | 1 | 1 | NA | DFS | Design port access features to fit smoke box port connection and pencil connections 1 st article measurements on Mold part features | NA | Test according to Protocol ENG-PRT-238 force to insert and shear force tests and for 1 st article measurements | See Report ENG-RPT-340 | See footnote 2 | 1 | 1 | 1 | NA | NA | NA |
| 54-d | ULPA Filter | Primary smoke filter and fluid trap | Wrong filter media | Reduced or no smoke capture (little to no suction) | Possible user exposure to potential carcinogens and infectious by-products | 5 | 5 | 25 | Yes | DFS | Design using correct ULPA filter media. Certificate of filter media performance | Not practical | Test according to Protocol ENG-PRT-238 for flow. | See Report ENG-RPT-340 | See footnote 2 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |
| 55-d | ULPA Filter | Primary smoke filter and fluid trap | Inadequate material used for housing | Unable to connect box or smoke pencil to filter (connection features are broken off) | User dissatisfaction | 1 | 1 | 1 | NA | DFS | Design and specify appropriate material | NA | Test according to Protocol ENG-PRT-238 for flow Specify correct material | See Report ENG-RPT-340 | See footnote 2 | 1 | 1 | 1 | NA | NA | NA |
| 56-d | ULPA Filter | Primary smoke filter and fluid trap PRS 1314 PRS 1315 | Inadequate material and/or weld joint design | Filter leaks fluid | Bio-hazard | 10 | 5 | 50 | Yes | DFS | Design and specify appropriate material and weld joint design | Not practical | Test according to Protocol ENG-PRT-238 for Leak test. Specify correct material and weld joint design | See Report ENG-RPT-340 | See footnote 2 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 56.5-d | ULPA Filter | Primary smoke filter and fluid trap PRS 1314 PRS 1315 | Inadequate material and/or weld joint design | Filter leaks smoke | Possible user exposure to potential carcinogens and infectious by-products | 5 | 5 | 25 | Yes | DFS | Design and specify appropriate material and weld joint design | Not practical | Test according to Protocol ENG-PRT-238 for Leak test. Specify correct material and weld joint design | See Report ENG-RPT-340 | See footnote 2 | 5 | 1 | 5 | No | N/A | See footnote 3 and 4 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 47 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation Risk reduction considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|---|--|---|---|--|-----------------|---------------------------|-----|--|--------------|--|--------------------|---|---|--|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification/ Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 57-d | Charcoal Filter | Capture odors and residual chemicals PRS 1114 | Threads on muffler damaged | Charcoal Filter will not thread to box port | Increased noise reduced filtration, User dissatisfaction | 1 | 1 | 1 | NA | DFS | Design packaging to protect filter threads | NA | Test according to Protocol ENG-PRT-238 for ship test | See Report ENG-RPT-340 | See footnote 2 | 1 | 1 | 1 | NA | NA | NA |
| 58-d | ULPA and Charcoal Filters | Filter air for surgical smoke evacuation PRS 1605 | Shipping Damage | Filter Failure | Biohazard | 10 | 5 | 50 | Yes | DFS | Design adequate packaging | Not practical | Test according to Protocol ENG-PRT-238 for ship test | See Report ENG-RPT-340 | See footnote 2 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 59-d | Filter Shipping Box | Protect product during shipping, handling, and storage PRS 1801 PRS 1802 PRS 1803 | Box sized incorrectly – does not protect product adequately | Damaged filter fails in use. | Biohazard | 10 | 5 | 50 | Yes | DFS | Design, material selection, and product validation | Not Practical | Test according to protocol ENG-PRT-238 See drawings 2211, 2220 | Shipping test report ENG-RPT-340 | See test report ENG-RPT-340 See drawings 2211, 2220 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 60-d | Filter Packaged Product | Protect filters during shipping, handling, and storage PRS 1201 PRS 1202 PRS 1803 | Transport and storage out of accepted limits | Filter material breakdown chemically / physically from exposure to extreme shipping / storage conditions (Temp and/or Humidity) | Biohazard | 10 | 5 | 50 | Yes | DFS | Design, material selection, and product validation | Not Practical | Ensure material requirements are defined. Test according to protocol ENG-PRT-238 | Material requirements are defined in drawing 2211 and 2220 Test Report ENG-RPT-340 | See test report ENG-RPT-340 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 61-d | Applied and non-applied part intended to supply heat to the patient | Deliver RF energy PRS 1319 PRS 1320 | RF Energy | Heat exceeding 41 C degrees | Patient Burn | 10 | 5 | 50 | Yes | DFS | Design, material selection, and product validation | Not Practical | Ensure material requirements are defined. Testing to standards IEC 60601-1 | Material requirements are defined in drawing. | Test Report ENG-RPT-395 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |
| 62-d | Instructions for use | Provide IFU with device | Do not provide IFU with device | Customer unable to use the device | Customer dissatisfaction | 1 | 1 | 1 | NA | IFS | NA | NA | NA | NA | NA | 1 | 1 | 1 | NA | NA | NA |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 48 of 52 |

Zip-Pen PROCESS FMEA resides with NDT (vendor)

Zip-Pen Megadyne EO Sterilization Process FMEA (applies to ME725M1E and ME725M1C)

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|---|--|--|---|--|-----------------|------------------------------|-----|---|---|--|--|---|---------------------------------|---------------------------------|-------------------------------------|------------------------------|-----|---|---|---------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification / Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 1-p | Receive Gaylord requiring sterilization | Sterilize Product | Product not marked as awaiting EO sterilization. | Not identifiable as product awaiting EO Sterilization | Use of non-sterile product (intended to be sterile). | 10 | 2 | 20 | Yes | Training and Traveler. | Traveler OPER-FRM-134 | Non-Sterile Product awaiting sterilization tape found on gaylords. | Validation of process traveler OPER-FRM-134 per protocol ENG-PRT-474 | See ENG-RPT-586 | See ENG-RPT-586 | 10 | 1 | 10 | No | N/A | See footnote 3 |
| 2-p | Transfer Gaylord to back building. | Product is not damaged. | Product dropped/damaged in transit. | Damaged Product – no sterility breach. | Damaged product not identified and used in surgery. Increased thermal damage to tissue, especially during plastic surgery. | 5 | 2 | 10 | Yes | DFS, Training and SOPs | Design, material selection, and process validation | QA-SOP-008 | Test according to protocol ENG-PRT-441 and validate process per ENG-PRT-474. | See ENG-RPT-535 and ENG-RPT-586 | See ENG-PRT-535 and ENG-RPT-586 | 5 | 1 | 5 | No | N/A | See footnote 3 |
| 3-p | Transfer Gaylord to back building. | Product is not damaged. | Product dropped/damaged in transit. | Damaged Product – sterility breach. | Use of non-sterile product (intended to be sterile). | 10 | 2 | 20 | Yes | DFS, Training and SOPs | Design, material selection, and process validation | QA-SOP-008 | Test according to protocol ENG-PRT-441 and validate process per ENG-PRT-474. | See ENG-RPT-535 and ENG-RPT-586 | See ENG-RPT-535 and ENG-RPT-586 | 10 | 1 | 10 | No | N/A | See footnote 3 |
| 4-p | Re-Palletize Product for sterilization | Sterilize Product | Product not properly palletized for sterilization. | Product not sterilized per validated process/cycle | Use of non-sterile product (intended to be sterile). | 10 | 2 | 20 | Yes | Training. Work instruction, and Traveler. | Traveler OPER-FRM-134 and OPER-WI-048 | Not Practical | Validation of process traveler OPER-FRM-134 per ENG-PRT-474 and update to OPER-WI-048 | See ENG-RPT-586 | See ENG-RPT-586 | 10 | 1 | 10 | No | N/A | See footnote 3 |

| | | | |
|---------------------------------|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | | Revision: D |
| | | | Page 49 of 52 |

| Risk Analysis | | | | | | Risk Estimation | | | Risk Evaluation considered (Yes/No) | Risk Control | | | | | | Risk Estimation of Residual Risk | | | Risk Evaluation of Residual Risk | | |
|---------------|---|---|--|---|--|-----------------|---------------------------|-----|---|---|-----------------------|--------------------|---|-----------------|-----------------|----------------------------------|---------------------------|-----|--|---|------------------------------|
| Item# | Application | Intended Use or Characteristic (Function) | Initiating Events / Circumstances / Causes | Potential Failure Mode (Hazard) | Potential Effects (Harm) | Severity | Probability of Occurrence | RPN | | Options | Preventive Measures | Detection Measures | Verification / Validation Plan | Verification | Validation | Severity | Probability of Occurrence | RPN | Additional Risk after Risk Control? (Yes/No) | Additional Risk Control (if applicable) | Risk Benefit Analysis/ Notes |
| 5-p | Re-Palletize Product for Sterilization | Prepare product for Sterilization. Product not damaged. | Product dropped/damaged during palletization | Damage to product packaging – no sterility breach | Cosmetic damage to packaging | 1 | 2 | 2 | Yes | Training and Work Instructions | OPER-WI-048 | Not Practical | Validation of process traveler OPER-FRM-134 per ENG-PRT-474 | See ENG-RPT-586 | See ENG-RPT-586 | 1 | 1 | 1 | No | N/A | See footnote 3 and 4 |
| 6-p | Final Release of Product for Distribution | Release product for distribution to customers | Incorrect handling | Final Release occurs without sterilization | Use of non-sterile product (intended to be sterile). | 10 | 2 | 20 | Yes | Training. Work instruction, and Traveler. | Traveler OPER-FRM-134 | Not Practical | Validation of process traveler OPER-FRM-134 per ENG-PRT-474 | See ENG-RPT-586 | See ENG-RPT-586 | 10 | 1 | 10 | No | N/A | See footnote 3 and 4 |

The risk assessment was reviewed and updated by Kiran Kumar (Lifecycle Quality Engineer) and Scot Harris (Lifecycle Quality Engineer Lead) on September 29, 2020.

Smoke Evacuation Pencil & Accessories IFU’s

| Product Code | MOD | IFU P/N |
|---|-----|------------|
| 252510* 2525-10 252510EC* 2525-10EC, 252515* 2525-15 252515EC* 2525-15EC | 00 | 3000312-01 |
| | 01 | 3000343-01 |
| ME7251C ME7251E ME725M1C ME725M1E | 01 | 3000317-01 |
| | 02 | 3000336-01 |
| | 00 | 3000313-01 |

| | | |
|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | Revision: D |
| | | Page 50 of 52 |

| | | |
|-------------------------------|----|------------|
| 2540J*, 2540, 2560, 2560J* | 01 | 3000344-01 |
| 2211, 2211J* | 00 | 3000190-01 |
| | 01 | 3000346-01 |
| 2220J*, 2220 | 00 | 3000191-01 |
| | 01 | 3000347-01 |

DRAFT

| | | |
|---------------------------------|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | Revision: D |
| | | Page 51 of 52 |

Table 1. Severity Risk Ranking

| Description | Impact | Rank |
|-------------|---|------|
| Critical | The potential problem could result in death or serious injury to the patient or user, such as alternate current site injury, patient burn, or infection. | 10 |
| Important | The potential problem could result in a non-serious injury to the patient or user that requires medical intervention, such as reduced or no surgical effect, pressure sore, or allergic reaction. | 5 |
| Minor | The potential problem could result in a non-serious injury to the patient or user that does not require medical intervention, such as a delay in the procedure. | 3 |
| Negligible | There is no risk of injury to the patient or user. | 1 |

Table 2. Probability of Occurrence Risk Ranking

| Qualitative Description | Quantitative Description | Rank |
|-------------------------------|--------------------------|------|
| Continual (expected to occur) | ≥ 10% | 5 |
| Frequent (likely to occur) | > 1% and < 10% | 3 |
| Occasional (can occur) | > 0.1% and ≤ 1% | 2 |
| Rare (unlikely to occur) | ≤ 0.1% | 1 |

Probability of Occurrence Risk Ranking

Determination of RPN

| | | | Severity Risk Ranking | | | |
|--|------------|---------|-----------------------|-------|-----------|----------|
| | | | Negligible | Minor | Important | Critical |
| | | Ranking | 1 | 3 | 5 | 10 |
| Probability of Occurrence Risk Ranking | Continual | 5 | 5 | 15 | 25 | 50 |
| | Frequent | 3 | 3 | 9 | 15 | 30 |
| | Occasional | 2 | 2 | 6 | 10 | 20 |
| | Rare | 1 | 1 | 3 | 5 | 10 |

FOOTNOTES

| | |
|----|--|
| #1 | This medical device is restricted to use by a physician. |
| #2 | A review of clinical experience (as noted in complaint analysis section) reveals that the Probability of Occurrence levels actually fall within the range listed in the Post Risk Control section. |

| | | |
|--|--|---------------------------------------|
| Megadyne Medical Products, Inc. | RISK ANALYSIS | Document Number <i>ENG-RMF-045</i> |
| | Smoke Evacuation Accessories Risk Analysis | Revision: D |
| | | Page 52 of 52 |

| | |
|----|---|
| #3 | Residual risk is deemed acceptable through verification / validation testing. |
| #4 | The risk is reduced ALAP (As Low As Possible). The risk is associated with the use of electrosurgery regardless of the medical device used. The benefit derived from the use of the medical device is greater than the risks associated with its use. The risks associated with electrosurgery are accepted by the patient through informed consent. The surgeons are trained in the proper application of electrosurgery and the prevention of injuries. |

DRAFT