



The Electrosurgical Authority®

DOCUMENT NUMBER: ENG-RPT-413

DOCUMENT TITLE: Zip Pen Shipping Test Report

DOCUMENT NOTES:

Document Information

Revision: 001

Vault: MEG-Rel

Status: Release

Document Type: ENG-RPT

Date Information

Effective Date: 16 Jan 2015 3:00 PM

Expiration Date:

Release Date: 16 Jan 2015 3:00 PM

Next Review Date:

Control Information

Author: MGLASSETT

Owner: MGLASSETT

Previous Number:

Change Number: 2015-ENG-DCO-006

Signature Manifest

Document Number: ENG-RPT-413

Revision: 001

Title: Zip Pen Shipping Test Report

All dates and times are in Mountain Standard Time.

Zip Pen Ship Test Validation

Change Request

Name/Signature	Title	Date	Meaning/Reason
Paul Borgmeier (PBORGMEIER)		14 Jan 2015, 08:58:43 AM	Approved

Collaboration

Name/Signature	Title	Date	Meaning/Reason
Balaji Sudabattula (BSUDABATTULA)			
Dave Shimkus (DSHIMKUS)			
Julie Barber (JBARBER)			
Jill Skoczen (JTRASAMAR)			
Paul Borgmeier (PBORGMEIER)			
Tara Rogers (TROGERS)			
Ronda Magneson (RMAGNESON)	RA Director	16 Jan 2015, 01:25:13 PM	Complete
Mark Glassett (MGLASSETT)		16 Jan 2015, 01:28:01 PM	Complete

Document Review

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		16 Jan 2015, 02:08:42 PM	Complete

RA-Approval

Name/Signature	Title	Date	Meaning/Reason
Ronda Magneson (RMAGNESON)	RA Director	16 Jan 2015, 02:15:45 PM	Approved

QA-Approval

Name/Signature	Title	Date	Meaning/Reason
Balaji Sudabattula (BSUDABATTULA)		16 Jan 2015, 02:23:15 PM	Approved

ENG-Approval

Name/Signature	Title	Date	Meaning/Reason
----------------	-------	------	----------------

Mark Glassett (MGLASSETT)

Paul Borgmeier
(PBORGMEIER)

16 Jan 2015, 02:58:59 PM Approved

Training Review

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		16 Jan 2015, 03:00:05 PM	Approved

Final Release

Name/Signature	Title	Date	Meaning/Reason
Lucy Richards (LRICHARDS)		16 Jan 2015, 03:00:17 PM	Approved

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 1 of 23

Authored By: Mark Glassett

TABLE OF CONTENTS	Page
1. ABSTRACT	1
2. OBJECTIVE	1
3. RESULTS	1
4. DISCUSSION	3
5. CONCLUSIONS	5
6. RECOMMENDATIONS	5

1. ABSTRACT

Zip Pen Catalog #2525-10 samples were tested per protocol XENG-PRT-229 Rev A to show compliance to ISO 11607-1:2009. The samples passed the requirements of the protocol after one revision for seal width. These products were exposed to three year accelerated aging, shipping temperature extremes and the ASTM D4169-5 performance test. The Zip Pen Catalog #2525-10 passed the testing.

2. OBJECTIVE

The objective of this test report is to document the ship testing that was done on the Zip Pen catalog item 2525-10 after accelerated aging to simulate three years expiration life, exposure to transport and storage conditions and package performance testing per ASTM D4169-5.

3. RESULTS

3.1. Accelerated Aging

The Zip Pen samples for testing were subjected to accelerated aging at 55°C to simulate 3 years of shelf life. The samples that were exposed to this aging are as follows:

Catalog Number	Lot Number	Quantity
2525-10	S140303	40

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 2 of 23

Documentation for aging is shown in Appendix IA. Note that this is a repeat of the shipping test for this product. Supporting documentation for package material compliance to ISO 11607-1 is located in ENG-RPT-330. The certificate of conformance and the sterilization certificate for the test product are also shown in appendix IA.

3.2. Transport and Storage Cycle and Pre-Conditioning.

The Zip Pen samples for testing were subjected to the transport and storage cycle required by the protocol. This cycle includes temperatures from -40°C to 55°C and humidity's from 15% to 95%. This treatment is shown in Appendix I. Following the transport and storage cycle, the samples were pre-conditioned per protocol ENG-PRT-229 Rev A and ASTM D4169-05. The pre-conditioning parameters were to 23°C and 50% RH for a minimum of 72 hours. The cycling and preconditioning were performed by the Lab Technician. The technician sign off for this treatment is shown in Appendix I.

3.3. Shipping Test

The shipping test was performed by the Lab Technician. The shipping test followed protocol XENG-PRT-229 Rev A. None of the boxes broke open during the test. There were minor indentations on the corners and edges of the boxes which is typical and acceptable per the protocol. This testing is documented on the log sheet in Appendix I.

3.4. Print and Lot Number Clear and Legible

Forty samples (two boxes) of 2525-10 were inspected per protocol XENG-PRT-229 Rev A. The package print was clear and legible and passed the protocol acceptance criteria. The Lot numbers were also clear and legible and passed the protocol acceptance criteria. This testing is documented on the log sheet in Appendix II.

3.5. Bubble Leak Test

The Bubble Leak Test was performed on 40 samples (two boxes) of 2525-10. The Bubble Leak test followed protocol XENG-PRT-229 Rev A. The 2525-10 samples passed the bubble leak test criteria of the protocol. The data is shown in appendix III.

3.6. Dye Penetration Test

The dye penetration test was performed per protocol XENG-PRT-229 rev A. The procedure follows ASTM F1929. Forty Zip Pen packages were dye tested and all 40 passed the protocol. See data in appendix IV.

3.7. Burst Test

The package Burst Test was performed on 40 samples of 2525-10. The package burst test followed protocol XENG-PRT-229 Rev A. The minimum burst for the protocol is 19 in. H₂O. The 40 samples passed the burst test with an average burst value of 27.9 in. H₂O

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 3 of 23

and a minimum value of 24.8 in. H₂O. See data in appendix V and summary statistics in Appendix VI.

3.8. Minimum Seal Width

Numerous tests of the Zip Pen package have resulted in passing all requirements except the seal width. The original seal width specification was set at 0.30" as a starting point without any data. The aged samples for this and previous protocols did not meet this minimum. As a result of all testing that has taken place on the Zip Pen package, the seal width specification has been re-evaluated. All of the data and the new seal width specification is summarized in Test report ENG-RPT-414. The new seal width specification is 0.20" minimum. The minimum seal width measured for this testing was 0.22". The average minimum seal width is 0.30" (average of the lowest measurement of each package). The Zip Pen packages for this protocol meet this new requirement. Seal width data is shown in Appendix VI.

3.9. Electrode Cap and Electrode Damage

Electrode cap and electrode damage was inspected and documented in a previous report, see test report ENG-RPT-330.

4. DISCUSSION

The samples used for the ship testing are shown in the table below:

Catalog Number	Lot Number	Quantity
2525-10	S140303	40

The Zip Pen samples were sterilized with a gamma irradiation dose equal to or greater than the minimum 50 kGy required by drawing X2525-10 rev 03. Sterilization certificate and C of C are shown in appendix IA.

4.1. Shipping Test

Two cases of 2525-10 were subjected to the shipping test. The testing followed protocol XENG-PRT-229 Rev A. This testing consisted of temperature extremes, pre-conditioning, drop test, compression test, vibration and a second drop test. The shipping test was completed with both boxes meeting the acceptance criteria.

4.2. Print and Lot Number Clear and Legible

The print and lot number inspections were performed by the Engineering Technician. Two cases of each 2525-10 were inspected per the protocol. The package print was clear and legible and passed the protocol acceptance criteria. The Lot numbers were clear and

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 4 of 23

legible and passed the protocol acceptance criteria. This inspection is documented on the log sheet in Appendix II.

4.3. Bubble Leak Test

The bubble leak test was performed per the requirements of the protocol and ASTM F2096 on two boxes of 2525-10. There are 20 pouches in each box for a total of 40 tests. All of the 2525-10 samples passed the bubble leak test criteria of the protocol. The data is shown in appendix III.

4.4. Dye Penetration Test

In an effort to better understand the seal width failures from previous tests, the dye penetration test was added to the protocol. The dye penetration test was conducted per ASTM F1929. The test creates a visible line where the seal exists. It allows for visual confirmation of a seal and an approximation of seal width before the burst test. The requirement of the protocol is that none of the seals are breached as evident from the dye test. The test samples passed this test. The data is shown in appendix IV.

4.5. Package Burst Test

The package burst test was performed on 40 each of 2525-10 per the protocol. All 40 samples passed the minimum requirement of 19 in. H₂O. The average burst value was 27.9 in. H₂O and the minimum value was 24.8 in. H₂O. See raw data in appendix V and a summary in appendix VI. Note that previous protocols have tested a smaller sample size for burst test. In this revision of the protocol the number of burst test samples was increased. This was done to show that a lower seal width specification does not correspond to lower burst values, see test report ENG-RPT-414.

4.6. Minimum Seal Width

The nominal seal width is designed to be .37" wide and the original minimum specification on drawing X2525-10 Rev A was 0.30". This minimum specification was chosen prior to any data collection. After testing many samples, both aged and non-aged, a justification has been written to lower the minimum seal width requirement to 0.20". Refer to test report ENG-RPT-414. A summary for the minimum seal widths of the packages in this protocol is as follows:

Minimum measured seal width	0.22"
Average measured Minimum	0.30"
Standard Deviation of minimum measurements	0.032

The seal width for the product passes this revised specification. See data in appendix VI.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 5 of 23

5. CONCLUSIONS

This testing demonstrates that the Zip Pencil 2525-10 (10 foot) complies with the requirements of the revised protocol ENG-PRT-229 Rev 002 and ISO 11607-2:2006. This testing included three year accelerated aging, Temperature Storage extremes, and package performance testing.

6. RECOMMENDATIONS

This testing was performed to demonstrate compliance of the Zip Pen packaging to ISO 11607-1:2009 and ISO 11607-2:2006 after three year accelerated aging to support the three year expiration life. The package met requirements for Package Performance. It is recommended that real time age samples be put aside for testing per Megadyne Protocol ENG-PRT-057. The product can be released to market upon completion of regulatory approvals and completion of design control.

The testing also establishes shipping and storage conditions for labeling per IEC 60601-1:2005 clause 7.2.17. The shipping box labels of the products will show the international symbols for shipping and storage with temperatures of 5°C to 50°C and relative humidity of 15% to 95%. The IFU will include the note "Normal storage conditions are assumed. Brief excursions to limits allowed".

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 6 of 23

APPENDIX IA
Aging, C of C and Sterilization Certificate

Accelerated Aging In Process

Product: 2525-10 E-Z Clear Pencil
Lot Numbers S140303
4 Boxes
Temperature 55°
Relative Humidity Ambient

Required Time 15 weeks 6 days
(111 days)

Thermotron ID Number 01268

Last Calibration Date 5-8-14

Calibration Due Date 5-31-15

Start	Time	Initials	Stop	Time	Total	Initials
4 Sep 2014	14:00	MSJ	9-13-2014	6:30	9:	Power off
9-16-2014	13:45	MSJ	9-21-2014	16:00	14	Power out
9-22-2014	8:40	MSJ	1-2-2015	8:00		

If this aging needs to be interrupted for any reason,
contact Mark Glassett at ext. 845

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 7 of 23

APPENDIX IA Aging, C of C and Sterilization Certificate



NEW DEANTRONICS TAIWAN LTD. WITH
QUALITY SYSTEM
CERTIFIED BY DNV
~ISO 13485:2003~

Certificate of Compliance

<p>Taiwan</p> <p>大德企業股份有限公司 New Deantronics Taiwan, Ltd. 新城市系統中心大樓 4 樓 51 號 12 樓 12F, No.51, Sec.4, Chong Yang Rd., Tu Cheng Dist., New Taipei City 23675, Taiwan R.O.C. Tel: +886 2 2268-1726 Fax: +886 2 2268-3800</p>	<p>Customer Name: MEGADYNE MEDICAL PRODUCTS, INC. Invoice Number : MD-493/14 P. O. Number: 24443 Customer P/N: 6020190-01 (2525-10) (With New Electrodes) Drawing Number: X2525-10 Rev.04 Redline New Deantronics P/N: PB352SM1 (With New Electrodes) Lot Number: S140303 Expiration Date: 2017-07 Quantity: 80 Pieces (4 Boxes) Carton Number: #1 ~#2 (2 Cartons)</p> <p>Signed: <i>D. Y. Chen</i> Printed Name: Da-Yu Chen Title: Q.A. Director</p>	<p>Date: 08-22-2014</p>
---	--	--------------------------------

<p>USA</p> <p>New Deantronics Ltd. 1690 North California Blvd. Suite 1040 Walnut Creek, CA 94596 Tel: +1 (925) 260-6388 Fax: +1 (925) 260-1788</p>	<table> <tr> <th><u>Materials</u></th> <th><u>ND Mat'l Lot#</u></th> <th><u>Megadyne Lot#</u></th> </tr> <tr> <td>Cable</td> <td>C302600 Conductor : 7/0.16*3BC(26 AWG) Bare copper Insulation: PPE, one red, one blue & one white Outer jacket: PVC, Gray</td> <td>140311Z</td> </tr> <tr> <td>Plug material</td> <td>F506100, TOP, ABS PA707 F506200, Bottom, ABS PA707</td> <td>140408-140423 140408-140423</td> </tr> <tr> <td>Terminal</td> <td>T101702, Nickel plated brass</td> <td>140310</td> </tr> <tr> <td>Overmold material</td> <td>R900801, TPR</td> <td>130801</td> </tr> <tr> <td>Swivel, male</td> <td>F916600, HDPE</td> <td>S140064</td> </tr> <tr> <td>Swivel, Female</td> <td>F917000, HDPE</td> <td>S130276-S140064</td> </tr> <tr> <td>Connector</td> <td>F916800, HDPE</td> <td>S140067</td> </tr> <tr> <td>Collet Terminal</td> <td>T202502, Phosphor Bronze contact plated nickel</td> <td>140408-140414</td> </tr> <tr> <td>Button</td> <td>F302500, NYLON 66</td> <td>140219-140317</td> </tr> <tr> <td>Tape</td> <td>A100401, PTFE</td> <td>140303-140521</td> </tr> <tr> <td>PCB</td> <td>H201602, PCB FR-4</td> <td>140311</td> </tr> </table>	<u>Materials</u>	<u>ND Mat'l Lot#</u>	<u>Megadyne Lot#</u>	Cable	C302600 Conductor : 7/0.16*3BC(26 AWG) Bare copper Insulation: PPE, one red, one blue & one white Outer jacket: PVC, Gray	140311Z	Plug material	F506100, TOP, ABS PA707 F506200, Bottom, ABS PA707	140408-140423 140408-140423	Terminal	T101702, Nickel plated brass	140310	Overmold material	R900801, TPR	130801	Swivel, male	F916600, HDPE	S140064	Swivel, Female	F917000, HDPE	S130276-S140064	Connector	F916800, HDPE	S140067	Collet Terminal	T202502, Phosphor Bronze contact plated nickel	140408-140414	Button	F302500, NYLON 66	140219-140317	Tape	A100401, PTFE	140303-140521	PCB	H201602, PCB FR-4	140311
<u>Materials</u>	<u>ND Mat'l Lot#</u>	<u>Megadyne Lot#</u>																																			
Cable	C302600 Conductor : 7/0.16*3BC(26 AWG) Bare copper Insulation: PPE, one red, one blue & one white Outer jacket: PVC, Gray	140311Z																																			
Plug material	F506100, TOP, ABS PA707 F506200, Bottom, ABS PA707	140408-140423 140408-140423																																			
Terminal	T101702, Nickel plated brass	140310																																			
Overmold material	R900801, TPR	130801																																			
Swivel, male	F916600, HDPE	S140064																																			
Swivel, Female	F917000, HDPE	S130276-S140064																																			
Connector	F916800, HDPE	S140067																																			
Collet Terminal	T202502, Phosphor Bronze contact plated nickel	140408-140414																																			
Button	F302500, NYLON 66	140219-140317																																			
Tape	A100401, PTFE	140303-140521																																			
PCB	H201602, PCB FR-4	140311																																			

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 8 of 23


APPENDIX IA

Aging, C of C and Sterilization Certificate



	Dome	M100600, 6mm ROUND W/DIMPLE & FEET	130820	
	Pen Body	F105400 ABS+TPR	140219-140321	
Taiwan 大悅企業股份有限公司 New Deantronics Taiwan, Ltd. 台北市中山路4段112號 12F, No.51, Sec.4, Chong Yang Rd., Tu Cheng Dist., New Taipei City 23675, Taiwan R.O.C. Tel: +886 2 2258-1726 Fax: +886 2 2258-3800	Nozzle	F916300, PC	140310-140326	
	Collet Holder	F918100, PC	140326-140407	
	Carriage	F916500, ABS	140326	
	Snap Swivel, Male	F916400, HDPE	S140065	
	Snap Swivel, Female	F918300, HDPE	S140065	
	Tubing, Connector	F916700, HDPE	S130318-S140066	
	Tubing, Convoluted	P305600, EVA	140303-140326	
	Tubing, Convoluted	P305700, EVA	140303-140326	
	Holster	F916900, HDPE	140318-140407	
	Blade	G102700, Coated Megadyne P/N: 0012BN5 (Provided by Megadyne)	140610B (80pcs)	142008
USA New Deantronics Ltd. 1930 North California Blvd. Suite 1040 Walnut Creek, CA 94596 Tel: +1 (925) 280-8388 Fax: +1 (925) 280-1788	Paper Band	A900300	140207	
	Tyvek	A000400	140425	
	Nylon Film	A000500 Nylon	140211	
	Glue	S400900, Loctite 4061	140210	
	Ink	S102000, Green, PMS356C	121127	
	IFU	I010100,P/N:3000185-01 Rev.A	140512	

APPENDIX IA
Aging, C of C and Sterilization Certificate



行政院原能會核准設立照射廠執照証號 IRRADIATION PLANT NO: 物字第 1100223 號

客戶名稱 CUSTOMER NAME : 大瓏企業(股)公司
NEW DEANTRONICS TAIWAN LTD.
照射日期 IRRADIATION RUN DATE : 2014/06/23

照射批號 IRRADIATION RUN NUMBER : NEW14372-01

客戶產品已照射 MATERIALS PROCESSED :

總數 4 箱 數
TOTAL BOXES

中國生化科技股份有限公司證明上述產品經本公司劑量偵測系統判讀，吸收劑量如下：
China Biotech Corporation certifies that the material listed above (has described by its manufacturer)
received the following doses within the precision limits of the dosimetry system employed

最低劑量 66.9 kGy ; 最高劑量 78.7 kGy
MINIMUM DOSAGE MAXIMUM DOSAGE

使用放射性同位素 ISOTOPE UTILIZED : 鈷 60 COBALT-60


客戶劑量要求 DOSE REQUIREMENT : 最低劑量 MIN 64.0 kGy ; 最高劑量 MAX 80.0 kGy

確 認 者: [Signature]
CERTIFIED BY
品保部主管
QUALITY ASSURANCE

確保人類健康，珍惜自然環境

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 10 of 23

APPENDIX I Transport and Storage Cycle, Preconditioning and Ship Testing



Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 11 of 23

APPENDIX I Transport and Storage Cycle, Preconditioning and Ship Testing

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-229
	Shipping Test – Zip Pen 2525-10	Revision: A
		Page 17 of 23

APPENDIX I Shipping Test Log Sheet

Preconditioning:

Start Date: 1-5-2015 Chamber Number: 01095
Completion Date: 1-8-2015 Last Calibration: 5-29-14
Signature/Date: Paul Valpreda 1-8-2015 Calibration due: 5-31-15

Drop Test:

Catalog 2525-10 Weight 13 Drop Height: 15
lot 8140303

Drop	Orientation	Specific face, edge or	Initials/Date
1	Top	Face 1	PV 1-8-2015
2	Edge	Edge 5-3	PV 1-8-2015
3	Edge	Edge 6-3	PV 1-8-2015
4	Corner	Corner 2-3-5	PV 1-8-2015
5	Corner	Corner 4-3-6	PV 1-8-2015
6	Bottom	Face 3	PV 1-8-2015

Comments: Pass

Signature: Paul Valpreda Date: 1-8-2015

Compression Test:

Catalog 2525-10 Pounds Force 275 lbs

Comments: Pass

Signature: Paul Valpreda Date: 1-8-2015

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:57:54 pm; Printed by: .

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 12 of 23

APPENDIX I
Transport and Storage Cycle, Preconditioning and Ship Testing

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-229
	Shipping Test – Zip Pen 2525-10	Revision: A
		Page 18 of 23

Appendix I Continued
Shipping Test Log Sheet

Vibration:

Low Frequency, 40 minutes, Initials PV High frequency 10 minutes, Initials PV

Completion Date: 1-8-2015

Signature: Paul Valprede Date: 1-8-2015

Second Drop Test:

Catalog 2525-10 Weight 13 Drop Height: 15/30
Lot 540303

Drop	Orientation	Specific face, edge or	Initials/Date
1	Edge	Edge 4-6	PV 1-8-2015
2	Face	Face 4	PV 1-8-2015
3	Face	Face 6	PV 1-8-2015
4	Corner	Corner 2-1-5	PV 1-8-2015
5	Edge	Edge 2-1	PV 1-8-2015
6	Bottom	Face 3, Increase height to 30 inches.	PV 1-8-2015

Comments: Pass

Signature: Paul Valprede Date: 1-8-2015

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:57:54 pm; Printed by: .

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 13 of 23

APPENDIX I
Transport and Storage Cycle, Preconditioning and Ship Testing

Document: ENG-PRT-229 Title: Shipping Test - Zip Pencil Revision: 001 Effective Date: 27 Mar 2014 12:00 AM Copy expires on: 02 Nov 2014 at 01:22:03 pm		
Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number 1150720-10
	Shipping Test - Zip Pencil MASTER DOCUMENT	Revision: A
		Effective Date: 2014 MAR 27
		Page 15 of 19

Appendix I Continued
Shipping Test Log Sheet

Vibration:

Low Frequency, 40 minutes, Initials PV High frequency 10 minutes, Initials PV

Completion Date: 11-7-2014

Signature: Paul Valpreda Date: 11-10-2014

Second Drop Test:

Catalog 2525-10 Weight 13 lbs. Drop Height: 15" and 30"

Drop	Orientation	Specific face, edge or	Initials/Date
1	Edge	Edge 4-6	PV 11-7-2014
2	Face	Face 4	PV 11-7-2014
3	Face	Face 6	PV 11-7-2014
4	Corner	Corner 2-1-5	PV 11-7-2014
5	Edge	Edge 2-1	PV 11-7-2014
6	Bottom	Face 3, Increase height to 30 inches.	PV 11-7-2014

Comments: _____

Signature: Paul Valpreda Date: 11-10-2014

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 31 Oct 2014, 02:22:33 pm; Printed by: MGLASSETT

Confidential - Controlled

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:57:54 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 14 of 23

APPENDIX II Print Legibility Log Sheet

Megadyne Medical Products, Inc.	TEST PROTOCOL	<u>Document Number</u> XENG-PRT-229
	Shipping Test – Zip Pen 2525-10	Revision: A
		Page 19 of 23

Appendix II Print Legibility Log Sheet

Inspect the product per the protocol and enter the number of units that pass or fail in the box below.

Catalog # <i>2525-10</i> <i>Lot 5140303</i>	Pass	Fail
Pouch Print	✓	
Lot Number Print	✓	

Comments: *Pass*

Paul Valprede
Inspected by: _____ Date completed: *1-8-2015*

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:57:54 pm; Printed by: .

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 15 of 23

APPENDIX III Bubble Leak Test Log Sheet

Megadyne Medical Products, Inc.	TEST PROTOCOL	<u>Document Number</u> XENG-PRT-229
	Shipping Test – Zip Pen 2525-10	Revision: A
		Page 20 of 23

Appendix III Bubble Leak Test Log Sheet

Catalog # 2525-10
Ref 5140303

Sample	Pass	Fail	Comment
1	X		
2	X		
3	X		
4	X		
5	X		
6	X		
7	X		
8	X		
9	X		
10	X		
11	X		
12	X		
13	X		
14	X		
15	X		
16	X		
17	X		
18	X		
19	X		
20	X		

Catalog # 2525-10
Ref 5140303

Sample	Pass	Fail	Comment
21	X		
22	X		
23	X		
24	X		
25	X		
26	X		
27	X		
28	X		
29	X		
30	X		
31	X		
32	X		
33	X		
34	X		
35	X		
36	X		
37	X		
38	X		
39	X		
40	X		

Signature: 

Date: 1-9-2015

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:57:54 pm; Printed by: .

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 16 of 23

APPENDIX IV
Dye Penetration Test Log Sheet

Megadyne Medical Products, Inc.	TEST PROTOCOL	<u>Document Number</u> XENG-PRT-229
	Shipping Test – Zip Pen 2525-10	Revision: A
		Page 21 of 23

Appendix IV
Dye Penetration Test Log Sheet

Catalog # 2525-10

Lot # 5140303

Sample	Pass	Fail	Comment
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
11			
12			
13			
14			
15			
16			
17			
18			
19			
20			

Sample	Pass	Fail	Comment
21			
22			
23			
24			
25			
26			
27			
28			
29			
30			
31			
32			
33			
34			
35			
36			
37			
38			
39			
40			

Signature: _____

Made Grant
Dan G. P.

Date: 1-12-15

1-12-15

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:57:54 pm; Printed by: .

Confidential - Controlled copy on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 17 of 23

APPENDIX V Burst Test Log Sheet

The user must ensure that they are using the correct/current revision of this document.
Document: OPER-FRM-004 Rev: 001 Effective: 30 Oct 2012 12:30 AM

Lucy Richards 2012 OCT 30
D.C. Verification:

Tim Kessinger 2012 OCT 25
Second Verification:

Box 1

Lot Number	Catalog Number	Description	Date Tested	Quantity Tested
3140303	2525-10	Zip Pen	1-9-2015	20

Equipment Serial #: 01397
Calibration Due Date: 9-30-2015
Equipment Identification #: 01397 1-9-15 01397

Unless otherwise specified, the sampling plan is C=0.
Burst / Creep Inspection method is specified by 1100016-10.

Specification, in. H ₂ O	Burst	Date	Time	Initials	Burst Data, Pass/Fail - List Lane #'s for Multivac Only	Sample Lane #	pouch 1	Lane #	pouch 2	Lane #	pouch 3
19	✓	10115				19-1	29.1	R			
						19-1	35.9	B			
						15-1	24.8	T			
						8-2	26.7	R			
						3-2	26.7	B			
						19-1	29.5	L			
						14-2	29.7	R			
						10-1	26.4	T			
						12-1	27.2	T			
						20-2	25.8	R			
						1-1	27.5	L			
						6-2	27.6	T			
						11-2	27.1	R			
						16-1	28.5	L			
						17-2	25.2	R			
						5-1	32.0	T			
						2-1	28.3	L			

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.
Printed on: 09 Jan 2015, 09:15:07 am; Printed by: MGLASSETT.

Confidential - Controlled copy on date printed

APPENDIX V

Burst Test Log Sheet

Confidential - Controlled copy on date printed

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 19 of 23

APPENDIX V Burst Test Log Sheet

The user must ensure that they are using the correct/current revision of this document.
Document: OPER-FRM-004 Rev: 001 Effective: 30 Oct 2012 12:00 AM

Lucy Richards 2012 OCT 30
D.C. Verification:

Tim Kessinger 2012 OCT 25
Second Verification:

Lot Number	Catalog Number	Description	Date Tested	Quantity Tested
5140303	2525-10	ZIP PEN	1-9-2015	20

Box 2

Equipment Serial #: 01397
Calibration Due Date: 9-30-2015
Equipment Identification #: 01397

Unless otherwise specified, the sampling plan is C=0.
Burst / Creep Inspection method is specified by 1100016-10.

Equipment Serial #: 01397
Calibration Due Date: 9-30-2015
Equipment Identification #: 01397

Initials: [Signature] Date: [Signature]

Specification, in. H ₂ O	Burst	Date	Time	Initials	Burst Data, Pass/Fail - List Lane #'s for Multivac Only	Lane #	pouch 1	Lane #	pouch 2	Lane #	pouch 3
19	✓	10:25		[Signature]	27-1	27.3	L				
					28-1	28.4	L				
					29-1	29.2	T				
					30-2	26.8	R				
					31-2	27.9	R				
					32-1	27.3	L				
					32-2	29.4	B				
					35-2	27.9	R				
					40-1	27.0	L				
					37-2	28.0	B				
					21-1	27.0	T				
					25-1	25.5	L				
					33-2	25.9	R				
					28-2	29.8	R				
					31-2	27.7	R				
					22-2	28.5	B				
					26-2	26.6	B				

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.
Printed on: 09 Jan 2015, 09:15:07 am; Printed by: MGLASSETT.

Confidential - Controlled copy on date printed

APPENDIX V

Burst Test Log Sheet

Confidential - Controlled copy on date printed

Printed on: 21 Jan 2020, 10:57:54 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 21 of 23

APPENDIX VI Package Seal Width Log Sheet

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-229 Rev: A Effective: 06 Jan 2015 1:35 PM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-229
	Shipping Test – Zip Pen 2525-10	Revision: A
		Page 22 of 23

Appendix V Seal Width Log Sheet

Paul Valprede 1-12-2015

Catalog # 2525-10 Lot #: 5140303

Sample	Cavity	Front	Back	Right	Left
1	1	.36	.28	.33	.34
2	1	.35	.27	.32	.33
3	2	.35	.31	.33	.36
4	1	.34	.27	.34	.32
5	1	.34	.29	.33	.33
6	2	.36	.34	.35	.37
7	2	.30	.33	.31	.26
8	2	.31	.28	.33	.33
9	1	.29	.35	.29	.33
10	1	.34	.30	.34	.33
11	2	.33	.31	.29	.35
12	1	.37	.31	.24	.32
13	1	.32	.32	.34	.32
14	2	.37	.29	.34	.33
15	1	.35	.22	.30	.32
16	1	.36	.34	.35	.31
17	2	.35	.32	.32	.32
18	1	.32	.35	.32	.31
19	1	.33	.32	.34	.34
20	2	.34	.27	.33	.37
21					
22					
23					
24					
25					
26					
27					
28					
29					
30					
31					
32					
33					
34					
35					
36					
37					
38					
39					
40					

Appendix VI

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 06 Jan 2015, 01:42:33 pm; Printed by: MGLASSETT.

started vib at 1:15 PM on 1-9-2015

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:57:54 pm; Printed by: .

Confidential - Controlled copy on date printed

Confidential - Pre-Production Controlled on date printed

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 22 of 23

APPENDIX VI Package Seal Width Log Sheet

The user must ensure that they are using the correct/current revision of this document.
Document: XENG-PRT-229 Rev: A Effective: 06 Jan 2015 1:35 PM

Megadyne Medical Products, Inc.	TEST PROTOCOL	Document Number XENG-PRT-229
	Shipping Test – Zip Pen 2525-10	Revision: A
		Page 22 of 23

Appendix V Seal Width Log Sheet

Catalog # 2525-10 Lot #: S140303

T. Corlyk 1/9/2015
Mitrukyo CD-8 (40110-4)
cut 439082 USA 01473
Cal. 4 max 13 / 0.05 31 max 15

Sample	Cavity	Front	Back	Right	Left
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20					
21	1	0.30	0.29	0.24	0.35
22	2	0.34	0.34	0.36	0.37
23	2	0.36	0.32	0.38	0.33
24	1	0.35	0.30	0.35	0.38
25	1	0.35	0.33	0.34	0.33
26	2	0.345	0.36	0.35	0.38
27	1	0.37	0.29	0.33	0.29
28	2	0.33	0.34	0.34	0.37
29	1	0.35	0.34	0.34	0.34
30	2	0.34	0.34	0.34	0.34
31	2	0.34	0.34	0.34	0.37
32	2	0.35	0.33	0.36	0.37
33	2	0.34	0.34	0.34	0.38
34	1	0.38	0.37	0.36	0.31
35	2	0.31	0.35	0.35	0.33
36	2	0.33	0.29	0.33	0.38
37	2	0.35	0.35	0.33	0.35
38	2	0.37	0.36	0.33	0.35
39	1	0.36	0.34	0.34	0.34
40	1	0.37	0.34	0.35	0.35

Appendix VI

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 06 Jan 2015, 01:42:33 pm; Printed by: MGLASSETT.

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.

Printed on: 21 Jan 2020, 10:57:54 pm; Printed by: .

Megadyne Medical Products, Inc.	TEST REPORT	Document Number ENG-RPT-413
	Zip Pen 2525-10 Shipping Test	Revision: 001
		Page 23 of 23

APPENDIX VI
Package Seal Width Summary

Burst Test Data Zip Pen 2525-10 Lot S140303, Box 1 and 2 Minimum Burst Specification is 19 in. H2O							
Sample-Lane/Box 1	Burst Value	Burst Loc	Minimum Seal/Loc	Sample-Lane/Box 2	Burst Value	Burst Loc	Minimum Seal/ Loc
1-1	27.5	L	.28/B	21-1	27.0	T	.24/R
2-1	28.3	L	.27/B	22-2	28.5	B	.34/F,B
3-2	26.7	B	.31/B	23-2	27.9	R	.32/B
4-1	28.3	L	.27/B	24-1	29.3	L	.30/B
5-1	32.0	T	.29/B	25-1	25.5	L	.33/B,L
6-2	27.6	T	.34/B	26-2	26.6	B	.34/F,B
7-2	29.1	R	.30/F	27-1	28.4	L	.28/L
8-2	26.7	R	.28/B	28-2	28.8	R	.33/F
9-1	35.9	B	.29/F,R	29-1	29.2	T	.34/B,R,L
10-1	26.4	T	.30/B	30-2	28.9	R	.34/ALL
11-2	27.1	R	.29/R	31-2	27.7	R	.34/F,B,R
12-1	27.2	T	.24/R	32-2	29.4	B	.33/B,L
13-1	28.1	L	.32/F,B,L	33-2	25.9	R	.34/B,R,L
14-2	29.7	R	.29/B	34-1	27.3	L	.31/L
15-1	24.8	T	.22/B	35-2	27.9	R	.31/F
16-1	28.5	L	.31/L	36-2	27.1	B	.29/B
17-2	25.2	R	.32/B,R,L	37-2	28.0	B	.25/B
18-1	28.6	T	.31/L	38-2	26.8	R	.33/R
19-1	29.5	L	.32/B,R,L	39-1	27.3	L	.34/B,R,L
20-2	25.8	R	.27/B	40-1	27.0	L	.34/B
					27.9	Average	0.30
					24.8	Minimum	0.22
						Std Dev	0.032
						Avg-3sigma	0.21

Possession of this document is an acknowledgment that the contents herein are the exclusive property of Megadyne Medical Products, Inc. This document may not be reproduced in any form whatsoever without written permission from Megadyne. The user of this document must ensure that they are using the most current revision of this document.